
Reprocessing Instructions

REF

900MR017	900MR024	900MR025	900MR026
900MR042	900MR042	900MR045	900MR047
900MR049	900MR050	900MR051	900MR061
900MR062	900MR067	900MR068	900MR070
900MR071	900MR072	900MR074	900MR075
900MR125	900MR126	900MR127	900MR128
900MR130	900MR131	900MR132	900MR135
900MR136	900MR138	900MR139	900MR143
900MR145	900MR146	900MR147	900MR148
900MR156	900MR157	900MR158	900MR159
900MR178	900MR405	900MR406	900MR510
900MR524	900MR526	900MR528	900MR529
900MR530	900MR531	900MR532	900MR533
900MR534	900MR535	900MR540	900MR711
900MR715	900MR749	900MR751	900MR754
900MR755	900MR761	900MR780	900MR781
900MR782	900MR784	900MR830	900MR923
MR340E	MR340F	MR340S	MR370


A**REF MR370****REF MR340****REF 900MR139**

English	4
български (Bulgarian).....	9
Hrvatski (Croatian).....	14
Česky (Czech)	19
Dansk (Danish).....	24
Nederlands (Dutch).....	29
Eesti keel (Estonian).....	34
Suomi (Finnish).....	39
Français (French)	44
Deutsch (German)	49
Ελληνικά (Greek)	54
Magyar (Hungarian).....	59
Italiano (Italian).....	64
Latviešu (Latvian)	69
Lietuvių (Lithuanian)	74
Norsk (Norwegian)	79
Polski (Polish)	84
Português (Portuguese)	89
Português (Brasil) (Brazilian Portuguese).....	94
Română (Romanian)	99
Slovenčina (Slovak).....	104
Slovenščina (Slovenian)	109
Español (Spanish).....	114
Svenska (Swedish)	119
Türkçe (Turkish)	124
繁體中文版 (Traditional Chinese).....	129
简体中文版 (Simplified Chinese)	134
Bahasa Indonesia (Indonesian)	139
ไทย (Thai).....	144
Tiếng Việt (Vietnamese).....	149
Русский (Russian)	154
한국어 (Korean).....	159
日本語 (Japanese)	164

Section 1: Component Reprocessing Overview

This instruction includes the manufacturer’s recommended reprocessing steps for each product but will be dependent on the reprocessing facilities available and hospital reprocessing guidelines.


1.1 Product Family Descriptions

1.2 Reprocessing Overview

Figure 1 below shows the recommended reprocessing steps for the respiratory products listed on the front page of this document. Automated washer disinfection and sterilization are strongly recommended for all respiratory products where possible.

Figure 1: Recommended Reprocessing Steps

For the products listed in Table 1 below only, use only manual cleaning (Step D) and high-level disinfection (chemical) (Step E). Do not use automated cleaning (Step F) or steam sterilization (Step I) on these products. Refer to Section 2 for detailed descriptions of each step.

Table 1

Heater-wires	900MR510, 900MR524, 900MR751, 900MR754, 900MR755
Hose Clips	900MR042, 900MR047, 900MR049
Draw wires	900MR070, 900MR071

Section 2: Specific Reprocessing Instructions

Clean the products prior to use and whenever the products are visibly soiled. Follow the cleaning instructions outlined in this document. Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.

Warnings:

- Before use on a patient all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process (all products except those listed in Table 1, Section 1.2).
- Products must be reprocessed according to the reprocessing instructions in this document. The use of unapproved cleaning methods may damage the product, reduce its usable life and lead to the risk of cross-infection between patients and potential serious harm. Fisher & Paykel Healthcare cannot be held responsible for correctly functioning components that are not reprocessed according to these instructions.
- Used products must be handled as contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and the correct disposal of these products.
- Follow the healthcare facility infection control procedures as well as local laws, standards and regulations. This applies in particular to regulations regarding the effective deactivation of prions.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These faults may cause gas leaks and a loss of ventilation or respiratory support.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Cautions:

- Before handling cleaning agents wear appropriate personal protective equipment, such as: disposable gloves, splash visor or safety glasses and a protective gown. Ensure disinfection is performed in a fume cupboard or with fume extraction and follow the disinfectant instructions for use.

Chemical Compatibility:

- Solutions containing ketones, formaldehyde, chlorinated hydrocarbons, hypochlorite, inorganic acids, aromatic hydrocarbons or phenol (>5 %) should be avoided as they may cause plastic components to crack, leading to compromised therapy.
- To avoid degradation of the aluminum chamber bases and black silicon rubber components, do not expose these components to hydrogen peroxide solutions.

Failure to comply with the above warnings and cautions may impair the performance of the device or compromise safety (including causing potential serious harm).

A Point of Use

It is recommended that reprocessing of the products should commence as soon as is reasonably practical following use to prevent gross contaminants from drying out on circuit components.

B Transportation

- Protect the product and components from mechanical damage during transport.
 - Package the product and components in accordance with hospital infection control protocols.
-

C Preparation Before Cleaning

- The products must be disconnected from both the power and gas source prior to cleaning.
 - Disconnect and separate all adaptors and tubes into individual items.
 - To avoid damage to the breathing tubes, attach or detach from equipment by handling end connectors only; do not pull or twist the tubing.
 - Disassemble all plugs and caps from ports and from the flexible retainers. Remove flexible retainers by stretching opposing retainer ends.
 - Disassemble the MR340 and MR370 reusable chambers as per diagram **A**.
 - Disassemble 900MR139 water trap as per diagram **A**.
 - Reprocess the components together.
-

D Manual Cleaning

Manual cleaning is required if the product has dried-on or gross soil or CANNOT be cleaned in an automated washer disinfectant.

1. Rinse in cold water (<43 °C, 109 °F).
2. Soak in enzymatic cleaning solution. Refer to the detergent manufacturer's instructions
Note: Ensure the product is fully submerged and absent of trapped air bubbles so the cleaning solution wets all surfaces.
3. If required, use a soft-bristle cleaning brush to remove gross soil.
4. Rinse thoroughly with potable water after removal from the cleaning solution.

 **Warning:** Do not expose 900MR510, 900MR524, 900MR751, 900MR754 and 900MR755 heater-wires to reprocessing temperatures >100 °C (212 °F).

Caution: It is recommended to use an automated washer disinfectant for the 900MR075 pressure tube due to the small inner bore and length of the tube.

E High-level Disinfection (Chemical)

The following chemical disinfectant is recommended: 0.55 % ortho-phthalaldehyde (e.g. Cidex OPA).

Note: Follow the disinfectant manufacturer's instructions for the recommended time to achieve high-level disinfection.

1. Ensure components are fully submerged and no air bubbles are present so the disinfectant solution wets all surfaces.
2. Rinse thoroughly with potable water after removal from the disinfection solution.

 **Warning:** Do not expose 900MR510, 900MR524, 900MR751, 900MR754 and 900MR755 heater-wires to reprocessing temperatures >100 °C (212 °F).

Caution: It is recommended to use an automated washer disinfectant for the 900MR075 pressure tube due to the small inner bore and length of the tube.

F Automated Cleaning and Disinfection (Washer Disinfector)

Caution: 900MR510, 900MR524, 900MR751, 900MR754 and 900MR755 heater-wires, hose clips and draw-wires are NOT suitable for automated cleaning and disinfection.

The washer disinfector must comply with the requirements of and be validated, maintained, and checked in accordance with, either ISO 15883-1 or for customers in the USA, ANSI or AAMI ST15883-1.

1. Breathing and pressure tubes must be loaded onto an auto-washer tray cart with anesthesia manifolds, where individual injectors or spigots circulate washing and rinsing solutions through the length of the tube.
2. Ensure the tube ends are firmly fixed to manifolds to prevent disconnection during the wash cycle.
3. Disassembled water traps, Y-pieces, 900MR711 and 900MR715 heater-wires, adaptors and connectors, etc. are to be placed in instrument trays with lids and be orientated so that all internal and external surfaces are accessible to cleaning agents and can drain freely.

Recommended washer disinfector cycle parameters:

- Cold rinse: 40 °C (104 °F) for at least 2 minutes.
- Wash with mildly alkaline cleaning agent: 60 °C (140 °F) for at least 3 minutes, or as per the detergent manufacturer's instructions.
- Rinse with a rinse aid at 80 °C (176 °F) for at least 2 minutes, or as per the rinse aid manufacturer's instructions.
- Disinfection phase: 90 °C (194 °F) for at least 5 minutes to ensure $A_0 = 3000$.

G Drying

Dry all products after manual disinfection or if additional drying is required. Ensure all products and components are completely dry before being sterilized.

- Air dry breathing and pressure tubes in a drying cabinet up to 90 °C (194 °F) for at least 18 minutes or until visibly dry.
- The tubes should hang vertically and be connected to a drying nozzle to allow warm air to circulate over all surfaces and prevent water pooling.
- Do not dry at temperatures >90 °C (194 °F).
- All other components are to be dried in perforated trays that allow warm air to circulate over all surfaces and prevent water pooling.
- Ensure all other components are oriented to allow water to drain away and to avoid water pooling.

Warning: Do not expose 900MR510, 900MR524, 900MR751, 900MR754 and 900MR755 heater-wires to reprocessing temperatures >100 °C (212 °F).

H Packaging

For high-level disinfected product:

- The product should be protected from contamination to maintain high-level disinfection.
- Package the product for usage as per hospital guidelines for high-level disinfected products.
- If the products will not be sterilized proceed to Storage (Step J in Figure 1).

For products intended to undergo steam sterilization:

- Package components individually in either non-woven wraps or steam sterilization pouches as per hospital guidelines.
-

I Steam Sterilization

Caution: 900MR510, 900MR524, 900MR751, 900MR754 and 900MR755 heater-wires, hose clips, and draw-wires are NOT suitable for steam sterilization.

- Autoclaves should comply with the requirements of and be validated, maintained and checked in accordance with EN 285 or EN 13060, EN ISO 17665, and ANSI/AAMI ST79. It is recommended that a pre-vacuum (forced air removal) cycle is used to ensure the steam accesses all surfaces. Follow hospital protocols for steam sterilization.
 - Autoclave cycle parameters: 134 °C (270 °F) at 207 kPa (30 psig) for at least 4 minutes, **or** 121 °C (250 °F) at 103 kPa (15 psig) for at least 15 minutes.
 - Drying stage: >30 minutes.
-

J Storage Post Reprocessing

- Follow hospital protocol for the storage of products that have been high-level disinfected or sterilized.
 - Reprocessed products should be stored to prevent any damage to the product packaging. Damaged packaging could affect the level of disinfection or sterilization.
 - Avoid storing under direct UV light.
-

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

The diagram illustrates the components of a respiratory system. It includes a Y-piece (1) connected to various tubes (3) and connectors (2). A water canister (6) is shown at the bottom left, and a pressure tube (7) is connected to a water cup (8). A coiled tube (9) is shown on the right. A tube clamp (4) is also depicted. The components are numbered 1 through 9 for identification.

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The flowchart shows the recommended reprocessing steps: A (Use location) → B (Transport) → C (Preparation for cleaning) → D (Manual cleaning) → F (Automatic cleaning and disinfection) → G (Drying) → H (Packaging) → I (Steam sterilization) → J (Storage after reprocessing). There is also a branch from D to E (High-level disinfection (chemical)) which then joins the main flow at J.

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