
Hose clip - 15 mm

USER INSTRUCTIONS REF 900MR047

| | |
|--|----|
| English | 3 |
| български (Bulgarian)..... | 5 |
| Hrvatski (Croatian)..... | 7 |
| Česky (Czech)..... | 9 |
| Dansk (Danish)..... | 11 |
| Nederlands (Dutch)..... | 13 |
| Eesti keel (Estonian)..... | 15 |
| Suomi (Finnish)..... | 17 |
| Français (French)..... | 19 |
| Deutsch (German)..... | 21 |
| Ελληνικά (Greek)..... | 23 |
| Magyar (Hungarian)..... | 25 |
| Italiano (Italian)..... | 27 |
| Latviešu (Latvian)..... | 29 |
| Lietuvių (Lithuanian)..... | 31 |
| Norsk (Norwegian)..... | 33 |
| Polski (Polish)..... | 35 |
| Português (Portuguese)..... | 37 |
| Português (Brasil) (Brazilian Portuguese)..... | 39 |
| Română (Romanian)..... | 41 |
| Slovenčina (Slovak)..... | 43 |
| Slovenščina (Slovenian)..... | 45 |
| Español (Spanish)..... | 47 |
| Svenska (Swedish)..... | 49 |
| Türkçe (Turkish)..... | 51 |
| 繁體中文版 (Traditional Chinese)..... | 53 |
| Bahasa Indonesia (Indonesian)..... | 55 |
| ไทย (Thai)..... | 57 |
| Tiếng Việt (Vietnamese)..... | 59 |
| Русский (Russian)..... | 61 |
| 한국어 (Korean)..... | 63 |
| 日本語 (Japanese)..... | 65 |

REF 900MR047**Hose clip – 15 mm****Intended Use**

The 900MR047 hose clip is used to secure temperature probe cables to the inspiratory limb. This accessory is intended for use with the Fisher & Paykel Healthcare 900MR749 pediatric reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

Warning, Cautions and Notes**WARNINGS**

- This product is reprocessible (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage.
- Do not modify this product.
- Do not use breathing circuits, chambers, accessories or combinations which are not approved by Fisher & Paykel Healthcare.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Before use on a patient, all components must be cleaned and disinfected.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.
- For use under the supervision of trained medical personnel.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

Reprocessing Instructions

For reprocessing instructions (UI-623716), please visit www.fphcare.com/fp-reusables

WARNINGS

- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

Symbol Definitions

| | | |
|--|---|---|
| <p>Product does not contain phthalates or pyrogens</p> | <p>For use on pediatric patients</p> | <p>European Union authorized representative</p> |
| <p>Transportation and storage temperature limits</p> | <p>Not made with natural rubber latex</p> | <p>Caution symbol for warnings</p> |
| <p>CE Mark</p> | <p>Use-by date</p> | <p>Prescription only</p> |
| <p>Date of manufacture</p> | <p>Instructions for Use</p> | <p>Lot number</p> |
| <p>Manufacturer</p> | <p>Medical Device</p> | <p>Reference number</p> |

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REF 623631 REV A 2022-05 © 2022 Fisher & Paykel Healthcare Limited