
Draw wire - 1.5 m

USER INSTRUCTIONS REF 900MR070

English	3
български (Bulgarian).....	5
Hrvatski (Croatian).....	7
Česky (Czech).....	9
Dansk (Danish).....	11
Nederlands (Dutch).....	13
Eesti keel (Estonian).....	15
Suomi (Finnish).....	17
Français (French).....	19
Deutsch (German).....	21
Ελληνικά (Greek).....	23
Magyar (Hungarian).....	25
Italiano (Italian).....	27
Latviešu (Latvian).....	29
Lietuvių (Lithuanian).....	31
Norsk (Norwegian).....	33
Polski (Polish).....	35
Português (Portuguese).....	37
Română (Romanian).....	39
Slovenčina (Slovak).....	41
Slovenščina (Slovenian).....	43
Español (Spanish).....	45
Svenska (Swedish).....	47
Türkçe (Turkish).....	49
繁體中文版 (Traditional Chinese).....	51
Bahasa Indonesia (Indonesian).....	53
ไทย (Thai).....	55
Tiếng Việt (Vietnamese).....	57
Русский (Russian).....	59
한국어 (Korean).....	61
日本語 (Japanese).....	63

REF 900MR070 Draw wire - 1.5 m**Intended Use**

The 900MR070 draw wire is used to assemble reusable heater-wires into reusable inspiratory limb tubes. This accessory is intended to assist with the assembly of Fisher & Paykel Healthcare's reusable breathing circuits in hospital or clinical environments.

Warning, Cautions and Notes**WARNINGS**

- This product is reprocessable (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect product before each use. Discard if faulty or if there is any sign of deterioration such as cracks or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Do not modify this product.
- Do not use breathing circuits, chambers, accessories or combinations which are not approved by Fisher & Paykel Healthcare.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Before use on a patient, the product must be cleaned and disinfected.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- Do not use the product during X-rays, MRIs or other radio-frequency procedures.
- For use under the supervision of trained medical personnel.
- If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

- Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

Reprocessing Instructions

For reprocessing instructions (UI-623716), please visit www.fphcare.com/fp-reusables

**WARNINGS**

- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

Symbol Definitions

 <p>Product does not contain phthalates or pyrogens</p>	 <p>Do not recycle</p>	 <p>European Union authorized representative</p>
 <p>Transportation and storage temperature limits</p>	 <p>Not made with natural rubber latex</p>	 <p>Caution symbol for warnings</p>
 <p>Use-by date</p>	 <p>Date of manufacture</p>	 <p>Manufacturer</p>
 <p>Lot number</p>	 <p>Reference number</p>	 <p>Instructions for Use</p>

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