
Reusable breathing tubing - 10 mm x 0.5 m

USER INSTRUCTIONS REF **900MR072**

English	3
български (Bulgarian).....	6
Hrvatski (Croatian).....	9
Česky (Czech).....	12
Dansk (Danish).....	15
Nederlands (Dutch).....	18
Eesti keel (Estonian).....	21
Suomi (Finnish).....	24
Français (French).....	27
Deutsch (German).....	30
Ελληνικά (Greek).....	33
Magyar (Hungarian).....	36
Italiano (Italian).....	39
Latviešu (Latvian).....	42
Lietuvių (Lithuanian).....	45
Norsk (Norwegian).....	48
Polski (Polish).....	51
Português (Portuguese).....	54
Português (Brasil) (Brazilian Portuguese).....	57
Română (Romanian).....	60
Slovenčina (Slovak).....	63
Slovenščina (Slovenian).....	66
Español (Spanish).....	69
Svenska (Swedish).....	72
Türkçe (Turkish).....	75
繁體中文版 (Traditional Chinese).....	78
Bahasa Indonesia (Indonesian).....	81
ไทย (Thai).....	84
Tiếng Việt (Vietnamese).....	87
Русский (Russian).....	90
한국어 (Korean).....	93
日本語 (Japanese).....	96

REF 900MR072**Reusable breathing tubing –
10 mm x 0.5 m****Intended Use**

The 900MR072 is a reusable tube used to convey breathing gases. This accessory is intended for use in neonatal reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

Product Specifications

Maximum Operating Pressure	80 hPa (cmH ₂ O) (8 kPa)
Gas leakage at 60 cmH₂O	<30 mL/min
Compliance at 60 cmH₂O	0.256 mL/cmH ₂ O
Resistance to flow at 2.5 L/min	0.004 cmH ₂ O/ L/min
Length	0.5 m (19.7 in)
Minimum internal diameter	10 mm (0.4 in)

Warning, Cautions and Notes**WARNINGS**

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing or accessories. Failure to comply may result in fire, burns and loss of system performance.
- This product is reprocessible (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect the product before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Ensure appropriate ventilator or flow source alarms are set before connecting the breathing circuit to the patient.
- Check all connections are tightly secured before use.
- Perform a pressure and leak test on the breathing circuit and check for occlusions before connecting to a patient.

- Use of non-compatible, non-ISO 5367, or non-ISO 80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Heated tubes should be kept out of direct contact with the patient's skin. Failure to comply may result in skin damage.
- Covering breathing tubes with blankets, pillows, towels or bed linen or heating them in an incubator or with an overhead heater can affect the quality of the therapy or injure the patient.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted, turn the humidifier off.
- Do not crush or stretch the breathing tube.
- Do not expose the breathing tube to excessive UV radiation.
- Do not modify this product.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Notes

- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- To avoid damage to the breathing tube, attach or detach from equipment by handling the end connectors only; do not pull or twist tubing.
- Keep the product away from heat or ignition sources, especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.

Monitoring Instructions**WARNING**

- Regularly monitor and drain condensate build-up in the circuit. Condensate build-up may result in loss of ventilation or respiratory support.

Note

- Avoid fans and air conditioner venting onto the breathing circuit during use as this may increase condensate in the circuit and potentially block a tube, resulting in loss of ventilation or respiratory support.

Reprocessing Instructions

For reprocessing instructions (UI-623716), visit <http://www.fphcare.com/fp-reusables>

**WARNINGS**

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.


Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

Symbol Definitions

 <p>Product does not contain phthalates or pyrogens</p>	 <p>Transportation and storage temperature limits</p>	 <p>For use on neonatal patients</p>
 <p>European Union authorized representative</p>	 <p>Not made with natural rubber latex</p>	 <p>Caution symbol for warnings</p>
 <p>Type BF Applied Part</p>	 <p>Do not recycle</p>	 <p>Instructions for Use</p>
<p>Rx only</p> <p>Prescription only</p>	 <p>Reference number</p>	 <p>Date of manufacture</p>
 <p>Manufacturer</p>	 <p>Lot number</p>	 <p>Medical Device</p>
<p>CE 0123</p> <p>CE Mark</p>	 <p>Use-by date</p>	

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REF 623671 REV A 2022-05 © 2022 Fisher & Paykel Healthcare Limited