
Y-piece - pediatric

USER INSTRUCTIONS

REF

900MR136

English	3
български (Bulgarian)	5
Hrvatski (Croatian)	7
Česky (Czech)	9
Dansk (Danish)	11
Nederlands (Dutch)	13
Eesti keel (Estonian)	15
Suomi (Finnish)	17
Français (French)	19
Deutsch (German)	21
Ελληνικά (Greek)	23
Magyar (Hungarian)	25
Italiano (Italian)	27
Latviešu (Latvian)	29
Lietuvių (Lithuanian)	31
Norsk (Norwegian)	33
Polski (Polish)	35
Português (Portuguese)	37
Português (Brasil) (Brazilian Portuguese)	39
Română (Romanian)	41
Slovenčina (Slovak)	43
Slovenčina (Slovenian)	45
Español (Spanish)	47
Svenska (Swedish)	49
Türkçe (Turkish)	51
繁體中文版 (Traditional Chinese)	53
Bahasa Indonesia (Indonesian)	55
ไทย (Thai)	57
Tiếng Việt (Vietnamese)	59
Русский (Russian)	61
한국어 (Korean)	63
日本語 (Japanese)	65

REF 900MR136**Y-piece – pediatric****Intended Use**

The 900MR136 Y-piece is used to connect the patient interface to the inspiratory and expiratory breathing tubes. This accessory is intended for use with pediatric reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

Product Specifications

Interface Connections	ISO 5356-1 conical connectors
Maximum Operating Pressure	80 hPa (cmH ₂ O) (8 kPa)

Warning, Cautions and Notes**WARNINGS**

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing, or accessories. Failure to comply may result in fire, burns and loss of system performance.
- This product is reprocessable (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Check all connections are tightly secured before use.
- The use of non-compatible, non-ISO-5367, or non-ISO-80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Do not modify this product.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

• California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, and suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

Reprocessing Instructions

For reprocessing instructions (UI-623716), please visit www.fphcare.com/fp-reusables

WARNINGS

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

Symbol Definitions

		
Product does not contain phthalates or pyrogens	Transportation and storage temperature limits	European Union authorized representative
		
Not made with natural rubber latex	For use on pediatric patients	Caution symbol for warnings
		
Type BF Applied Part	Do not recycle	Instructions for Use
Rx only	REF	
Prescription only	Reference number	Date of manufacture
	LOT	MD
Manufacturer	Lot number	Medical Device
CE 0123		
CE Mark	Use-by date	

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