
Hose heater 0.7 m 29 ohm 17 W inspiratory

USER INSTRUCTIONS **REF** **900MR524**


A


1


2

3

4

English	4
Български (Bulgarian).....	7
Hrvatski (Croatian).....	10
Česky (Czech).....	13
Dansk (Danish).....	16
Nederlands (Dutch).....	19
Eesti keel (Estonian).....	22
Suomi (Finnish).....	25
Français (French).....	28
Deutsch (German).....	31
Ελληνικά (Greek).....	34
Magyar (Hungarian).....	37
Italiano (Italian).....	40
Latviešu (Latvian).....	43
Lietuvių (Lithuanian).....	46
Norsk (Norwegian).....	49
Polski (Polish).....	52
Português (Portuguese).....	55
Português (Brasil) (Brazilian Portuguese).....	58
Română (Romanian).....	61
Slovenčina (Slovak).....	64
Slovenščina (Slovenian).....	67
Español (Spanish).....	70
Svenska (Swedish).....	73
Türkçe (Turkish).....	76
繁體中文版 (Traditional Chinese).....	79
Bahasa Indonesia (Indonesian).....	82
ไทย (Thai).....	85
Tiếng Việt (Vietnamese).....	88
Русский (Russian).....	91
한국어 (Korean).....	94
日本語 (Japanese).....	97

REF 900MR524

Hose heater 0.7 m 29 ohm 17 W inspiratory

Intended Use

The 900MR524 is a reusable heater-wire used to heat breathing gases within a 0.7 m (27.5 in) 13 mm diameter inspiratory breathing tube. This accessory is intended for use in neonatal or pediatric reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

Product Specifications

Maximum Operating Pressure 80 hPa (cmH₂O)
(8 kPa)

Warning, Cautions and Notes

WARNINGS

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
 - Do not lubricate connections, tubing, or accessories. Failure to comply may result in fire, burns, and loss of system performance.
 - This product is reprocessable (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
 - Ensure appropriate ventilator or flow source alarms are set before connecting the breathing circuit to patient.
 - Check all connections are tightly secured before use.
 - Perform a pressure and leak test on the breathing circuit and check for occlusions before connecting to a patient.
 - The circuit is powered through the heater-wire adaptor from the MR850. Refer to MR850 humidifier user instructions.
 - Use of non-compatible, non-ISO-5367, or non-ISO-80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
 - Heated tubes should be kept out of direct contact with the patient's skin. Failure to comply may result in skin damage.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted, turn the humidifier off.
 - Do not modify this product.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

Notes

- Before use on a patient all components must be cleaned and disinfected.
- Wear and degradation of electrical connectors may affect the product performance. Inspect the product components between uses.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- To avoid damage to breathing tubes, attach and detach from equipment by handling end connectors only; do not pull or twist tubing.
- Keep product away from heat or ignition sources, especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen, or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- Do not use the product during X-rays, MRIs, or other radio frequency procedures.
- Refer to MR850 Respiratory Humidifier user instructions for Electromagnetic Compatibility information.
- This product is safe for use on children.

Assembly Instructions

Refer to diagram **A** for heater-wire assembly.

**WARNINGS**

- Do not assemble heater-wire into single-use tubing.
- To avoid damaging heater-wire insulation, the draw-wire hook should only be attached to the heater-wire o-ring.
- Ensure heater-wire is extended along the tube length and is not bunched at any point.
- Do not use excessive force when drawing the heater-wire through tubes.
- Do not use the draw wire to tighten the tube and heater-wire elbow together.

Notes

- Ensure the chamber end probe is securely inserted into the heater-wire elbow port.

Monitoring Instructions**WARNING**

- Visually inspect the heater-wire before each use on a patient. Discard if faulty or if there is any sign of insulation deterioration such as cracks, tears, or damage. These may lead to electrical hazards and possible risk of serious patient harm.

Reprocessing Instructions

For reprocessing instructions (UI-623716), visit <http://www.fphcare.com/fp-reusables>.

WARNINGS

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever it is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.
- Do not expose heater-wires to reprocessing temperatures >100 °C (212 °F).

Notes

- It is recommended that reprocessing should commence as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state, and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories should be stored in a clean plastic bag or container and away from direct sunlight or heat if not reusing immediately.

Symbol Definitions

Product does not contain phthalates or pyrogens	Unsafe for use in MRI environments	European Union authorized representative
Transportation and storage temperature limits	Not made with natural rubber latex	Reference number
		Rx only
Caution symbol for warnings	Use-by date	Prescription only
Type BF Applied Part	Date of manufacture	Instructions for Use
Lot number	Manufacturer	Medical Device
Do not recycle	CE Mark	

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REF 623698 REV A 2022-05 © 2022 Fisher & Paykel Healthcare Limited