
Reusable breathing tubing 15 mm x 1.3 m

**USER INSTRUCTIONS** REF **900MR529**

---

English .....	3
Български (Bulgarian).....	6
Hrvatski (Croatian).....	9
Česky (Czech).....	12
Dansk (Danish).....	15
Nederlands (Dutch).....	18
Eesti keel (Estonian).....	21
Suomi (Finnish).....	24
Français (French).....	27
Deutsch (German).....	30
Ελληνικά (Greek).....	33
Magyar (Hungarian).....	36
Italiano (Italian).....	39
Latviešu (Latvian).....	42
Lietuvių (Lithuanian).....	45
Norsk (Norwegian).....	48
Polski (Polish).....	51
Português (Portuguese).....	54
Português (Brasil) (Brazilian Portuguese).....	57
Română (Romanian).....	60
Slovenčina (Slovak).....	63
Slovenščina (Slovenian).....	66
Español (Spanish).....	69
Svenska (Swedish).....	72
Türkçe (Turkish).....	75
繁體中文版 (Traditional Chinese).....	78
Bahasa Indonesia (Indonesian).....	81
ไทย (Thai).....	84
Tiếng Việt (Vietnamese).....	87
Русский (Russian).....	90
한국어 (Korean).....	93
日本語 (Japanese).....	96

**REF 900MR529****Reusable breathing tubing  
15 mm x 1.3 m****Intended Use**

The 900MR529 is a reusable tube used to convey breathing gases. This accessory is intended for use in Fisher & Paykel Healthcare 900MR749 pediatric reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

**Product Specifications**

<b>Maximum Operating Pressure</b>	80 hPa (cmH <sub>2</sub> O) (8 kPa)
<b>Gas leakage at 60 cmH<sub>2</sub>O</b>	<40 mL/min
<b>Compliance at 60 cmH<sub>2</sub>O</b>	0.460 mL/cmH <sub>2</sub> O
<b>Resistance to flow at 15 L/min</b>	0.004 cmH <sub>2</sub> O/ L/min
<b>Length</b>	1.3 m (51.2 in)
<b>Minimum internal diameter</b>	15 mm (0.59 in)

**Warning, Cautions and Notes****WARNINGS**

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing or accessories. Failure to comply may result in fire, burns and loss of system performance.
- This product is reprocessible (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect the product before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Ensure appropriate ventilator or flow source alarms are set before connecting the breathing circuit to the patient.
- Check all connections are tightly secured before use.
- Perform a pressure and leak test on the breathing circuit and check for occlusions before connecting to a patient.

- Use of non-compatible, non-ISO 5367, or non-ISO 80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Heated tubes should be kept out of direct contact with the patient's skin. Failure to comply may result in skin damage.
- Covering breathing tubes with blankets, pillows, towels or bed linen or heating them in an incubator or with an overhead heater can affect the quality of the therapy or injure the patient.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted, turn the humidifier off.
- Do not crush or stretch the breathing tube.
- Do not expose the breathing tube to excessive UV radiation.
- Do not modify this product.
- Do not use breathing circuits, chambers, accessories or combinations which are not approved by Fisher & Paykel Healthcare.

**Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).**

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)

**Notes**

- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- To avoid damage to the breathing tube, attach or detach from equipment by handling the end connectors only; do not pull or twist tubing.
- Keep the product away from heat or ignition sources, especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

**Monitoring Instructions****WARNING**

- Regularly monitor and drain condensate build-up in the circuit. Condensate build-up may result in loss of ventilation or respiratory support.

**Note**

- Avoid fans and air conditioner venting onto the breathing circuit during use as this may increase condensate in the circuit and potentially block a tube, resulting in loss of ventilation or respiratory support.

**Reprocessing Instructions**

For reprocessing instructions (UI-623716), please visit [www.fphcare.com/fp-reusables](http://www.fphcare.com/fp-reusables)

**WARNINGS**

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

**Notes**

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

**Storage**

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

## Symbol Definitions

 <p>Product does not contain phthalates or pyrogens</p>	 <p>Transportation and storage temperature limits</p>	 <p>European Union authorized representative</p>
 <p>For use on pediatric patients</p>	 <p>Not made with natural rubber latex</p>	 <p>Caution symbol for warnings</p>
 <p>Reference number</p>	<p><b>Rx only</b></p> <p>Prescription only</p>	 <p>Type BF Applied Part</p>
 <p>Date of manufacture</p>	 <p>Instructions for Use</p>	 <p>Do not recycle</p>
 <p>Lot number</p>	 <p>Use-by date</p>	<p>CE 0123</p> <p>CE Mark</p>
 <p>Manufacturer</p>	 <p>Medical Device</p>	

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