
Hose heater 20 ohm 9.4 ohm/m inspiratory

USER INSTRUCTIONS **REF** **900MR755**

Rx only

CE 0123

MD

50 °C
122 °F
-10 °C
14 °F

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REF 623713 REV A 2021-07 © 2021 Fisher & Paykel Healthcare Limited


Fisher & Paykel
HEALTHCARE

A


1


2

3

4

English	4
български (Bulgarian).....	7
Hrvatski (Croatian).....	10
Česky (Czech).....	13
Dansk (Danish).....	16
Nederlands (Dutch).....	19
Eesti keel (Estonian).....	22
Suomi (Finnish).....	25
Français (French).....	28
Deutsch (German).....	31
Ελληνικά (Greek).....	34
Magyar (Hungarian).....	37
Italiano (Italian).....	40
Latviešu (Latvian).....	43
Lietuvių (Lithuanian).....	46
Norsk (Norwegian).....	49
Polski (Polish).....	52
Português (Portuguese).....	55
Română (Romanian).....	58
Slovenčina (Slovak).....	61
Slovenščina (Slovenian).....	64
Español (Spanish).....	67
Svenska (Swedish).....	70
Türkçe (Turkish).....	73
繁體中文版 (Traditional Chinese).....	76
简体中文版 (Simplified Chinese).....	79
Bahasa Indonesia (Indonesian).....	82
ไทย (Thai).....	85
Tiếng Việt (Vietnamese).....	88
Русский (Russian).....	91
한국어 (Korean).....	94
日本語 (Japanese).....	97

Hose Heater 20 ohm 9.4 ohm/m inspiratory

Intended Use

The 900MR755 is a reusable heater-wire used to heat breathing gases within an inspiratory breathing tube. This accessory is intended for use in Fisher & Paykel Healthcare 900MR781 neonatal reusable breathing circuits within hospital or clinical environments and should only be prescribed by a physician. Addition of heat and humidity to the supply of cool and dry respiratory gases provided through ventilation is beneficial to prevent drying of the patient's airways.

Product Specifications

Operating Flow Range

Infant Warmer

Invasive:	4–30 L/min
Noninvasive:	2–30 L/min

Incubator

Invasive:	2–30 L/min
Noninvasive:	2–30 L/min

Maximum Operating Pressure	8 kPa
-----------------------------------	-------

Technical Specifications

WARNINGS

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing, or accessories. Failure to comply may result in fire, burns, and loss of system performance.
- This product is reprocessable (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Ensure appropriate ventilator or flow source alarms are set before connecting the breathing circuit to patient.
- Check all connections are tightly secured before use.
- Perform a pressure and leak test on the breathing circuit and check for occlusions before connecting to a patient.
- The circuit is powered through the heater-wire adaptor from the MR850. Refer to MR850 humidifier user instructions.

- Use of non-compatible, non-ISO-5367, or non-ISO-80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
 - Heated tubes should be kept out of direct contact with the patient's skin. Failure to comply may result in skin damage. At flow rates below 4 L/min, the tube surface temperature may exceed 44 °C (111 °F).
 - Do not use the heated breathing circuit without gas flow. If gas flow is interrupted, turn the humidifier off.
 - Do not modify this product.
 - Do not use breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare.
- Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).**
- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Before use on a patient all components must be cleaned and disinfected.
- Wear and degradation of electrical connectors may affect the product performance. Inspect the product components between uses.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- To avoid damage to breathing tubes, attach and detach from equipment by handling end connectors only; do not pull or twist tubing.
- Keep product away from heat or ignition sources, especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen, or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- Do not use the product during X-rays, MRIs, or other radio frequency procedures.
- Refer to MR850 Respiratory Humidifier user instructions for Electromagnetic Compatibility information.
- This product is safe for use on children.
- Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.
- If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Assembly Instructions

Refer to Diagram **A** for heater-wire assembly.

**WARNINGS**

- Do not assemble heater-wire into single-use tubing.
- To avoid damaging heater-wire insulation, the draw-wire hook should only be attached to the heater-wire O-ring.
- Ensure heater-wire is extended along the tube length and is not bunched at any point.
- Do not use excessive force when drawing heater-wire through tubes.
- Do not use the draw wire to tighten the tube and heater-wire elbow together.

Notes

- Ensure the chamber end probe is securely inserted into the heater-wire elbow port.

Monitoring Instructions**WARNING**

- Visually inspect the heater-wire before each use on a patient. Discard if faulty or if there is any sign of insulation deterioration such as cracks, tears, or damage. These may lead to electrical hazards and possible risk of serious patient harm.

Reprocessing Instructions

For reprocessing instructions (UI-623716), visit <http://www.fphcare.com/fp-reusables>

WARNINGS

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever it is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.
- Do not expose heater-wires to reprocessing temperatures >100 °C (212 °F).




Notes

- It is recommended that reprocessing should commence as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state, and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories should be stored in a clean plastic bag or container and away from direct sunlight or heat if not reusing immediately.

Symbol Definitions

		
Product does not contain phthalates or pyrogens	Unsafe for use in MRI environments	European Union authorized representative
		
Transportation and storage temperature limits	Not made with natural rubber latex	Caution symbol for warnings
		
CE Mark	Use-by date	Prescription only
		
Type BF Applied Part	Date of manufacture	Instructions for Use
		
Lot number	Manufacturer	Medical Device
		
Reference number	Do not recycle	

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