

Neonatal reusable low-flow circuit with watertrap









USER INSTRUCTIONS

REF **900MR781**

Rx only  0123     

Manufacturer  Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

C**1****2****3****4****5****D****E****F**

English	2
български (Bulgarian).....	7
Hrvatski (Croatian).....	12
Česky (Czech)	17
Dansk (Danish).....	22
Nederlands (Dutch).....	27
Eesti keel (Estonian).....	32
Suomi (Finnish).....	37
Français (French)	42
Deutsch (German)	47
Ελληνικά (Greek)	52
Magyar (Hungarian).....	57
Italiano (Italian).....	62
Latviešu (Latvian)	67
Lietuvių (Lithuanian)	72
Norsk (Norwegian).....	77
Polski (Polish)	82
Português (Portuguese)	87
Português (Brasil) (Brazilian Portuguese)	92
Română (Romanian)	97
Slovenčina (Slovak).....	102
Slovenščina (Slovenian)	107
Español (Spanish).....	112
Svenska (Swedish).....	117
Türkçe (Turkish)	122
繁體中文版 (Traditional Chinese).....	127
Bahasa Indonesia (Indonesian)	132
ไทย (Thai).....	137
Tiếng Việt (Vietnamese).....	142
Русский (Russian)	147
한국어 (Korean).....	152
日本語 (Japanese)	157

REF 900MR781 Neonatal reusable low-flow circuit with watertrap**Intended Use**

The 900MR781 is a reusable breathing circuit for delivery of heated and humidified respiratory gases to neonatal patients requiring respiratory support. It is intended for use with Fisher & Paykel Healthcare MR850 humidifiers in hospital or clinical environments and should only be prescribed by a physician. Addition of heat and humidity to the supply of cold and dry respiratory gases provided through ventilation is beneficial to prevent drying of the patient's airways.

Product Specifications

The following product specifications are based on testing done with the MR850 humidifier and the MR340 neonatal reusable chamber.

Tidal volume	≤50 mL
Interface connections	ISO 5356-1 Conical Connectors
Resistance to flow at 2.5 L/min (maximum water level)	
Inspiratory	0.12±0.02 cmH ₂ O
Expiratory	0.06±0.02 cmH ₂ O
Resistance to flow at 2.5 L/min (minimum water level)	
Inspiratory	0.13±0.02 cmH ₂ O
Expiratory	0.07±0.02 cmH ₂ O
Compliance at 60 cmH₂O (maximum water level)	0.67±0.02 mL/cmH ₂ O
Compliance at 60 cmH₂O (minimum water level)	0.83±0.02 mL/cmH ₂ O
Maximum operating pressure	80 hPa (cmH ₂ O) (8 kPa)
Operating flow range with Infant Warmer	
Invasive (>33 mg/L)	4–30 L/min
Noninvasive (>12 mg/L)	2–30 L/min
Operating flow range with Incubators	
Invasive (>33 mg/L)	2–30 L/min
Noninvasive (>12 mg/L)	2–30 L/min
Operating ambient temperatures	20–26 °C (68–79 °F)
Inlet gas temperature range	Same as ambient temperature
Gas leakage rate	<30 mL/min (at 60 cmH ₂ O)
Minimum tube internal diameter	10 mm (0.4 in)
Shelf life	5 years from date of manufacture

Warnings, Cautions and Notes

Warnings

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing or accessories. Failure to comply may result in fire, burns and loss of system performance.
- Heated tubes should be kept out of direct contact with the patient's skin. Failure to comply may result in skin damage. At flow rates below 2 L/min, the tube surface temperature may exceed 44 °C (111 °F).
- Covering breathing tubes with blankets, pillows, towels or bed linen or heating them in an incubator or with an overhead heater can affect the quality of the therapy or injure the patient.
- Do not exceed maximum peak flow rate of 70 L/min. Excessive flow rate can cause liquid to splash into the breathing circuit.
- These products are reprocessable (reusable). The products must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Use of nebulizers may affect the performance of the humidifier and temperatures of delivered respiratory gases. It is recommended to place nebulizers downstream of all airway temperature probes.
- Do not use the product above altitudes of 3000 m (9842 ft) or outside the recommended ambient temperature or flow range. Using the humidification system outside of specified ranges or above this altitude can affect the quality of the therapy or injure the patient.
- Use of non-compatible, non-ISO 5367, or non-ISO 80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Do not modify these products.
- Do not use breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare.

Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Keep product away from heat or ignition sources especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- Do not use the product during X-rays, MRIs or other radio-frequency procedures.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.

If a serious incident has occurred while using this device please contact your local Fisher & Paykel Healthcare representative and Competent Authority.

Assembly Instructions

Warnings

- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damaged insulation. These may cause gas leaks, electrical hazards, loss of ventilation or respiratory support and risk of serious patient harm.
- Before connecting to a patient, ensure the complete circuit functions correctly with required ventilator settings and appropriate ventilator or flow source alarms are set.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted turn the humidifier off.
- Perform a pressure and leak test on the breathing circuit and check for occlusions before connecting to a patient.
- Check all connections are tightly secured before use.
- To avoid damage to breathing tubes, attach and detach from equipment by handling end connectors only. Do not stretch, crush or twist tubing.
- Ensure the water trap and chamber are positioned below the patients head so that condensate flows away from the patient interface.
- Do not fill the humidification chamber above the maximum fill-level line. Liquid could enter the breathing circuit if the chamber is overfilled.
- Do not fill the chamber with water with a temperature in excess of 37 °C (99 °F).
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Ensure the MR340 chamber o-ring is not twisted when inserted into the chamber base groove and that the chamber dome is pushed fully onto the chamber base as per the chamber assembly instructions.
- Ensure the humidifier and chamber are level and positioned lower than the patient. Operating the chamber at an angle in excess of 10° may result in water entering the breathing circuit.
- The circuit is powered through the heater-wire adaptor from the MR850. Refer to MR850 humidifier user instructions for further set up information including Electromagnetic Compatibility.

- 1 Refer to diagram **A** for heater-wire assembly.
 -  Ensure heater-wire is extended along the tube length and is not bunched at any point.
 -  To avoid damaging heater-wire insulation the draw wire hook should only be attached to the heater-wire o-ring.
 -  Do not use excessive force when drawing the heater-wire through tubes.
 -  Do not use the draw wire to tighten the tube and heater-wire elbow together.
- 2 Assemble reusable circuit components as per diagram **B**. Refer to the 900MR139 water trap user instructions for assembly instructions and additional information.
 - Do not use 0.3 m (12 in) inspiratory extension tube for infant radiant warmer patients or in ambient environments. Use under these environments may affect the quality of the therapy. For incubator patients, ensure the patient-end temperature probe remains outside the incubator.
- 3 Assemble the chamber onto the heater base as per diagram **C**. Refer to the MR340 chamber user instructions for chamber assembly instructions and additional information.
- 4 Connect the heater-wire adaptor, chamber-end probe and patient-end probe to the circuit.
 - Refer to diagram **D** for circuit setup.
 - Ensure the chamber-end probe is securely inserted into the heater-wire elbow port.
 - Ensure the pressure port is on the top of the patient Y-piece. Do not have the pressure elbow positioned so condensate can block the elbow hole. This may cause loss of ventilation or respiratory support.

Notes

- Refer to diagram **E** for individual component REF numbers.
- Wear and degradation of electrical connectors may affect the product performance. Inspect the product components between uses.

Monitoring Instructions

Warnings

- Appropriate patient monitoring (e.g., oxygen saturation) must be used at all times. Failure to monitor the patient may result in serious harm or death (e.g., in the event of an interruption to gas flow).
- Regularly monitor and drain condensate build-up in the circuit. Condensate build-up may result in loss of ventilation or respiratory support.
- If a water feed set is used, ensure the water bag is placed below the level of the humidifier after the chamber has been filled.
- Do not touch the heater-plate or chamber base. Surfaces may exceed 85 °C (185 °F). Failure to comply may result in a skin burn.

Chamber

-  Ensure the chamber contains water at all times. Check the water level in the humidification chamber regularly as per diagram **E** and refill as required. The absence of water will result in low levels of humidity.

Water Trap

-  Monitor the water trap's liquid level regularly and empty as required. The water container should be reassembled as soon as possible after emptying, to prevent gas leakage. Refer to diagram **F**.

Temperature Probes

-  Ensure the patient-end probe is securely inserted into the inspiratory limb probe port (as per diagram **D**) and positioned vertically downwards to prevent condensate pooling around the sensor and affecting the airway gas temperature measurement.

Ensure both temperature probes remain fully inserted into the inspiratory limb during use.

Notes

- Avoid fans and air conditioner venting onto the breathing circuit during use as this may increase condensate in the circuit and potentially block a tube, resulting in loss of ventilation or respiratory support.

Reprocessing Instructions

For reprocessing instructions please refer to (UI-623716), or go to www.fphcare.com/fp-reusables

Warnings

- The products must be disconnected from both the power and gas source prior to cleaning.
- Clean product prior to use and whenever product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

Warning

- Do not expose breathing tubes to excessive UV radiation.

Note

- As a minimum, reprocessed breathing circuits and accessories should be stored in a clean plastic bag or container and away from direct sunlight or heat if not reusing immediately.

Symbol Definitions

 Product does not contain phthalates or pyrogens	 For use on neonatal patients	 Transportation and storage temperature limits	 Type BF Applied Part	 European Union authorized representative
 Not made with natural rubber latex	 Instructions for Use	 Unsafe for use in MRI environments	 Caution symbol for warnings	 Date of manufacture
 Manufacturer	 Incorrect water level	 Correct water level	 No water	 CE Mark
 Do not recycle	 Use-by date	Rx only Prescription only	 Medical Device	 Lot number
 Reference number				

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