

Reusable Adaptor – Flow Probe

**USER INSTRUCTIONS** REF **900MR830**

---

**A**

English .....	4
български (Bulgarian).....	6
Hrvatski (Croatian).....	8
Česky (Czech).....	10
Dansk (Danish).....	12
Nederlands (Dutch).....	14
Eesti keel (Estonian).....	16
Suomi (Finnish).....	18
Français (French).....	20
Deutsch (German).....	22
Ελληνικά (Greek).....	24
Magyar (Hungarian).....	26
Italiano (Italian).....	28
Latviešu (Latvian).....	30
Lietuvių (Lithuanian).....	32
Norsk (Norwegian).....	34
Polski (Polish).....	36
Português (Portuguese).....	38
Português (Brasil) (Brazilian Portuguese).....	40
Română (Romanian).....	42
Slovenčina (Slovak).....	44
Slovenščina (Slovenian).....	46
Español (Spanish).....	48
Svenska (Swedish).....	50
Türkçe (Turkish).....	52
繁體中文版 (Traditional Chinese).....	54
Bahasa Indonesia (Indonesian).....	56
ไทย (Thai).....	58
Tiếng Việt (Vietnamese).....	60
Русский (Russian).....	62
한국어 (Korean).....	64
日本語 (Japanese).....	66

**REF 900MR830****Reusable Adaptor – Flow Probe****Intended Use**

The 900MR830 adaptor allows chambers and circuits compatible with MR600 and MR700 series humidifiers to be used with MR850 humidifiers. This accessory is intended for use with Fisher & Paykel Healthcare reusable breathing circuits in hospital or clinical environments and should only be prescribed by a physician.

**Product Specifications**

**Maximum Operating Pressure** 80 hPa (cmH<sub>2</sub>O)  
(8 kPa)

**Warning, Cautions and Notes****WARNINGS**

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing or accessories. Failure to comply may result in fire, burns and loss of system performance.
- This product is reprocessible (reusable). It must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Check all connections are tightly secured before use.
- Use of non-compatible, non-ISO 5367, or non-ISO 80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Do not modify this product.

**Failure to comply with the above warnings and cautions may impair the performance of the product or compromise safety (including causing potential serious harm).**

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)

Contact your local Fisher & Paykel Healthcare representative for a printed copy of these user instructions.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

**Notes**

- Refer to diagram **A** for assembly instructions.
- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended if products can withstand the process.
- Ensure the complete circuit functions correctly, with the required ventilator settings, before connecting to a patient.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on pregnant women and breastfeeding women.
- This product is safe for use on children.

## Reprocessing Instructions

For reprocessing instructions (UI-623716), visit <http://www.fphcare.com/fp-reusables>

### WARNINGS

- The product must be disconnected from both the power and gas source prior to cleaning.
- Clean the product prior to use and whenever the product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

### Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

### Storage

- As a minimum, reprocessed breathing circuits and accessories are to be stored in a clean plastic bag or container and away from direct sunlight or heat if not being reused immediately.

## Symbol Definitions

Product does not contain phthalates or pyrogens	European Union authorized representative	Transportation and storage temperature limits
Not made with natural rubber latex	Reference number	Caution symbol for warnings
Use-by date	Prescription only	Type BF Applied Part
Date of manufacture	Instructions for Use	Lot number
Manufacturer	Medical Device	Do not recycle
CE Mark		

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