
Temperature Probe

USER INSTRUCTIONS

Rx only








-10°C
14°F 50°C
122°F

10 95
%
%


Fisher & Paykel
HEALTHCARE

Reprocessing Diagram


C**c1**

Equipment	Setting	Reprocessing Cycles
STERIS V-PRO 1	Standard Cycle	
STERIS V-PRO 1Plus		78
STERIS V-PRO 60	Non-Lumen Cycle	
STERIS V-PRO maX		
Sterrad 100NX	Standard Cycle	*78
VP4 STERIZONE	Cycle 1	78

D**d1****d2**

Setup Diagram

English	6	Ελληνικά	32
Français	7	Português	34
Deutsch	9	ไทย	35
Nederlands	11	Hrvatski	37
Español	13	Česky	38
Italiano	15	Magyar	39
Português (Brasil)	16	Latviešu	40
Svenska	17	Lietuvių	41
Русский	18	Română	42
Türkçe	20	Slovenčina	43
日本語	21	Slovenština	44
한국어	22	български	45
Dansk	23	Eesti keel	47
Norsk	24	Tiếng Việt	48
Suomi	25	Español (latinoamérica)	49
Polski	26	Українська	50
简体中文版	28	العربية	52
繁體中文版	30	Srpski	53
Bahasa Indonesia	31	Македонски	54

Temperature Probe User Instructions

Intended Use:

The Temperature Probe is designed for measuring gas temperature and flow with Fisher & Paykel Healthcare MR850 and HC550 Humidifiers and Fisher & Paykel Healthcare Breathing Circuits.

Refer to the MR850 Respiratory Humidifier instructions for use for additional safety information, clinical claims and instructions for setup and use.

Specification:

MODEL OVERALL LENGTH

900MR860 1820 mm (71.5")

900MR861 2400 mm (94.4")

900MR863 1240 mm (48.7")

900MR868 1620 mm (63.7")

900MR869 2100 mm (82.6")

Expected Service Life:

3 years; if used in accordance with this User Instruction. Note: Service life may vary depending on usage pattern.

Usage: Multiple patient multiple use.

Setup:

Refer to the setup diagram.

1. Before use remove and discard **Protective Sleeve** (D), if present, from **Chamber Probe** (E).
2. Visually inspect entire product for damage before use. Discard if damaged.
3. Push the **Airway Probe** (A) and **Chamber Probe** (E) into the Breathing Circuit making sure they are correctly located and pushed into place.
4. The probe lead can be restrained using Breathing Circuit Clips.
- If patient is in an incubator, refer to Breathing Circuit User Instructions for probe positioning information.
5. Insert the **Temperature Probe Connector** (C) into the blue socket on the Humidifier.
6. Perform ventilator leak test on the Breathing Circuit before use.

Maintenance:

Probe maintenance and accuracy checks should be routinely performed. For further information, refer to the MR850 or HC550 Technical Manual (available from your local Fisher & Paykel Healthcare representative).

Reprocessing:

The following parts are in contact with the gas path and can become contaminated with body fluids or expired gases:

Airway Probe (A), Chamber Probe (E)

Recommended cleaning procedures between uses (refer to reprocessing diagram):

1. Remove any visible soil using a soft bristle brush, water and a mild non-abrasive detergent.
2. Do not immerse the **Temperature Probe Connector** (C) in any liquid.

3. Rinse in water after cleaning.
4. Dry with a non-shedding cloth.
5. Check for visible signs of damage. Discard if damaged. Repeat the cleaning steps if parts are not visually clean.
6. Proceed to either step A, B, C or D depending on hospital guidelines.

A. MANUAL DISINFECTION - WIPES

- a1. Follow the wiping instructions supplied by the manufacturer.
- a2. Check for visible signs of damage. Discard if damaged.
- a3. If applicable, product can be packaged for reuse.

B. MANUAL DISINFECTION - CHEMICAL

- b1. Follow the soaking instructions supplied by the manufacturer of the disinfection solution.
- Do not immerse the **Temperature Probe Connector** (C) in any liquid. Cover container while soaking.
- b2. Rinse in water after chemical soaking.
- b3. Dry with a non-shedding cloth.
- b4. Check for visible signs of damage. Discard if damaged.
- b5. If applicable, product can be packaged for reuse.

C. VAPORISED HYDROGEN PEROXIDE (VHP) STERILIZATION

- c1. Follow hospital guidelines for using the equipment (**Equipment**) and settings (**Setting**).

*Note: Some tackiness may be observed after 25 cycles; probe performance is unaffected by this condition.

D. ETHYLENE OXIDE (ETO) STERILIZATION

- d1. Follow hospital guidelines for process. d2. Allow at least 12 hours aeration.

6. Store as per hospital guidelines.

Notes:

- Fisher & Paykel Healthcare do not recommend exceeding the number of reprocessing cycles stated (**Reprocessing Cycles**).
- The product is not designed for autoclave sterilization processes.
- Follow these recommended cleaning procedures and the user instructions of the cleaning agents or disinfection/sterilization equipment. It is the user's responsibility to qualify any deviations from these procedures, both for disinfecting efficacy and physical effect.

Warnings:

- The use of respiratory humidifier, breathing circuits, chambers, other accessories or parts which are not approved by Fisher & Paykel Healthcare may impair performance or compromise safety.
- Failure to perform routine probe accuracy and visible damage checks may impair performance or compromise safety.
- Ensure that both Temperature Probe sensors are correctly and securely fitted. Failure to do so may result in gas temperatures in excess of 41 °C being delivered to the patient.
- Do not touch the glass tip of the **Chamber Probe** (E) during use; it may cause a skin burn.
- Gas mixes, such as helium-oxygen mixtures, that have different physical or thermal properties than air or air-oxygen mixture, may impair system performance or compromise safety.
- No automated method for cleaning and disinfection is recommended.
- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

Notes

- For use under the supervision of trained medical personnel.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

Disposal:

Dispose according to standard hospital guidelines for electrical and electronic equipment.

Symbol Definition:

	Operating Instructions
	Date and Country of Manufacture NZ: New Zealand
	Manufacturer
	Dispose of product in correct manner
	Transportation and storage humidity limitation
	Transportation and storage temperature limitation
	Lot number
	Reference number
	European union authorised representative
	Prescription only
	European Conformity
	Medical Device
	Switzerland authorised representative
	UK responsible person
	Importer
	Distributor

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