


Temperature Probe

USER INSTRUCTIONS



Rx only

Fisher & Paykel
HEALTHCARE

Reprocessing Diagram

C

c1

STERILE VHP

Equipment	Setting	Reprocessing Cycles
STERIS V-PRO 1	Standard Cycle	78
STERIS V-PRO 1 Plus	Non-Lumen Cycle	
STERIS V-PRO 60		
STERIS V-PRO maX		
Sterrad 100NX	Standard Cycle	*78
VP4 STERIZONE	Cycle 1	78


D

d1


STERILE EO

d2

Setup Diagram

English	6	繁體中文版	24
Français	7	Bahasa Indonesia	25
Deutsch	8	Ελληνικά	26
Nederlands	9	Português	27
Español	10	ไทย	28
Italiano	11	Hrvatski	29
Português (Brasil)	12	Česky	30
Svenska	13	Magyar	31
Русский	14	Latviešu	32
Türkçe	15	Lietuvių	33
日本語	16	Română	34
한국어	17	Slovenčina	35
Dansk	18	Slovenščina	36
Norsk	19	български	37
Suomi	20	Eesti keel	38
Polski	21	Tiếng Việt	39
简体中文版	22	Español (latinoamérica)	40

Intended Use:

The Temperature Probe is designed for measuring gas temperature and flow with Fisher & Paykel Healthcare MR850 and HC550 Humidifiers and Fisher & Paykel Healthcare Breathing Circuits.

Specification:


MODEL	OVERALL LENGTH
900MR860	1820 mm (71.5")
900MR861	2400 mm (94.4")
900MR863	1240 mm (48.7")
900MR868	1620 mm (63.7")
900MR869	2100 mm (82.6")

Expected Service Life:

3 years, if used in accordance with this User Instruction. Note: Service life may vary depending on usage pattern.

Setup:

Refer to the setup diagram.

- Before use remove and discard **Protective Sleeve** , if present, from **Chamber Probe** .
- Visually inspect entire product for damage before use. Discard if damaged.
- Push the **Airway Probe**  and **Chamber Probe**  into the Breathing Circuit making sure they are correctly located and pushed into place.
- The probe lead can be restrained using Breathing Circuit Clips.
 - If patient is in an incubator, refer to Breathing Circuit User Instructions for probe positioning information.
- Insert the **Temperature Probe Connector**  into the blue socket on the Humidifier.
- Perform ventilator leak test on the Breathing Circuit before use.

Maintenance:

Probe maintenance and accuracy checks should be routinely performed. For further information, refer to the MR850 or HC550 Technical Manual (available from your local Fisher & Paykel Healthcare representative).

Reprocessing:

The following parts are in contact with the gas path and can become contaminated with body fluids or expired gases:

Airway Probe , **Chamber Probe** .

Recommended cleaning procedures between

uses (refer to reprocessing diagram):

- Remove any visible soil using a soft bristle brush, water and a mild non-abrasive detergent.
Do not immerse the **Temperature Probe Connector**  in any liquid.
- Rinse in water after cleaning.
- Dry with a non-shedding cloth.
- Check for visible signs of damage. Discard if damaged. Repeat the cleaning steps if parts are not visually clean.
- Proceed to either step A, B, C or D depending on hospital guidelines.

A. MANUAL DISINFECTION – WIPES

- Follow the wiping instructions supplied by the manufacturer.
- Check for visible signs of damage. Discard if damaged.
- If applicable, product can be packaged for reuse.

B. MANUAL DISINFECTION – CHEMICAL

- Follow the soaking instructions supplied by the manufacturer of the disinfection solution.
Do not immerse the **Temperature Probe Connector**  in any liquid. Cover container while soaking.
- Rinse in water after chemical soaking.
- Dry with a non-shedding cloth.
- Check for visible signs of damage. Discard if damaged.
- If applicable, product can be packaged for reuse.

C. VAPORISED HYDROGEN PEROXIDE (VHP) STERILIZATION

- Follow hospital guidelines for using the equipment (**Equipment**) and settings (**Setting**).

*Note: Some tackiness may be observed after 25 cycles; probe performance is unaffected by this condition.

D. ETHYLENE OXIDE (ETO) STERILIZATION

- Follow hospital guidelines for process.
- Allow at least 12 hours aeration.


6. Store as per hospital guidelines.

Notes:

- Fisher & Paykel Healthcare do not recommend exceeding the number of reprocessing cycles stated (**Reprocessing Cycles**).

- The product is not designed for autoclave sterilization processes.
- Follow these recommended cleaning procedures and the user instructions of the cleaning agents or disinfection/sterilization equipment. It is the user's responsibility to qualify any deviations from these procedures, both for disinfecting efficacy and physical effect.





Warnings:

- The use of breathing circuits, chambers, other accessories or parts which are not approved by Fisher & Paykel Healthcare may impair performance or compromise safety.
- Failure to perform routine probe accuracy and visible damage checks may impair performance or compromise safety.
- Ensure that both Temperature Probe sensors are correctly and securely fitted. Failure to do so may result in gas temperatures in excess of 41 °C being delivered to the patient.
- Do not touch the glass tip of the **Chamber Probe**  during use; it may cause a skin burn.
- Gas mixes, such as helium-oxygen mixtures, that have different physical or thermal properties than air or air-oxygen mixture, may impair system performance or compromise safety.
- No automated method for cleaning and disinfection is recommended.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

Disposal:

Dispose according to standard hospital guidelines for electrical and electronic equipment.

Symbol Definition:

	Operating Instructions
	Date of manufacture
	Manufacturer
	Dispose of product in correct manner
	Transportation and storage humidity limitation
	Transportation and storage temperature limitation
	Lot number
	Reference number
	European union authorised representative
	Prescription only
	European Conformity
	MD symbol

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