



Single Use **Rx only**

For patent information, see www.fphcare.com/ip

Fisher & Paykel HEALTHCARE

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

Indications for use - Airvo nebulizer adapter:**When used with a nebulizer:**

The nebulizer adapter is a medical device accessory for single-patient use to facilitate aerosolization of Albuterol sulfate for inhalation to adult patients receiving high-flow humidified breathing gases via tracheostomy patient interface. Intended for use by healthcare professionals in hospitals or long-term care facilities.

When used without a nebulizer:

For use in hospitals and long-term care facilities, for the delivery of humidified respiratory gases to patients via nasal, tracheostomy and mask interfaces.

System Specifications:

Circuit Length: 1.8 m (6 ft) + interface.

Use WITHOUT a nebulizer

- If never used with a nebulizer, the 900PT563 kit must be disposed of after a maximum of 14 days' use.
- Do not soak, wash or sterilize.

Setup:

- Airvo 2 humidifier ("Airvo") on a flat surface.
- For use at flows from 2 to 60 L/min depending on the patient interface.
- OPT316/OPT318/OPT942/OPT944/OPT946/OPT946/OPT1042/OPT1044/OPT1046 nasal interfaces OR OPT970 tracheostomy interface OR OPT980 mask adapter interface

WARNING: Use of a non-approved accessory could impair performance or compromise safety.

Steps 1 - 7:

Assemble and connect water chamber.

Step 8:

Ensure water level is rising but not going over the maximum water level.



Incorrect water level
Replace MR290 chamber



Correct water level in the
MR290 chamber

Use USP sterile/distilled water for inhalation, or equivalent.

Step 9:

Connect breathing tube to Airvo.

Step 10:

Make sure blue connector is fully located into place.

Step 11:

Connect breathing tube to patient interface once the system has warmed up. Position breathing tube below patient interface so that condensate flows away from the patient.

Step 12:

Connect breathing tube clip to patient clothing or bedding. Note: connecting to patient clothing may not be suitable for all patient groups.

Step 13:

In the absence of a nebulizer, make sure the nebulizer port plug is firmly in place.

- **WARNING:** Loss of therapy may occur if the nebulizer port plug is not in place.

Use WITH a nebulizer

- If used with a nebulizer at any time, the 900PT563 kit must be disposed of after a maximum of 7 days' use, whether nebulization is continuous or intermittent during this time.
- Do not soak, wash or sterilize.

CAUTION: Using the kit with a nebulizer for longer than 7 days may cause damage to the kit, particularly the water chamber, and eventually the Airvo itself.

Setup:

WARNING: Nasal delivery of nebulized drugs for lung deposition is not approved by FDA.

WARNING: Use of a non-approved accessory could impair performance or compromise safety.

- Airvo 2 humidifier ("Airvo") on a flat surface.
- For use at flows from 10 to 30 L/min using the OPT970 patient interface.
- OPT970 tracheostomy interface
- Aerogen Solo nebulizer

Compatible medications:

- **This kit is compatible with the following medications only:**
 - **Albuterol sulfate**
- Refer to drug manufacturers' instructions for correct usage and restrictions.
- **WARNING:** Do not nebulize the following medications, as they are incompatible with Airvo and its consumables and may cause patient harm:
 - Alcohol-based medications
 - Significantly acidic or basic medications eg. epoprostenol (Flolan)
 - Saline >0.9%

Step 14 (Nebulizer use):

- Follow Steps 1-13.
- Remove the nebulizer port plug from the nebulizer port and connect the nebulizer in its place.
 - CAUTIONS:**
 - Refer to the nebulizer manufacturer's instructions for setup and usage information.
 - Do not add oxygen or other gases through the nebulizer port or use a flow-generating nebulizer (eg. jet nebulizer) as this will cause discrepancies in humidity, flow and oxygen outputs.
 - Ensure the Airvo is switched on with humidified gas flowing before the nebulizer is used.
 - The medication delivery efficiency decreases as the Airvo flow increases. Ensure that the patient receives adequate flow at all times when adjusting flow during nebulization. Titrate medication delivery as required.
 - Nebulized medications can form residue and cause cosmetic discoloration of the breathing tube and chamber during normal operation.
 - Nebulization may cause increased levels of condensate in the breathing tube and chamber.
 - Avoid skin contact with condensate as it may contain nebulized medication.

Performance:

Albuterol sulphate (1 mg/mL, 2.5 mL)	10 L/min	20 L/min	30 L/min
Total delivered dose (µg)*	1362.9 - 2087.7	548.6 - 1938	383.3 - 1461.8
Total respirable dose (µg, 1-5µm)*	1035.9 - 1550.3	470.7 - 1428.2	387.9 - 837.3
Respirable fraction (%)	63 - 87.1	59.3 - 93.8	31.3 - 100
Coarse particle dose (µg, >4.7µm)*	206.6 - 528.2	0 - 507.8	0 - 531.9
Coarse particle fraction (>4.7µm)	14.9 - 27.6	4.7 - 32.4	0 - 47.2
Fine particle dose (µg, <4.7µm)*	1111.3 - 1604.5	519 - 1500.8	448.5 - 987.4
Fine particle fraction (<4.7µm)	72.4 - 85.1	67.6 - 95.4	52.8 - 100
Mass-Median Aerosol Diameter (MMAD) (µm)*	2.67 - 2.94	2.41 - 3.05	2.08 - 3.00
Geometric Standard Deviation (GSD) (µm)*	1.56 - 2.30	1.18 - 2.58	0.1 - 4.91

*** 95% confidence intervals**

Testing carried out with one Aerogen Solo nebulizer, three sets of Airvo 2, 900PT563 kit and OPT970 interface and three tests per set. Cascade impactor testing conducted at 15 L/min using NGI.

WARNINGS:

- For single patient use only. Reuse may result in transmission of infectious substances. Attempting to reprocess will result in degradation of materials and render the product defective.
- Using the breathing tube for longer than the specified time can result in serious injury, including infection.
- Do not use the auto-fill MR290 chamber if it has been dropped or been allowed to run dry as this could lead to the chamber over-filling.
- Do not use the MR290 chamber if the water level rises above the maximum water level line as this may lead to water entering the patient's airway.
- Never operate the unit if the breathing tube has been damaged with holes, tears or kinks.
- Do not block the flow of air through the unit and breathing tube.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

To avoid burns:

- Do not allow the breathing tube to remain in direct contact with skin for prolonged periods of time.
- Do not use in the presence of a naked flame, to avoid fires.
- Do not add heat to any part of the breathing tube e.g. covering with a blanket, or heating it in an incubator or overhead heater for a neonate, as this could result in serious injury.

Cautions:

- Avoid product contact with chemicals, cleaning agents, or hand sanitizers.
- Use only USP sterile/distilled water for inhalation, or equivalent.
- Nebulized medications may be emitted into the surrounding area.
- Condensate containing nebulized medication may be present after use. Wear appropriate personal protective equipment such as gloves.

To prevent condensation:

- Use in a room warmer than 18 °C (64 °F).
- Remove/minimize the impact of anything that may cool the heated breathing tube, e.g. a fan, air conditioning, open window.

To manage excessive condensation:

- Place the Airvo below patient head height.
- Drain condensate back into the water chamber. At higher target flow settings it may be necessary to first reduce the target flow setting to 30 L/min, to ensure that the condensate drains safely and effectively.
- Disconnect the patient interface from the heated breathing tube.
- Lift the patient end of the heated breathing tube, allowing the condensate to run into the water chamber.