CE 0123 Rx only  fphcare.com/
neopuff-ifu   

F&P and Neopuff are trademarks of Fisher & Paykel Healthcare Limited.
REF 185043544 REV K 2023-05 © 2023 Fisher & Paykel Healthcare Limited

For patent information,
see www.fphcare.com/ip

 Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box
14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz
Web: www.fphcare.com

 **Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street,
Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel
do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil
Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有
限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486
France (FR)  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5,
BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201
Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49,
73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032
India (IN) Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)**
Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel
Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464
Russia (RU) Tel: +7 495 782 21 50 **Switzerland (CH)**  Fisher & Paykel Healthcare
GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)**
Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret
Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara,
Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit
16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189
USA (US)/Canada (CA) Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)**
Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)**
Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346
Sweden (SE) Tel: +46 8 564 76 680


www.fphcare.com

Fisher & Paykel
HEALTHCARE

REF 900RD110


F&P Neopuff™ Humidified Resuscitation Circuit

Instructions for use


www.fphcare.com



Fisher & Paykel
HEALTHCARE

1

2

English	6	Português (Brasil)	72
العربية (Arabic)	9	(Brazilian Portuguese)	
български (Bulgarian)	12	Română (Romanian)	75
Čeština (Czech)	15	Русский (Russian)	78
Dansk (Danish)	18	Slovenčina (Slovak)	81
Deutsch (German)	21	Slovenščina (Slovenian)	84
Ελληνικά (Greek)	24	Svenska (Swedish)	87
Español (Spanish)	27	ภาษาไทย (Thai)	90
Eesti keel (Estonian)	30	Türkçe (Turkish)	93
Suomi (Finnish)	33	اردو (Urdu)	96
Français (French)	36	Tiếng Việt (Vietnamese)	99
Hrvatski (Croatian)	39	繁體中文 (Traditional Chinese)	102
Magyar (Hungarian)	42	简体中文 (Simplified Chinese)	105
Bahasa Indonesia (Indonesian)	45		
Italiano (Italian)	48		
한국어 (Korean)	51		
Lietuvių (Lithuanian)	54		
Latviešu (Latvian)	57		
Nederlands (Dutch)	60		
Norsk (Norwegian)	63		
Polski (Polish)	66		
Português (Portuguese)	69		

900RD110 F&P Neopuff Humidified Resuscitation Circuit

English (en)

Intended purpose

The F&P Neopuff Humidified T-piece Circuit (900RD110) is a single use breathing circuit intended to provide heated, humidified ventilatory support to neonates and infants with respiratory insufficiency. The device is designed to connect to the F&P Neopuff Infant T-piece Resuscitator and F&P MR850 humidifier. This device is designed for use in the hospital environment and must be prescribed by a physician. It is intended to be used by medical professionals including physicians, nurses, midwives and respiratory therapists. The intended population for use is neonates and infants.

Instructions for use

Consult the Labelling for your Neopuff T-piece Resuscitation Device, MR850/MR810* Humidifier and MR290/MR225* chambers for additional set-up & use instructions and additional warnings.

*Not available in all markets

An alternative non-gas powered means of resuscitation (such as a self-inflating bag) must be available.

Set-up

Do not use if package is damaged

Use Sterile water in the humidifier's Chamber

① Fit a chamber to the humidifier.

Slide a MR290 chamber onto the humidifier base. Remove the blue caps. Hang the water bag, unwind the waterfeed set and spike water bag. The bag should be at least 50 cm above the chamber. Use USP sterile water for inhalation or equivalent.

Or

Slide a MR225 chamber onto the humidifier base. Remove the blue caps. Fill with 30 mL of sterile water or equivalent. This will allow for 60 minutes of running time (@ 8 L/min flow rate) and then refill as necessary during resuscitation.

② Connect Circuit

- Connect the dry-line between the chamber inlet port and the Neopuff gas outlet.
- Connect the blue inspiratory limb to the remaining chamber port.
- Ensure the temperature probes and heater wire adaptors are fully inserted. If temperature probes are not used (i.e. using the MR810), ensure caps (BC410) are placed over the temperature probe ports.

③ Connect to flow source

- Connect a blended oxygen/air supply or oxygen to the gas inlet port using gas supply line.
 - The Humidified Resuscitation System is to be used with medical grade oxygen/air**
- Adjust the flow rate to the prescribed input flow rate between 8 and 10 L/min.
 - DO NOT attempt to use product outside this flow range. Recommended flow rate 8 L/min. For infants weighing ≥ 8 kg, a flow rate of 10 L/min is recommended.**

④ Set humidifier

Turn on the humidifier base.

MR850 – Check the humidifier is set to invasive mode. If not, hold the mode button until the invasive mode light is lit.

MR810 – Set the humidifier to high mode. Ensure the green heater wire light is lit.  Caution is required when delivering gases that are not body temperature saturated, to patients with bypassed airways.

⑤ Set PIP

Ensure the test lung is inserted into the patient T-piece. Occlude PEEP (Positive End Expiratory Pressure) cap and adjust blue PIP (Peak Inspiratory Pressure) knob on the Neopuff (or similar on other T-piece Resuscitators) to set PIP.

Connect the test lung to the T-piece circuit and test resuscitator function before connecting to the patient.

⑥ Set PEEP

While the test lung is still attached to the patient T-piece, adjust the PEEP cap to set PEEP.

Resuscitate by placing and removing finger or thumb over the PEEP cap to allow inspiration and expiration.

Ensure pressures are monitored throughout resuscitation.

Technical specifications

Input Flow Range	8–10 L/min
Positive End-Expiratory (PEEP)	@8 L/min 1 to 10 cmH ₂ O [mbar] @10 L/min 2 to 13 cmH ₂ O [mbar]
T-piece Connection	15 mm female (conical connector)
Circuit Length	1.3 m
Patient Weight Range	Up to 10 kg (22 lb)
Dead Space	6 mL
Intended T-piece Connections	Standard neonatal interface
Intended Duration of Use	Maximum of 7 days
Operating Temperature	18 °C to 26 °C (+64 to +79 °F), 30 to 75% humidity
Storage Temperature	Storage: -10 to 50 °C (+14 to +122 °F), up to 95% humidity
Humidity	MR850 >33 mg/L MR810 >10 mg/L

Contraindications

There are no contraindications for using this device. Please refer to general warnings.

GENERAL WARNINGS

-  The humidified resuscitation system is designed for SINGLE USE ONLY. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- Do not smoke, use naked flames and other sources of ignition near the unit when it is in use.
 - The F&P Neopuff Humidified T-piece Circuit and F&P Neopuff T-piece Resuscitator are intended for use by medical professionals trained in infant resuscitation.
 - DO NOT soak, wash, sterilize or re-use this product.
 - DO NOT occlude the resuscitation tubing.
 - DO NOT operate the humidifier or the circuit with no gas flow.
 - Ensure the humidifier is mounted below the level of the patient.
 - DO NOT cover the circuit with materials such as blankets, towels or bed linen.
 - Avoid prolonged contact with patient's skin.
 - DO NOT stretch or milk tubing.
 - DO NOT fill the chamber above the maximum fill level line.
 - California residents please be advised of the following, pursuant to Proposition 65: This

product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

- As this medical device uses an alternative small-bore connector design different from those specified in the ISO 80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need be taken by the user to mitigate these reasonable foreseeable risks.

Risks

- There are no reported side effects associated with the use of this device.
- The use of this product is not without risk even when used as intended. Following all instructions and warnings provided, risks of hypoxia and barotrauma remain. These risks may result in serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority

Disposal

Discard according to hospital protocols.

Symbol definitions

	Consult Instructions for use www.fphcare.com/neopuff-ifu		CE Mark
	Latex free		Do not reuse
	Authorized representative in the European Community		Not made with phthalates (DEHP, DBP, BBP)
	Date of manufacture and country of manufacture		Use by date
	Storage temperature limit		Medical device
Rx only	Prescription only		Catalogue number
	Manufacturer		Batch code
	Distributor		Importer
	UKCA mark		Do not use if package is damaged

For patent information, see www.fphcare.com/ip

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