
**F&P 950**

# Bi-Level / CPAP Single Limb Circuit Kit

## USER INSTRUCTIONS

**REF** 950A60J  
950A61J

Single Use 14 Days Maximum Use

Rx only

**Fisher & Paykel  
HEALTHCARE**

English

### Indications for use

The breathing set is an accessory to the F&P 950™ Respiratory Humidifier for use in hospitals and long-term care facilities. It is intended for delivery of heated humidified respiratory gases to adult, adolescent and child patients, within the limits of its stated technical specifications.

### Contraindications

There are no contraindications for heated respiratory humidification.

Do not use the exhalation port or noninvasive therapy on patients who are non-spontaneously breathing; uncooperative or unresponsive; have copious secretions; are at risk of nausea or vomiting; or are at high risk of aspiration of vomitus.

### Side effects

There are no known side effects of heated respiratory humidification.

### Technical specifications

Compatible with the F&P 950 Respiratory Humidifier. Refer to humidifier, accessories, and patient interface user instructions for additional warnings, cautions, contraindications, side effects, and system information.

**Patient weight:** >10 kg  
**Interface connections:** ISO 5356-1 Conical Connectors

**Breathing tube length:**  
Inspiratory: 1.60 m

**Gas leakage (60 cmH<sub>2</sub>O, BTPS):**  
Single limb: <40 mL/min

**Resistance to flow (30 L/min, BTPS):**  
Inspiratory: 0.31 ± 0.03 cmH<sub>2</sub>O

*With accessories:*  
Inspiratory: 0.66 ± 0.09 cmH<sub>2</sub>O

includes measurement uncertainty of 0.02 cmH<sub>2</sub>O

**Compliance (60 cmH<sub>2</sub>O, BTPS):**  
Single limb: 1.51 ± 0.21 mL/cmH<sub>2</sub>O

*With accessories:*  
Single limb: 1.39 ± 0.29 mL/cmH<sub>2</sub>O

includes measurement uncertainty of 0.08 mL/cmH<sub>2</sub>O

**Tube minimum internal diameter:** 17 mm

**Maximum tube surface temperature:** 44 °C

**Minimum operating pressure:** 4 cmH<sub>2</sub>O

**Maximum operating pressure:** 80 cmH<sub>2</sub>O

**Operating flow range**

Mask mode: 5–120 L/min

@ 31 °C humidifier set point

Optiflow mode: 5–70 L/min

**Room temperature:** 18–26 °C

**Chamber maximum volume of water:** 160 mL

**Filter:**  
Viral: > 99.94%  
Organism: ΦX174 bacteriophage  
Bacterial: > 99.996%  
Organism: staphylococcus aureus  
Mean particle size: 3 μm  
NaCl: 96.73%

### Warnings, cautions, and notes

#### ⚠ WARNINGS

Failure to comply with the following warnings may impair the performance of the device or compromise safety (including potentially causing serious harm):

- Do not use beyond the 14-day maximum duration of use.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- Single use. Do not reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm, or death.
- Exhalation port must be used with a non-vented interface in a single limb system. Failure to comply may lead to patient inhaling excess carbon dioxide resulting in hypercapnia.
- Do not use product near any ignition source. Exposure to oxygen increases the risk of fire that may cause serious injury or death.
- Do not use this product in or near a magnetic resonance imaging (MRI) scanner.
- Do not cover the circuit with materials such as towels, pillows, or bed linen. Failure to comply may result in skin burn.
- Avoid leaving tubing in prolonged contact with patient's skin. Failure to comply may result in a skin burn.
- The use (or modification) of breathing circuit/chamber combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification, ventilator malfunction, or serious harm.
- Do not crush, stretch, or milk the tubing.
- Do not clean or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Set appropriate ventilator or flow source alarms to monitor therapy delivery.
- Ensure there is gas flow through the tubing before connecting to a patient.
- Perform a pressure and leak test on the breathing system before connecting to a patient.
- Check all connections are tight before use.
- Monitor circuit condensate to prevent occlusion or build-up of fluid. Drain as required.
- The humidifier should always be level and positioned lower than the patient.
- Do not use the chamber if it has been visibly damaged or dropped.
- Do not use the chamber if any of the sensor covers are dislodged or missing from the chamber.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not fill the chamber with water in excess of 37 °C.
- Do not block or seal the vent holes on the exhalation port.
- Replace the filter every 24 hours following standard hospital procedure. Failure to do so may result in reduced filtration efficacy.
- Monitor the filter when in use. Replace if blockage occurs due to condensation, patient secretions or nebulized drugs. Failure to do so may result in excessive carbon dioxide rebreathing.
- Do not use for ventilation via tracheostomy tube. Copious secretions produced by tracheostomy patients may block the disposable exhalation port and lead to excessive carbon dioxide rebreathing.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)


#### ⚠ CAUTIONS

- Do not touch the hot surface of the chamber base. Failure to comply may result in a skin burn.
- Do not operate the chamber without water or gas flow. Failure to comply may result in low levels of humidity.
- The water source must be at least 50 cm higher than the chamber. Failure to comply may result in low levels of humidity.

#### NOTES

- Use USP Sterile Water for Inhalation, or equivalent. Adding other substances may have adverse effects.
- User may be exposed to breathing tract fluids during disposal. Breathing circuits must be treated as contaminated. Dispose of according to local, state or country specific environmental laws and regulations.
- This product is designed for the delivery of air and/or oxygen. This product is not intended for use with Heliox gas. Use with Heliox gas may result in false No Water alarms.

### Kit Contents

### Additional contents in 950A61J:


#### RT016 Filter

1 Press down on the blue finger guard and insert the humidification chamber.

2 Pull ring to remove the yellow float retainer from the humidification chamber. Recycle or dispose of this yellow part.


3 Unwind the feed set and spike the water bag. Recycle or dispose of the yellow cover.

4 Connect the dryline to the humidification chamber.

CONTINUED OVER PAGE

5 Connect the dryline to the ventilator. If required, an inspiratory filter may be used.

6 Attach the inspiratory limb to the humidification chamber.

**7 For non-vented masks:**


7 Connect the inspiratory limb to the Disposable Exhalation Port (DEP) and connect to the non-vented mask.

8 Connect the pressure line to the Disposable Exhalation Port (DEP).


**7 For vented masks:**

7 Connect the inspiratory limb to the vented mask.

8 Connect the pressure line to the vented mask.

9 Connect pressure line to the ventilator.

10 Ensure the water level in the humidification chamber is below the maximum line.

**Symbol definitions**

	Follow instructions for use – safety.
	Consult instructions for use. <a href="http://fphcare.com/950IFU">fphcare.com/950IFU</a>
	Single use.
	Magnetic resonance (MR) unsafe.
	14 days maximum use.
	Correct water level in water chamber.
	Incorrect water level; do not use water chamber.
	Component can be disposed of in setup.
	Type BF applied part.
	This product is not made with natural rubber latex.
	Temperature limitations.
	Manufacturer.
	Date of manufacture.
	Date of expiration.

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