Neonatal / Pediatric Ventilator Circuit Kit

USER INSTRUCTIONS

REF 950N80J
950N81J

Single Use 14 Days Maximum Use

Rx only Fisher & Paykel HEALTHCARE

English en

Indications for use

The breathing set is an accessory to the F&P 950™ Respiratory Humidifier. It is intended for delivery of heated humidified respiratory gases to neonatal, infant and child patients, within the limits of its stated technical specifications.

Contraindications

There are no contraindications for heated respiratory humidification.

Side effects

There are no known side effects of heated respiratory humidification.

Technical specifications

Compatible with the F&P 950 Respiratory Humidifier. Refer to humidifier, accessories, and patient interface user instructions for additional warnings, cautions, contraindications, side effects, and system information.

Tidal volume:	<300 mL
Interface connections:	ISO 5356-1 Conical Connectors
Breathing tube length:	
Inspiratory:	1.75 m
Expiratory:	1.60 m
Gas leakage (60 cmH₂O, BTPS):	
Dual limb:	<30 mL/min
Resistance to flow (15 L/min, BTPS):	
Inspiratory:	1.07 ± 0.07 cmH ₂ O
Expiratory:	1.13 ± 0.10 cmH ₂ O
With accessories:	
Inspiratory:	1.59 ± 0.07 cmH ₂ O
Expiratory:	1.24 ± 0.08 cmH ₂ O
includes measurement uncertainty of 0.02 cmH ₂ O	
Compliance (60 cmH₂O, BTPS):	
Dual limb:	1.01 ± 0.13 mL/cmH ₂ O
With accessories:	
Dual limb:	1.02 ± 0.13 mL/cmH ₂ O
includes measurement uncertainty of 0.08 mL/cmH ₂ O	
Tube minimum internal diameter:	11.2 mm
Maximum tube surface temperature:	44 °C (<25 cm from patient port)
	80 cmH ₂ O
Maximum operating pressure:	
Operating flow range:	
Invasive ventilation:	0.5–40 L/min
Ventilator CPAP or NIV:	1.5–40 L/min
Optiflow:	1.5–40 L/min
Room temperature:	20–26 °C
Maximum incubator temperature:	37 °C
Chamber maximum volume of water:	160 mL

Warnings, cautions, and notes

⚠ WARNINGS

Failure to comply with the following warnings may impair the performance of the device or compromise safety (including potentially causing serious harm):

- ⌚ Do not use beyond the 14-day maximum duration of use.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- ⊗ Single use. Do not reuse this product. Reuse may result in infection of the airway.
- Do not use product near any ignition source. Exposure to oxygen increases the risk of fire that may cause serious injury or death.
- Ⓜ Do not use this product in or near a magnetic resonance imaging (MRI) scanner.
- Do not cover the circuit with materials such as towels, pillows, or bed linen. Failure to comply may result in skin burn.
- Avoid leaving tubing in prolonged contact with patient's skin. Failure to comply may result in a skin burn.
- Do not crush, stretch, or milk the tubing.
- Do not clean or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Set appropriate ventilator or flow source alarms to monitor therapy delivery.
- Perform a pressure and leak test on the breathing system before connecting to a patient.
- Check all connections are tight before use.
- Monitor circuit condensate to prevent occlusion or build-up of fluid. Drain as required.
- The humidifier should always be level and positioned lower than the patient.
- Do not use the chamber if it has been visibly damaged or dropped.
- Do not use the chamber if any of the sensor covers are dislodged or missing from the chamber.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not fill the chamber with water in excess of 37 °C.
- Do not use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- The use (or modification) of breathing circuit/chamber combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification, ventilator malfunction, or serious harm.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65


⚠ CAUTIONS

- Do not touch the hot surface of the chamber base. Failure to comply may result in a skin burn.
- Do not operate the chamber without water or gas flow. Failure to comply may result in low levels of humidity.
- The water source must be at least 50 cm higher than the chamber. Failure to comply may result in low levels of humidity.
- If using a Puritan Bennett™ 840 or Puritan Bennett 980 ventilator, use the yellow SST adaptor to complete the short self-test (SST). Refer to ventilator manufacturer for consequences of failure to comply with SST.


NOTES

- Use USP Sterile Water for Inhalation, or equivalent. Adding other substances may have adverse effects.
- User may be exposed to breathing tract fluids during disposal. Breathing circuits must be treated as contaminated. Dispose of according to local, state or country specific environmental laws and regulations.
- This product is designed for the delivery of air and/or oxygen. This product is not intended for use with Heliox gas. Use with Heliox gas may result in false No Water alarms.

Kit Contents

Additional contents in 950N81J:

Press down on the blue finger guard and insert the humidification chamber.


Pull ring to remove the yellow float retainer from the humidification chamber. Recycle or dispose of this yellow part.

Unwind the feed set and spike the water bag. Recycle or dispose of the yellow cover.

Connect the dryline to the humidification chamber.

Connect the inspiratory limb to the humidification chamber.


Connect the dryline to the ventilator inspiratory port and the expiratory limb to the ventilator expiratory port.

Connect the expiratory heater wire adapter to the sensor cartridge.

Connect the expiratory heater wire adapter to the expiratory limb.


Ensure the water level in the humidification chamber is below the maximum line.

Symbol definitions

	Follow instructions for use — safety.
	Consult instructions for use. fphcare.com/950IFU
	Single Use.
	Magnetic resonance (MR) unsafe.
	14 days maximum use.
	Correct water level in water chamber.
	Incorrect water level; do not use water chamber.
	Component can be disposed of in setup.
	Type BF applied part.
INSP	Ventilator inspiratory port (flow to patient).
EXP	Ventilator expiratory port (flow from patient).
	This product is not made with natural rubber latex.
SST	Indicates component required for ventilator short self-test.
	Temperature limitations.
	Manufacturer.
	Date of manufacture.
	Date of expiration.

For Optional Accessories:

For ventilators with 22 mm ports only:

Connect the dryline and expiratory limb to the ventilator using the 22F-15F adapters.

When a pressure line is required:


Connect pressure line to pressure port on the expiratory limb.

For Puritan Bennett Ventilators:

To complete the short shelf test (SST), attach the SST adapter between the inspiratory limb and the inspiratory port of the ventilator.

When an 8 mm CPAP interface is required:

Remove swivel and connect offset adapters to the inspiratory and expiratory limbs before connecting limbs to the interface.

When an Optiflow Junior interface is required:

Remove swivel and expiratory limb. Connect the Optiflow adapter to the inspiratory limb before connecting to the interface.

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