F&P 950

Adult Bi-Level/CPAP Heated Circuit Kit

USER INSTRUCTIONS





REF **950A60**
950A61

CE 0123

Single
Use

REF: ZB5002 REV H 2021-12 © 2021 Fisher & Paykel Healthcare Limited.

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraisópolis, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/ 售后服务机构 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât. F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Östım Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

F&P and F&P 950 are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

www.fphcare.com

Fisher & Paykel
HEALTHCARE

F&P 950™ Adult Bi-Level/CPAP Heated Circuit Kit

Indications for use

The breathing set is an accessory to the F&P 950 Respiratory Humidifier. It is intended for delivery of heated humidified respiratory gases to adult and paediatric patients, within the limits of its stated technical specifications.

Contraindications

There are no contraindications for heated respiratory humidification.

Do not use the exhalation port or noninvasive therapy on patients who are non-spontaneously breathing; uncooperative or unresponsive; have copious secretions; at risk of nausea or vomiting; or at high risk of aspiration of vomitus.

Side effects

There are no known side effects of heated respiratory humidification.

Technical specifications

Compatible with the F&P 950 Respiratory Humidifier. Refer to humidifier, accessory, and patient interface user instructions for additional warnings, cautions, contraindications, side effects, and system information.

Patient weight:	>10 kg
Interface connections:	ISO 5356-1 Conical Connectors

Breathing tube length:	
Inspiratory:	1.60 m

Gas leakage (60 cmH₂O, BTPS):	
Single limb:	<40 mL/min

Resistance to flow (30 L/min, BTPS):	
Single limb:	0.31 ± 0.03 cmH ₂ O

<i>With accessories:</i>	
Single limb:	1.12 ± 0.10 cmH ₂ O

(includes measurement uncertainty of 0.02 cmH₂O)

Compliance (60 cmH₂O, BTPS):	
Single limb:	1.51 ± 0.21 mL/cmH ₂ O

<i>With accessories:</i>	
Single limb:	1.47 ± 0.20 mL/cmH ₂ O

(includes measurement uncertainty of 0.02 mL/cmH₂O)

Tube minimum internal diameter:	17 mm
Maximum tube surface temperature:	44 °C
Maximum operating pressure:	80 cmH ₂ O
Operating flow range	
Mask mode:	5–120 L/min
@ 31 °C humidifier set point	
Optiflow mode:	5–70 L/min
Room temperature:	18–26 °C
Chamber maximum volume of water:	160 mL

Filter:	
Viral:	>99.99%
Organism:	ΦX174 bacteriophage
Bacterial:	>99.9997%
Organism:	bacillus subtilis
Mean particle size:	3 μm
NaCl:	98.17%

Symbol definitions

	Follow instructions for use - safety.
	Consult instructions for use. fphcare.com/950IFU
	Single use.
	Correct water level in water chamber.
	Incorrect water level; do not use water chamber.
	Component can be disposed of in setup.
	Type BF applied part.
	Ventilator inspiratory port (flow to patient).
	This product is not made with natural rubber latex.
	Temperature limitations.

	Manufacturer.
	Date of manufacture.
	Date of expiration.

Warnings, cautions, and notes

-  **WARNINGS**
 - Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm, or death.
 - Single use. Do not reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm, or death.
 - Exhalation port must be used with a non-vented interface in a single limb system. Failure to comply may lead to patient inhaling excess carbon dioxide resulting in hypercapnia.
 - The use (or modification) of breathing circuit/chamber combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification, ventilator malfunction, or serious harm.
 - Do not use product near any ignition source. Exposure to oxygen increases the risk of fire that may cause serious injury or death.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- This product is designed for the delivery of air and/or oxygen. It is not suitable for the delivery of flammable anaesthetic gas mixes or Heliox gas.
- Do not crush, stretch, or milk the tubing.
- Do not clean or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Set appropriate ventilator or flow source alarms to monitor therapy delivery.
- Perform a pressure and leak test on the breathing system before connecting to a patient.
- Check all connections are tight before use.
- Monitor circuit condensate to prevent occlusion or build-up of fluid. Drain as required.
- The humidifier should always be level and positioned lower than the patient.
- Do not use the chamber if it has been visibly damaged or dropped.
- Do not use the chamber if any of the sensor covers are dislodged or missing from the chamber.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not fill the chamber with water in excess of 37 °C.
- Do not block or seal the vent holes on the exhalation port.
-  Do not use beyond the 14-day maximum duration of use.
- When nebulized drugs are used resistance to flow should be monitored and the filter replaced, following standard hospital procedure.
- Change filter provided every 24 hours or sooner if notable deterioration occurs, following standard hospital procedure.

-  **CAUTIONS**
 - Do not touch the hot surface of the chamber base. Failure to comply may result in a skin burn.
 - Do not cover the circuit with materials such as towels, pillows, or bed linen. Failure to comply may result in a skin burn.
 - Avoid leaving tubing in prolonged contact with patient's skin. Failure to comply may result in a skin burn.
 - Do not operate the chamber without water or gas flow. Failure to comply may result in low levels of humidity.
 - The water source must be at least 50 cm higher than the chamber. Failure to comply may result in low levels of humidity.

NOTES

- Use USP Sterile Water for Irrigation, or equivalent. Adding other substances may have adverse effects.
- User may be exposed to breathing tract fluids during disposal. Breathing circuits must be treated as contaminated. Dispose of according to local, state or country specific environmental laws and regulations.
- If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and, for European Union member countries, the Competent Authority in your country.

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