F&P 950

Bubble CPAP Circuit Kit (for Hudson Interface)

USER INSTRUCTIONS

REF 950N62J

Single Use **14** 14 Days Maximum Use

Rx only

Fisher & Paykel
HEALTHCARE

English ^{en}

Indications for use

The breathing set is an accessory to the F&P 950™ Respiratory Humidifier. It is intended to provide heated humidified CPAP to spontaneously breathing premature and full term neonates and infant patients (up to a weight of 10 kg) who require breathing support, due to conditions associated with prematurity (such as Respiratory Distress Syndrome) or other conditions where CPAP is required or desired. It is intended for use within a hospital clinical environment such as the NICU (Neonatal Intensive Care Unit) and PICU (Pediatric Intensive Care Unit), within the limits of its stated technical specification.

Contraindications

There are no contraindications for heated respiratory humidification. Bubble CPAP therapy must not be used for non-spontaneously breathing patients, patients with congenital abnormalities or malformations (e.g., diaphragmatic hernia, choanal atresia, tracheoesophageal fistula), nasal trauma/severe deformity that might be exacerbated by use of nasal prongs or nasal mask.

Side effects

There are no known side effects of heated respiratory humidification.

Technical specifications

Compatible with the F&P 950 Respiratory Humidifier. Refer to humidifier, accessories, and patient interface user instructions for additional warnings, cautions, contraindications, side effects, and system information.

Patient weight: <10 kg

Patient interface: Hudson RCI™ Infant Nasal CPAP

Breathing tube length:

Inspiratory: 1.75 m

Expiratory: 1.60 m

Gas leakage (60 cmH₂O, BTPS):

Dual limb: <30 mL/min

Resistance to flow (2.5 L/min, BTPS):

Inspiratory: 0.07 ± 0.03 cmH₂O

Expiratory: 0.07 ± 0.03 cmH₂O

includes measurement uncertainty of 0.02 cmH₂O

Compliance (60 cmH₂O, BTPS):

Dual limb: 0.89 ± 0.18 mL/cmH₂O

includes measurement uncertainty of 0.08 mL/cmH₂O

Tube minimum internal diameter: 11.2 mm

Maximum tube surface temperature: 44 °C (<25 cm from patient port)

BCPAP Pressure manifold maximum pressure limit:

17 cmH₂O
@ Test Flow: 8 L/min

Operating flow range (BCPAP): 4–15 L/min

Chamber maximum volume of water: 160 mL

Room temperature: 20–26 °C

Maximum incubator temperature: 37 °C

Warnings, cautions, and notes

⚠ WARNINGS

Failure to comply with the following warnings may impair the performance of the device or compromise safety (including potentially causing serious harm):

- ⌚ Do not use this breathing circuit kit beyond the 14-day maximum duration of use.
- 👁 Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death. Ensure alternative means for providing breathing support is available such as a bag-valve and mask (BVM).
- ⊗ Single use. Do not reuse this product. Reuse may result in infection of the airway.
- 🔥 Do not use product near any ignition source. Exposure to oxygen increases the risk of fire that may cause serious injury or death.
- 🧲 Do not use this product in or near a magnetic resonance imaging (MRI) scanner.
- 👉 The BCPAP pressure relief manifold (white cap) must be used with Bubble CPAP therapy. Failure to do so may result in harmful pressures and lead to barotrauma.
- 🚫 Do not use if pressure relief manifold cap has been removed or is missing.
- 📏 Ensure the Bubble CPAP Generator is mounted below the level of the patient.
- 🛏 Do not cover the circuit with materials such as towels, pillows, or bed linen. Failure to comply may result in skin burn.
- 🧴 Avoid leaving tubing in prolonged contact with patient's skin. Failure to comply may result in a skin burn.
- 👉 Do not crush, stretch, or milk the tubing.
- 🧼 Do not clean or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- 🔊 Set appropriate ventilator or flow source alarms to monitor therapy delivery.
- 🔧 Perform a pressure and leak test on the breathing system before connecting to a patient.
- 🔧 Check all connections are tight before use.
- 🚰 Monitor circuit condensate to prevent occlusion or build-up of fluid. Drain as required.
- 📏 The humidifier should always be level and positioned lower than the patient.
- 🚫 Do not use the chamber or Bubble CPAP Generator if it has been visibly damaged or dropped.
- 🚫 Do not use the chamber if any of the sensor covers are dislodged or missing from the chamber.
- 🚫 Do not use the chamber if the water level rises above the maximum water level line.
- 🚫 Do not fill the chamber with water in excess of 37 °C.
- 🚫 Do not use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- 🚫 The use (or modification) of breathing circuit/chamber combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification, ventilator malfunction, or serious harm.
- 🇺🇸 California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65


⚠ CAUTIONS

- 🔥 Do not touch the hot surface of the chamber base. Failure to comply may result in a skin burn.
- 🚫 Do not operate the chamber without water or gas flow. Failure to comply may result in low levels of humidity.
- 📏 The water source must be at least 50 cm higher than the chamber. Failure to comply may result in low levels of humidity.
- 🔧 Check Bubble CPAP generator regularly. Failure to comply may result in compromised pressure delivered to the patient.

NOTES

- 👉 For use with nasal CPAP interfaces with nasal prongs or nasal mask. Ensure interface is adjusted for best fit and minimal leak as per the interface user instructions.
- 👉 The delivery of nebulized medication via the nasal route has not been cleared by the FDA.
- 👉 Use USP Sterile Water for Inhalation, or equivalent, in the humidification chamber. Adding other substances may have adverse effects.
- 👉 Use USP Sterile Water for Inhalation, or equivalent, to fill the Bubble CPAP Generator. A solution of 0.25% acetic acid may also be used.
- 👉 User may be exposed to breathing tract fluids during disposal. Breathing circuits must be treated as contaminated. Dispose of according to local, state or country specific environmental laws and regulations.
- 👉 This product is designed for the delivery of air and/or oxygen. This product is not intended for use with Heliox gas. Use with Heliox gas may result in false No Water alarms.


Kit Contents

1 Press down on the blue finger guard and insert the humidification chamber.

2 Pull ring to remove the yellow float retainer from the humidification chamber. Recycle or dispose of this yellow part.

3 Unwind the feed set and spike the water bag. Recycle or dispose of the yellow cover.

4 Connect the BCPAP pressure relief manifold to the humidification chamber.

CONTINUED OVER PAGE


Connect the gas supply line to the BCPAP pressure relief manifold.

Connect the inspiratory limb to the humidification chamber.

Mount the Bubble CPAP generator below the level of the patient.

Fill the Bubble CPAP generator with water using the funnel. Ensure water is above the fill line.


Connect the expiratory limb to the Bubble CPAP generator.

Connect the expiratory heaterwire adapter to the sensor cartridge.


Connect the expiratory heaterwire adapter to the expiratory limb.


Connect the flow test elbow to the inspiratory and expiratory limbs. Set flow rate to 1 L/min. Check the Bubble CPAP generator has gentle bubbling.

Set the Bubble CPAP pressure to the prescribed level. The number above the lid indicates the CPAP pressure in cmH_2O .


Adjust the flow rate to the prescribed input flow rate.

Connect the inspiratory and expiratory limbs to the Hudson interface.

Regularly check the Bubble CPAP generator has gentle bubbling during use.

Ensure the water level in the humidification chamber is below the maximum line.

Symbol definitions

	Follow instructions for use — safety.
	Consult instructions for use. fphcare.com/950IFU
	Single use.
	Magnetic resonance (MR) unsafe.
	14 days maximum use.
	Correct water level in water chamber.
	Incorrect water level; do not use water chamber.
	Fill Bubble CPAP generator with water above line.
	Maximum overflow water level in Bubble CPAP generator.
	Component can be disposed of in setup.
	Type BF applied part.
	This product is not made with natural rubber latex.
	Temperature limitations.
	Manufacturer.
	Date of manufacture.
	Date of expiration.

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