# User Manual

**LANGUAGE**

**SECTION**

English ..... A  
Español..... B  
Français..... C  
العربية..... D  
繁體中文..... E

## Before you start

- The Airvo 2 and this User Manual is intended for use by appropriately qualified healthcare professionals.
- Read this User Manual, including all warnings. Failure to do so may result in injury. Keep in a safe place for future reference.
- Training to use the Airvo 2 is necessary, where instructions in this User Manual are provided as training material.
- Before the Airvo™ 2 is used for the first time, it must be set up according to the instructions in the Airvo 2 Technical Manual. The Airvo 2 needs special precautions regarding electromagnetic compliance (EMC) therefore must be installed and put into service according to the EMC information provided in this User Manual and the Technical Manual.
- Some accessories may not be available in certain countries. Please contact your local Fisher & Paykel Healthcare representative for more information.

### Other references

- For Disinfection instructions, refer to the Disinfection Kit Manual provided with the Disinfection Kit (900PT600).
- Refer to user instructions supplied with individual accessories for correct use and additional safety information.
- Refer to the Airvo 2 Technical Manual for initial setup, maintenance, servicing and additional troubleshooting instructions.
- Download the Airvo 2 Simulator App to learn how to use the Airvo 2. You can change settings, simulate faults and test your skills. Available for download from [Apple](#), [Google Play](#), and the [Windows](#) App Store. Please note that not all applications are available for all countries.
- For assistance from your Fisher & Paykel Healthcare representative, contact us at: [www.fphcare.com/contact-us](http://www.fphcare.com/contact-us).

### Conventions used in this manual

**WARNINGS** A warning alerts the user to a potential hazard with use or misuse of the device which, if not avoided, could result in death or serious injury.

**CAUTIONS** A caution alerts the user to a potential hazard with use or misuse of the device which, if not avoided, could result in minor or moderate injury.

## Table of contents

<b>1. Overview</b>	<b>2</b>
Intended use.....	2
Contraindications.....	2
Side-effects.....	2
Warnings.....	2
Airvo 2 and accessories.....	4
<b>2. Setting up Airvo 2</b>	<b>6</b>
<b>3. Using Airvo 2</b>	<b>10</b>
Target dew-point temperature.....	11
Target flow.....	11
Oxygen.....	11
Alarms.....	13
<b>4. Reprocessing</b>	<b>17</b>
Schedule for changing accessories.....	18
Filter replacement.....	18
Servicing.....	18
<b>5. Technical information</b>	<b>19</b>
Symbol definitions.....	19
Specifications.....	20

# 1. Overview

The Airvo 2 is a humidifier with integrated flow generator that delivers high flow warmed and humidified respiratory gases to spontaneously breathing patients through a variety of patient interfaces.

A blower inside the Airvo 2 entrains flows of room air of 2 - 60 L/min, which may be blended with oxygen from low-pressure sources (such as flowmeters). The air-oxygen mixture is then warmed and humidified in the water chamber, before being transported through the heated breathing tube to a nasal, tracheostomy or mask patient interface.

## Intended use

The Airvo 2 is for the treatment of spontaneously breathing infant, child, adolescent and adult patients who would benefit from receiving high flow warmed and humidified respiratory gases. This includes patients who have had upper airways bypassed. The flow may be from 2-60 L/min depending on the patient interface. The Airvo 2 is for use in hospitals and long-term care facilities by appropriately qualified healthcare professionals.

## Contraindications

Contraindications are therapy-specific. Refer to instructions of patient interfaces and/or tube and chamber kits for therapy-specific contraindications.

## Side-effects

Side-effects are therapy-specific. Refer to instructions of patient interfaces and/or tube and chamber kits for therapy-specific side-effects.

## WARNINGS

- The unit is not intended for life support. Do not use the Airvo 2 on patients who cannot tolerate a brief interruption of therapy.
- Appropriate patient monitoring must be used at all times. Loss of power means loss of therapy.
- Nasal delivery of respiratory gases may generate flow-dependent dynamic positive airway pressure. This must be taken into account where positive airway pressure could have adverse effects on a patient.

### To avoid burns:

- Use only interfaces, water chambers and breathing tubes specified in this user manual.
- Do not use accessories beyond the maximum periods of use specified in this manual.
- Before using oxygen with the unit, read all warnings in the “Oxygen” section of this manual.
- Never operate the unit if:
  - the heated breathing tube has been damaged with holes, tears or kinks,
  - it is not working properly,
  - the case screws have ever been loosened.
- Do not block the flow of the air through the unit and breathing tube.
- Locate the unit in a position where ventilation around the unit is not restricted.
- Never block the air openings of the unit or place it on a soft surface, such as a bed or couch/sofa, where the filter area may be blocked. Keep the air openings free of lint, hair etc.
- Only use lotions and/or salves that are labeled as oxygen-compatible to avoid the risk of fire and burns.

### To avoid electric shock:

- Do not store or use the unit where it can fall or be pulled into water. If water has entered the unit enclosure, disconnect the power cord and discontinue use.
- Never operate the unit if:
  - it has been dropped or damaged,
  - it has a damaged power cord or plug,
  - it has been dropped into water.
- Avoid unnecessary removal of the power cord from the rear of the device. If removal is necessary, hold the connector during removal. Avoid pulling on the power cord.
- Return the unit to an authorized service center for examination and repair, except as outlined in this manual.

### To avoid choking, or inhalation of a foreign object:



- Ensure an air filter is fitted when operating your unit.
- Never drop or insert any object into any opening or tube.

**Miscellaneous:**

- Prior to each patient use, ensure that the auditory alarm signal is audible by conducting the alarm system functionality check described in the Alarms section.
- Do not use the Airvo 2 at an altitude or temperature outside the rated range listed in the specifications section of the manual. Using outside of these ranges can compromise the equipment performance which consequently can result in degradation of the health of the patient.
- Do not use Airvo 2 when outside the operating conditions listed in the specifications section. Humidity output will be compromised below 18 °C (64 °F) and above 28 °C (82 °F).
- To prevent disconnection during use, especially during ambulatory use, use only heated breathing tubes specified in this manual.
- Do not use the Airvo 2 system in a magnetic resonance imaging (MRI) environment.
- The unit is not suitable for use in the presence of a flammable, anesthetic mixture with air or oxygen or nitrous oxide.
- Do not use the Airvo 2, or accessories, during defibrillation.
- Do not use the Airvo 2, or accessories, during electrocautery.
- Explosion hazard: Do not use the Airvo 2 in the presence of flammable anesthetics or other flammable substances in combination with air, oxygen-enriched environments, or nitrous oxide.
- Do not use the Airvo 2 in a hyperbaric chamber.
- Avoid using the Airvo 2, or accessories, adjacent to, or stacked with, other equipment, which could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- The Airvo 2 is not a sealed system. Follow hospital infection control guidelines to reduce risk of cross-contamination
- Use of accessories or power cables not specified by Fisher & Paykel Healthcare could result in increased electromagnetic emissions, decreased electromagnetic immunity and/or improper operation.
- Do not use any patient consumables, accessories or spare parts that are not listed in this user manual, or the Airvo 2 Technical Manual. Incompatible consumables, parts or accessories could affect the quality of therapy, injure the patient, decrease electromagnetic immunity or increase electromagnetic emissions.
- As the low-pressure oxygen (LPO) inlet port uses an alternative small-bore connector design different from those specified in the ISO-80369 series, there is a possibility that a misconnection can occur with a medical device using an alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need to be taken by the user to mitigate these reasonably foreseeable risks.
- The therapy delivered to the patient can be impacted by the use of a pneumatic/jet nebuliser. Refer to compatible accessories and drug manufacturer instructions for use for correct usage.
- Do not connect the Airvo 2 to the battery of a battery-powered wheelchair, which may compromise device performance and therapy delivery.
- Use only patient interfaces, heated breathing tubes, water chambers and filters specified in this manual to prevent disconnection during use, especially when moving the Airvo 2.
- Do not lubricate fittings, connections, tubing, or other accessories of the equipment to avoid the risk of fire and burns.
- Do not use sealed patient interfaces with Optiflow high flow therapy, to avoid the risk of suffocation or barotrauma.
- Ensure a sufficient intended leakage between the breathing system and the patient to allow the patient to exhale.
- It is the responsibility of the responsible organization to ensure that the oxygen source is compatible with the rated range of pressure, flowrate and oxygen concentration as marked on the equipment and indicated in the instructions for use as this can affect the performance of the equipment or pipeline system that can consequently result in serious deterioration of health.
- Open flames during Optiflow high flow therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m of the equipment or any oxygen-carrying accessories.
- Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions, if the device is turned on, but not in use; the oxygen will make the materials more flammable. Turn the device off when not in use to prevent oxygen enrichment.
- Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking or open flames within the same room as the equipment or any oxygen-carrying accessories. If the patient intends to smoke, always turn the equipment off, remove the cannula and leave the room where the equipment is located. If unable to leave the room, wait 10 minutes after the equipment has been turned off.
- There is a risk of fire associated with oxygen enrichment during oxygen therapy. Do not use the equipment or accessories near sparks or open flames.

**Notice to user:** If a serious incident has occurred while using this device, please contact your local Fisher & Paykel Healthcare representative and Competent Authority.

# Airvo 2 and accessories

**Airvo 2**  
PT101xx

Oxygen inlet port



Heated breathing tube connection port

Measurement point of displayed dew-point temperature

Chamber ports

Heater plate

Finger guard

**Water chamber**

Auto-fill water chamber (MR290) (with adapter fitted)

## Parts and accessories

Description	Part number	Size	Pack size
Optiflow Junior nasal interface	OPT316	Infant	20
	OPT318	Pediatric	20
Optiflow Junior 2 nasal interface	OJR416	L	20
	OJR418	XL	20
Optiflow Junior 2+ nasal interface	OJR520	XXL	10
Optiflow+ nasal interface	OPT942	Small	20
	OPT944	Medium	20
	OPT946	Large	20
Optiflow+ Duet nasal interface	OPT962	Small	20
	OPT964	Medium	20
	OPT966	Large	20
Optiflow 3S nasal interface	OPT1042	Small	20
	OPT1044	Medium	20
	OPT1046	Large	20
Optiflow+ tracheostomy interface	OPT970	15 mm	20
Optiflow+ mask interface adapter	OPT980	22 mm	20
AirSpiral tube and chamber kit	900PT561		10
AirSpiral tube and chamber kit (with nebulizer adapter)	900PT562/900PT563*		10

\*For USA only


**Note:** The Airvo 2 is compatible with delivering specified nebulized medications. See the user instructions for the 900PT562/900PT563 for more information (including warnings and cautions).

Accessories		Cleaning and disinfection	
900PT405	Pole mounting tray	900PT600	Disinfection kit
900PT411	UPS mounting kit	900PT601	Disinfection filter (2-pack)
900PT420	Mobile pole stand (extendable)	900PT602	Cleaning sponge-stick (20-pack)
900PT421	Mobile pole stand	900PT603	Clean storage cover (20-pack)
900PT426	Plastic basket		
900PT427	Oxygen bottle holder		
900PT427L	Oxygen bottle holder (large)		
900PT428	Pole clamp		
900PT445	Adjustable handle for hospital stand Airvo series		
OPT012/WJR112	Wigglepads for Optiflow Junior (20-pack)		
		Spare parts	
		900PT422	Oxygen inlet extension kit
		900PT912	Filter holder
		900PT913	Air filter (2-pack)

**Note:** Some accessories may not be available in all countries. Please contact your local Fisher & Paykel Healthcare representative for more information.

The responsible organisation is required to ensure that all patient consumables are compliant to ISO 18190. Fisher & Paykel Healthcare has assessed the following items as compliant.


## 2. Setting up Airvo 2

### 1. Before you begin


Check that the Airvo 2 is attached securely to a pole mounting tray (900PT405) and mobile pole stand (900PT421) and is below the patient's head height.

Position the Airvo 2 so that the power cord connection to the wall power supply is easily accessible and can be disconnected if necessary.


### 2. Assemble the water chamber


Open the tube and chamber kit and remove the MR290 auto-fill water chamber and chamber adapter. Remove the blue port caps from the chamber by pulling the tear tab upwards then remove the bracket holding the water supply tube.

Fit the supplied adapter over the two vertical ports on the chamber and push on fully then clip the water supply tube into position.

Check that the chamber adapter sits on the two vertical chamber ports and there are no gaps.

Fit the water chamber to the unit by pressing down the finger guard and sliding the chamber on, carefully aligning with the blue chamber port ends.

Push the chamber on firmly until the finger guard clicks into place.

## WARNINGS

### To avoid burns:


- Do not start the unit without the water chamber in place.
- Do not touch the heater plate, water chamber or chamber base during use.
- The water in the chamber becomes hot during use. Exercise caution when removing and emptying the chamber.

**To avoid electric shock:**

- When handling the unit with the water chamber in place, avoid tilting the machine to prevent any chance of water entering the unit enclosure.
- Empty all the water from the water chamber before transporting the unit.

**CAUTIONS****To ensure optimal therapy (MR290 only):**

- Do not use the auto-fill MR290 chamber if it has been dropped or been allowed to run dry. This could lead to the chamber over filling.

**3. Connect water bag**


Attach the sterile water bag to the hanging bracket 20 cm (8") above the unit, and push the bag spike into the fitting at the bottom of the bag. Open the vent cap on the side of the bag spike. The chamber will now automatically fill to the required level and maintain that level until the water bag is empty.

To ensure continual humidification, always ensure that the water chamber and/or water bag are not allowed to run out of water.

**CAUTIONS**

Adding substances other than water can adversely affect the humidifier and delivered therapy.

Check that water flows into the chamber and is maintained below the maximum water level line. If the water level rises above the maximum water level line, replace the chamber immediately.


**MR290: Flow setting vs. usage time**

(1 litre sterile water bag, at 37 °C target temperature)

L/min	2	5	10	15	20	25	30	35	40	45	50	55	60
Hours	129	60	31	21	16	12	10.5	9	8	7	6.5	5.5	5

**CAUTIONS****To ensure optimal therapy (MR290 only):**

- Do not use the MR290 chamber if the water level rises above the maximum water level line as this may lead to water entering the patient's airway.

**4. Install heated breathing tube**

Connect the breathing tube to the Airvo 2. Make sure the blue connector is fully locked into place.

**WARNINGS****To avoid burns:**

- Do not modify the breathing tube or interface in any way.
- Do not allow the breathing tube to remain in direct contact with skin for prolonged periods of time. The healthcare professional shall assess the conditions for safe contact, such as duration and skin condition.
- Do not add heat above ambient levels to any part of the breathing tube or interface, e.g. by covering with a blanket or by heating with infrared radiation, an overhead heater, or an incubator. This can affect the quality of the therapy or injure the patient.
- Do not use an insulating sleeve or any similar accessories that are not recommended by Fisher & Paykel Healthcare.


**CAUTIONS**

- Position the heated breathing tube away from any electrical monitoring leads (EEG, ECG/EKG, EMG, etc.) to minimize any possible interference with the monitored signal.

## 5. Visual check of inlet and outlet ports

Check that inlet and outlet ports are not damaged and shows no signs of deformation.

Ensure the water chamber adapters are fully inserted into the ports and that there are no gaps.

## 6. Select patient interface

The Airvo 2 can be used with a variety of patient interfaces. Read the separate user instructions for the patient interface that will be used, including all warnings.

### Nasal interface


**Optiflow Junior 2**  
**Optiflow Junior 2+**


**Optiflow, Optiflow+**  
**Optiflow+ Duet, Optiflow 3S**

### Tracheostomy interface

**Optiflow+ tracheostomy interface**

### Mask interface

**Optiflow+ mask interface adapter**

All patient interfaces are Type BF applied parts.

The following table shows the target flow settings able to be used with these interfaces.

	L/min									
Patient interface	2	5	10	15	20	25	...	50	55	60
<b>Junior Mode</b>										
OPT316/OJR416	2		20							
OPT318/OJR418	2		25							
<b>Default Mode</b>										
OJR520			10			50				
OPT942/OPT962/OPT1042	<b>S</b>			10			50			
OPT944/OPT964/OPT1044	<b>M</b>			10			60			
OPT946/OPT966/OPT1046	<b>L</b>			10			60			
OPT970			10			60				
OPT980			10			60				

\*The rated flow rate depends on the device, breathing circuit and patient interface. These flow rates are an indication of the achievable range. However, the rated flow rates in the breathing circuit and patient interface user instructions need to also be considered. The smallest range across these instructions is applicable. For any inconsistencies, please refer to your FPH representative to confirm the rated flow rate of your system configuration.


Low temperature ambient conditions may prevent the unit from reaching a 37 °C target temperature setting at high target flow settings. In these cases, consider decreasing the target flow setting. At altitude, the maximum flow rates achievable may be lower than those in the above table, by approximately 5 L/min per 1000 m (3000 ft).

## WARNINGS

### To avoid burns:

- Do not modify the breathing tube or interface in any way.
- Do not use any patient interfaces not listed here.

## 3. Using Airvo 2


### 1. Switch on unit

Plug the unit's power cord into the mains/utility power socket. The connector at the other end of the power cord should be well secured to the rear of the unit

#### WARNINGS

- Do not cover or position next to materials that can block the gas intake, thereby interfering the patient therapy.

#### To avoid electric shock:

- Ensure that the unit is dry before plugging into the mains/utility power socket.

Switch on the unit by pressing the On/Off button for 5 seconds.

### 2. Check disinfection status


The unit will show you whether it is safe for use on a new patient.

This Airvo 2 is safe for use on a new patient.

This Airvo 2 has not been cleaned and disinfected since last use.  
This Airvo 2 is NOT safe for use on a new patient.

### 3. Warmup

The unit will begin to warm up. You will see numbers showing the current output dew-point temperature, flow and oxygen values. These numbers will pulse until they approach their target settings.

This screen is called the "Summary screen".

### 4. Junior Mode

If the patient will be using an Optiflow Junior nasal interface (OPT316/OJR416/ OPT318/OJR418), you must activate Junior Mode. Do not use Junior Mode for other patient interfaces.


Junior mode limits the target settings to: 34 °C and 2 – 25 L/min, in increments of 1 L/min.

#### To activate Junior Mode:


Hold the Mode button for 5 seconds.


#### New target settings

The target settings for dew-point temperature and flow will be changed automatically. The colorful icons in the corners of the screen indicate that this unit is in Junior Mode.

To deactivate Junior Mode, follow the same procedure: hold the Mode button for 5 seconds.

## 5. Configure target settings

Press the Mode button to view target settings.

 These settings are locked by default.

### Target dew-point temperature

You can set the Airvo 2 to three target dew-point temperature settings:

- 37 °C (98.6 °F) [always use for patients with bypassed airways (tracheostomies)]
- 34 °C (93 °F) [if compliance at 37 °C is a problem]
- 31 °C (88 °F) [for face masks only]

You may not have access to all settings, if:


- the unit is in Junior Mode (limited to 34 °C),
- the unit was initially set up with tighter limits.

The Airvo 2 will return to its default setting (37 °C) after every disinfection cycle.

#### ➔ To change the target dew-point temperature setting:


Hold the Up and Down buttons for 3 seconds to “unlock” the setting.

The lock will disappear and be replaced by an arrow showing the minimum and maximum accessible settings. Press the Up and Down buttons to choose the new setting.

When you have finished, press the Mode button to ‘lock’ the setting again.

The lock will reappear.

Press the Mode button to move on to the next screen.

### Target flow

You can set the Airvo 2 to flows between 10 L/min and 60 L/min, in increments of 1 L/min (10 – 25 L/min) and 5 L/min (25 – 60 L/min).


You may not have access to all settings, if:

- the unit is in Junior Mode (limited to 2 – 25 L/min, in increments of 1 L/min),
- the unit was initially set up with tighter limits.

The Airvo 2 will remember its target flow setting when you switch it off.

#### ➔ To change the target flow setting:

Follow the same sequence of steps as above in “To change the target dew-point temperature setting”.

Press the Mode button to move on to the next screen.

## Oxygen

You can connect up to 60 L/min of supplementary oxygen from a regulated supply to the Airvo 2. The Airvo 2 contains an oxygen analyzer to help you determine the oxygen fraction you are delivering to the patient. Your unit may have been initially set up with tighter limits.

Use continuous oxygen monitoring on patients who would desaturate significantly in the event of disruption to their oxygen supply.

### WARNINGS


Before using the Airvo 2 with oxygen, read all of the following warnings:

- The use of oxygen requires that special care be taken to reduce the risk of fire. Accordingly, for safety it is necessary that all sources of ignition, e.g. electrocautery or electrosurgery be kept away from the unit and preferably out of the room in which it is being used. Oxygen should not be used while smoking or in the presence of an open flame. The unit should be located in a position where ventilation around the unit is not restricted.
- A spontaneous and violent ignition may occur if oil, grease or greasy substances come in contact with oxygen under pressure. These substances must be kept away from all oxygen equipment.
- Ensure that the Airvo 2 is switched on before connecting oxygen.
- Oxygen must only be added through the special oxygen inlet port on the back of the unit. To ensure that oxygen enters the unit correctly, the oxygen inlet port must be fitted properly to the filter holder and the filter holder must be fitted properly to the unit. The power cord connector should also be well secured.
- Do not connect supplementary oxygen to the Airvo 2 at flow rates higher than the Airvo 2 target flow rate, as excess oxygen will be vented into the surroundings, or 60 L/min.
- The oxygen concentration delivered to the patient can be affected by changes to the flow setting, oxygen setting, patient interface or if the airpath is obstructed.
- When finished, turn off the oxygen source. Remove the output of the oxygen source from the oxygen inlet port on the back of the unit. The oxygen flow must be turned off when the unit is not operating so that oxygen does not build up inside the device.
- The oxygen analyzer within the Airvo 2 uses ultrasonic measurement technology. It does not require in-field calibration or barometric pressure compensation. It is designed for use with pure or concentrator oxygen - connecting any other gases or mixtures of gases will cause it to function incorrectly.
- The 'Oxygen concentration' setting must be configured appropriately to maintain oxygen measurement accuracy when using pure or concentrator oxygen. Refer to the Airvo 2 Technical Manual for further information.

### Connect oxygen

Connect the output from the oxygen source to the oxygen inlet port on the side of the unit. Make sure you push the oxygen tube firmly onto this connection port.

### Adjust oxygen

Adjust the level of oxygen from the oxygen source until the desired oxygen fraction is displayed onscreen. It may take the reading several minutes to settle. You can set the oxygen fraction between the maximum and minimum values displayed above and below the arrow.


Real-time O<sub>2</sub> measurement is displayed when O<sub>2</sub> >25% and O<sub>2</sub> <95%. However, note that oxygen fractions below 25% and above 95% will be displayed as 21% and 100% respectively. If the oxygen fraction exceeds 95%, the oxygen reading will pulse red and the device will beep.

### WARNINGS

- Note that if the patient's peak inspiratory demand exceeds the flow delivered by the unit, the fraction of oxygen inspired by the patient will be lower than the value shown onscreen due to the additional entrainment of ambient air.
- Check that suitable blood saturation levels are achieved at the prescribed flow.

Press the Mode button to return to the Summary screen.


### 6. Connect your patient

Wait until the “Ready for use” symbol is displayed on the Summary screen.

“Ready for use” symbol

Connect the patient interface to the heated breathing tube.

Monitor the flow and oxygen values displayed on the Summary screen. Adjust the level of oxygen from the oxygen source as necessary.

When the patient first uses the unit, the air will feel warm. This is normal. The patient should continue to breathe normally through the nose and/or mouth, or tracheostomy.

### 7. During use

If the “Ready for use” symbol has been displayed for 2 minutes and no button has been pushed in this time, a screensaver will be launched.

#### Condensate management

The unit must be placed below head height and flat. This allows condensate to drain towards the water chamber and away from the patient.

If excess condensate accumulates in the heated breathing tube, disconnect the patient interface from the heated breathing tube. Drain the condensate by lifting the patient end of the tube, allowing the condensate to run into the water chamber.


At higher target flow rates, it may be necessary to first reduce the target flow rate to 30 L/min or below to ensure the condensate drains into the water chamber.

Minimize local sources of cooling acting on the heated breathing tube, such as a fan to cool the patient or an air-conditioning unit/vent.

If condensate persists, consider turning the target temperature down.

Note: a lower target temperature will decrease the humidity output of the unit, decreasing the level of condensation.

Note: The temperature and humidity level delivered to the patient will also be reduced.


### 8. After use

Switch off the unit by pressing the On/Off button for 5 seconds.

## Alarms

The Airvo 2 has visual and auditory alarms to warn you about interruptions to your patient's treatment. These alarms are generated by an intelligent alarm system, which processes information from the sensors and target settings of the unit and compares this information to pre-programmed limits.

### Alarm signals

Alarm signal	Symbols	Meaning
<b>Visual alarm signal</b>  <p>All medium priority alarm displays will feature a “flashing” effect whereby the pixels on the top part of the screen (and the screensaver) will alternate between yellow and black.</p>	 	<p>Alarm condition</p> <p>Audio paused</p>
<b>Auditory alarm signal</b> 3 beeps in 3 seconds. Repeated every 5 seconds.		<p>Press this button to mute the auditory alarm for 115 seconds.            The auditory alarm can be reactivated by pressing this button again.</p>

### Alarm conditions

All of the alarms listed below have been assessed as “Medium Priority”. These priorities have been allocated for an operator's position within 1 meter of the device. The unit also uses an internal priority-ranking system. If multiple alarm conditions occur simultaneously, the unit will display the highest-priority alarm.

The following table lists all of the alarm conditions from highest to lowest priority, their causes, possible solutions and delays. Alarm conditions that affect oxygen delivery require an immediate response to assess the patient's saturation levels. Alarm conditions that affect humidity delivery require a prompt response to assess potential drying of mucus and associated blockages.

The following alarm delays assume operation in ‘Ready for use’ mode.

Message	Meaning	Affects delivery of:	Delays
<b>Fault (E###)</b>	<p><b>The unit has detected an internal fault.</b></p> <p>Switch the unit off and then restart. If the problem persists, note the fault code and contact your Fisher &amp; Paykel Healthcare representative.</p>	Oxygen, humidity	<5 seconds
<b>Check tube</b>	<p><b>The unit cannot detect the heated breathing tube.</b></p> <p>Check that the heated breathing tube is not damaged and that it is plugged in correctly. If the problem persists, then change the heated breathing tube.</p>	Oxygen, humidity	<5 seconds
<b>Check for leaks</b>	<p><b>The unit has detected a leak in the system.</b></p> <p>The check for leaks alarm is designed to alert the user if the water chamber has been completely disconnected from the Airvo's chamber ports and heater plate, or, if an Optiflow Junior interface is disconnected from the heated breathing tube while operating in Junior mode. It is not designed to reliably detect other possible sources of leak.</p> <p>If the Check for leaks alarm is activated, check that the:</p> <ul style="list-style-type: none"> <li>• Water chamber is correctly inserted onto the Airvo 2 chamber ports.</li> <li>• If in Junior mode, check the interface is correctly connected to the heated breathing tube.</li> <li>• Check for any other possible sources of leak.</li> </ul>	Oxygen, humidity	<120 seconds

Message	Meaning	Affects delivery of:	Delays
Check for blockages	<p><b>The unit has detected a blockage in the system.</b></p> <p>Check the heated breathing tube or patient interface for blockage. Check the air filter and filter holder for blockage. Check whether the unit should be in Junior Mode. If the patient will be using an Optiflow Junior nasal interface (OPT316/OJR416/OPT318/ OJR418), you must activate Junior Mode.</p>	Oxygen, humidity	<15* seconds
O <sub>2</sub> too low	<p><b>The measured oxygen level has fallen below the allowed limit.</b></p> <p>Check that the oxygen source is still operational and is correctly connected. Adjust the level of oxygen from the oxygen source as necessary.</p>	Oxygen	<20 seconds
O <sub>2</sub> too high	<p><b>The measured oxygen level has exceeded the allowed limit.</b></p> <p>Check that the Airvo flow rate has been set correctly. Adjust the level of oxygen from the oxygen source as necessary.</p>	Oxygen	<20 seconds
Cannot reach target flow	<p><b>The unit cannot reach the target flow setting.</b></p> <p>Check the heated breathing tube or patient interface for blockage. Check whether the target flow setting is too high for the patient interface being used (refer to “Setting up Airvo 2” – “Select Patient Interface”). You will be prompted for acknowledgement.</p> <p><b>WARNINGS</b></p> <ul style="list-style-type: none"> <li>The oxygen concentration delivered to the patient can be affected by changes to the flow setting. Adjust the level of oxygen from the oxygen source as necessary.</li> </ul>	Oxygen	<120 seconds
Check water	<p><b>The chamber has run out of water.</b></p> <p>When a chamber runs dry, the chamber float may be damaged. Replace the chamber and water bag. To ensure continual humidification, always ensure that the water chamber and/or water bag do not run out of water.</p>	Humidity	<30 minutes
Cannot reach target temperature	<p><b>The unit cannot reach the target temperature setting.</b></p> <p>You will be prompted for acknowledgement. The most likely cause for this is that the unit is operating at a high flow rate in low ambient conditions. Consider decreasing the target flow setting.</p> <p><b>WARNINGS</b></p> <ul style="list-style-type: none"> <li>The oxygen concentration delivered to the patient can be affected by changes to the flow setting. Adjust the level of oxygen from the oxygen source as necessary.</li> </ul>	Humidity	30 +/- 3 minutes
Check operating conditions	<p><b>The unit has detected that it is operating in unsuitable ambient conditions.</b></p> <p>This alarm may be caused by a sudden change in ambient conditions. Leave the unit running for 30 minutes. Switch the unit off and then restart.</p>	Humidity	60 +/- 6 seconds
Power out	<p><b>The unit has been disconnected from the mains/utility power socket.</b></p> <p>No visual alarm. The auditory alarm will sound for at least 120 seconds. If power is reconnected in this time, the unit will automatically restart.</p> <p><b>WARNINGS</b></p> <ul style="list-style-type: none"> <li>Appropriate patient monitoring must be used at all times. Loss of power means loss of therapy.</li> </ul>	Oxygen, humidity	<5 seconds

\*Delays measured from unit detection.

## Alarm limits

Most alarm limits are pre-programmed. The exceptions are listed below. These alarm limits may be changed to other values by authorized personnel. Changes will be preserved during or after any power loss.

Alarm condition	Factory-set alarm limit	Possible preset values
O <sub>2</sub> too low	21% O <sub>2</sub>	21 or 25% O <sub>2</sub>
O <sub>2</sub> too high	95% O <sub>2</sub>	30 - 100% O <sub>2</sub> , in 5% increments

## WARNINGS

- A hazard can exist if different alarm presets are used on different units within any single area, e.g., an intensive care unit.
- Alarm limits set to extreme values can render the alarm system useless.

## Checking alarm system functionality

The functionality of the alarm system can be checked at any time when the unit is turned on.

Remove the heated breathing tube. You should see the “Check tube” visual alarm signal and hear the auditory alarm signal. If either alarm signal is absent, do not use the unit and refer to the Airvo 2 Technical Manual for a guide on troubleshooting. If problems persist, contact your Fisher & Paykel Healthcare representative.

## Auditory information signals

In addition to auditory alarm signals, auditory information signals are provided. These are described below.

Melody	Meaning
Ascending sequence of 5 tones	The “Ready for use” symbol has appeared
Ascending sequence of 3 tones	Activation/deactivation of Junior mode
Single tone every 5 seconds	Measured oxygen level $\geq 33\%$ at turn-off
Single tone every 30 seconds	Measured oxygen level $>95\%$

## 4. Reprocessing

Standard aseptic techniques should be followed to minimize contamination when handling the Airvo 2 and accessories. The patient interface, heated breathing tube, water chamber and outlet elbow may become contaminated during use.

As soon as possible after using the Airvo 2:

1. Remove the single-use accessories from the Airvo 2 and dispose of them in accordance with local laws, regulations and hospital protocols for disposing of contaminated products.
  - Squeeze the sides of the breathing tube connector and lift to remove it from the Airvo 2.
  - Grip the port adapter and pull the water chamber away from the Airvo 2 to remove it.
2. Reprocess the Airvo 2 exterior by following the instructions in section 4.1.
3. Clean and high-level disinfect the Outlet Elbow by following the instructions in the Disinfection Kit Manual (900PT600).
4. Replace accessories within the maximum use period shown in the “Schedule for changing accessories” table.

### 4.1 Device exterior reprocessing

#### 4.1.1. Airvo 2 device exterior reprocessing

##### Equipment

- Mild detergent and clean water
- Clean lint-free cloths
- Protective gloves

##### Instructions:

1. Mix a solution of warm clean water and mild detergent (refer to the detergent manufacturer’s instructions for use).
2. Dampen a clean cloth with the cleaning solution.
3. Thoroughly wipe all outside surfaces of the device (including the Outlet Elbow) for at least one minute to remove any visible soil. Use the corner/edge of the cloth to clean all crevices of the device.
4. Dampen a clean cloth with clean water.
5. Thoroughly wipe all outside surfaces of the device with the damp cloth to rinse and remove any detergent residue.
6. Thoroughly wipe all outside surfaces of the device with a dry cloth until it is visibly dry.
7. Allow to air dry until completely dry.

#### 4.1.2 Device exterior disinfection

**Perform disinfection only after all cleaning steps are complete**

##### Equipment

- Disinfectant wipes
- Clean lint-free cloths
- Clean water
- Protective gloves

##### Instructions

1. Use pre-soaked disinfectant wipes, thoroughly wipe all outside surfaces of the device (including the Outlet Elbow).
2. Ensure that these surfaces remain visibly wet as directed by the manufacturer of the disinfectant wipes. Use additional wipes as needed to ensure that these surfaces remain wet for the required length of time.
3. Dampen a clean cloth with clean water.
4. Thoroughly wipe all outside surfaces of the device with the damp cloth to remove any disinfectant residue.
5. Thoroughly wipe all outside surfaces of the device with a dry cloth until it is visibly dry.
6. Allow to air dry until completely dry.

### 4.2 Outlet elbow reprocessing

The Outlet Elbow requires cleaning and high-level disinfection. Clean and high-level disinfect the Outlet Elbow by following the instructions in the Disinfection Kit Manual (900PT600).

**WARNING**

Other cleaning agents may be used if they are: non-abrasive, non-toxic and non-corrosive. Do not use any cleaning agents that are not compatible with polycarbonate plastic. Cleaning agents that are not suitable for use with the Airvo 2 include: ammonia, ammonium hydroxide, caustic soda, iodine, methanol, methylated spirits, turpentine and alkaline bleaches such as sodium hypochlorite. Use of any of these products will damage the Airvo 2.

**Schedule for changing accessories**


The accessories for the unit must be changed frequently to avoid the risk of infection. Parts should be replaced immediately if they are damaged or discolored; otherwise they must be replaced within the periods shown in the following table:

Accessory	Part Number	Maximum use
Optiflow Junior nasal interfaces	OPT316/OPT318	1 week (single-patient use) *replace wigglepads as required
Optiflow Junior 2 nasal interfaces	OJR416/OJR418	2 weeks (single-patient use) *replace wigglepads as required
Optiflow Junior 2+ nasal interface	OJR520	
Optiflow+ tracheostomy interface	OPT970	2 weeks (single-patient use)
Optiflow +mask Interface adapter	OPT980	
Optiflow+ nasal Interfaces	OPT942/OPT944/OPT946	
Optiflow+ Duet nasal interfaces	OPT962/OPT964/OPT966	
Optiflow 3S nasal interfaces	OPT1042/OPT1044/OPT1046	3 months or 1000 hours
All tube and chamber kits	900PT561/900PT562/900PT563	
Air filter	900PT913	


Some products may not be available in your country. Please contact your local Fisher & Paykel Healthcare representative.

**Filter replacement**

After the Airvo 2 has been switched on for 1000 hours, a prompt will appear at the start of the next disinfection cycle indicating that an air filter change is due. Follow the steps below if a filter change is due:

1. Take the filter holder from the back of the unit and remove the filter.
2. Replace the old filter with a new filter (900PT913).
3. Reattach the filter holder to the unit (clip the bottom of the filter holder in first, then rotate it upwards until the top clips into place).

4. Press the Mode button to move on to the "Replace now" screen.
5. Press the Up button to select "Now".
6. Press the Mode button to confirm. The hours counter will be reset to zero.

If you choose the "Later" option, the prompt will continue to appear at the start of subsequent disinfection cycles.

**Servicing**

Airvo 2 humidifiers do **NOT** require routine servicing and calibration. Refer to the Airvo 2 Technical Manual for a list of external spare parts.

## 5. Technical information

### Symbol definitions

	For safety reasons, refer to the instructions for use		Class II equipment (double insulated)
	Caution		Catalogue number
<b>LPO</b>	Low-pressure oxygen		Serial number
	Warning, hot surface		Batch code
	Manufacturer		Humidity range
 YYYY-MM-DD	Date of manufacture		Temperature range
	Date of shelf life expiry	<b>IP22</b>	Protected against ingress of small objects and water drops
	Type BF applied part (body floating)		Magnetic Resonance (MR) unsafe
<b>Rx only</b>	(USA) Federal Law restricts this device to sale by, or on the order of a physician.*		EU Authorised representative*
	Alarm symbol	<b>CE 0123</b>	European conformity - TÜV SÜD*
	Alarm pause		Power on/off (standby)
	Regulatory Compliance Mark (RCM)*		Medical Device*
	Importer		Distributor
	Contains hazardous substances Pb Lead, Co Cobalt*		

\*symbol seen on selected models

## Specifications

### General

<b>Dimensions</b>	295 mm x 170 mm x 175 mm (11.6" x 6.7" x 6.9")
<b>Weight</b>	2.2 kg (4.8 lb) unit only, 3.4 kg (7.5 lb) packaged in bag incl. accessories
<b>Supply frequency</b>	50 – 60 Hz
<b>Supply voltage/current</b>	100 – 115 V 2.2 A (2.4 A max <sup>1</sup> ) 220 – 240 V 1.8 A (2.0 A max <sup>1</sup> )
<b>Auditory alarm</b>	
<b>Sound pressure level</b>	>45 dBA @ 1 m
<b>Audio pause duration</b>	115 seconds
<b>Sound power level</b>	<60 dBA
<b>Sound pressure level</b>	<50 dBA @ 1m
<b>Ingress protection</b>	IP22
<b>Expected service life</b>	5 years <sup>7</sup>
<b>Serial port</b>	The serial port is used for downloading product data, using F&P Infosmart™ software.
<b>Warmup time</b>	10 minutes to 31 °C (88 °F), 30 minutes to 37 °C (98.6 °F) using a MR290 chamber with flow rate of 35 L/min <sup>5</sup> and starting temperature 23 ± 2 °C (73 ± 3 °F)
<b>Target temperature settings</b>	37, 34, 31 °C
<b>Humidity performance<sup>2,5,6</sup></b>	>33 mg/L at 37 °C target humidity, 10-60 L/min target flow <sup>4</sup> >16 mg/L for all other settings
<b>Maximum temperature of delivered gas<sup>2</sup></b>	43 °C (109 °F)
<b>Maximum surface temperature of applied parts<sup>2</sup></b>	44 °C (111 °F)
<b>Static temperature stability<sup>6</sup></b>	± 2 °C
<b>Flow range (default)<sup>5</sup></b>	10 – 60 L/min
<b>Flow range (Junior Mode)<sup>5</sup></b>	2 – 25 L/min
<b>Maximum oxygen input<sup>3</sup></b>	60 L/min
<b>Maximum limited pressure</b>	60 cmH <sub>2</sub> O
<b>Maximum operating pressure</b>	30 cmH <sub>2</sub> O
<b>Oxygen response time</b>	<60 s
<b>Low-pressure oxygen (LPO) inlet port</b>	
<b>Line pressure</b>	0 – 50 kPa
<b>Flow rate<sup>3</sup></b>	0 – 60 L/min
<b>% Concentration</b>	93% , >99%

### Operating Conditions

<b>Ambient temperature</b>	18 – 28 °C (64 - 82 °F)
<b>Humidity</b>	10 – 95% RH
<b>Altitude</b>	0 – 2000 m (6000 ft)
<b>Mode of operation</b>	Continuous operation

### Storage and transport conditions

<b>Ambient temperature</b>	- 10 – 60 °C (14 - 140 °F)
<b>Humidity</b>	10 – 95% RH, non-condensing

The unit may require up to 24 hours to warm up or cool down from the minimum or maximum storage temperature before it is ready for use.

1. Inrush current may reach 50 A.
2. In accordance with ISO 80601-2-74. Tested to an accuracy of ± 1 °C or ± 1 mg/L, as appropriate.
3. Flow rate is expressed in STPD (standard temperature and pressure, dry) as per ISO 80601-2-74.
4. Applies to use with bypassed airway patient interfaces only.
5. Flow rates are measured in Body Temperature/Pressure, Saturated (BTPS).
6. To meet safety and performance standards, the system has been tested using the worst case breathing system setup. This setup includes 900PT561 with OPT970 interface for bypassed airway, and nasal interfaces for all other cases.
7. Fisher & Paykel Healthcare recommends that this device be replaced at the end of its expected service life."

## Range and accuracy of delivered parameters

Parameter	Range	Accuracy
Flow** <sup>§</sup>	2 – 60 L/min	± (1 + 8% of the target flow rate) L/min
Humidity** <sup>§</sup>	37 °C 31 – 37 °C	>33 mg/L at 10 – 60 L/min target flow rate <sup>#</sup> >16 mg/L

\* Test equipment and methods are selected and controlled to ensure the uncertainty coverage is no more than 30% of the disclosed tolerance.

\* Values are expressed in body temperature, pressure, saturated (BTPS), in accordance with ISO 80601-2-74.

<sup>#</sup> Applies to use with bypassed airway patient interfaces only.

<sup>§</sup> To meet safety and performance standards, the system has been tested using the worst case breathing system setup. This setup includes 900PT561 with OPT970 interface for bypassed airway, and nasal interfaces for all other cases.

## Range and accuracy of measured parameters

Parameter	Range	Accuracy
Flow** <sup>§</sup>	2 – 60 L/min	± (1 + 8% of reading) L/min
Humidity	31 – 37 °C	Not specified
Oxygen concentration** <sup>§</sup>	25 – 95%	Lower of: ± 4%, or ± (2.5% + 2.5% of reading) - provided “Oxygen concentration” setting is correct Operating conditions: 18 – 28 °C (64 – 82 °F) 30 – 70% RH

\* Test equipment and methods are selected and controlled to ensure the uncertainty coverage is no more than 30% of the disclosed tolerance.

\* Values are expressed in body temperature, pressure, saturated (BTPS), in accordance with ISO 80601-2-74.

<sup>§</sup> To meet safety and performance standards, the system has been tested using the worst case breathing system setup. This setup includes 900PT561 with OPT970 interface for bypassed airway, and nasal interfaces for all other cases.

## WARNINGS

Do not use the unit at an altitude above 2000 m (6000 ft) or outside a temperature range of 18 – 28 °C (64 – 82 °F). Doing so may affect the quality of the therapy or injure the patient.

Designed to conform to the requirements of:

IEC 60601-1:2005+AMD1:2012 +AMD2:2020

IEC 60601-1-2:2014+AMD1:2020

ANSI/AAMI ES 60601-1:2005 and A1:2012 and A2:2021

CAN/CSA-C22.2 NO.60601-1:14+A2:2022 (R2022)

IEC 60601-1-6:2010+AMD1:2013+AMD2:2020

IEC 60601-1-8:2006+AMD1:2012+AMD2:2020

ISO 80601-2-74:2021

ISO-80601-2-90:2021

The unit complies with the electromagnetic compatibility requirements of IEC 60601-1-2. In certain circumstances, the unit may affect or be affected by nearby equipment due to the effects of electromagnetic interference. Excessive electromagnetic interference may affect the therapy delivered by the unit. If this should happen, try moving the unit or the location of the unit causing interference, or alternatively consult your healthcare provider. To avoid potential interference, do not place any part of the device or accessories within 30 cm (12”) of any portable or mobile radio frequency communication equipment.

Accessory equipment connected to the serial port of the device must be certified to either IEC 60601-1 or IEC 60950-1. Furthermore, all configurations shall comply with the system standard IEC 60601-1-1. Anyone who connects additional equipment to the signal input part or signal output part configures a medical system and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical services department or your local representative.

## Disposal instructions

### Unit disposal instructions

This unit contains electronics. Please do not discard with regular waste. At the end of the unit’s service life, return to Fisher & Paykel Healthcare or dispose according to local guidelines for disposing of electronics. Dispose according to Waste Electrical and Electronic Equipment (WEEE) directive in European Union.

### Consumables disposal instructions

Dispose of accessories, spare parts and packaging according to local guidelines. Place the breathing tube and water chamber in a waste bag at the end of use and discard with normal waste. Hospitals should discard according to their standard method for disposing of contaminated product.

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REF 185045495 REV 1 2026-05 © 2026 Fisher & Paykel Healthcare Limited
F&P, Airvo, Optiflow and AirSpiral are trademarks of Fisher & Paykel Healthcare Limited.