
F&P Bubble CPAP System - REF BC151/BC161 User Instructions

F&P and Flexitrunk are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, see www.fphcare.com/ip

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **Colombia (CO)** Tel: +57 3142852934 **France (FR)** **EC REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115

Korea (KR) Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 2150 **Sweden (SE)** Fisher & Paykel Healthcare AB, Svetsarvägen 15, 2tr, 17141 Solna, Sweden Tel: +46 8 564 76 680 **Switzerland (CH)** **CH REP** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **UK REP** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000











Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 45 25 70 85 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346

Rx only CE 0123

REF 185042847 REV V 2024-05 © 2024 Fisher & Paykel Healthcare Limited

Fisher & Paykel
HEALTHCARE

Pressure-Flow Table:

F&P Bubble CPAP Generator used with F&P FlexiTrunk™ Interface, Typical mean pressure values (cmH₂O) with no leak

		Probe setting (cmH ₂ O)															
Flow (L/min)	3	4	5	6	7	8	9	10	Flow (L/min)	3	4	5	6	7	8	9	10
4	3.1	4.1	5.1	6.1	7.1	8.1	9.1	10.0	10	4.3	5.3	6.2	7.2	8.1	9.0	9.9	10.8
5	3.2	4.2	5.3	6.3	7.3	8.3	9.2	10.1	11	4.6	5.6	6.5	7.4	8.4	9.2	10.1	10.9
6	3.4	4.4	5.4	6.4	7.4	8.4	9.3	10.2	12	4.9	5.8	6.8	7.7	8.5	9.4	10.2	11.1
7	3.6	4.6	5.6	6.6	7.6	8.5	9.5	10.4	13	5.2	6.1	7.0	7.9	8.8	9.6	10.4	11.3
8	3.8	4.8	5.8	6.8	7.8	8.7	9.6	10.5	14	5.5	6.4	7.3	8.2	9.0	9.8	10.6	11.4
9	4.0	5.0	6.0	7.0	7.9	8.9	9.8	10.6	15	5.8	6.7	7.5	8.4	9.2	10.0	10.8	11.6

F&P Bubble CPAP System – BC151/BC161 User Instructions

English en

Intended Use

The Fisher & Paykel Healthcare Bubble CPAP System is intended to provide CPAP to spontaneously breathing neonates and infants who require breathing support due to conditions associated with prematurity (such as Respiratory Distress Syndrome) or other conditions where CPAP is required or desired and is prescribed by a physician.

The Bubble CPAP System is for use in the hospital clinical environment such as the NICU (Neonatal Intensive Care Unit) and PICU (Pediatric Intensive Care Unit).

The intended patient population is premature and full term neonates up to a weight of 10 kg.

All components of the Bubble CPAP System are single use only.

Product Specifications

- CPAP range (mean): 3 – 10 cmH₂O
- Operational input flow range is 4 – 15 L/min.
- An initial range of 6 – 8 L/min is recommended.
- Gas Leakage @ 60 cmH₂O: <30 mL/min
- Humidification Output (Invasive Mode) >33 mg/L from 4 to 15L/min over temp range 20 °C to 26 °C
- Resistance to flow @ 2.5 L/min
Inspiratory Limb: 0.09 ± 0.032 cmH₂O
Expiratory Limb: 0.05 ± 0.022 cmH₂O (including 0.02 cmH₂O measurement uncertainty)
- Compliance @ 60 cmH₂O:
0.8 ± 0.16 mL/cmH₂O (including 0.07 mL/cmH₂O measurement uncertainty).
- Tube minimum internal diameter:
Inspiratory: 10.9 mm
Expiratory: 10.0 mm

- The **BC151** is for use with third party **short binasal prongs** (e.g. Hudson prongs).
- The **BC161** is for use with the **F&P FlexiTrunk Interface**.
- For use with F&P MR850 and MR730 humidifiers only.
- Dispose of product according to Hospital protocol. User may be exposed to breathing tract fluids during disposal.

Symbol Definitions

	Chamber water level incorrect		Chamber water level correct
	Consult Instructions for use		Type BF applied part
	Phthalate Free (PHT, DEHP, BBP, DBP Free)		Latex free
	Do not reuse		Catalogue number
	EC REP		RX only
	Date of manufacture		Manufacturer
	Use-by		7 Day use
	Temperature limitation		Do not use if package is damaged
	Batch code		CE Mark
	Caution		Overflow water level

Contraindications

- Non-spontaneous breathing
- Congenital abnormalities or malformations where bi-nasal prongs or nasal mask are contraindicated (e.g. choanal atresia).

- Congenital abnormalities or malformations where positive pressure therapies are contraindicated (e.g. diaphragmatic hernia and tracheo-oesophageal fistula).
- Nasal trauma/severe deformity that might be exacerbated by use of nasal prongs or nasal mask

WARNING

- DO NOT reuse any part of the **BC151** or **BC161**. It is manufactured for single patient use only. **Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.**
- **The Fisher & Paykel Bubble CPAP System is intended for use by adequately trained medical practitioners. Please contact a Fisher & Paykel Healthcare representative for suggested training material.**
- This product is intended to be used for a maximum of 7 days.
- DO NOT use if product or packaging is damaged.
- US Federal law restricts this device to sale by, or on the order of, a physician.
- Use patient oxygen monitoring.
- **Always use pressure monitoring to verify that the patient is receiving the prescribed CPAP level.**
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.
- Remove any sources of ignition, such as cigarettes or open flames
- DO NOT add or modify any attachments or accessories to the humidifier that are not specified or approved for use with the humidifier or accessory, or the humidifier

might not function correctly, affecting the quality of the therapy or injuring the patient.

- The water source must be at least 50 cm higher than the chamber.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- To prevent disconnection of the tubing or tubing system during use, especially during ambulatory use, only tubes in compliance with ISO 5367 or ISO 80601-2-74 should be used.
- Ensure any unused ports have their caps and/or plugs in place before use.
- DO NOT use an in-line nebulizer between the pressure relief manifold and the dry side of the humidification chamber.
- DO NOT bypass the pressure relief valve as this may cause damage/injury.
- DO NOT use if pressure relief manifold shroud has been removed or is missing.
- DO NOT soak, wash, sterilize or reuse this product. Avoid contact with non-approved chemicals.

Bubble CPAP Circuit

This Bubble CPAP circuit is designed for use with Fisher & Paykel Healthcare humidifiers only.

Product Specifications

Circuit Length	BC151	BC161
	Inspiratory 1.4 m Expiratory 1.2 m	Inspiratory 1.1 m Expiratory 1.2 m

Length includes the end adaptor and any assembled ends.

- Interface Connections: ISO 5356-1 Conical Connectors.
- For use with F&P MR850 and MR730 humidifiers only.

-  Latex free.

WARNING

- Check that the heater wire is evenly distributed along the blue inspiratory tube and is not bunched or kinked.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted turn the humidifier off.
- Check all connections are tight before use and after any adjustment.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT crush or pinch tubing.
- **This circuit is designed for use with the F&P Bubble CPAP System only.**
- Avoid prolonged contact with patient skin.
- Test circuit for occlusions and pressure leaks using the flow elbow provided before connection to the infant.

MR290 Humidification Chamber

A single use auto-fill chamber for use with Fisher & Paykel Healthcare humidifiers.

The chamber is to be supplied with sterile water, (USP or equivalent), from a flexible bag or rigid water container. DO NOT add other substances to the water. The water level is controlled by a dual float mechanism, one float acting as a backup to prevent the chamber from flooding should damage occur to the primary float mechanism.

WARNING

- An MR290 chamber MUST be used with the bubble CPAP system.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10 degrees.

Water Level Symbol

-  denotes correct water level
-  denotes water level too high


- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT spike the water source until the blue caps have been removed.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- Ensure the humidifier is mounted below the level of the patient.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.

Bubble CPAP Pressure Manifold

The pressure manifold incorporates a pressure relief valve which activates automatically in the event of a downstream circuit occlusion and which resets upon restoration of normal circuit pressure. The pressure manifold also provides additional ports to allow connection of external pressure monitoring devices and air/oxygen analysers.

Product Specifications

Maximum Pressure Limit:	17 cmH ₂ O at 8 L/min 24 cmH ₂ O at 15 L/min
① Manifold Pressure Port:	Female luer
② Manifold O ₂ Ports:	22 mm male or 15 mm female white
③ Pressure relief shroud:	Air/O ₂ inlet adaptor
④ Manifold Inlet Port:	22 mm/15 mm female
⑤ Manifold Outlet Port:	
⑥ Flow direction arrow	

WARNING

- To reduce the risk of unsafe circuit pressure, the pressure manifold must be used.
- DO NOT connect the pressure manifold to the outlet port of the humidification chamber as incorrect operation of the pressure relief mechanism may occur.
- Remove the correct port cap from the pressure manifold before connecting monitoring devices.
- USE ONLY the pressure manifold supplied with the Bubble CPAP System (White Cover with text "DO NOT REMOVE").

Bubble CPAP Generator


A single use water-bubble CPAP generator that is capable of producing CPAP levels from 3–10 cmH₂O. This unique CPAP generator produces CPAP with pressure oscillations. The CPAP generator is made up of three major components; the Bubbler, an Overflow container and a CPAP probe. The mean CPAP level remains constant due to the generator's "Auto-level" mechanism.

Product Specifications

Gas flow inlet port:	15 mm female
Gas flow exit port:	22 mm male
CPAP range (mean):	3 – 10 cmH ₂ O

Setting the CPAP Level

Using the fill funnel provided, fill the Bubble CPAP generator until water flows into the overflow container. Setting the CPAP level is as simple as raising or lowering the CPAP probe to the desired level. The number directly above the lid of the CPAP generator indicates the CPAP pressure set. (Refer to illustration 7 "Set CPAP Level").

CPAP Generator Water Level Symbol

-  denotes correct water level
-  denotes water level too low

Refill the CPAP generator if the water level falls below the minimum water level line.

Overflow Water Level Symbol Check and empty the overflow container once every 8 hours.

To remove, lift and then slide Overflow container out.

WARNING

- Ensure audible continuous bubbling is maintained.
- To ensure set CPAP is maintained, DO NOT operate this CPAP Generator at an angle of more than 5 degrees.
- If the CPAP generator isn't bubbling the CPAP pressure will be less than the amount set. Check for and minimize air leaks in the system and at the patient.
- If continuous bubbling is still not observed after minimizing air leaks, input flow rate may be increased up to a maximum of 15 L/min.
- DO NOT use this CPAP generator if it has been dropped or damaged in any way.
- Ensure the CPAP generator is mounted below the level of the patient.

Set up Instructions

Refer to illustrations on page 2.

1. Insert Humidification Chamber

-  Discard the chamber if the seals are not intact when received.
- Slide the chamber onto the humidifier base.
- Remove the blue caps.

2. Connect Waterbag

- Hang the water bag, unwind the water feed set and spike the water bag.
- Open the vent on the vented spike.
- Check the humidification chamber for water flow from the bag.
- If no water is visible in the chamber or water consumption is low, check that the bag is spiked properly and that the feed tube is not kinked or blocked. Try gently squeezing the bag to promote water flow. Ensure that the height of the water bag is at least **50 cm** higher than the humidification chamber. If in doubt replace the chamber.
-  Discard the chamber if the water level exceeds the maximum water level line.

3. Fill Bubble CPAP Generator

- Using the fill-funnel provided, fill the CPAP generator with sterile water until water flows into the overflow container.
- Set the CPAP probe to 10 ready for the leak test.
-  Ensure that the humidifier and Bubble CPAP generator are mounted below the patient.

4. Connect Bubble CPAP Pressure Manifold and Breathing Circuit

- Connect the oxygen tubing between the pressure manifold and the flow source.
- Connect the pressure manifold to the chamber inlet port.
- Connect the blue inspiratory tube to the remaining chamber port.

- Remove the blue caps from the inspiratory tube and connect the temperature probes and heater wire adapter.
-  Ensure that the heater wire is evenly distributed along the circuit and is not bunched or kinked.
- Connect test elbow to the expiratory and inspiratory tubes.
- Connect the expiratory tube to the CPAP probe.

5. Leak Test

-  Check all connections are tight before use.
- Set input flow rate to 1 L/min.
- Observe the CPAP generator. Gentle, audible bubbling is acceptable, no bubbling means unacceptable leakage.
- If no bubbling is observed, check the entire system.
-  Remove the correct port cap from the pressure manifold before connecting any monitoring device.

6. Set Flow Rate

- Adjust the flow rate to the prescribed input flow rate.
- Recommended 6 – 8 L/min.
- Allowable flow range: 4 – 15 L/min.

Note: Refer to the Pressure-Flow table on page 2

7. Set CPAP Level

- The number on the probe above the lid indicates the CPAP pressure in cmH₂O.
- Set the CPAP probe at the prescribed level (3 – 10 cmH₂O).

8. Connect Bubble CPAP Circuit to Infant Interface

-  Remove the flow test elbow and connect the circuits to the Infant Interface using the instructions provided with the interface.

- BC161 connects to F&P FlexiTrunk Interface.
- BC151 connects to short binasal prongs. (e.g. Hudson prongs)

9. Set the Humidifier

-  Ensure there is air flow present before turning on the humidifier.
- Turn on the humidifier.
 - If using the MR850, ensure it is on invasive mode.
 - If using the MR730, set the temperature to 40°, -3.
- Refer to MR850 or MR730 user instructions for further information.

10. Attach Interface to Infant

- Connect the interface to the infant using the instructions provided with the interface.

Checks During Operation

- Regularly observe that the water is feeding into the humidification chamber.
- Should the water level exceed the maximum level marked on the humidification chamber, the chamber should be replaced.
- Check all connections are tight before use and after any adjustment.
- Ensure air flow is present at all times. If air flow is interrupted, turn off the humidifier.
- Regularly observe the circuit for condensate. Drain as required in order to maintain set CPAP.
- Regularly observe the CPAP generator for bubbling. If bubbling is not observed, check for and minimize air leaks in the system and at the patient. If air leaks have been minimized, air flow may be increased to achieve continuous bubbling.
- Regularly observe the water level in the CPAP generator and overflow container. Refill the CPAP generator if the water level drops below the minimum water level line.
Check and empty the overflow container once every 8 hours.

- Monitor Patient Oxygen levels.
- **Always use pressure monitoring to verify that the patient is receiving the prescribed CPAP level.**

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