CPAP and NHF in neonatal intensive care

An overview of current neonatal literature and practices globally 1)



www.fphcare.com

Mechanisms of action and key benefits of CPAP and NHF



Magnenant, E. et al. Pediatr. Pulmonol. 37, 485-491 (2004). Bhutani, V. K. Manual of Neonatal Respiratory Care 3-15 (Springer US, 2012). Gupta, S. & Donn, S. M. Semin. Fetal Neonatal Med. 21, 204-211 (2016). Courtney, S. E. & Barrington, K. J. Clin. Perinatol. 34, 73-92, vi (2007). Pollett H.F., Reid W.D. Can Anaesth Soc. J. 24(5): 615-617 (1997).

6. de Klerk, A. In: Physiology of Humidification in

Critically III Neonates. Springer Berlin Heidelberg (2012).

8 Hough LL et al Pediatr Crit Care Med 13 A46 451 (2013

www.fphcare.com

NHF



While there is an overlap in the key mechanisms delivered by pressurebased and flow-based therapies, the primary mechanism of each therapy is different. This is important to consider when deciding choice of therapy:

NHF

When adequate flow is delivered on NHF (meeting or exceeding peak inspiratory demand), reducing anatomical dead space is the key mechanism.

CPAP

When delivering CPAP, the continuous distending pressure helps to establish and maintain functional residual capacity (FRC).

ten Brink, F. et al. Pediatri Crit Care Med. 14, e326-31 (2013).
 Woodhead, D. D. et al. J Perinatol. 26, 481-5 (2006).
 de Klerk, A. Adv. Neonatal Care. 8, 98-106 (2008).
 Saslow, J. et al. J Perinatol. 26, 476-80 (2006).
 Milési, C. et al. Intensive Care Med. 39, 1088-1094 (2013).

Hough, J. et al. Pediatr Crit Care Med. 15, e214-9 (2014).
 Sinha, I. et al. Chest. 148, 810-823 (2015).
 S. Collins, C. L. et al. Eur J. Pediatr. 173, 181-186 (2014).
 Roberts, C. et al. N Engl J Med. 375, 1142-51 (2016).



Evidence based guidance supporting the use of CPAP and NHF

The following information collates data from published literature. The body of evidence suggests that CPAP continues to be the standard of care in neonates <28 weeks gestational age (GA). There is also evidence to support the introduction of NHF (with rescue CPAP available) in the care of neonates ≥ 28 weeks GA.

< 28 weeks GA

	POSTEXTUBATION SUPPORT ¹⁶	CPAP CPAP continues to be the standard of care, however, NHF may be considered once infants are stable.
	PROLONGED CPAP ¹⁷⁸	CPAP CPAP continues to be the standard of care, however, NHF may be considered once infants are stable.
	PRIMARY RESPIRATORY SUPPORT ^{1,942}	CPAP CPAP continues to be the standard care and is used routinely as an alternative to invasive mechanical ventilation. There is emerging data describing the use of NHF in this population.

Disclaimer: Clinical judgement is necessary to assess the appropriate treatment for an individual patient. The use of NHF therapy is not typically supported for infants with extreme prematurity, severe respiratory distress syndrome, or untreated surfactant deficiency.

www.fphcare.com

1. Wilkinson, D. et al. Cochrane database Syst. Rev. 2, CD006405 (2016). 5. Collins, C. L. et al. J Pediatr. 162, 949–54.e1 (2013). 2. Campbell, D. M. et al. J of Peds. 26, 546 (2006). 6. Liu et al. J of Peds. 52, 271-6 (2014) 4. Manley, B. et al. N Engl J Med. 369, 1425–33 (2013).

7. Yoder, B.A. et al. J of Perinat. 37, 809 (2017) 8. Roehr. C.C. et al. Cli in Perinat. 43. 693-705 (2016)

≥ 28 weeks GA	Evidence
NHF + Rescue CPAP The use of NHF as an alternative to CPAP is associated with no difference in the rate of treatment failure, reintubation, and adverse outcomes (such as death, BPD and pneumothorax) and significantly less nasal trauma.	STRONG SUPPORT e.g. Cochrane Review
CPAP, then NHF once stable at the clinician's discretion ⁴	 STRONG SUPPORT e.g. consensus of published
The use of NHF as an alternative to prolonged CPAP may be considered once infants are stable. The benefits of NHF for older preterm infants are demonstrated in the body of literature.	expert opinion
CPAP or NHF + Rescue CPAP	✓ GENERAL SUPPORT
Emerging evidence suggests that the use of NHF (with rescue CPAP available) may be	e.g. emerging RCT data
considered once infants are stable, with no significant difference to intubation rates.	and consensus of published
The benefits of NHF for older preterm infants are demonstrated in the body of literature.	expert opinion



CPAP and NHF for postextubation support in neonates ≥ 28 weeks GA

High flow nasal cannula for respiratory support in preterm infants.

Cochrane Review¹ Wilkinson et al. 2016.

Data from six post-extubation RCTs was analyzed to assess the efficacy and safety of NHF:

Manley et al. 2013 (Australia)²

N Engl J Med.

- ightarrow 303 infants < 32 weeks GA
- → Primary outcome: Treatment failure within 7 days
- → Result: NHF was non-inferior to CPAP (Risk difference: 8.5%, margin of non-inferiority: 20%)

Campbell et al. 2006 (USA)³

J Perinatol.

- \rightarrow 40 infants \leq 1250 g at birth
- → Primary outcome: Need for intubation within 7 days
- → Result: Statistically significant difference favoring CPAP compared with NHF.

Liu et al. 2016 (China)⁴

Chinese J Pediatr.

- \rightarrow 256 infants 150 preterm < 7 days old
- → Primary outcome: Treatment failure within 7 days.
- → Result: No statistically significant difference between CPAP and NHF.

Collins et al. 2013 (Australia)⁵

J Pediatr.

- ightarrow 132 infants < 32 weeks GA
- → Primary outcome: Treatment failure within 7 days
- → Result: No statistically significant difference between CPAP and NHF.

Mostafa - Gharehbahgi et al. 2015 (Iran)

Zahedan J Res Med Sci.

- \rightarrow 85 infants 1250-2000 g at birth
- → Primary outcome: Treatment failure within 3 days
- → Result: No statistically significant difference between CPAP and NHF.

Yoder et al. 2013 (USA & China)⁷

Pediatrics.

- \rightarrow 432 infants (226 in post-extubation arm)
- → Primary outcome: Need for intubation within 3 days
- → Result: No statistically significant difference between CPAP and NHF.

www.fphcare.com

Wilkinson, D. et al. Cochrane database Syst. Rev. 2, CD006405 (2016).
 Manley, B. et al. N Engl J Med. 369, 1425–33 (2013).
 Campbell, D. M. et al. J of Peds. 26, 546 (2006).
 Liu et al. J of Peds. 52, 271-6 (2014)

Collins, C. L. et al. J Pediatr. 162, 949–54.e1 (2013).
 Mostafa-Gharehbaghi et al. Zahedan J Res Med Sci. 17 (2015).
 Yoder, B. et al. Pediatrics. 131, e1482-90 (2013).

The Cochrane Review ¹ found that compared with CPAP, the	No statistically significant difference in rate of treatment failure	Typical relative risk: 1.21, 95% CI 0.95 to 1.55 Data from 5 studies, 786 neonates
use of NHF for postextubation is associated with:	No statistically significant difference in rate of reintubation	Typical relative risk: 0.91, 95% CI 0.68 to 1.20 Data from 6 studies, 934 neonates
	No statistically significant difference in adverse outcomes i.e. pneumothorax	Typical relative risk: 0.35, 95% CI 0.11 to 1.06 Data from 5 studies, 896 neonates
	Significant reduction in rates of nasal trauma	Typical relative risk: 0.64, 95% CI 0.51 to 0.79 Typical risk difference: -0.14, 95% CI -0.20 to -0.08 Data from 4 studies, 645 neonates
SA		

1. Wilkinson, D. et al. Cochrane database Syst. Rev. 2, CD006405 (2016).

Cochrane Reviews are internationally recognized as the highest standard in evidence-based health care. All the existing primary research on a topic is collated to establish whether or not there is a conclusive evidence about a specific treatment. (Cochrane Collaboration, 2076)



An overview of key evidence: Flow rates on NHF and pressure settings



To achieve similar clinical outcomes as described in the Cochrane Review¹, it is important to consider how much flow is sufficient to provide the physiological benefits associated with NHF.

Setting Flow (NHF)

Results from the Cochrane Review demonstrate that adequate starting flow rates provide physiological and clinical benefits. Data from RCTs and guidance from leading experts suggests that NHF can be initiated between 4-6 L/min:

RCTs investigating NHF & CPAP: Comparison of flow rates and pressure settings



RCT: Randomized Controlled Tria

www.fphcare.com

Wilkinson, D. et al. Cochrane database Syst. Rev. 2, CD006405 (2016).
 Collins, C. L. et al. J Pediatr. 162, 949–54.e1 (2013).
 Manley, B. et al. N Engl J Med. 369, 1425–33 (2013).

Yoder, B. et al. Pediatrics. 131, e1482-90 (2013)
 Roberts, C.T. et al. N Engl J Med. 375, 1142-1151 (2016).
 Lavizzari, A. et al. JAMA Pediatrics. (2016)

on CPAP

Setting Pressure (CPAP)

Data from RCTs demonstrates that CPAP is typically initiated between 5-7 cmH₂O:







CPAP and NHF: A comparison of setup and interface design



CPAP

CPAP therapy is typically delivered using a dual limb circuit and, often a bubble generator.

Key characteristics of an interface designed to deliver pressure are:



Sealed interface



Prescribe pressure



Larger tubes lower resistance to flow

www.fphcare.com



Key characteristics of an interface designed to deliver flow are:



Unsealed interface



Prescribe flow



Narrower tubes (compared to CPAP) increase resistance to flow









For more information about CPAP or NHF, please contact your local representative or scan the QR Code to request a sample.



F&P, Optiflow and FlexiTrunk are trademarks of Fisher & Paykel Healthcare. For patent information, refer to www.fphcare.com/ip



www.fphcare.com

PM-619022 REV A © 2020 Fisher & Paykel Healthcare Limited