
Nasal Cannula

F&P Optiflow + Duet

- OPT962 (S)
- OPT964 (M)
- OPT966 (L)

Single Use **CE 0123 Rx Only**

Fisher & Paykel
HEALTHCARE

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** [EC REP] Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 2150 **Switzerland (CH)** [CH REP] Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** [UK REP] Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

REF 901629 REV D 2023-04

© 2022 Fisher & Paykel Healthcare Limited

Nasal Cannula

INDICATIONS FOR USE

Nasal cannula patient interface for delivery of humidified respiratory gases.

SYSTEM SPECIFICATIONS

- Permissible transport and storage conditions: -10°C to 50°C
- Refer to Airvo™ 2/myAirvo™ 2 User Manual for permissible operating conditions.

SETUP

- Airvo 2 or myAirvo 2 Humidifier with an AirSpiral™ heated breathing tube and chamber kit (e.g. 900PT561).

Flow Range:

| | |
|--------|-------------|
| OPT962 | 10-50 L/min |
| OPT964 | 10-60 L/min |
| OPT966 | 10-60 L/min |

or

- MR850 Respiratory Humidifier in invasive mode, RT series kit with 22 mm heated inspiratory tube and chamber (e.g. RT232).

Flow Range: 5-60 L/min

or

- F&P 950™ Respiratory Humidifier in Optiflow™ mode with F&P 950 Adult Heated Circuit kit (e.g. 950A40)

Flow Range: 5-70 L/min

Note: F&P 950 Respiratory Humidifier may not be available in all countries.

SETUP INSTRUCTIONS

Step 1

Select appropriate size. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.

Step 2

Adjust headstrap to fit. Do not over-tighten.

Step 3

- Ensure headstrap clip is attached to the headstrap, to prevent cannula from being pulled out of the nares.
- Cannula can become unattached if not used with the headstrap clip.

Step 4

Attach breathing tube clip to a secure location (e.g. clothing or bedding) to prevent cannula from pulling off face.

Step 5

If using an MR850 Respiratory Humidifier, attach the breathing tube clip to the heated breathing tube but ensure probe cable is not crushed by tubing clip.

⚠ WARNINGS

- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in an event of an interruption to gas flow) may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption of therapy, to avoid serious injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Do not use with an air entrainer, to prevent loss of therapy.
- Do not use system near a naked flame or any ignition source, including electrosurgery, electrocautery, or laser surgery instruments. Exposure to oxygen increases the risk of fire that may result in patient injury or death.
- Do not crush or stretch tube, to prevent loss of therapy.
- For single patient use only. Reuse may result in transmission of infectious substances. Attempting to reprocess will result in degradation of materials and render the product defective.

⚠ CAUTIONS

- Failure to use the setup described above can compromise performance and affect patient safety.
- Achievable flow rates may depend on individual patient and/or flow source.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Do not allow the prong(s) to seal in the nares. Occlusion may result in nasal injury.
- Check for condensate regularly. Drain as required.
- Do not use if packaging is not sealed.
- Ensure that patients in an awake prone position (i.e. on their stomach) are not lying on the interface tubing or heated breathing tube. Reposition the breathing tube clip if required. Failure to comply may result in damage to the device, causing loss of therapy.

DURING USE

As appropriate for the patient's condition, monitor for any disruptions to the patient receiving flow.

Notice to User: If a serious incident has occurred while using this device please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

HOSPITAL USE

- This product is intended to be used for a maximum of 14 days.
- Do not soak, wash or sterilize.

HOME AND LONG-TERM CARE FACILITY USE

This product is intended to be used for a maximum of 30 days providing daily and weekly cleaning instructions are followed (see the myAirvo 2 User Manual).


DISPOSAL INFORMATION

Dispose of products and packaging according to local guidelines. Hospitals and long-term care facilities should discard according to their standard method for disposing of contaminated product.

INTELLECTUAL PROPERTY INFORMATION

F&P, Airvo, myAirvo, 950, Optiflow and Duet are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

SYMBOL DEFINITIONS

| | | | | | |
|---|------------------------------------|---|--|---|--------------------------------------|
|  REF | Catalogue number |  | Warning or Caution |  | Consult instructions for use |
|  | Manufacturer |  | Conforms with medical device directive 93/42/EEC |  | Small nasal cannula size |
|  | Authorised European Representative |  | Storage temperature range -10°C min. (+14°F) to +50°C max. (+122°F) |  | Medium nasal cannula size |
|  | Batch code |  | Rx Only Prescription only |  | Large nasal cannula size |
|  | Date of manufacture |  | Single use |  | Correct FITTING INSTRUCTION method |
|  | Use-by date |  | Polyethylene Terephthalate |  | Incorrect FITTING INSTRUCTION method |

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F&P, Airvo, myAirvo, 950, Optiflow, dan Duet adalah merek dagang Fisher & Paykel Healthcare Limited.