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Ventilation alvéolaire  
7,9 L/min  
**+23 %**  
(↑ 1,5 L/min)

Ventilation alvéolaire  
8,8 L/min  
**+37 %**  
(↑ 2,4 L/min)

1. Morán I, Bellapart J, Vari A, Mancebo J. Heat and moisture exchangers and heated humidifiers in acute lung injury/acute respiratory distress syndrome patients. Effects on respiratory mechanics and gas exchange. *Intensive Care Med.* 2006 Apr;32(4):524–31.
2. Pitoni S, D'Arrigo S, Grieco DL, Idone FA, Santantonio MT, Di Giannatale P, et al. Tidal Volume Lowering by Instrumental Dead Space Reduction in Brain-Injured ARDS Patients: Effects on Respiratory Mechanics, Gas Exchange, and Cerebral Hemodynamics. *Neurocrit Care.* 2021;34(1):21–30.
3. Shimoda T, Sekino M, Higashijima U, Matsumoto S, Sato S, Yano R, et al. Removal of a catheter mount and heat-and-moisture exchanger improves hypercapnia in patients with acute respiratory distress syndrome: A retrospective observational study. *Medicine (Baltimore).* 2021 Sep 10;100(36):e27199.
4. Lellouche F, Taillé S, Lefrançois F, Deye N, Maggiore SM, Jouvét P, et al. Humidification performance of 48 passive airway humidifiers: comparison with manufacturer data. *Chest.* 2009 Feb;135(2):276–86.
5. Lellouche F, Delorme M, Brochard L. Impact of Respiratory Rate and Dead Space in the Current Era of Lung Protective Mechanical Ventilation. *Chest.* 2020 Jul;158(1):45–7.
6. Griffiths MJD, McAuley DF, Perkins GD, Barrett N, Blackwood B, Boyle A, et al. Guidelines on the management of acute respiratory distress syndrome. *BMJ Open Respir Res.* 2019;6(1):e000420.
7. Papazian L, Aubron C, Brochard L, Chiche JD, Combes A, Dreyfuss D, et al. Formal guidelines: management of acute respiratory distress syndrome. *Ann Intensive Care.* 2019 Jun 13;9(1):69.