


F&P 950

Adult Ventilator
Dual Heated
Circuit Kit

USER INSTRUCTIONS

REF 950A80
950A81
950A82

CE 0123

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraisópolis, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/ 售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** [EC REP] Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

F&P and F&P 950 are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

www.fphcare.com

Fisher & Paykel
HEALTHCARE

REF Z85001 REV J 2021.12 © 2021 Fisher & Paykel Healthcare Limited.

F&P 950™ Adult Ventilator Dual Heated Circuit Kit

Indications for use

The breathing set is an accessory to the F&P 950 Respiratory Humidifier. It is intended for delivery of heated humidified respiratory gases to adult and paediatric patients, within the limits of its stated technical specifications.

Contraindications

There are no contraindications for heated respiratory humidification.

Side effects

There are no known side effects of heated respiratory humidification.

Technical specifications

Compatible with the F&P 950 Respiratory Humidifier. Refer to humidifier, accessory, and patient interface user instructions for additional warnings, cautions, contraindications, side effects, and system information.

Tidal volume:	>200 mL
Interface connections:	ISO 5356-1 Conical Connectors

Breathing Tube Lengths:	
Inspiratory	1.60 m
Expiratory	1.60 m
Gas leakage (60 cmH₂O, BTPS):	
Dual limb	<40 mL/min
Resistance to flow (30 L/min, BTPS):	
Inspiratory	0.35 ± 0.03 cmH ₂ O
Expiratory	0.20 ± 0.03 cmH ₂ O
<i>With accessories:</i>	
Inspiratory	0.35 ± 0.03 cmH ₂ O
Expiratory	1.00 ± 0.11 cmH ₂ O
(includes measurement uncertainty of 0.02 cmH ₂ O)	

Compliance (60 cmH₂O, BTPS):	
Dual limb	2.27 ± 0.25 mL/cmH ₂ O

<i>With accessories:</i>	
Dual limb	2.36 ± 0.28 mL/cmH ₂ O
(includes measurement uncertainty of 0.02 mL/cmH ₂ O)	

Tube minimum internal diameter:	17 mm
Maximum tube surface temperature:	44 °C
Maximum operating pressure:	80 cmH ₂ O
Operating flow range	
Invasive mode:	5–60 L/min
Mask mode:	5–120 L/min
Optiflow mode:	5–70 L/min
Room temperature:	18–26 °C
Chamber maximum volume of water:	160 mL
Filter:	
Viral:	>99.99%
Organism:	ΦX174 bacteriophage
Bacterial:	>99.9997%
Organism:	bacillus subtilis
Mean particle size:	3 µm
NaCl:	98.04%

Symbol definitions

	Follow instructions for use - safety.
	Consult instructions for use. fphcare.com/950IFU
	Single use.
	Correct water level in water chamber.
	Incorrect water level; do not use water chamber.
	Indicates port for delivery of Metered Dose Inhaler medication

SST

Indicates component required for ventilator Short Self-Test

	Component can be disposed of in setup.
	Type BF applied part.
	Ventilator inspiratory port (flow to patient).
	Ventilator expiratory port (flow from patient)
	This product is not made with natural rubber latex.

	Temperature limitations.
	Manufacturer.
	Date of manufacture.
	Date of expiration.

Warnings, cautions, and notes

⚠ WARNINGS

- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- Single use. Do not reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm, or death.
- The use (or modification) of breathing circuit/chamber combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification, ventilator malfunction, or serious harm.
- Do not use product near any ignition source. Exposure to oxygen increases the risk of fire that may cause serious injury or death.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- This product is designed for the delivery of air and/or oxygen. It is not suitable for the delivery of flammable anaesthetic gas mixes or Heliox gas.
- Do not crush, stretch, or milk the tubing.
- Do not clean or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Set appropriate ventilator or flow source alarms to monitor therapy delivery.
- Perform a pressure and leak test on the breathing system before connecting to a patient.
- Check all connections are tight before use.
- Monitor circuit condensate to prevent occlusion or build-up of fluid. Drain as required.
- The humidifier should always be level and positioned lower than the patient.
- Do not use the chamber if it has been visibly damaged or dropped.
- Do not use the chamber if any of the sensor covers are dislodged or missing from the chamber.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not fill the chamber with water in excess of 37 °C.
- Do not use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- 14 Do not use beyond the 14-day maximum duration of use.
- When nebulized drugs are used resistance to flow should be monitored and the filter replaced, following standard hospital procedure.
- Change filter provided every 24 hours or sooner if notable deterioration occurs, following standard hospital procedure.

⚠ CAUTIONS

- Do not touch the hot surface of the chamber base. Failure to comply may result in a skin burn.
- Do not cover the circuit with materials such as towels, pillows, or bed linen. Failure to comply may result in a skin burn.
- Avoid leaving tubing in prolonged contact with patient’s skin. Failure to comply may result in a skin burn.
- Do not operate the chamber without water or gas flow. Failure to comply may result in low levels of humidity.
- If using a Puritan Bennett™™ 840 or Puritan Bennett™ 980 ventilator, use the yellow SST adapter to complete the short self-test (SST). Refer to ventilator manufacturer for consequences of failure to comply with SST.
- The water source must be at least 50 cm higher than the chamber. Failure to comply may result in low levels of humidity.

NOTES

- Use USP Sterile Water for Irrigation, or equivalent. Adding other substances may have adverse effects.
- User may be exposed to breathing tract fluids during disposal. Breathing circuits must be treated as contaminated. Dispose of according to local, state or country specific environmental laws and regulations.
- If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and, for European Union member countries, the Competent Authority in your country.

*Puritan Bennett is a trademark of Medtronic, Inc.

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*Puritan Bennett là thương hiệu của Medtronic, Inc.