**F&P** **950**

## F&P 950™ Respiratory Humidifier

### USER INSTRUCTIONS

## Languages

Language	Section
<b>English</b> .....	A
<b>Dansk</b> (Danish).....	B
<b>Nederlands</b> (Dutch).....	C
<b>Suomi</b> (Finnish).....	D
<b>Français</b> (French).....	E
<b>Deutsch</b> (German).....	F
<b>Italiano</b> (Italian) .....	G
<b>Norsk</b> (Norwegian) .....	H
<b>Polski</b> (Polish) .....	I
<b>Português</b> (Portuguese) .....	J
<b>EspaÑol</b> (Spanish) .....	K
<b>Svenska</b> (Swedish) .....	L

## Contents

---

Indications for use .....	A-4
Operating principle .....	A-4
Package contents .....	A-4
F&P 950 Respiratory Humidifier set up .....	A-5
User interface.....	A-7
Alarms.....	A-9
Information and service menus.....	A-12
Cleaning and maintenance .....	A-13
Warnings, cautions, and notes.....	A-14
Symbol definitions.....	A-16
Technical specifications .....	A-17
<i>Product specifications.</i> .....	A-17
<i>Operating conditions.</i> .....	A-18
<i>Storage conditions.</i> .....	A-18

## Indications for use

The F&P 950 Respiratory Humidifier is intended to provide heat and humidity to respiratory gases delivered to patients. It is for use in a professional healthcare facility by a health professional.

## Operating principle

The F&P 950 Respiratory Humidifier provides heat and humidity to medical gases by passing the gas through a heated water chamber and heated breathing tubes.

The amount of heating is controlled based on the gas temperature measured at different parts of the humidifier.

## Package contents

**F&P 950 Heaterbase**

(e.g. 950ANZ)

### ACCESSORIES TO COMPLETE THE F&P 950 RESPIRATORY HUMIDIFIER

**Equipment mount**

(e.g. 900MR030)


**F&P 950 Sensor Cartridge**

(e.g. 950S01)

**Power cord**

(e.g. 950XPE)

**F&P 950 Expiratory Heater Wire Adapter**

(e.g. 950X00)


**F&P 950 Breathing Circuit Kit**

(e.g. 950A81, 950N80)

## F&P 950 Respiratory Humidifier set up

Attach the power cord and power cord retainer to the heaterbase.

Attach the sensor cartridge to the heaterbase.

### **WARNING**

The heaterbase must be mounted on an equipment mount capable of supporting 4 kg. Failure to comply may result in damage to the equipment mount and heaterbase and potentially cause serious patient harm.

**NOTE:** Ensure the heaterbase does not block access to the power supply outlet.

## F&P 950 Respiratory Humidifier set up


Set up the breathing circuit as described in the user instructions provided with the selected circuit kit.

When turning on the humidifier, an audible single beep sound should be heard.

## User interface

### SCREEN NAVIGATION

### MODES

The modes available will depend on the type of breathing circuit connected. The availability and operating principles for each mode are shown below.

#### Breathing Circuit Kit

Adult Breathing Circuit Kits

	Invasive	Mask	Optiflow
Adult Breathing Circuit Kits	<b>Invasive mode</b> is intended for patients whose upper airways have been bypassed by either a tracheostomy or endotracheal tube.	<b>Mask mode</b> is intended for patients whose upper airways have not been bypassed but are receiving gas via a face mask or similar.	<b>Optiflow mode</b> is intended for patients who require respiratory therapy through an Optiflow interface.

Neonatal Breathing Circuit Kit

<b>Neonatal</b>
-----------------

**Neonatal mode**  
is intended for neonatal patients who require respiratory support.

Optiflow Oxygen Kit

<b>Optiflow</b>
-----------------

**Optiflow mode**  
is intended for patients who require respiratory therapy through an Optiflow interface.

## User interface

When multiple modes exist for a type of breathing circuit kit, selection can be accessed via the drop-down menu button.

### COMFORT SETTINGS

In Mask and Optiflow modes it is possible to change the target temperature to provide conditions which may encourage patient comfort.

The available comfort settings are:

Mode	Default	Medium	Low
Invasive	37 °C	-	-
Mask	31 °C	29 °C	27 °C
Optiflow	37 °C	35 °C	33 °C
Neonatal	37 °C	-	-

The humidifier will reset to the default set-point if the mode is changed or the humidifier is turned off and back on. It is possible for service personnel to change the default set-point for Mask and Optiflow modes in the service menu.

## Alarms

### ALARM SIGNALS

The F&P 950 Respiratory Humidifier has visual and audible alarms to warn about interruptions to treatment. These alarms are generated by an intelligent alarm system, which processes information from the sensors and target settings of the unit and compares this information to pre-programmed limits.

### ALARM CONDITIONS

All possible alarm conditions are listed on the following pages, and all are classified as medium priority.

As the F&P 950 Respiratory Humidifier does not include patient monitoring, these alarms are considered technical indicators of humidifier performance. It is possible to have multiple alarm conditions occur simultaneously; under these conditions the humidifier uses an internal ranking system to display the highest-ranked alarm.

All alarms have been designed to be detectable within one meter of the humidifier, with the alarm signal being three beeps repeated every five seconds.

### CHECKING ALARM SYSTEM FUNCTIONALITY

To check alarm functionality, remove the heated breathing tube at any time while the humidifier is powered on **but not connected to a patient**. This action should activate the "Disconnection" visual and audible alarms. If either signal is absent, do not use the humidifier. Contact your servicing department for assistance.

In the event of an unexpected shutdown, the humidifier shall resume the operating mode and alarm settings (except algorithm-based alarms) prior to the reset if the interruption is for less than or equal to 30 seconds.

## Alarms

Alarm conditions	Required Action
<p><b>The Disconnection alarm</b> activates when the humidifier detects a disconnection of the inspiratory circuit.</p>	<p>Connect inspiratory circuit and fully insert the chamber for complete connection. Delay: &lt;10 seconds</p>
<p><b>The No Water alarm</b> activates when the humidifier detects that the chamber is empty or almost empty of water.</p> <p>The time-to-alarm signal generation is dependent on operating mode set-point and flow rates as these determine the water evaporation rate.</p>	<p>Replace the empty water bag. Delay: &lt;60 minutes*</p>
<p><b>The Low Temperature alarm</b> activates when the humidifier detects a low temperature condition at the patient end or chamber outlet for a continuous period of time.</p> <p>The alarm threshold is 2 °C below the set-point temperature.</p> <p>The time-to-alarm signal generation is dependent on the flow rates.</p>	<p>Check the humidifier is receiving flow within the range stated in this user instruction. Delay: &gt;10 minutes</p>
<p><b>The High Temperature alarm</b> activates when the humidifier detects a high temperature condition at the patient end or chamber outlet.</p> <p>The alarm threshold is a patient end temperature of &gt;43 °C.</p>	<p>Check the humidifier is receiving flow within the range stated in this user instruction. Check connections to the flow source. Delay: &lt;10 seconds</p>
<p><b>The Cartridge Disconnection alarm</b> activates when the humidifier detects that the sensor cartridge is not electrically connected.</p>	<p>Connect the sensor cartridge. Delay: &lt;10 seconds</p>
<p><b>The Tube Fault alarm</b> activates when the humidifier detects a potential fault in the breathing circuit.</p>	<p>Replace the breathing circuit when safe to do so. Delay: &lt;10 seconds</p>


\*For flow rates <1 L/min the alarm will be delayed by >90 minutes upon start-up.

## Alarms

ALARM CONDITIONS	REQUIRED ACTION
<b>The Service Required alarm</b> activates when the humidifier detects a potential fault that requires the humidifier to be serviced.	Turn off the humidifier as soon as appropriate, remove from service, and contact a technician. Delay: 10 seconds to 5 minutes
<b>The Caution Indicator LED light</b> illuminates when the humidifier detects that there is a potential fault with the humidifier and the screen is not operational.	Turn off the humidifier as soon as appropriate, remove from service, and contact a technician. Delay: <10 seconds
<b>The Cartridge Authentication alarm</b> activates when the humidifier does not recognize the sensor cartridge.  If this occurs, the user may choose to press "Accept" to acknowledge that the sensor cartridge is not Fisher & Paykel Healthcare approved. A sensor cartridge authentication failure icon will appear at the bottom right of the display.	To remove the sensor cartridge authentication failure icon, contact technician to replace sensor cartridge as soon as appropriate.
<b>The Cartridge Service Life alarm</b> activates when the humidifier detects the sensor cartridge has exceeded the recommended service life.  The sensor cartridge should be replaced at the next opportunity that it is safe to do so (when not in use by a patient).	Press "Pause Alarm" button to dismiss the alarm screen. Contact technician to replace sensor cartridge as soon as appropriate. Delay: 7 years from the date of manufacture or after 15,000 hours of use. If the alarm is paused, it will reappear 4 hours later.

## Alarms

### INFORMATION SIGNALS




INFORMATION SIGNALS	POSSIBLE ACTIONS
<p><b>The Cartridge Service Life warning</b> activates when the humidifier detects the sensor cartridge is approaching the end of its recommended service life. At this point the sensor cartridge has one month of service life remaining and a sensor cartridge should be made available for replacement.</p>	<p>Press "Remind me later" button to dismiss the warning screen. Contact technician to replace sensor cartridge as soon as appropriate. Delay: 30 days prior to expiry and will reappear every 24 hours, or every 8 hours if less than 7 days remaining.</p>

## Information and service menus

### OPTIONS SCREEN

The "Options" screen contains additional information about the humidifier and can be accessed by pressing the "Information" button. Tapping on each option enables navigation through the screens.

## Information and service menus

---

The servicing functions are password protected and should only be accessed by technical personnel. Refer to the Product Technical Manual for more information.

**NOTE:** The readings displayed in the Operating Functions page under the Information directory are additional information for troubleshooting purposes only. These values are not intended to be used to specify patient treatment or for patient diagnosis.

## Cleaning and maintenance

---

### CLEANING

Clean the heaterbase, sensor cartridge, or expiratory heater wire adapter using a cloth dampened with either isopropyl alcohol or neutral detergent. Always disconnect the humidifier from the power supply before cleaning.

#### NOTES:

- Do not immerse or autoclave the heaterbase, sensor cartridge, or expiratory heater wire adapter.
- Do not spray liquid into the vents or onto electrical connectors. Failure to comply may result in irreparable damage to the humidifier.

### ROUTINE MAINTENANCE

A full technical description, including routine maintenance and service data, is contained in the Product Technical Manual available from your supplier or Fisher & Paykel Healthcare.

**WARNING:** The Product Technical Manual must be followed for all servicing and maintenance of the humidifier. Failure to comply may impair performance of the humidifier or compromise safety (including potentially causing serious harm).

## Warnings, cautions, and notes

---

### **WARNINGS**

- This product is only designed and verified for use with accessories and spare parts approved by Fisher & Paykel Healthcare. Unauthorized accessories or spare parts which are used with the humidifier may impair performance of the humidifier, or compromise safety (including potentially causing serious patient harm), or result in increased electromagnetic emissions, or decreased electromagnetic immunity, resulting in improper operation.
- This product is designed for the delivery of air and/or oxygen. It is not suitable for the delivery of flammable anesthetic gas mixes or Heliox gas. Failure to comply may impair performance of the humidifier or compromise safety (including potentially causing patient harm).
- The humidifier should always be level and positioned lower than the patient. Failure to comply may impair performance of the humidifier or compromise safety (including potentially causing serious patient harm).
- Visually inspect components and accessories for damage before use and replace if damaged. Use of damaged components or accessories may impair performance of the humidifier or compromise safety (including potentially causing serious harm).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- Do not touch the electrical connectors and the patient simultaneously. Failure to comply may result in serious harm.
- Operation of the humidifier outside of the recommended operating conditions (as described in these user instructions) may impair performance of the humidifier or compromise safety (including potentially causing patient harm).
- Monitor circuit condensate every six hours to prevent occlusion or build-up of fluid. Drain as required. Failure to comply may impair performance of the humidifier or compromise safety (including potentially causing serious patient harm).
- Remove any sources of ignition, such as: cigarettes, an open flame, or materials which ignite easily at high oxygen concentrations.
- Follow the instructions of the oxygen device provider; keep oxygen regulators, cylinder valves, tubing, connections, and all other oxygen equipment away from oil, grease, or greasy substances. Spontaneous and violent ignition may occur if these substances come into contact with oxygen under pressure.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)
- The operation of high frequency surgical apparatus, or shortwave, or microwave equipment in the vicinity of the humidifier may adversely affect its performance. If this occurs, the humidifier should be removed from the vicinity of such devices.
- Do not use this product in or near a magnetic resonance imaging (MRI) scanner.

## Warnings, cautions, and notes

---


### CAUTIONS

- Ensure that Invasive mode is set for patients who have bypassed airways. Prolonged exposure to reduced humidity will result in patient harm including decreased mucociliary clearance, atelectasis, or pneumonia.
- Do not touch the hot surface of the heater plate, chamber base, or probes. Failure to comply may result in a skin burn.

### NOTES

- Use USP sterile water for irrigation, or equivalent. Adding other substances may have adverse effects.
- The F&P 950 Respiratory Humidifier contains an embedded software system licensed to Fisher & Paykel Healthcare by Microsoft. The license contains certain restrictions that are relevant to the use of the F&P 950 Respiratory Humidifier. Visit [www.fphcare.com/microsoftlicensing](http://www.fphcare.com/microsoftlicensing) for more information about such restrictions.

## Symbol definitions

				
Follow instructions for use – safety	Consult instructions for use	Manufacturer	Date of manufacture	Catalogue reference number
				
Batch code	Serial number	Type BF Applied Part	Class II equipment	Alternating current
	IP21			
Standby (On/Off)	IP Classification	Temperature limitations	Humidity limitations	USB 2.0
				
WEEE (Waste Electrical and Electronic Equipment)*	European representative*	CE Marking 93/42/EEC*	Regulatory Compliance Mark*	Raise finger guard
				
For USA: prescription only	Fragile, handle with care	Keep dry	Recyclable	Caution
				
Alarm	Warning: hot surface	Information screen	Alarm audible pause	Alarm audible paused
				
Invasive mode	Mask mode	Optiflow mode	Neonatal mode	Sensor cartridge authentication failure
				
Sensor cartridge service life warning	Accept	Cancel	Back arrow	Date of expiration

\* symbol displayed on select models only

## Technical specifications

### PRODUCT SPECIFICATIONS

	ADULT	NEONATAL
SPECIFICATION	VALUE	VALUE
<b>Dimensions (heaterbase only)</b>	240 mm (D) x 154 mm (W) x 253 mm (H)	240 mm (D) x 154 mm (W) x 253 mm (H)
<b>Weight (heaterbase and power cord only)</b>	3.45 kg	3.45 kg
<b>Supply frequency</b>	50/60 Hz	50/60 Hz
<b>Supply voltage</b>	<small>[REF] 950AXX<sup>1</sup></small> 230 V <small>[REF] 950JXX<sup>1</sup></small> 115 V <small>[REF] 950GXX<sup>1</sup></small> 100 V	<small>[REF] 950AXX<sup>1</sup></small> 230 V <small>[REF] 950JXX<sup>1</sup></small> 115 V <small>[REF] 950GXX<sup>1</sup></small> 100 V
<b>Power rating</b>	350 VA	350 VA
<b>Sound pressure level</b>	Alarms exceed 45 dBA @ 1 m	Alarms exceed 45 dBA @ 1 m
<b>Auditory alarm pause</b>	120 seconds	120 seconds
<b>Humidity performance (except in the event of a humidifier alarm or power failure or electromagnetic disturbance)</b>	<b>Invasive mode (37 °C):</b> >33 mg/L <b>Mask mode (31 °C):</b> >12 mg/L <b>Optiflow mode (37 °C):</b> >12 mg/L	<b>Neonatal mode (37 °C):</b> >33 mg/L
<b>Maximum temperature of delivered gas</b>	43 °C	43 °C
<b>Operating flow range (at 22 °C room temperature)</b>	<b>Invasive mode:</b> 5–60 L/min <b>Mask mode:</b> 5–120 L/min <b>Optiflow mode:</b> 5–70 L/min	<b>Neonatal mode:</b> 0.5–40 L/min See breathing circuit user instructions for specific flow ranges
<b>Time to reach set temperature (gas flow is required)</b>	<30 minutes	<30 minutes
<b>Maximum surface temperature of the breathing circuit (applied part section)</b>	44 °C	44 °C
<b>Temperature measurement accuracy</b>	± 2 °C	± 2 °C
<b>Component service life</b>	<b>Heaterbase:</b> 7 years	<b>Heaterbase:</b> 7 years

<sup>1</sup>XX represents the country code

## Technical specifications

### OPERATING CONDITIONS

	ADULT	NEONATAL
SPECIFICATION	VALUE	VALUE
<b>Room temperature</b>	18–26 °C	20–26 °C
<b>Incoming gas temperature</b>	Minimum = Room temperature Maximum = 10 °C above room temperature (at 30% relative humidity)	Minimum = Room temperature Maximum = 10 °C above room temperature (at 30% relative humidity)
<b>Operator position</b>	<1 m from heaterbase	<1 m from heaterbase
<b>Atmospheric pressure</b>	Minimum of 70 kPa (equivalent to a maximum altitude of 3000 m) Maximum of 106 kPa	Minimum of 70 kPa (equivalent to a maximum altitude of 3000 m) Maximum of 106 kPa

### STORAGE CONDITIONS

	ADULT	NEONATAL
SPECIFICATION	VALUE	VALUE
<b>Temperature</b>	-20–60 °C	-20–60 °C
<b>Humidity</b>	10–95% relative humidity non-condensing	10–95% relative humidity non-condensing

Dispose of according to standard hospital procedure and national regulations for electrical and electronic equipment.

Note that a lithium cell battery is included in the sensor cartridge.

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REF 618942 REV B 2021-02 © 2021 Fisher & Paykel Healthcare Ltd
F&P, F&P 950, Optiflow, AirSpiral, Evaqua, SensiDomes, and Thermadapt are trademarks of Fisher & Paykel Healthcare Limited.