
Chinstrap

Instructions for use


REF BC351-10, BC353-10, BC355-10, BC357-10
Pack of 10

BC351
20-26 cm


BC353
26-32 cm

BC355
32-38 cm

BC357
38-44 cm

English	4	Latviešu (Latvian)	46
العربية (Arabic)	6	македонски (Macedonian)	48
български (Bulgarian)	8	Bahasa Melayu (Malay)	50
Česky (Czech)	10	Nederlands (Dutch)	52
Dansk (Danish)	12	Norsk (Norwegian)	54
Deutsch (German)	14	Polski (Polish)	56
Ελληνικά (Greek)	16	Português (Portuguese)	58
Español (Spanish)	18	Português (Brasil) (Brazilian Portuguese)	60
Eesti keel (Estonian)	20	Română (Romanian)	62
Suomi (Finnish)	22	Русский (Russian)	64
Français (French)	24	Slovenčina (Slovak)	66
Français Canadien (French Canadian)	26	Slovenščina (Slovenian)	68
עברית (Hebrew)	28	Srpski (Serbian)	70
Hrvatski (Croatian)	30	Svenska (Swedish)	72
Magyar (Hungarian)	32	ไทย (Thai)	74
Bahasa Indonesia (Indonesian)	34	Türkçe (Turkish)	76
Íslenska (Icelandic)	36	Urdu (Urdu)	78
Italiano (Italian)	38	Tiếng Việt (Vietnamese)	80
日本語 (Japanese)	40	繁體中文版 (Traditional Chinese)	82
한국어 (Korean)	42	简体中文版 (Simplified Chinese)	84
Lietuvių (Lithuanian)	44		

FlexiTrunk™ Chinstrap User Instructions

English (en)

Intended use

Chinstraps are intended to be used to control mouth leak during infant and neonatal nasal CPAP by gently holding the mouth closed.

Contraindications

The following are contraindications related to nasal CPAP and the use of nasal CPAP interfaces:

- non-spontaneous breathing
- congenital abnormalities or malformations where bi-nasal prongs or nasal mask are contraindicated (e.g. choanal atresia)
- congenital abnormalities or malformations where positive pressure therapies are contraindicated (e.g. diaphragmatic hernia and tracheoesophageal fistula)
- nasal trauma/severe deformity that might be exacerbated by the use of nasal prongs or nasal mask.

Side effects

The following are the known side effects related to nasal CPAP and the use of CPAP interfaces:

- Gastric distension and feeding intolerance of minor severity are undesirable side effects of using CPAP therapy.
- Nasal trauma and skin injury of minor to moderate severity are undesirable side effects of using CPAP interfaces.

General warnings

- The components of the patient interface are designed for single patient use only. Reuse may result in transmission of infectious substances or interruption to treatment leading to serious harm or death.
- This product is intended to be used for a maximum of 7 days.
- DO NOT soak, wash or sterilize this product.
- Strangulation/suffocation hazard – ensure that the chinstrap is positioned as shown, at all times.
- DO NOT overtighten the chinstrap. Overtightening can cause skin damage.
- When applied, the chinstrap should not restrict the patient from making deliberate jaw movements such as yawning.
- Ensure the label is visible. This makes sure that the Velcro® is facing away from the patient.
- If using with the F&P FlexiTrunk nasal CPAP interface, use only the recommended Fisher & Paykel Healthcare patient interface components. Ensure that the patient interface instructions provided are followed.

- DO NOT use if the product or packaging is damaged.
- DO NOT modify this product.
- Patient interface components are supplied in a clean state. Employ suitable methods to prevent bio-contamination during use and disposal.
- The Fisher & Paykel Healthcare chinstrap may only be used in a hospital or clinical setting where the patient is adequately monitored by trained medical staff.
- The use of this device is not without risk. Even if used as intended following all instructions and warnings provided, the following risks remain: infection and adverse reactions. These risks may result in serious injury.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

General cautions






- US federal law restricts this device to sale by, or on the order of, a physician.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

Checks during operation

Periodically:

- Ensure that the chinstrap is positioned as shown, at all times. Failing to do so may cause a strangulation or suffocation hazard.
- Inspect and adjust to ensure that the chinstrap is not overtightened by checking the skin under the strap.
- Discontinue use if any skin irritation is observed.

Symbol definitions

	Consult instructions for use		Catalog number
	Medical device		Prescription only
	Do not reuse		Manufacturer
	Authorized representative in the European Community		Use by date
	Date of manufacture and country of manufacture		CE mark
	Storage temperature limit		Batch code
	Importer		Distributor
	Authorized representative for Switzerland		UK responsible person

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REF 185043545 REV I 2023-08 © 2023 Fisher & Paykel Healthcare Limited
F&P and FlexiTrunk are trademarks of Fisher & Paykel Healthcare Limited.