


**F&P** FlexiTrunk





## Nasal CPAP Interface





### Instructions for use


---



13

15

14

<b>English</b>	6	<b>Norsk</b> (Norwegian)	143
العربية (Arabic)	11	<b>Polski</b> (Polish)	148
<b>български</b> (Bulgarian)	16	<b>Português</b> (Portuguese)	154
<b>Česky</b> (Czech)	22	<b>Português (Brasil)</b> (Brazilian Portuguese)	160
<b>Dansk</b> (Danish)	27	<b>Română</b> (Romanian)	166
<b>Deutsch</b> (German)	32	<b>Русский</b> (Russian)	172
<b>Ελληνικά</b> (Greek)	38	<b>Slovenčina</b> (Slovak)	178
<b>Español</b> (Spanish)	44	<b>Slovenščina</b> (Slovenian)	183
<b>Eesti keel</b> (Estonian)	50	<b>Srpski</b> (Serbian)	188
<b>Suomi</b> (Finnish)	55	<b>Svenska</b> (Swedish)	193
<b>Français</b> (French)	60	ไทย (Thai)	198
<b>Français Canadien</b> (French Canadian)	66	<b>Türkçe</b> (Turkish)	204
עברית (Hebrew)	72	<b>Urdu</b> (Urdu)	209
<b>Hrvatski</b> (Croatian)	77	<b>Tiếng Việt</b> (Vietnamese)	214
<b>Magyar</b> (Hungarian)	82	繁體中文版 (Traditional Chinese)	219
<b>Bahasa Indonesia</b> (Indonesian)	88	简体中文版 (Simplified Chinese)	224
<b>Íslenska</b> (Icelandic)	94		
<b>Italiano</b> (Italian)	99		
日本語 (Japanese)	105		
한국어 (Korean)	110		
<b>Lietuvių</b> (Lithuanian)	115		
<b>Latviešu</b> (Latvian)	120		
<b>македонски</b> (Macedonian)	125		
<b>Bahasa Melayu</b> (Malay)	131		
<b>Nederlands</b> (Dutch)	137		

# **FlexiTrunk™ Nasal CPAP Interface User Instructions**

English [en](#)

## **Intended use**

The FlexiTrunk nasal CPAP interface is designed to deliver noninvasive positive pressure therapy to a spontaneously breathing patient weighing up to 10 kg in a hospital or clinical setting where the patient is adequately monitored by trained medical staff.

## **Product specifications**

BC190 BC191 BC192	The FlexiTrunk nasal CPAP interface is designed to connect to the Fisher & Paykel Healthcare Bubble CPAP System. The adapters allow connection to other manufacturers' CPAP circuits with 7.8 or 10 mm connectors.
Maximum input flow	15 L/min
Maximum allowable operating pressure	15 cmH <sub>2</sub> O

Please refer to the component section for product codes and dimensions.

Not made with natural rubber latex

Phthalate free (PHT, DEHP, BBP, DBP free)

## **Contraindications**

The following are contraindications related to nasal CPAP and the use of nasal CPAP interfaces:

- non-spontaneous breathing
- congenital abnormalities or malformations where bi-nasal prongs or nasal mask are contraindicated (e.g. choanal atresia)
- congenital abnormalities or malformations where positive pressure therapies are contraindicated (e.g. diaphragmatic hernia and tracheoesophageal fistula).
- nasal trauma/severe deformity that might be exacerbated by the use of nasal prongs or nasal mask.

## **Side effects**

The following are the known side effects related to nasal CPAP and the use of nasal CPAP interfaces:

- Gastric distension and feeding intolerance of minor severity are undesirable side effects of using CPAP therapy.

- Nasal trauma and skin injury of minor to moderate severity are undesirable side effects of using CPAP interfaces.

## **Warnings**

- The components of the patient interface are designed for single patient use only. Reuse may result in transmission of infectious substances or interruption to treatment leading to serious harm or death.
- This product is intended to be used for a maximum of 7 days.
- DO NOT soak, wash or sterilize this product.
- Handle with care. Use caution when positioning or disconnecting the interface. Avoid excessive pull forces, sharp objects and tubing holders. Damage to the tubing may cause loss of pressure and require immediate replacement.
- DO NOT use if the product or packaging is damaged.
- US federal law restricts this device to sale by, or on the order of, a physician.
- DO NOT modify this product.
- Use patient oxygen monitoring.
- Connection to other manufacturers' humidification systems increases the risk of condensate and/or high gas temperature delivery to the patient.
- Only air and oxygen mixtures are intended to be used with the FlexiTrunk interface. The materials used may not be compatible with anesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.
- Connect only pressure monitoring devices to the pressure monitoring port on the nasal tubing.
- The FlexiTrunk interface components are supplied in a clean state. Employ suitable methods to prevent bio-contamination during use and on disposal.
- DO NOT overtighten the bonnet straps.
- DO NOT cover the patient's face with the bonnet as suffocation may occur.
- USE ONLY the recommended FlexiTrunk interface components. Impaired positioning of the nasal prongs and mask can occur from non-approved substitutes.
- Septal necrosis can occur if the correct sizing or placement of the prongs and mask is not observed.
- Always use pressure monitoring to verify that the patient is receiving the prescribed CPAP level.
- DO NOT use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to loss of CPAP pressure.

- The use of this device is not without risk. Even if used as intended following all instructions and warnings provided, the following risks remain: infection and adverse reactions. These risks may result in serious injury.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

## **Set-up**

- Refer to the F&P Bubble CPAP System instructions prior to set up of the FlexiTrunk interface.
- If using the adapters, refer to the set-up instructions provided by the circuit manufacturer prior to set-up of the FlexiTrunk Interface.
- Refer to illustrations 1-15 for correct set up.
- For correct product sizing, refer to the components section.
- A maximum pressure relief device MUST be used. Pressure relief settings shall be suitable for infant and neonatal CPAP.

## **Checks before operation**

- When pressure monitor port on nasal tubing is not in use, the cap must be closed.
- Check that all circuit connections are tight before use and after any adjustment.

## **Set-up instructions**

Refer to the illustrations on page 1.

- Measure the baby's head circumference in centimeters.
- Choose the correct-sized bonnet. The bonnet should be snug fitting. Stretch the bonnet with your hands initially for ease of placement.
- Slip the bonnet onto the baby's head, completely covering the ears and with the back edge of the bonnet at the base of the neck. The front edge of the bonnet should be just above or on the eyebrows.

④ Use the sizing guide to choose suitable prongs or mask. Prongs should fill the nares completely without stretching the skin. Use the biggest possible size. The mask should not touch the edge of the nose, septum or the eyes.

- Choose the correct length of the nasal tubing. Use the shortest length possible. The clear tubing should not extend over the patient's forehead. The following may be used as a guide: 50 mm  $\leq$ 1.5 kg, 70 mm  $\leq$ 2.5 kg, 100 mm  $>$ 2.5 kg.
- Connect the prongs or mask to the nasal tubing ensuring that they are inserted fully. If using prongs: squeeze the sides of the prongs firmly to expose the grooves. Start from one end and insert prong grooves into the nasal tubing. Push the end in firmly.

- Remove the tear-off strips as required to adjust the nasal tubing angle to optimise the seal. DO NOT discard the tear-off strips as they may be needed when alternating between nasal prongs and masks.
- Attach the circuit to the nasal tubing in either orientation. Refer to the instructions provided by the circuit manufacturer for more detail.  
Note: Disconnect the FlexiTrunk interface gently by twisting the connectors. DO NOT pull forcefully, as this may damage the tubing. Handle with care.
- Set the gas flow to the prescribed level.

 Turn the humidifier on. Refer to the manufacturer's instructions.

 Place your hand close to the nasal prongs or mask to ensure there is gas flow present.

- If using prongs: Clear any nasal secretions before inserting the nasal prongs.
-  Ensure the nasal prongs are positioned at least 2 mm (0.08 inches) from the septum to avoid pressure necrosis. Adjust as necessary. Hourly checks of the septum's integrity are recommended.

- (11) If using a mask: Connect to the patient by placing the mask around the nose. The mask should sit comfortably around the patient's nose. It must not occlude the nostrils or touch the septum and should not be over the lip or eyes.
- (12) Set the height of the foam block by adding or removing the tear-off strips. Set-ups with masks may require adding strips for proper positioning. Place the strap over the foam block and tighten to hold the tubing firmly in place.
- (13) Hook the clips from the side straps of the bonnet to the glider. Pull both straps at the same time for central positioning. Affix the round Velcro® tabs to the blue strap of the bonnet. Use the least tension possible to maintain the prescribed CPAP level and stability of the interface.
- (14) DO NOT use angel frames, Christmas trees or similar devices as they can puncture the FlexiTrunk tubing and cause a loss of pressure.
- (15) Ensure the weight of the circuit is gently supported to reduce tension on the nasal tubing. A pillow or folded towel is recommended. The FlexiTrunk tubing must have free movement to allow the baby's head to move without restriction. The circuit should be kept below the level of the nasal tubing to minimize condensate at the baby.


#### **During use:**



- Check the patient frequently for signs of redness, pressure sores or irritation, which may result from extended prong or mask use. Hourly checks are recommended.
- Alternating between prongs and mask can reduce the risk of irritation. Alternating every 4 to 6 hours is recommended.
- Positioning of the mask should be adjusted for best fit and minimal leak. If excessive leak is still present, the flow may be increased to compensate. The maximum allowable flow rate is 15 L/min.
- Check for condensate regularly and drain as required.
- Replace the bonnet if it becomes loose during usage.

#### **Symbol definitions**


	Consult instructions for use		CE mark
	Latex free		Do not reuse
	Authorized representative in the European Community		Not made with phthalates (DEHP, DBP, BBP)
	Date of manufacture and country of manufacture YYYY-MM-DD		Use by date
	CE mark		Medical device
	Storage temperature limit		Caution
	Size Small		Size Medium
	Size Extra Large		Size Large
	Prescription only		Catalog number
	Manufacturer		Polyethylene terephthalate recycling
	Batch code		Handle with care
	Importer		Distributor
	Authorized representative for Switzerland		UK responsible person

## Components

FlexiTrunk Nasal Prongs	Product Code <b>[REF]</b>	A	B	Quantity
	BC3020-10	3.0 mm	2.0 mm	Pack of 10
	BC3520-10	3.5 mm	2.0 mm	
	BC4030-10	4.0 mm	3.0 mm	
	BC4540-10	4.5 mm	4.0 mm	
	BC5040-10	5.0 mm	4.0 mm	
	BC5050-10	5.0 mm	5.0 mm	
	BC5550-10	5.5 mm	5.0 mm	
	BC5560-10	5.5 mm	6.0 mm	
	BC6060-10	6.0 mm	6.0 mm	
	BC6070-10	6.0 mm	7.0 mm	
	BC6570-10	6.5 mm	7.0 mm	

FlexiTrunk Nasal Tubing		
 		
BC19x series		
Product Code <b>[REF]</b>	Size	Quantity
BC190-05	50 mm	Pack of 5
BC191-05	70 mm	
BC192-05	100 mm	

Chinstrap	Product Code <b>[REF]</b>	Size	Quantity
	BC351-10	20-26 cm	Pack of 10
	BC353-10	26-32 cm	
	BC355-10	32-38 cm	
	BC357-10	38-44 cm	

FlexiTrunk Nasal Mask	Product Code <a href="#">REF</a>	Size	Quantity
	BC800-10	S	Pack of 10
	BC801-10	M	
	BC802-10	L	
	BC803-10	XL	

FlexiTrunk Bonnet	Product Code <a href="#">REF</a>	Size	Quantity
	BC300-05	17-22 cm	Pack of 5
	BC303-05	22-25 cm	
	BC306-05	25-29 cm	
	BC309-05	29-36 cm	

FlexiTrunk Headgear	Product Code <a href="#">REF</a>	Size	Quantity
	BC325-05	29-36 cm	Pack of 5
	BC328-05	35-40 cm	
	BC331-05	40-45 cm	

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REF 185042952 REV Y 2024-05 © 2024 Fisher & Paykel Healthcare Limited
F&P and FlexiTrunk are trademarks of Fisher & Paykel Healthcare Limited.