

Flow Matters



Focus: First published clinical protocol detailing the application of nasal high flow (NHF) therapy for adult patients with acute hypoxemic respiratory failure (AHRF). Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. Ischaki E, Pantazopoulos I, Zakynthinos S. Eur Respir Rev. 2017

In this issue of *Flow Matters,* we introduce the first published clinical protocol (algorithm, as referred to) detailing the application of NHF therapy for adult patients with AHRF.

The development and application of evidencebased protocols within the ever-evolving global healthcare landscape serve to maximize efficiency and quality of care. Robust clinical protocols are supported by the best-available evidence. These protocols can also be regarded as clinical decision trees and/or algorithms.

The authors of this publication provide a narrative review of 99 NHF studies which were located through the use of a robust search strategy. They appraise and differentiate research findings that suggest a benefit from those that do not. The included studies compare NHF to conventional oxygen devices, and/or noninvasive ventilation (NIV). The narrative review then provides the foundation for a proposed algorithm for the application of NHF for adult patients with AHRF regardless of cause.

What is unique about this protocol?

- It is the first peer-reviewed, published clinical protocol for this patient population.
- This protocol is currently in active use in clinical practice.
- It was developed using robust clinical evidence (including Frat et al. 2015 and Hernández et al. 2016), which we have outlined in previous editions of this newsletter.

Click here for the newsletter featuring Frat et al. 2015.

Click here for the newsletter featuring Hernández et al. 2016.

Who does this protocol apply to?

- Adult patients with AHRF from almost any cause, defined as PaO₂/FiO₂ ratio < 300.
- Immunocompromised patients are included.
- Hypercapnic patients defined as PaCO₂
 > 45 mmHg and pH < 7.35 were excluded.

Why were hypercapnic patients excluded?

 The algorithm was developed to reflect its foundation trials in which, in the majority of cases, hypercapnic patients were excluded.

Refresher:

what's the PaO₂/FiO₂ ratio?

- It is the ratio of arterial oxygen partial pressure to the fraction of inspired oxygen.
- This may help to describe respiratory efficiency.
- For ARDS patients, a decreasing PaO₂/FiO₂ ratio may be associated with increased mortality.
- A normal ratio is > 500.

ARDS Severity	PaO ₂ /FiO ₂	Mortality Risk
Mild	200-300	27%
Moderate	100-200	32%
Severe	< 100	45%



 It is to be noted that NHF can be used effectively in patients with hypercapnia. NHF has been shown to decrease PaCO₂ in some populations, when compared to conventional low-flow oxygen, CPAP, and even BiPAP (Fraser et al. 2016, Braunlich et al. 2016, Biselli et al. 2016, and Fricke et al. 2016).

What are the key points for this protocol?

- This protocol is designed for use on adults with AHRF.
- This protocol describes strategies for:
- Initiation of NHF at 60 L/min, 37 °C, FiO₂ 100%
- Use of NHF for pre-oxygenation and peri-laryngoscopy oxygenation
- Monitoring whilst on NHF
- Titration of NHF
- Weaning from NHF

Ischaki E, Pantazopoulos I, Zakynthinos S. Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. Eur Respir Rev. 2017.

ARDS Definition Task Force, Ranieri VM, Rubenfeld GD, Thompson BT, Ferguson ND, Caldwell E, Fan E, Camporota L, Slutsky AS. Acute respiratory distress syndrome: the Berlin Definition. JAMA. 2012 Jun 20; 307(23):2526-2533.

Frat J, Thille A, Mercat A, Girault C, Ragot S, Perbet S et al. High-Flow Oxygen through Nasal Cannula in Acute Hypoxemic Respiratory Failure. N Engl J Med. 2015;372(23):2185-2196.

Hernández G, Vaquero C, González P, Subira C, Frutos-Vivar F, Rialp G et al. Effect of Postextubation High-Flow Nasal Cannula vs Conventional Oxygen Therapy on Reintubation in Low-Risk Patients. JAMA. 2016;315(13):1354-1361.


Fraser J, Spooner A, Dunster K, Anstey C, & Corley A. Nasal high flow oxygen therapy in patients with COPD reduces respiratory rate and tissue carbon dioxide while increasing tidal and end-expiratory lung volumes: A randomised crossover trial. Thorax. 2016;71(8), 759–761.


Bräunlich J, Köhler M, Wirtz H. Nasal highflow improves ventilation in patients with COPD. International Journal of Chronic Obstructive Pulmonary Disease. 2016;1077-1085.

Biselli PJ, Kirkness JP, Grote L, Fricke K, Schwartz AR, Smith P, Schneider H. Nasal high-flow therapy reduces work of breathing compared with oxygen during sleep in COPD and smoking controls: a prospective observational study. J Appl Physiol (1985). 2017 Jan 1; 122(1):82-88.

Fricke K, Tatkov S, Domanski U, Franke KJ, Nilius G, Schneider H. Nasal high flow reduces hypercapnia by clearance of anatomical dead space in a COPD patient. Respir Med Case Rep. 2016 Aug 26;19:115-117.


Algorithm used under Creative Commons license 4.0





Monitoring

Presence of one of the following: respiratory rate >35 breaths·min⁻¹, S_{DO2} <88–90%, thoraco-abdominal asynchrony and/or persistent auxiliary muscle use, respiratory acidosis (P_{aCO_2} >45 mmHg with pH <7.35)



Presence of one of the following within hours (maximum 48 h), besides optimum NHF titration: respiratory rate >35 breaths·min ⁻¹, SpO₂ <88–90%, thoraco-abdominal asynchrony and/or persistent auxiliary muscle use, respiratory acidosis (PacO₂ >45 mmHg with pH <7.35), haemodynamic instability [§]



FIGURE 1 Recommended algorithm for high-flow nasal cannula use in acute hypoxaemic respiratory failure in immunocompetent or immunocompromised patients. #: arterial oxygen tension (P_{aO_2}) /inspiratory oxygen fraction $(F_{iO_2}) <300$ (patients with arterial carbon dioxide tension $(P_{aCO_2}) >45$ mmHg and pH <7.35 are excluded); *: systolic arterial blood pressure <90 mmHg despite adequate fluid administration; *: the rationale for change in nasal high flow (NHF) settings are as follows. 1) Flow rate could be adjusted downwards by $5-10 \text{ L}\cdot\text{min}^{-1}$ per 1–2 h if none of the negative prognostic factors are present. However, if targets of arterial oxygen saturation measured by pulse oximetry (S_{PO_2}) and respiratory rate are not achieved, while the flow rate is <60 L·min⁻¹, increase of flow rate by $5-10 \text{ L}\cdot\text{min}^{-1}$ is preferred to raising Fio₂; 2) increase in Fio₂ causes increases in P_{aO_2} and S_{PO_2} ; 3) temperature can be set at 37°C or lower (31–34°C), based on the patient's comfort; [§]: haemodynamic instability is defined by heart rate >140 beats·min⁻¹ or change >20% from baseline. MV: mechanical ventilation; SOT: standard oxygen treatment.

Click here for more information about Ischaki et al. 2017. Click here to contact us.

Disclaimer: Any clinical opinions in this newsletter are the opinions of the contributing authors and are given for information purposes only. The clinical opinions are not intended as and do not substitute medical advice. There are no known financial interests by any study authors in the product or manufacturer. Fisher & Paykel Healthcare provided equipment for this study. Fisher & Paykel Healthcare is not aware of any significant risks or safety concerns related specifically to this use and patient population, not discussed in this article. F&P and Optiflow are trademarks of Fisher & Paykel Healthcare Ltd.



614370 REV A © 2019 Fisher & Paykel Healthcare Limitec