## Édition COVID-19

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- Wilson NM, Marks GB, Eckhardt A, et al. The effect of respiratory activity, non-invasive respiratory support and facemasks on aerosol generation and its relevance to COVID-19. *Anaesthesia*. 2021 Mar 30;10.1111/anae.15475.
- World Health Organization COVID-19 Clinical management: living guidance. 25 January 2021. Available from: <https://www.who.int/publications/i/item/WHO-2019-nCoV-clinical-2021-1/> [Accessed 5 Jul 2021].
- National Institutes of Health. Coronavirus Disease 2019 (COVID-19) Treatment Guidelines. Available from: <https://www.covid19treatmentguidelines.nih.gov/> [Accessed 5 Jul 2021].
- National Health Commission of the People Republic of China. Diagnosis and Treatment Protocol for COVID-19 Patients (Tentative 8th Edition). Available from: [http://en.nhc.gov.cn/2020-09/07/c\\_81565.htm](http://en.nhc.gov.cn/2020-09/07/c_81565.htm) [Accessed 5 Jul 2021].
- Alhazzani W, Evans L, Alshamsi F, et al. Surviving Sepsis Campaign Guidelines on the Management of Adults With Coronavirus Disease 2019 (COVID-19) in the ICU: First Update. *Crit Care Med*. 2021 Mar 1;49(3):e219-e234.
- National COVID-19 Clinical Evidence Taskforce. Australian guidelines for the clinical care of people with COVID-19. Available from: <https://covid19evidence.net.au/> [Accessed 5 Jul 2021].
- Chalmers JD, Crichton ML, Goeminne PC, et al. Management of hospitalised adults with coronavirus disease 2019 (COVID-19): a European Respiratory Society living guideline. *Eur Respir J*. 2021 Apr 15;57(4):2100048.
- Nasa P, Azoulay E, Khanna AK, et al. Expert consensus statements for the management of COVID-19-related acute respiratory failure using a Delphi method. *Crit Care*. 2021 Mar 16;25(1):106.
- French Intensive Care Societies guidelines. Expert recommendations on the intensive care of patients during a SARS-CoV-2 epidemic. Available from: <https://www.srlf.org/article/coronavirus/> [Accessed 5 Jul 2021].
- Patel M, Gangemi A, Marron R, et al. Retrospective analysis of high flow nasal therapy in COVID-19-related moderate-to-severe hypoxaemic respiratory failure. *BMJ Open Respir Res*. 2020;7(1):e000650.
- Bonnet N, Martin O, Boubaya M, et al. High flow nasal oxygen therapy to avoid invasive mechanical ventilation in SARS-CoV-2 pneumonia: a retrospective study. *Ann Intensive Care*. 2021 Feb 27;11(1):37.
- Agarwal A, Basmaji J, Muttalib F, et al. High-flow nasal cannula for acute hypoxemic respiratory failure in patients with COVID-19: systematic reviews of effectiveness and its risks of aerosolization, dispersion, and infection transmission. *Can J Anaesth*. 2020 Sep;67(9):1217-1248.
- Demoule A, Vieillard Baron A, Darmon M, et al. High-Flow Nasal Cannula in Critically Ill Patients with Severe COVID-19. *Am J Respir Crit Care Med*. 2020 Oct 1;202(7):1039-1042.
- Anesi GL, Jablonski J, Harhay MO, et al. Characteristics, Outcomes, and Trends of Patients With COVID-19-Related Critical Illness at a Learning Health System in the United States. *Ann Intern Med*. 2021 May;174(5):613-621.
- Wendel Garcia PD, Aguirre-Bermeo H, Buehler PK, et al. Implications of early respiratory support strategies on disease progression in critical COVID-19: a matched subanalysis of the prospective RISC-19-ICU cohort. *Crit Care*. 2021 May 25;25(1):175.
- Calligaro GL, Lalla U, Audley G, et al. The utility of high-flow nasal oxygen for severe COVID-19 pneumonia in a resource-constrained setting: A multi-centre prospective observational study. *EClinicalMedicine*. 2020 Nov;28:100570.
- Franco C, Facciolo N, Tonelli R, et al. Feasibility and clinical impact of out-of-ICU noninvasive respiratory support in patients with COVID-19-related pneumonia. *Eur Respir J*. 2020 Nov 5;56(5):2002130.
- Guy T, Créac'hcadec A, Ricordel C, et al. High-flow nasal oxygen: a safe, efficient treatment for COVID-19 patients not in an ICU. *Eur Respir J*. 2020 Nov 12;56(5):2001154.
- Mellado-Artigas R, Ferreyro BL, Angriman F, et al. High-flow nasal oxygen in patients with COVID-19-associated acute respiratory failure. *Crit Care*. 2021 Feb 11;25(1):58.
- Deng L, Lei S, Jiang F, et al. (2020). The Outcome Impact of Early vs Late HFNC Oxygen Therapy in Elderly Patients with COVID-19 and ARDS. *medRxiv*. <https://doi.org/10.1101/2020.05.23.20111450>
- Vianello A, Arcaro G, Molena B, et al. High-flow nasal cannula oxygen therapy to treat patients with hypoxemic acute respiratory failure consequent to SARS-CoV-2 infection. *Thorax*. 2020 Nov;75(11):998-1000.
- Duan J, Chen B, Liu X, et al. Use of high-flow nasal cannula and noninvasive ventilation in patients with COVID-19: A multicenter observational study. *Am J Emerg Med*. 2020 Jul 29;S0735-6757(20)30666-5.
- Westafer LM, Soares WE 3rd, Salvador D, et al. No evidence of increasing COVID-19 in health care workers after implementation of high flow nasal cannula: A safety evaluation. *Am J Emerg Med*. 2021 Jan;39:158-161.
- Gershengorn HB, Hu Y, Chen JT, et al. The Impact of High-Flow Nasal Cannula Use on Patient Mortality and the Availability of Mechanical Ventilators in COVID-19. *Ann Am Thorac Soc*. 2021 Apr;18(4):623-631.
- Gaekle NT, Lee J, Park Y, et al. Aerosol Generation from the Respiratory Tract with Various Modes of Oxygen Delivery. *Am J Respir Crit Care Med*. 2020 Oct 15;202(8):1115-1124.
- Bem RA, van Mourik N, Klein-Blommert R, et al. Risk of Aerosol Formation During High-Flow Nasal Cannula Treatment in Critically Ill Subjects. *Respir Care*. 2021 Jun;66(6):891-896.
- Li J, Fink JB, Elshafei AA, et al. Placing a mask on COVID-19 patients during high-flow nasal cannula therapy reduces aerosol particle dispersion. *ERJ Open Res*. 2021 Jan 25;7(1):00519-2020.
- Iwashyna T, Boehman A, Capeceletro J, et al. (2020). Variation in Aerosol Production Across Oxygen Delivery Devices in Spontaneously Breathing Human Subjects. *medRxiv*. <https://doi.org/10.1101/2020.04.15.20066688>
- Jermy MC, Spence CJT, Kirton R, et al. Assessment of dispersion of airborne particles of oral/nasal fluid by high flow nasal cannula therapy. *PLoS One*. 2021 Feb 12;16(2):e0246123.
- Kotoda M, Hishiyama S, Mitsui K, et al. Assessment of the potential for pathogen dispersal during high-flow nasal therapy. *J Hosp Infect*. 2020 Apr;104(4):534-537.
- Leung CCH, Joynt GM, Gomersall CD, et al. Comparison of high-flow nasal cannula versus oxygen face mask for environmental bacterial contamination in critically ill pneumonia patients: a randomized controlled crossover trial. *J Hosp Infect*. 2019 Jan;101(1):84-87.
- Miller DC, Beamer P, Billheimer D, et al. Aerosol Risk with Noninvasive Respiratory Support in Patients with COVID-19. *J Am Coll Emerg Physicians Open*. 2020 May 21;1(4):521-6.
- McGain F, Humphries RS, Lee JH, et al. Aerosol generation related to respiratory interventions and the effectiveness of a personal ventilation hood. *Crit Care Resusc*. 2020 May 26.
- Gall ET, Laguerre A, Noelck M, et al. (2020). Aerosol generation in children undergoing high flow nasal cannula therapy. *medRxiv*. <https://doi.org/10.1101/2020.12.10.20245662>
- Hui DS, Chow BK, Lo T, et al. Exhaled air dispersion during high-flow nasal cannula therapy versus CPAP via different masks. *Eur Respir J*. 2019 Apr 11;53(4):1802339.
- Loh NW, Tan Y, Taculod J, et al. The impact of high-flow nasal cannula (HFNC) on coughing distance: implications on its use during the novel coronavirus disease outbreak. *Can J Anaesth*. 2020 Jul;67(7):893-894.
- Hamada S, Tanabe N, Inoue H, et al. Wearing of medical mask over the high-flow nasal cannula for safer oxygen therapy in the COVID-19 era. *Pulmonology*. 2021 Mar-Apr;27(2):171-173.
- Dellweg D, Kerl J, Gena AW, et al. Exhalation Spreading During Nasal High-Flow Therapy at Different Flow Rates. *Crit Care Med*. 2021 Jul 1;49(7):e693-e700.
- Kaur R, Weiss TT, Perez A, et al. Practical strategies to reduce nosocomial transmission to healthcare professionals providing respiratory care to patients with COVID-19. *Crit Care*. 2020 Sep 23;24(1):571.
- Takazono T, Yamamoto K, Okamoto R, et al. Effects of surgical masks on droplet dispersion under various oxygen delivery modalities. *Crit Care*. 2021 Feb 27;25(1):89.
- Hamilton F, Arnold D, Bzdek BR, et al. Aerosol generating procedures: are they of relevance for transmission of SARS-CoV-2? *Lancet Respir Med*. 2021 Jul;9(7):687-689.
- Hamilton F, Gregson F, Arnold D, et al. (2021). Aerosol emission from the respiratory tract: an analysis of relative risks from oxygen delivery systems. *medRxiv*. <https://doi.org/10.1101/2021.01.29.21250552>
- Lyons C, Callaghan M. The use of high-flow nasal oxygen in COVID-19. *Anaesthesia*. 2020 Jul;75(7):843-847.
- Li J, Fink JB, Ehrmann S. High-flow nasal cannula for COVID-19 patients: low risk of bio-aerosol dispersion. *Eur Respir J*. 2020 May 14;55(5):2000892.
- Li J, Fink JB, Ehrmann S. High-flow nasal cannula for COVID-19 patients: risk of bio-aerosol dispersion. *Eur Respir J*. 2020 Oct 8;56(4):2003136.
- Dhand R, Li J. Coughs and Sneezes: Their Role in Transmission of Respiratory Viral Infections, Including SARS-CoV-2. *Am J Respir Crit Care Med*. 2020 Sep 1;202(5):651-659.
- Cook TM, El-Boghdady K, Brown J, et al. The safety of anaesthetists and intensivists during the first COVID-19 surge supports extension of use of airborne protection PPE to ward staff. *Clin Med (Lond)*. 2021 Mar;21(2):e137-e139.
- Morgensten J. 5 Dec 2020. "Aerosol Generating Medical Procedure" is a faulty paradigm. <https://first0em.com/aerosol-generating-medical-procedure-is-a-faulty-paradigm/> [Accessed 5 Jul 2021].
- Torjesen I. Covid-19: Risk of aerosol transmission to staff outside of intensive care is likely to be higher than predicted. *BMJ*. 2021 Feb 5;372:n354.
- Addleman S, Leung V, Asadi L, et al. Mitigating airborne transmission of SARS-CoV-2. *CMAJ*. 2021 Jul 5;193(26):E1010-E1011.

Fabricant : Fisher & Paykel Healthcare Ltd

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Organisme notifié : TÜV SÜD Product Service GmbH