COVID-19 Edition

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1. Gershengorn H, Hu Y, Chen J, et al. The Impact of High-Flow Nasal Cannula Use on Patient Mortality and the Availability of Mechanical Ventilators in COVID-19 [published online ahead of print, 2020 Oct 13]. *Ann Am Thorac Soc*. 2020.
2. Duan J, Chen B, Liu X, et al. Use of high-flow nasal cannula and noninvasive ventilation in patients with COVID-19: A multicenter observational study [Journal pre-proof]. *Am J Emerg Med*. 2020.
3. Guy T, Créac'hacdec A, Ricordel C, et al. High-flow nasal oxygen: a safe, efficient treatment for COVID-19 patients not in an ICU [published online ahead of print, 2020 Sep 9]. *Eur Respir J*. 2020;2001154.
4. Patel M, Gangemi A, Marron R, et al. Retrospective analysis of high flow nasal therapy in COVID-19-related moderate-to-severe hypoxaemic respiratory failure. *BMJ Open Respir Res*. 2020;7(1):e000650.
5. Vianello A, Arcaro G, Molena B, et al. High-flow nasal cannula oxygen therapy to treat patients with hypoxemic acute respiratory failure consequent to SARS-CoV-2 infection [published online ahead of print, 2020 Jul 23]. *Thorax*. 2020;thoraxjnl-2020-214993.
6. Westafer L, Soares W, Salvador D, et al. No evidence of increasing COVID-19 in health care workers after implementation of high flow nasal cannula: A safety evaluation [published online ahead of print, 2020 Oct 7]. *Am J Emerg Med*. 2020;S0735-6757(20)30881-0.
7. World Health Organization interim guidance for the Clinical management of COVID-19. 27 May 2020. Available from: [https://www.who.int/publications/i/item/clinical-management-of-severe-acute-respiratory-infection-when-novel-coronavirus-\(ncov\)-infection-is-suspected/](https://www.who.int/publications/i/item/clinical-management-of-severe-acute-respiratory-infection-when-novel-coronavirus-(ncov)-infection-is-suspected/) [Accessed 18 Sep 2020].
8. NIH. Coronavirus Disease 2019 (COVID-19) Treatment Guidelines. Available from: <https://www.covid19treatmentguidelines.nih.gov/> [Accessed 18 Sep 2020].
9. Alhazzani W, Møller MH, Arabi YM, et al. Surviving Sepsis Campaign: guidelines on the management of critically ill adults with Coronavirus Disease 2019 (COVID-19). *Intensive Care Med*. 2020;46(5):854-887.
10. National COVID-19 Clinical Evidence Taskforce. Australian guidelines for the clinical care of people with COVID-19. Available from: <https://covid19evidence.net.au/> [Accessed 18 Sep 2020].
11. Gaekle N, Lee J, Park Y, et al. Aerosol Generation from the Respiratory Tract with Various Modes of Oxygen Delivery. *Am J Respir Crit Care Med*. 2020;202(8):1115-1124.
12. Iwashyna T, Boehman A, Capecehatro J, et al. (2020). *Variation in Aerosol Production Across Oxygen Delivery Devices in Spontaneously Breathing Human Subjects*. medRxiv. <https://doi.org/10.1101/2020.04.15.20066688>
13. Jermy M, Spence C, Kirton R, et al. (2020). *Assessment of dispersion of airborne particles of oral/nasal fluid by high flow nasal cannula therapy*. medRxiv. <https://doi.org/10.1101/2020.05.14.20102517>
14. Kotoda M, Hishiyama S, Mitsui K, et al. Assessment of the potential for pathogen dispersal during high-flow nasal therapy [published online ahead of print, 2020 Apr]. *J Hosp Infect*. 2020;104(4):534-537.
15. Leung CCH, Joynt GM, Gomersall CD, et al. Comparison of high-flow nasal cannula versus oxygen face mask for environmental bacterial contamination in critically ill pneumonia patients: a randomized controlled crossover trial. *J Hosp Infect*. 2019;101(1):84-8.
16. Kaur R, Weiss TT, Perez A, Fink JB, Chen R, Luo F, Liang Z, Mirza S, Li J. Practical strategies to reduce nosocomial transmission to healthcare professionals providing respiratory care to patients with COVID-19. *Crit Care*. 2020 Sep 23;24(1):571.
17. Hui DS, Chow BK, Lo T, et al. Exhaled air dispersion during high-flow nasal cannula therapy versus CPAP via different masks. *Eur Respir J*. 2019;53(4):1802339.
18. Hui DS, Chan MTV, Chow B. Aerosol dispersion during various respiratory therapies: A risk assessment model of nosocomial infection to health care workers. *Hong Kong Med J*. 2014;20(suppl 4):S9-13.
19. Lyons C, Callaghan M. The use of high-flow nasal oxygen in COVID-19 [published online ahead of print, 2020 Apr 20]. *Anaesthesia*. 2020;75(7):843-847.
20. Li J, Fink JB, Ehrmann S. High-flow nasal cannula for COVID-19 patients: low risk of bio-aerosol dispersion. *European Respiratory Journal*. 2020 May 1;55(5).
21. Li J, Fink JB, Ehrmann S. High-flow nasal cannula for COVID-19 patients: risk of bio-aerosol dispersion [published online ahead of print, 2020 Oct 28]. *Eur Respir J*. 2020;56(4):2003136.

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620236 REV B © 2021 Fisher & Paykel Healthcare Limited 

www.fphcare.com

Fisher & Paykel
HEALTHCARE