F&P **Optiflow Junior 2**

Next generation care

Instructions for use
Home Nasal Cannula

OJR416HM

OJR418HM

REF


www.fphcare.com

Fisher & Paykel
HEALTHCARE

Sizing information

1. Prongs must not create a seal in the nares.
A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.

		WEIGHT GUIDE (KG)														
		1	2	3	4	5	6	7	8	9	10	15	20	25	30	SPARES
 OJR416HM L																
 OJR418HM XL																

**F&P
Wigglepads™ 2
WJR112**

3.  Expected to fit patient  May fit patient




Technical specifications

Maximum flow rates (L/min)

Compatible Platform	 myAirvo™ 2
Compatible Circuit	MYAIRVOKITI 900PT560 900PT560E MYAIRSPIRAL
 OJR416HM L	20
 OJR418HM XL	25

- Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
- Some products may not be available in your country. Please contact your local Fisher & Paykel Healthcare representative.
- myAirvo 2 flow rates are expressed in BTPS.
- Refer to circuit user instructions for minimum flow rates.



1 Apply Cannula



II

Remove Cannula

III

Replace F&P Wigglepads 2

F&P Optiflow Junior 2

Next generation care

English	6	Norsk (Norwegian)	76
العربية (Arabic)	9	Polski (Polish)	79
български (Bulgarian)	12	Português (Portuguese)	82
Česky (Czech)	15	Português (Brasil) (Brazilian Portuguese)	85
Dansk (Danish)	18	Română (Romanian)	88
Deutsch (German)	21	Русский (Russian)	91
Ελληνικά (Greek)	24	Slovenčina (Slovak)	95
Español (Spanish)	28	Slovenščina (Slovenian)	98
Eesti keel (Estonian)	31	Српски (Serbian)	101
Suomi (Finnish)	34	Svenska (Swedish)	104
Français (French)	37	ไทย (Thai)	107
Français Canadien (French Canadian)	40	Türkçe (Turkish)	110
עברית (Hebrew)	43	اردو (Urdu)	113
Hrvatski (Croatian)	46	Tiếng Việt (Vietnamese)	116
Magyar (Hungarian)	49	简体中文版 (Simplified Chinese)	119
Íslenska (Icelandic)	52	繁體中文版 (Traditional Chinese)	122
Italiano (Italian)	55		
日本語 (Japanese)	58		
한국어 (Korean)	61		
Lietuvių (Lithuanian)	64		
Latviešu (Latvian)	67		
Bahasa Melayu (Malay)	70		
Nederlands (Dutch)	73		

F&P Optiflow™ Junior 2 home nasal cannula

English 

Intended use

The Fisher & Paykel Healthcare Optiflow Junior 2 nasal cannula is a single use nasal cannula intended for use with a nasal high flow therapy system to deliver heated and humidified nasal high flow therapy and nasal low flow therapy to spontaneously breathing patients. This product is designed for use in homes and long term care facilities and must be prescribed by a physician.

The intended pediatric subpopulation targeted for use on the F&P Optiflow Junior 2 Home Nasal Cannula includes:

- Infants, 1 month up to 2 years of age
- Children, 2 years up to 12 years of age

Contraindications

Suspected or confirmed abnormalities or trauma to the skull or airway involving any abnormal connection may allow pressure to be transmitted to unintended anatomical structures or tissues. Use of NHF in this clinical scenario can lead to serious injury or death.

Side effects

Nasal trauma and skin injury are side effects of using NHF interfaces.

GENERAL WARNINGS

- This product is not intended for patients who cannot endure a disruption to therapy.
- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Using unauthorized equipment, accessories or spare parts with this product may impair performance or potentially cause serious patient injury.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient may result in serious injury or death. For example, in the event of an interruption to gas flow.
- If using supplemental oxygen, keep ignition sources away from the patient. For example, keep away from candles, cigarettes and open flames.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate pre-existing, untreated serious air leak syndrome which may lead to further serious injury or death.
- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of barotrauma, hypoxic injury and skin damage remain. These risks may result in serious injury or death.

GENERAL CAUTIONS

- This cannula is intended to be used for a maximum of 14 days. Replace Wigglepads 2 as required.
- Repeated removal and replacement of the Wigglepads 2 may cause skin injury.
- Wigglepads 2 may cause skin irritation. If symptoms continue, please contact your healthcare professional.
- Ensure the patient does not lie on the tubing as this may result in skin injury to the ears or face.
- Do not add heat above room temperature to the cannula. Directly heating any part of the cannula with an external heat source could result in serious injury.
- Insulating any part of the cannula could result in patient injury. For example, do not cover the cannula with bedding.
- Minimize direct cooling on the cannula, such as a fan to cool the patient or an air-conditioning unit/vent.
- When the cannula is not in use, store securely out of the reach of children and pets.
- Regularly monitor the patient to ensure skin is intact and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to minimize irritation.
- DO NOT soak or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers. Secretions on the cannula and the prongs can be removed by gently wiping with a damp cloth.
- DO NOT use this cannula on more than one patient.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious harm or death.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- Tubing may pose a risk of strangulation or airway restriction.
- DO NOT use if the product or its packaging has been tampered with.
- If using with supplemental oxygen, only use medical grade gas supplies. The gas supply used with this device may unexpectedly fail to deliver oxygen or flow.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or potentially cause serious patient harm.
- Notice to user: Any serious incident that has occurred in relation to the use of the device must be reported to Fisher & Paykel Healthcare.
- Only air and oxygen mixtures are intended to be used with the Optiflow Junior 2 Home Interface. The materials used may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.

Fitting instructions

Apply Cannula I

1. Select appropriate cannula size; a clear gap must be visible around each prong.
 - Ensure the patient's cheeks are clean and dry.
 - Connect the system to the myAIRVO 2 and ensure there is gas flow through the prongs.
2. Remove the first backing tabs from the F&P Wigglepads 2 and avoid touching the adhesive.
3. Insert the cannula into the nostrils. Ensure the cannula rests close to the nose without touching the septum. Do not stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.
4. Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.

WARNINGS

- DO NOT allow the prongs to seal in the nostrils. Sealing may result in skin or pressure injury.

CAUTIONS

- DO NOT place the Wigglepads 2 on the patient's eyes, ears or injured skin.
- Ensure the cannula is placed directly on the Wigglepads 2. Incorrect placement may result in skin injury.

Remove Cannula II

Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2. Starting from the outside, peel towards the nose.

Replace F&P Wigglepads 2 III

- Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.
1. Attach the replacement Wigglepads 2 to the cannula, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.
 2. Remove the second backing tabs and stick onto the patient's cheeks.

CAUTIONS

- Failure to follow the Wigglepads 2 removal/replacement instructions will increase the risk of skin injury.

If you require help with setting up, using or maintaining the Optiflow Junior 2 Home Nasal Cannula, or wish to report a problem or injury, please contact your healthcare professional or your Fisher & Paykel Healthcare representative.

Checks during operation

- Regular monitoring of the patient is necessary to ensure a slight gap between the cannula and septum is maintained, as well as the correct placement of the prongs in the nostrils. Reposition the cannula on the Wigglepads 2 if required.
- Regularly monitor the patient to ensure skin is intact and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to minimize irritation.
- To prevent nostril sealing, clean secretions from the cannula and the patient's face as required.
- Check the cannula remains secure. Replace the Wigglepads 2 if required.
- Ensure all connections are secure during use. Check the cannula is undamaged and that the tubes are free from kinks or blockages. Replace if necessary.
- Under excessive load, the cannula may disconnect.

CAUTIONS

- DO NOT wrap, insulate, stretch or crush the tubing as this may impair the performance of this product or potentially cause patient injury.
- Monitor condensate to prevent blockage or build-up of fluid. Drain away from the patient as required.

Disposal instructions

The product may be disposed of with general waste.

CAUTION

Using this product beyond 14 days may impair the performance of this product or compromise safety (including potentially causing serious patient injury).

Operating conditions

- Ambient temperature range: 18 to 26 °C.

Approved compatible equipment

- myAIRVO 2 Humidifier

Approved compatible heated breathing tube kits

- MYAIRVOKIT1 - AirSpiral Tube and Auto-fill Chamber kit (1 pack)
- 900PT560 - AirSpiral Tube (10 pack)
- 900PT560E - AirSpiral Tube (1 pack)
- MYAIRSPIRAL - AirSpiral Tube (1 pack)

Operating flow rates

The operating flow rates for each F&P Optiflow Junior 2 cannula size are dependent on the breathing circuit and humidifier in use. Refer to the operating flow rates for each cannula.



Sizing notes

1. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.
3.  Expected to fit patient.  May fit patient.

Technical specification notes

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. Some products may not be available in your country. Please contact your local Fisher & Paykel Healthcare representative.
6. myAirvo 2 flow rates are expressed in BTPS.
7. Refer to circuit user instructions for minimum flow rates.

Symbol definitions

	Single patient multiple use		European Conformity
	Not made with phthalates (DEHP, DBP, BBP).		Storage temperature range
	Batch code		Catalogue number
	Consult instructions for use fphcare.com/oj2-ifu		This product is not made with natural rubber latex.
	Authorized representative in the European Community		Switzerland authorized representative
 YYYY-MM-DD	Use by date		UK responsible person
 YYYY-MM-DD NZ: New Zealand	Date and Country of Manufacture NZ: New Zealand		Manufacturer
 >PET<	Recyclable PET		Maximum 14 days use
	Importer		Distributor
	Unique device identifier		Prescription only
	Medical device		

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型號：OJR416HM, OJR418HM

F&P, myAirvo, AirSpiral, Optiflow, Wigglepads, and the Seahorse, Crab, Starfish, Octopus and Turtle images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

CE 0123 Rx only  

REF 900673 REV C 2024-08 © 2024 Fisher & Paykel Healthcare Limited

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)**  Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

www.fphcare.com

Fisher & Paykel
HEALTHCARE