

Evidence-based Humidification

Patient Indications

ISSUE

1

2022



There are no contraindications to heating and humidifying inspired gases to physiological conditions or using a heated humidifier to achieve this. There are instances where heat and moisture exchangers (HMEs) are contraindicated and contexts where significant concerns over their use have been expressed in the literature, precluding recommendations for their general use.¹

These contraindications to HME use can be classified into two main categories of concerns:



Sub-optimal humidity

- Passive devices result in a net loss of moisture from the respiratory tract, which increases the risk of mucosal dysfunction and airway occlusion.
- Any situation that reduces the delivery of expired gas to the device or reduces the heat and humidity in expired gas further impacts performance.



Sub-optimal ventilation

- The position of a passive device in the circuit increases dead space, resistance to flow and work of breathing (WOB).
 - Without compensation, this means reduced alveolar ventilation and increased PaCO₂.
- Any condition which increases the risk of fluids contacting the HME risks occlusion and further elevation of both resistance to flow and WOB.

Using a heated humidifier assists natural defenses in the airway, promotes efficient gas exchange and ventilation, supports lung-protective strategies and can improve comfort with noninvasive respiratory support.

Absolute contraindications to HME use

- Secretions that are thick, copious, bloody or tenacious¹⁻⁷
- Bronchopleural fistula^{1,4,7,8}
- Large mask leak^{1,9}
- Incompetent or absent endotracheal tube cuff^{1,4,7,10,11}
- Expired tidal volume (Vt) < 70% of delivered Vt¹
- Neonates and infants¹²
- Hypothermia (body temperature < 32 °C)^{1,3,5-8}
- Nebulized medications, including during aerosol treatments^{1,3,5}
- High minute volume (> 10 L/min)^{1,8,15}

Relative contraindications* to HME use

- Noninvasive ventilation^{1,7,9,13}
- Low tidal volumes, including lung-protective ventilation^{1,6,7,14}
- Acute respiratory distress syndrome (ARDS)^{6,8,9,11,16,17}
- Hypercapnia management^{6,13}
- Long-term ventilation^{2,5,8,9}
- Airway burns^{6,8,18}
- Dehydration
- Acute respiratory failure^{1,9,13,19,20}
- Chronic respiratory failure²¹
- Asthma¹⁸
- Atelectasis⁸
- Chronic obstructive pulmonary disease (COPD)^{6,8}
- Difficult to wean patients^{7,8,21}
- Airway hemorrhage or trauma^{8,18}
- Pulmonary edema^{3,7}
- Immunocompromise⁸

1. Restrepo RD, Walsh BK. AARC Clinical Practice Guideline. Humidification during invasive and noninvasive mechanical ventilation: O12. *Respir Care*. 2012 May;57(5):782–8.
2. Branson RD. Secretion management in the mechanically ventilated patient. *Respir Care*. 2007 Oct;52(10):1327–8.
3. Rathgeber J. Devices used to humidify respired gases. *Respir Care Clin N Am*. 2006 Jun;12(2):165–82.
4. Branson RD. Preventing Moisture Loss from Intubated Patients. *Clin Pulm Med*. 2000;7(4):187–98.
5. Branson RD, Davis Jr K, Brown R, Rashkin M. Comparison of three humidification techniques during mechanical ventilation: patient selection, cost, and infection considerations. *Respir Care*. 1996;41(9):809–816.
6. Plotnikow GA, Accoce M, Navarro E, Tiribelli N. Humidification and heating of inhaled gas in patients with artificial airway. A narrative review. *Rev Bras Ter intensiva*. 2018 Mar;30(1):86–97.
7. Branson RD. The ventilator circuit and ventilator-associated pneumonia. *Respir Care*. 2005 Jun;50(6):774–7.
8. Ricard J-D, Cook D, Griffith L, Brochard L, Dreyfuss D. Physicians' attitude to use heat and moisture exchangers or heated humidifiers: a Franco-Canadian survey. *Intensive Care Med*. 2002 Jun;28(6):719–25.
9. Branson RD. Humidification of respired gases during mechanical ventilation: mechanical considerations. *Respir Care Clin N Am*. 2006 Jun;12(2):253–61.
10. Campbell RS, Davis KJ, Johannigman JA, Branson RD. The effects of passive humidifier dead space on respiratory variables in paralyzed and spontaneously breathing patients. *Respir Care*. 2000 Mar;45(3):306–12.
11. Morán I, Bellapart J, Vari A, Mancebo J. Heat and moisture exchangers and heated humidifiers in acute lung injury/acute respiratory distress syndrome patients. Effects on respiratory mechanics and gas exchange. *Intensive Care Med*. 2006 Apr;32(4):524–31.
12. Doctor TN, Foster JP, Stewart A, Tan K, Todd DA, McGrory L. Heated and humidified inspired gas through heated humidifiers in comparison to non-heated and non-humidified gas in hospitalised neonates receiving respiratory support. *Cochrane Database Syst Rev*. 2017 Feb 24;2017(2):CD012549.
13. Lellouche F, Maggiore SM, Deye N, Taillé S, Pigeot J, Harf A, et al. Effect of the humidification device on the work of breathing during noninvasive ventilation. *Intensive Care Med*. 2002 Nov;28(11):1582–9.
14. Hinkson CR, Benson MS, Stephens LM, Deem S. The effects of apparatus dead space on P(aCO₂) in patients receiving lung-protective ventilation. *Respir Care*. 2006 Oct;51(10):1140–4.
15. Cerpa F, Cáceres D, Romero-Dapuerto C, Giugliano-Jaramillo C, Pérez R, Budini H, et al. Humidification on Ventilated Patients: Heated Humidifications or Heat and Moisture Exchangers? *Open Respir Med J*. 2015 Jun 26;9:104–11.
16. Prat G, Renault A, Tonnelier J-M, Goetghebeur D, Oger E, Boles J-M, et al. Influence of the humidification device during acute respiratory distress syndrome. *Intensive Care Med*. 2003 Dec;29(12):2211–5.
17. Papazian L, Aubron C, Brochard L, Chiche J-D, Combes A, Dreyfuss D, et al. Formal guidelines: management of acute respiratory distress syndrome. *Ann Intensive Care*. 2019 Jun;9(1):69.
18. Boots RJ, George N, Faoagali JL, Druery J, Dean K, Heller RF. Double-heater-wire circuits and heat-and-moisture exchangers and the risk of ventilator-associated pneumonia. *Crit Care Med*. 2006 Mar;34(3):687–93.
19. Pelosi P, Solca M, Ravagnan I, Tubiolo D, Ferrario L, Gattinoni L. Effects of heat and moisture exchangers on minute ventilation, ventilatory drive, and work of breathing during pressure-support ventilation in acute respiratory failure. *Crit Care Med*. 1996 Jul;24(7):1184–8.
20. Jaber S, Chanques G, Matecki S, Ramonatxo M, Souche B, Perrigault P-F, et al. Comparison of the effects of heat and moisture exchangers and heated humidifiers on ventilation and gas exchange during non-invasive ventilation. *Intensive Care Med*. 2002 Nov;28(11):1590–4.
21. Girault C, Breton L, Richard J-C, Tamion F, Vandelet P, Aboab J, et al. Mechanical effects of airway humidification devices in difficult to wean patients*. *Crit Care Med*. 2003;31(5):1306–11.