




Humidification Chamber

USER INSTRUCTIONS REF MR225

CE 0123 Rx only -10 °C 50 °C MD

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构-费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** **EC REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s.@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** **CH REP** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **UK REP** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

English (en)

HUMIDIFICATION CHAMBER

The chamber is designed for humidifying gas with Fisher & Paykel Healthcare breathing circuits.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to the MR850 Respiratory Humidifier instructions for use for additional safety information, clinical claims and instructions for setup and use.

Unless otherwise stated, flow rates are expressed at BTPS (body temperature and pressure, saturated) conditions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
GAS PORTS	Inlet: 22 mm Male Outlet: 22 mm Male
MAXIMUM WATER CAPACITY	210 mL
COMPRESSIBLE VOLUME	Full: 90 mL Empty: 300 mL
COMPLIANCE	Full: 0.12 mL/cmH ₂ O Empty: 0.3 mL/cmH ₂ O
MAXIMUM OPERATING PRESSURE	20 kPa
MAXIMUM PEAK FLOW	80 L/min
GAS LEAKAGE @ 60 cmH₂O	<10 mL/min

SYMBOL DEFINITIONS

	Water level too high		European Conformity
	Maximum water level		Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Consult operating instructions		Manufacturer
	Medical device		Use-by date
	Distributor		Caution/Consult instructions for use
	Single use		7 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative
Rx only	Prescription only		Importer
	UK Responsible person		Switzerland Authorised representative

WARNINGS, CAUTIONS, AND NOTES

⚠️ WARNINGS

- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/user.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.

- Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):**
- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
 -  DO NOT use beyond 7 days maximum duration of use
 - DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
 - DO NOT use the chamber if the seals are not intact when received.
 - Ensure the humidifier is not tilted. Tilting the humidifier may result in water entering the breathing circuit.
 - Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
 - Check all connections are tight before use.
 - Ensure that the water level in the chamber is periodically monitored. Refill when necessary.
 - DO NOT fill the chamber with water in excess of 37 °C.
 - Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
 - DO NOT fill the chamber above the maximum fill level line. Liquid could enter the breathing circuit if the chamber is overfilled.
 - If disconnecting the breathing circuit from the chamber, the humidifier **MUST** be put into stand-by or turned off to prevent excessive temperatures.
 - Remove water feed inlet plug only prior to inserting water feed set. Never reinsert water feed inlet plug as this may cause leakage.

- Ensure the water bag is placed below the level of the humidifier after the chamber has been filled.
- In the event of the chamber leaking, switch the humidifier off and replace the chamber.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

NOTES

- For use under the supervision of trained medical personnel.
- Refer to Hospital protocol for disposal methods. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

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