

Reusable humidification chamber

USER INSTRUCTIONS










REF **MR340**

Rx only  0123     

Manufacturer  Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English	3
български (Bulgarian).....	6
Hrvatski (Croatian).....	9
Česky (Czech)	12
Dansk (Danish).....	15
Nederlands (Dutch).....	18
Eesti keel (Estonian).....	21
Suomi (Finnish).....	24
Français (French).....	27
Deutsch (German).....	30
Ελληνικά (Greek)	33
Magyar (Hungarian).....	36
Italiano (Italian).....	39
Latviešu (Latvian)	42
Lietuvių (Lithuanian)	45
Norsk (Norwegian).....	48
Polski (Polish).....	51
Português (Portuguese)	54
Português (Brasil) (Brazilian Portuguese).....	57
Română (Romanian)	60
Slovenčina (Slovak).....	63
Slovenščina (Slovenian)	66
Español (Spanish).....	69
Svenska (Swedish).....	72
Türkçe (Turkish)	75
繁體中文版 (Traditional Chinese).....	78
Bahasa Indonesia (Indonesian).....	81
ไทย (Thai).....	84
Tiếng Việt (Vietnamese).....	87
Русский (Russian).....	90
한국어 (Korean).....	93
日本語 (Japanese)	96

A**1****2****3****4****B****1****2****3****4****5****C**

REF MR340

Intended Use

The MR340 is a reusable chamber for heating and humidifying respiratory gases for pediatric and neonatal patients requiring respiratory support. It is intended for use with Fisher & Paykel Healthcare MR850 humidifiers in hospital or clinical environments and should only be prescribed by a physician. Addition of heat and humidity to the supply of cold and dry respiratory gases provided through ventilation is beneficial to prevent drying of the patient's airways.

The following product specifications are based on testing done with the MR850 humidifier and the MR340 reusable chamber.

Product Specifications

Maximum Flow Rate	70 L/min
Maximum Operating Pressure	80 hPa (cmH ₂ O) (8 kPa)

Warnings, Cautions, and Notes

Warnings

- Breathing circuits and accessories should be prescribed and initially set up by a healthcare professional qualified in respiratory care and used under the supervision of trained medical personnel.
- Do not lubricate connections, tubing or accessories. Failure to comply may result in fire, burns and loss of system performance.
- Ensure the MR340 o-ring is not twisted when inserted into the chamber base groove and that the chamber dome is pushed fully onto the chamber base as per the chamber assembly instructions.
- Do not exceed maximum peak flow rate of 70 L/min. Excessive flow rate can cause liquid to splash into the breathing circuit.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Do not fill the humidification chamber above the maximum fill-level line. Liquid could enter the breathing circuit if the chamber is overfilled.
- Do not touch the heater-plate or chamber base. Surfaces may exceed 85 °C (185 °F). Failure to comply may result in a skin burn.
- This product is reprocessable (reusable). The product must be reprocessed according to the reprocessing instructions for this product. Incorrect reprocessing may lead to the risk of cross-infection between patients and potential serious harm.
- Visually inspect products before each use on a patient. Discard if faulty or if there is any sign of deterioration such as cracks, tears or damage. These may cause gas leaks and loss of ventilation or respiratory support.
- Do not use the product above altitudes of 3000 m (9842 ft) or outside the recommended ambient temperature or flow range. Using the humidification system outside of specified ranges or above this altitude can affect the quality of the therapy or injure the patient.
- Use of non-compatible, non-ISO 5367, or non-ISO 80601-2-74 compliant circuit components may cause unintentional tube disconnection and loss of ventilation or respiratory support.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted, turn the humidifier off.
- Do not modify this product.
- Do not use breathing circuits, chambers, accessories or combinations which are not approved by Fisher & Paykel Healthcare.

Failure to comply with the above warnings and cautions may impair the performance of the device or compromise safety (including causing potential serious harm).

- California residents, please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Notes

- Keep product away from heat or ignition sources, especially when operating with oxygen equipment.
- This product is designed for the delivery of air, oxygen or a mix of both gases. It is not suitable for the delivery of flammable anesthetic gas mixes, Heliox or gas solutions, suspensions or emulsions that have not been evaluated.
- Do not use the product during X-rays, MRIs or other radio-frequency procedures.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used.
- This product is safe for use on children.

If a serious incident has occurred while using this device, please contact your Fisher & Paykel Healthcare representative and Competent Authority.

Assembly Instructions**Warnings**

- Before use on a patient, all components must be cleaned and disinfected. Sterilization is highly recommended.
- Before connecting to a patient, ensure the complete circuit functions correctly with required ventilator settings and appropriate ventilator or flow source alarms are set.
- Ensure that the chamber is always positioned lower than the patient so that condensate flows away from the patient interface. Refer to humidifier user instructions for further set up and correct use information.
- Do not fill the chamber with water with a temperature in excess of 37 °C (99 °F).
- Check all connections are tightly secured before use.
- If a water feed set is used, ensure the water bag is placed below the level of the humidifier after the chamber has been filled.

Refer to diagram **A** to assemble the chamber after reprocessing.

Refer to diagram **B**

- Assemble chamber onto the heater base as shown in **1**.
- Connect the tubes to the correct ports following the breathing circuit user instructions as shown in **2**.
- Remove the water fill port plug as shown in **2**.
- Connect the water line to the chamber and water bag as shown in **3**.
- Fill the chamber with water to the maximum water level **4**.
- Place the water bag below the level of the humidifier after filling the chamber **5**.

Refer to the breathing circuit user instructions for additional information and warnings.

Monitoring Instructions**Warnings:**

- Ensure the chamber contains water at all times. Check the water level in the humidification chamber regularly (as per diagram **C**) and top up as required. Absence of water will result in low levels of humidity.
- Ensure the humidifier and chamber are level and positioned lower than the patient. Operating the chamber at an angle in excess of 10° may result in water entering the breathing circuit.

Refer to the monitoring diagram **C** for instructions on monitoring and refilling the chamber.

Reprocessing Instructions

For reprocessing instructions (UI-623716), visit <http://www.fphcare.com/fp-reusables>

Warnings:

- The products must be disconnected from both the power and gas source prior to cleaning.
- Clean product prior to use and whenever product is visibly soiled. Follow the reprocessing instructions (UI-623716). Alternative cleaning methods may deteriorate or damage the product, reducing the product life and compromising therapy.
- Discard the product after 50 reprocessing cycles, or five years from date of manufacture, whichever occurs earlier.

Notes

- It is recommended that reprocessing commences as soon as is reasonably practical following use to prevent gross contaminants from drying on the product.
- Used products must be treated as being contaminated. The user may be exposed to breathing tract fluids. Follow all local, state and federal regulations with respect to environmental protection and correct disposal of these products.

Storage

- As a minimum, reprocessed breathing circuits and accessories should be stored in a clean plastic bag or container and away from direct sunlight or heat if not reusing immediately.

Symbol Definitions

Product does not contain phthalates or pyrogens	Unsafe for use in MRI environments	European Union authorized representative	Transportation and storage temperature limits	Not made with natural rubber latex
Caution symbol for warnings	Use-by date	Type BF Applied Part	Date of manufacture	Manufacturer
Rx only				
Prescription only	Instructions for Use	Lot number	Reference number	Medical Device
CE Mark	Do not recycle	Incorrect water level	Correct water level	No water

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