# Optiflow™+ Duet Nasal Cannula

REF MYDUET **SMALL**

REF MYDUET **MEDIUM**

REF MYDUET **LARGE**




Thank you for choosing F&P myAirvo™

Please read this entire guide before you first use your nasal cannula kit. Following the steps in this guide and referring to the images and diagrams will assist you with daily use of your F&P myAirvo device and its consumables.

## LANGUAGES

PAGE	LANGUAGE	PAGE	LANGUAGE
1	English (en)	43	Italiano (it)
4	(ar) العربية	46	日本語 (ja)
7	简体中文 (zh)	49	Bahasa Melayu (ms)
10	繁體中文 (zht)	52	Norsk (no)
13	Česky (cs)	55	Polski (pl)
16	Dansk (da)	58	Português (pt)
19	Nederlands (nl)	61	Português (Brasil) (ptbr)
22	Suomi (fi)	64	Română (ro)
25	Français (fr)	67	Español (es)
28	Français (frca)	70	Svenska (sv)
31	Deutsch (de)	73	ภาษาไทย (th)
34	Ελληνικά (el)	76	Türkçe (tr)
37	(he) עברית	79	Tiếng Việt (vi)
40	Bahasa Indonesia (id)		

**i****iii****iv**

 Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: [info@fphcare.co.nz](mailto:info@fphcare.co.nz) Web: [www.fphcare.com](http://www.fphcare.com)

 **Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号电话: +86 20 32053486 **Colombia (CO)** Tel: +57 3142852934 **France (FR)**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: [c.s@fphcare.fr](mailto:c.s@fphcare.fr) **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Indonesia (ID)** Tel: +62 21 509 13476 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Sweden (SE)** Fisher & Paykel Healthcare AB, Svetsarvägen 15, 2tr, 17141 Solna, Sweden Tel: +46 8 564 76 680 **Switzerland (CH)**  Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

 **Austria (AT)** Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 45 25 70 85 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 910 38 81 18

# myAirvo Optiflow+ Duet Nasal Cannula

**REF** MYDUETSMALL  
myAirvo Small Duet Nasal Cannula 2-Pack

**REF** MYDUETMEDIUM  
myAirvo Medium Duet Nasal Cannula 2-Pack

**REF** MYDUETLARGE  
myAirvo Large Duet Nasal Cannula 2-Pack

## INDICATIONS FOR USE

The F&P Optiflow+ Duet is a nasal cannula patient interface for use with specified respiratory gas humidifiers to treat spontaneously breathing pediatric (10 years and older) and adult patients who would benefit from receiving high-flow warmed and humidified respiratory gases to the upper airway.

This device is designed to be used in a long-term care facility by appropriately qualified healthcare professionals, or in the home by lay users operating the device as prescribed by a healthcare professional.

## SPECIFICATIONS

- myAirvo series humidifiers with MYAIRVOKITI/MYAIRSPIRAL/900PT560/900PT560E Series tubes or kits
- Flow range: MYDUETSMALL 10–60\* L/min  
MYDUETMEDIUM 10–60\* L/min  
MYDUETLARGE 10–60\* L/min

\*Maximum flow rate depends on cannula size and humidifier model. Do not adjust flow rates without first consulting your healthcare provider. Refer to your myAirvo Series User Manual for more information.

*Note: Flow rates are expressed in body temperature, pressure, saturated (BTPS).*

### Additional information:

- This product is intended to be used for a maximum of 30 days provided that daily and weekly care instructions are followed.
- Replace the interface every 30 days, or earlier if damage or deterioration occurs.
- Ambient operating temperature: 18 °C to 28 °C (64 °F to 82 °F)
- Permissible transport and storage conditions: -10 °C to 50 °C (14 °F to 122 °F)
- Refer to the myAirvo Series User Manual for permissible operating conditions.

## SET-UP INSTRUCTIONS

### Step 1

Select the appropriate size nasal cannula. The prongs must not create a seal in the nares. A clear gap must be visible around each prong.

### Step 2

Adjust the headstrap to fit. Do not over-tighten.

### Step 3

- Ensure the headstrap clip is attached to the headstrap, to prevent the cannula from being pulled out of the nares.
- The cannula can become unattached if not used with the headstrap clip.

### Step 4

Attach the breathing tube clip to a secure location (e.g., clothing or bedding) to prevent the cannula from pulling on the ears or off the face.

## ii DURING USE

- As appropriate for the patient's condition, monitor for any disruptions to the patient receiving flow.
- To prevent condensate, remove or minimize the impact of anything that may cool the tubing (e.g., a fan, air conditioning).
- Notice to user:** If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

## iii DAILY CARE

### Step 1

Disconnect the patient interface from the heated breathing tube.

### Step 2

Rinse the interface in drinking-quality water then reconnect to the breathing tube while still in Drying Mode (refer to the myAirvo Series User Manual) to dry the interface.

## iv WEEKLY CARE

### Step 1

Handwash the interface in a tub of warm water with a mild dishwashing detergent. Ensure the interface is visually clean.

### Step 2

Rinse the interface thoroughly in drinking-quality water. Ensure that all soap residue has been removed. Leave the interface to dry away from direct sunlight or heat, then reconnect it to the breathing tube.

## SIDE EFFECTS

- Skin reddening may occur during nasal high flow therapy.

## DISPOSAL

Dispose of products and packaging according to local guidelines. Packaging and User Manual should follow recycling guidelines as applicable. A used product should be disposed of as contaminated waste. In the home, place the product in the general waste collection. In long-term care facilities, follow the applicable waste disposal procedure.

## WARNINGS

- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption of therapy, to avoid serious injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Failure to use the set-up described above can compromise performance and affect patient safety.
- If you are experiencing severe symptoms, contact your healthcare provider directly.

### Exposure to oxygen increases the risk of fire.

#### To avoid risks of burns, patient injury or death:

- Do not use system near a naked flame or any ignition source, including while smoking.
- Only use lotions and/or salves that are labeled as oxygen compatible.
- Do not lubricate connections or tubing.
- Do not leave the interface on bed coverings or cushions if the device is turned on but not in use. Oxygen can make the materials more flammable.

## CAUTIONS

- For single-patient use only. Reuse between patients may result in transmission of infectious substances.
- Do not use if individual interface packaging is not sealed to avoid contamination.
- Do not allow the prong(s) to seal in the nares. Occlusion may result in nasal injury.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Deviations from the recommended methods of daily and weekly care instructions may result in inadequate cleaning which may lead to adverse health effects.

- Do not sterilize or disinfect this product, or exceed the maximum duration of use to avoid degradation of materials and patient harm.
- Using the cannula may cause skin reddening. Monitor for signs of skin damage or irritation and adjust the cannula if necessary.

### Failure to comply with the following may result in loss of therapy:

- Check for condensate regularly. Drain condensate away from patient as required, to prevent blockages or condensate from reaching the patient.
- Do not crush or stretch tube. Do not obstruct or modify any part of the device.
- Do not use if tubing is damaged with holes, tears or kinks.

## SYMBOL DEFINITIONS

	Catalog number
	Manufacturer
	Importer
	Distributor
	Authorized European Representative
	Authorized Swiss Representative
	Authorized UK Representative
	Batch code
	Date and Country of Manufacture CC NZ: New Zealand MX: Mexico
	Use-by date
	Storage and transport temperature range
<b>Rx Only</b>	Prescription only
	Single patient multiple use
	Polyethylene Terephthalate

Polyethylene Terephthalate,  
Low Density Polyethylene

Refer to instructions

Correct FITTING INSTRUCTION  
method

Incorrect FITTING INSTRUCTION  
method

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