# Optiflow™+ Tracheostomy Interface

REF MYOPT9TRACHE


Thank you for choosing F&P myAirvo™

Please read this entire guide before you first use your tracheostomy interface kit. Following the steps in this guide, and referring to the images and diagrams will assist you with daily use of your F&P myAirvo device and its consumables.

## LANGUAGES

PAGE	LANGUAGE	PAGE	LANGUAGE
1	English (en)	14	Türkçe (tr)
2	Italiano (it)	15	العربية (ar)
3	Deutsch (de)	16	עברית (he)
4	Français (fr)	17	繁體中文 (zht)
5	Español (es)	18	Polski (pl)
6	ภาษาไทย (th)	19	Română (ro)
7	Português brasileiro (pt-br)	20	Ελληνικά (el)
8	Português (pt)	21	Slovenčina (sk)
9	Nederlands (nl)	22	Čeština (cs)
10	Dansk (da)	23	한국어 (ko)
11	Suomi (fi)	24	Русский (ru)
12	Norsk (no)	25	日本語 (ja)
13	Svenska (sv)		

Single patient multiple use

CE 0123 Conforms with medical device directive 93/42/EEC

Rx only Prescription only

F&P, Optiflow and myAirvo are trademarks of Fisher & Paykel Healthcare Limited.

For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip)

#### Manufacturer

Fisher & Paykel Healthcare Ltd,  
15 Maurice Paykel Place,  
East Tamaki, Auckland 2013  
PO Box 14 348 Panmure,  
Auckland 1741,  
New Zealand

Tel: +64 9 574 0100  
Fax: +64 9 574 0158  
Email: [info@fphcare.co.nz](mailto:info@fphcare.co.nz)  
Web: [www.fphcare.com](http://www.fphcare.com)

Importer/Distributor

**Australia** (Sponsor)  
Fisher & Paykel Healthcare Pty Ltd,  
19-31 King Street, Nunawading,  
Melbourne, Victoria 3131  
Tel: +61 3 9871 4900  
Fax: +61 3 9871 4998

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 MYOPT9TRACHE

## INDICATIONS FOR USE

Tracheostomy patient interface for delivery of humidified respiratory gases.

For use with the myAirvo Series humidifiers.

## SETUP

- myAirvo Series humidifiers with MYAIRVOKIT1/MYAIRSPIRAL/900PT560/900PT560E series tubes or kits.
- Flow range: 10-60 L/min.

## SYSTEM SPECIFICATIONS

### Tracheostomy tube connection:

15 mm female conical connector

### Spare part:

Sputum Guard Spare (OPT971)

## HOME AND LONG-TERM CARE FACILITY USE

This product is intended to be used for a maximum of 30 days providing daily and weekly cleaning instructions are followed (see the myAirvo Series User Manual).

## SETUP INSTRUCTIONS

### Step 1

- To ensure loading and movement on tracheostomy tube is kept to a minimum, make sure the lanyard is secured properly.
- Lanyard position as shown in 1a is not suitable for patients with damaged skin.

### Step 2

- Connect to the humidification system. Make sure all connections are secure.
- Attach the breathing tube clip to a secure location (e.g. clothing or bedding) to support the interface.

### Step 3

Orientate the sputum guard so that sputum drains downwards. The sputum guard can be removed if required.

### Warnings:

- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption to therapy, to avoid serious injury or death.
- Do not obstruct and or modify any part of the interface to avoid hypoxia from loss of therapy.
- Do not use near a naked flame, to avoid fires.
- Do not crush or stretch tube, to prevent loss of therapy.
- For single patient use only. Reuse may result in transmission of infectious substances. Attempting to reprocess will result in degradation of materials and render the product defective.
- The lanyard poses a potential risk of strangulation to vulnerable patients (e.g. children).
- The lanyard is designed to minimize loading and movement of the tracheostomy tube. Secure the lanyard properly to avoid accidental decannulation or airway damage.
- The sputum guard is intended to be detachable which may pose a choking risk.

### Cautions:

- Failure to use the set-up described above can compromise performance and affect patient safety.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Check for condensate regularly. Drain as required.
- Do not use if individual interface packaging is not sealed.

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