

Nasal High Flow Therapy: Neonates

Clinical Paper Summaries



Table of Contents



Reviews: Overview of Nasal High Flow Literature

R1	Mechanisms of action (Dysart et al., Manley et al., Saslow et al., Sivieri et al., Wilkinson et al., Lavizzari et al.)	2
R2	Post-extubation (Wilkinson et al., Colin et al., Yoder et al., Manley et al., Roberts et al.)	3
R3	Primary support (Roberts et al., Lavizzari et al.)	5
R4	Consensus on implementing NHF (Roehr et al., Yoder et al.)	7



Summaries: Key Randomized Controlled Trials

NHF POST-EXTUBATION

S1	Collins et al. 2013	10
S2	Manley et al. 2013	12
S3	Yoder et al. 2013	14

NHF AS PRIMARY SUPPORT
















S4	Lavizzari et al. 2016	16
S5	Roberts et al. 2016	18



All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom.



Neonatal Key References

1.  Collins CL, Holberton JR, Barfield C, Davis PG. A randomised control trial to compare heated humidified high flow nasal cannulae with nasal continuous positive airway pressure post-extubation. *Journal of Pediatrics*. 2013; 162(5): 949-954
2.  Dysart K, Miller TL, Wolfson MR, et al. Research in high flow therapy: mechanisms of action. *Respiratory Medicine*. 2009; 103(10): 1400-1405.
3.  Lavizzari A, Veneroni C, Colnaghi M, et al. Respiratory mechanics during NCPAP and HHHFNC at equal distending pressures. *Archives of Disease in Childhood-Fetal and Neonatal Edition*. 2014; 99(4): F315-320
4.  Lavizzari A, Colnaghi M, Ciuffini F, et al. Heated, Humidified High-Flow Nasal Cannula vs Nasal Continuous Positive Airway Pressure for Respiratory Distress Syndrome of Prematurity: A Randomized Clinical Noninferiority Trial. *JAMA Pediatrics*. 2016
5.  Manley BJ, Owen LS. High-flow nasal cannula: Mechanisms, evidence and recommendations. *Seminars in Fetal and Neonatal Medicine*. 2016; 21(3): 139-145
6.  Manley BJ, Owen LS, Doyle LW, et al. High-flow nasal cannulae in very preterm infants after extubation. *New England Journal of Medicine*. 2013; 369 (15): e1425-1433.
7.  Roberts CT, Hodgson KA. Nasal high flow treatment in preterm infants. *Maternal Health, Neonatology and Perinatology*. 2017 3(1):15.
8.  Roberts CT, Owen LS, Manley BJ, et al. Nasal high-flow therapy for primary respiratory support in preterm infants. *New England Journal of Medicine*. 2016;375(12):1142-51.
9.  Roehr CC, Yoder BA, Davis PG, et al. Evidence Support and Guidelines for Using Heated, Humidified, High-Flow Nasal Cannulae in Neonatology: Oxford Nasal High-Flow Therapy Meeting, 2015. *Clinics in Perinatology*. 2016; 43(4):693-705.
10.  Saslow JG, Aghai ZH, Nakhla TA, et al. Work of breathing using high-flow nasal cannula in preterm infants. *Journal of Perinatology*. 2006; 26(8): 476-480.
11.  Sivieri EM, Gerdes JS, Abbasi S. Effect of HFNC flow rate, cannula size, and nares diameter on generated airway pressure: an in vitro study. *Pediatric Pulmonology*. 2013; 48(5): 506-514
12.  Wilkinson DJ, Andersen CC, Smith K, et al. Pharyngeal pressure with high-flow nasal cannulae in premature infants. *Journal of Perinatology* 2008; 28(1): 42-47.
13.  Wilkinson D, Andersen C, O'Donnell CP, et al. High flow nasal cannula for respiratory support in preterm infants. *The Cochrane Library*. 2016
14.  Yoder BA, Stoddard RA, Li M, King J, et al. Heated, humidified high-flow nasal cannula versus nasal CPAP for respiratory support in neonates. *Pediatrics*. 2013; 131(5):1482-1490
15.  Yoder BA, Manley B, Collins C, et al. Consensus approach to nasal high-flow therapy in neonates. *Journal of Perinatology*. 2017.



All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom.



Dysart et al. 2009; Lavizzari et al. 2014; Manley et al. 2016; Saslow et al. 2006; Sivieri et al. 2013; Wilkinson et al. 2008

INTRODUCTION

- Nasal high flow (NHF) therapy is increasingly being used as an alternative form of respiratory support in neonatal intensive care units (NICU) in preterm infants with respiratory illnesses (Manley et al. 2016; Dysart et al. 2016). Continuous positive airway pressure (CPAP) has previously been regarded as the gold-standard for noninvasive support in this population. NHF is a form of noninvasive respiratory support that uses conditioned (heated and humidified), gas flow applied via small nasal cannulae. The following mechanisms of action are associated with NHF when used at effective flow rates

MECHANISMS OF ACTION

Washout of dead space

- NHF allows for a very effective flushing of the nasopharyngeal dead space, resulting in alveolar ventilation being a higher fraction of minute ventilation (the volume of gas inhaled or exhaled from the lungs per minute) (Dysart et al. 2009, Manley et al. 2013)
- The flow of gas during NHF flushes CO₂ from the nasopharynx, facilitating clearance of CO₂ and improving oxygenation by creating a reservoir of fresh gas in the airway

Improved mechanics

- NHF uses gas that is humidified (100% relative humidity) and adequately heated (maintained at 34-37° C) which reduces conductance and pulmonary compliance compared with the use of dry, cooler gas (Dysart et al. 2009)
- The use of gas that is already warmed and humidified presumably reduces the metabolic energy that might have been expended had the inspired gases been cooler or drier (Dysart et al. 2009)

Reduction in the work of breathing

- The distensibility of the nasopharynx results in an increased resistance on inspiration compared with expiration
- NHF provides gas flow rates that meet or exceeds the peak inspiratory flow, thereby leading to a reduction of inspiratory resistance. This translates to a decrease in the resistive work of breathing (Dysart et al. 2009)
- The work of breathing with NHF (air flow rates of 3-5 L/min) was comparable to that with CPAP set to 6 cm H₂O, according to a study in 18 preterm infants (Saslow et al. 2006)

Provision of distending airway pressure

- Providing distending pressure to the lungs improves ventilatory mechanics by optimizing lung compliance and assists with gas exchange by maintaining the patency of the alveoli (Dysart et al. 2009)
- Initial data from premature infants showed that increasing NHF flow rates resulted in increased pharyngeal pressures (Wilkinson et al. 2008)
 - Pressures generated in the nasopharynx with NHF were within the range of commonly used CPAP pressures.
- Appropriate prong-to-nares ratio is essential to achieve adequate, but not excessive, positive pressure support (Sivieri et al. 2013)
- When similar end-expiratory pressures were applied, there was no difference in breathing pattern, gas exchange, lung mechanics or work of breathing between NHF and CPAP, according to a study in 20 preterm infants with mild to moderate respiratory distress (Lavizzari et al. 2014)

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Collins et al. 2013; Manley et al. 2013; Roberts et al. 2017; Wilkinson et al. 2016; Yoder et al. 2013

INTRODUCTION

- Preterm infants are prone to respiratory failure and often require mechanical ventilation through an endotracheal tube after birth. Once these infants are extubated, traditionally Continuous Positive Airway Pressure (CPAP) has been used as the standard of care for noninvasive respiratory support. However, there is increasing evidence that Nasal High Flow (NHF) therapy is a suitable alternative form of noninvasive respiratory support in infants ≥ 28 weeks GA.
- To date there have been a number of Randomized Controlled Trials (RCTs) comparing the use of NHF to CPAP for post-extubation support in preterm infants. Of particular interest are the following three RCTs that have also been included in a Cochrane Review meta-analysis by Wilkinson et al (2016.):

	Manley et al. 2013	Collins et al. 2013	Yoder et al. 2013
Patient population	303 infants <32 weeks GA	132 infants <32 weeks GA	432* term and preterm Infants >28 GA or >1000 g
Study design	Multicentre RCT Non-inferiority trial (20% margin)	Single-centre RCT	Multicentre RCT
NHF flow rate	5-8 L/min	4-8 L/min	3-8 L/min
CPAP cmH₂O	5-8 cmH ₂ O	4-8 cmH ₂ O	5-8 cmH ₂ O
Primary outcome	Treatment failure within 7 days	Extubation failure within 7 days	Extubation failure within 72 hours

*226 preterm infants were randomized to post-extubation treatment. 125 preterm infants were randomized to primary treatment. Subgroup analysis are unavailable.

RESULTS

1. No difference to rate of treatment failure

Treatment failure rates were not significantly different between neonates randomized to NHF or CPAP across all three RCTs. However, Collins et al. (2013) found that across both treatment groups, extubation failure was significantly greater in infants with a GA <28 weeks (P<0.001)

Extubation Failure within 72 hours, %			
	NHF (n=212)	CPAP (n=220)	P-value
Yoder et al. (2013)	10.8	8.2	NS
Extubation Failure within 7 days, %			
	NHF (n=67)	CPAP (n=65)	P-value
Collins et al. (2013)	22	34	NS
Treatment Failure within 7 days, %			
	NHF (n=152)	CPAP (n=151)	P-value
Manley et al. (2013)	34.2	25.8	NS

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom.



2. No difference to rate of reintubation

There was no significant difference between NHF and CPAP with respect to reintubation within 7 days in all three studies conducted by Collins et al. (2013) Manley et al. (2013) and Yoder et al. (2013)

Reintubation within 7 days, %			
	NHF	CPAP	P-value
Collins et al. (2013)	(n=67)	(n=65)	
	10	12	NS
Manley et al. (2013)	NHF (n=152)	CPAP (n=151)	
	17.8	25.2	NS
Yoder et al. (2013)	NHF (n=212)	CPAP (n=220)	
	11	10	NS

3. Significant reduction in rates of nasal trauma

NHF was associated with significantly less nasal trauma than CPAP across all three studies.

Nasal Trauma, %			
	NHF	CPAP	P-value
Yoder et al. (2013)	(n=212)	(n=220)	
	9	16	0.047
Manley et al. (2013)	NHF (n=152)	CPAP (n=151)	
	39.5	54.3	0.01
Collins et al. (2013) [Nasal trauma score]	NHF (n=67)	CPAP (n=65)	
	[3.1]	[11.8]	<0.01

4. No difference in rates of other adverse outcomes i.e. death, pneumothorax, bronchopulmonary dysplasia (BPD)

Apart from significant differences in reported nasal trauma, there was no difference in the adverse events profile with NHF compared to CPAP.

CONCLUSION & KEY POINTS

- The evidence indicates that for preterm infants aged ≥ 28 weeks' GA, NHF and CPAP are associated with similar efficacy and safety.
- The use of NHF resulted in significantly lower rates of nasal trauma, with no additional risk of adverse events, compared to CPAP.
- For infants < 28 weeks, there are limited data and insufficient evidence to change the current practice of using CPAP for post-extubation respiratory support.

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Lavizzari et al. 2016; Roberts et al. 2016

BACKGROUND

- Traditionally, invasive ventilation and noninvasive ventilation in the form of CPAP have been the standard treatment for premature infants requiring respiratory support after birth and initial stabilization. There is emerging evidence comparing the efficacy of CPAP with that of Nasal High Flow (NHF) therapy.
- Two non-inferiority, Randomized Controlled Trials (RCTs) have investigated the use of NHF compared with CPAP for the primary treatment of preterm neonates.

	Lavizzari et al. 2016	Roberts et al. 2016
Patient population	316 infants ≥29 weeks' GA	564 infants ≥28 weeks GA
Study design	Single-center study in Italy	Multi center study in Australia and Norway
Non inferiority margin	10%	10%
NHF flow rate	4-6 L/min	6-8 L/min
CPAP cmH₂O	4-6 cmH ₂ O	6-8 cmH ₂ O
Type of CPAP	SIPAP	SIPAP, Bubble, Ventilator
Primary outcome	Intubation and mechanical ventilation	Failure of initial therapy
Surfactant use	NHF 44.3% CPAP 46.2%	NHF 14.4% CPAP 10.5%
	High usage/ low threshold	Low usage/threshold not specified
Rescue therapy	“Rescue CPAP” was used for infants in the NHF group at clinician discretion	“Rescue CPAP” was used for infants in the NHF group who met the failure criteria

A note on recruitment for Roberts et al.

Analysis of the primary outcome data from the first 515 recruited infants indicated a significant difference between groups (P<0.001). The trial was stopped since continued recruitment was extremely unlikely to indicate the non-inferiority of NHF to CPAP in terms of primary outcome.

RESULTS

Rate of treatment failure

Roberts et al. reported significantly higher rates of initial treatment failure in the NHF group, while Lavizzari et al. found that NHF was non-inferior to CPAP with regards to intubation and mechanical ventilation

Treatment Failure within 72 hours, %	NHF	CPAP	P-value	Risk difference
Roberts et al. (2016)	25.5	13.3	<0.001	12.3% (95% CI, 5.8 to 18.7)
Lavizzari et al. (2016)	Not reported			

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



The role of “Rescue CPAP” on the rate of intubation

The use of “Rescue CPAP” in both studies resulted in no significant difference to the rates of intubation between therapies. Lower rates of intubation were found among older preterm infants (>32 weeks’ GA).

Intubation rates within 72 hours, %	NHF	CPAP	P-value	Risk difference
Roberts et al. (2016)	15.5	11.5	0.17	3.9% (95% CI, -1.7 to 9.6)
Lavizzari et al. (2016)	10.8	9.5	0.71	1.3% (95% CI, -6.0 to 8.6)

Nasal trauma

Nasal trauma was significantly lower with NHF than with CPAP (Roberts et al.)

Rate of adverse events

- There were no significant differences in the secondary outcomes between the NHF and CPAP groups, including the duration of respiratory support, need for surfactant, duration of hospitalization, enteral feeding, weight or adverse events
- Both studies report surfactant use and neither observed a difference between CPAP and NHF. However, the use of surfactant appears to differ between studies.
- In Lavizzari et al., >44% of infants were briefly intubated for surfactant delivery. These infants continued in the study. In Roberts et al., the use of surfactant was much lower and it appears that all infants who received surfactant were classified as treatment failure.
- This difference between study populations is important to note when comparing outcomes between studies.

CONCLUSION & KEY POINTS

- Compared to NHF, CPAP may be associated with a lower rate of initial treatment failure when used for primary treatment.
- However, the need to intubate does not appear to differ significantly between the two groups. This suggests that NHF may be a feasible option when CPAP is available as rescue therapy in preterm infants, especially for those more mature preterm infants.
- While CPAP continues to be the gold standard for primary treatment of preterm infants requiring noninvasive respiratory support, further research into the use of NHF in this neonate group is still warranted



Roehr et al. 2016; Yoder et al. 2017

AIM

- To provide an overview of clinical data regarding Nasal High Flow (NHF) therapy in neonates
- To establish consensus on the mechanisms of action and the clinical indications for the use of NHF compared with Continuous Positive Airway Pressure (CPAP) in neonates

BACKGROUND

- CPAP is the current gold standard for providing noninvasive respiratory support to neonates with respiratory insufficiency
- NHF is an alternative form of respiratory support that uses conditioned (heated and humidified), high flows of gas applied via small nasal cannulae
 - NHF is typically defined as flows of 2L/min or higher
 - The high flows of gas leads to the effective washout of carbon dioxide, a reduction of inspiratory resistance in the upper airways and a reduction in the work of breathing (WOB)

METHODS

NHF therapy expert meeting (Roehr, et al.)

- A group of 24 international physicians with a particular interest in neonatal and pediatric NHF gathered in Oxford, United Kingdom in June 2015 to discuss the present state of research into the respiratory management of newborn infants and young infants
- Respiratory physiologists, epidemiologists and authors of publications and reviews on NHF also attended the meeting
- A summary of discussions from the meeting and treatment recommendations was published, based on the latest available evidence from RCTs and the collective experience of the attendees

Consensus questionnaire (Yoder, et al.)

- Seven international NHF clinical researchers were questioned regarding their approaches to initiation, escalation, weaning and discontinuing NHF
- A prospective, modified Delphi approach was taken to construct a series of tables that answered questions related to specific aspects of NHF practice
- Completed tables were reviewed independently by each investigator, the results were discussed and areas of consensus were determined



CONSENSUS AGREEMENT AND MEETING RECOMMENDATIONS

Humidification

There was strong consensus that NHF should be humidified (100% relative humidity) and adequately heated (maintained at 34-37° C)

Flow rate

- The maximum flow rate should be determined by the manufacturer's design or approval
 - There was general agreement that the initial gas flow should start at 4-6 L/min for preterm and term infants

Nares occlusion

- When selecting an NHF interface, there should be allowance for the generous egress of gas around the cannula tip and from the nares by ensuring that the prong diameter is approximately half that of the nostril

Appropriate population for NHF

- Current evidence supports the application of NHF for:
 - Post-extubation support of neonates \geq 28 weeks' in lieu of nasal CPAP
 - As an alternative to CPAP for stable infants who continue to require respiratory support above the level of standard Oxygen therapy
 - There was general agreement, but not consensus, that NHF can be used as the primary mode of support for neonatal respiratory distress at the clinician's discretion (based on the neonate's GA and the level of oxygen support required)

Contraindications for NHF

- Infants with signs of severe respiratory distress ($FiO_2 > 0.7$), severe apnea, significant active air leak, or craniofacial or airway anomalies should not be considered for NHF

Feeding

- Further study is required regarding the approaches to safe oral feeding in infants on NHF
- There is very limited evidence regarding the use of nasogastric or orogastric feeding tubes in infants on NHF

Post-extubation

- Post-extubation (support initiated after a period of intubation with resolving respiratory distress) evidence comes from randomized, clinical trials involving more than 1100 preterm infants and meta-analyses
- NHF was as effective and as safe as nasal CPAP
 - Failure rates were similar and there was no increase in adverse events, air leaks, or durations of oxygen use or hospital stay with NHF
 - Comfort scores of neonates were similar, but care of the infant was easier with NHF
 - NHF was associated with reduced rates of nasal trauma and potentially a reduced incidence of pneumothorax in preterm infants
 - The authors note that the majority of infants enrolled in trials to date have been >28 weeks GA

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom.



- The approach to either escalation or weaning NHF should be based on signs of distress or increase in WOB, as well as the FiO_2 needed to maintain targeted SpO_2
 - Typically, weaning is considered in infants who have been stable for 12-24 hours
 - FiO_2 should be weaned first to <0.3 , followed by flow rates in decrements of 0.5 to 1 L/min and every 12 to 24 h as tolerated, guided by the WOB of the infant
 - There was no consensus as to when to stop NHF, with discontinuation occurring at flow rates between 1-4 L/min
 - Escalation of gas flow is recommended for increasing FiO_2 and WOB, to a maximum recommended flow rate of 8 L/min
 - Alternative approaches to noninvasive support should be considered in neonates requiring elevated FiO_2 (typically >0.40) or with increased WOB or distress

Primary treatment

- There was general agreement for NHF as the primary mode of support for neonatal respiratory distress, but GA and the requirement for oxygen support must be considered
 - Only limited evidence is available from clinical trials for this indication, including two randomized clinical trials which suggested no difference in efficacy between NHF and either nasal ventilation or nasal CPAP (Roberts et al, 2016; Lavizzari et al. 2016) and one randomized trial which found higher treatment failure with NHF than with CPAP (Roberts et al. 2016)

AREAS FOR FUTURE RESEARCH

- Further studies are required to:
 - Determine the safety and efficacy of NHF in various groups of infants including extremely preterm infants, and infants with other neonatal lung disorders (i.e., meconium aspiration, diaphragmatic hernia)
 - Compare different initial NHF gas flow rates
 - Evaluate different approaches to discontinuing NHF
 - Compare different NHF devices and cannula interfaces (i.e., dual versus single nasal prong cannulas)
 - Analyze the economics of NHF versus other noninvasive modes
 - Assess the efficacy and safety of NHF for initial delivery room stabilization and during neonatal transport
 - Assess the efficacy and safety of NHF in resource-limited countries
 - Develop and assess the use of NHF for post-discharge noninvasive support
 - Assess approaches to oral feeding on NHF



A randomized controlled trial to compare heated humidified high-flow nasal cannulae with nasal continuous positive airway pressure postextubation in premature infants

AIM

- To determine whether Nasal High Flow (NHF) therapy increases the rate of successful extubation in premature infants following endotracheal positive pressure ventilation, compared with nasal Continuous Positive Airway Pressure (CPAP).

METHOD

Patient group

- 132 infants were randomized to either NHF (n=67) or nCPAP (n=65) and <32 weeks' GA, stratified into two groups (<28 weeks vs ≥28 weeks' GA)

Study design

- Single center RCT

Primary outcome

- Extubation failure within 7 days, defined as one or more of the following:
Apnea (respiratory pause >20 sec), >6 episodes in 6 hours or an episode requiring intermittent positive pressure ventilation; acidosis (pH <7.25 and PCO₂ >66 mmHg); sustained increase in FiO₂ of >15% from extubation
 - Patients were re-intubated at the treating physician's discretion

Secondary outcomes

- Nasal trauma, duration of respiratory support, supplemental oxygen requirement, bronchopulmonary dysplasia
 - Nasal trauma was assessed at the internal and external nares, philtrum and septum, using the sum of thrice-daily nasal trauma score recordings from 0 (normal) to 3 (skin tear)

Treatment regimen

- NHF (Vapotherm) at a starting flow of 8L/min. Flow rate was weaned to a minimum of 4L/min OR
- nCPAP (Hudson Respiratory Care) at a starting positive end-expiratory pressure (PEEP) of 8cm H₂O or 7cm H₂O for FiO₂ values of >0.3 and <0.3, respectively; PEEP was weaned to a minimum of 5cm H₂O

RESULTS

Primary outcome

- NHF and nCPAP were associated with similar extubation failure rates at 7 days (see table)
 - Stratification of patients by GA (<28 weeks vs. ≥28 weeks) also showed no significant difference between NHF and nCPAP therapy in extubation failure rates at 7 days
 - Overall extubation failure rates were higher among infants born at <28 weeks' GA than those born at 28–32 weeks (44% vs. 15%; P<0.001)



Secondary outcomes

- NHF was associated with significantly reduced nasal trauma compared with nCPAP
 - At 7 days after extubation, 20% of the infants randomized to nCPAP switched to NHF due to nasal trauma
- No significant differences were seen in bronchopulmonary dysplasia (BPD) rates or the duration of supplemental Oxygen and respiratory support:

Variable	HFNC (n=67)	nCPAP (n=65)	P-value
Extubation failure at 7 days (primary outcome), % pts.	22	34	NS
Apnoea, % pts.	21	26	NS
Acidosis, % pts.	0	5	NS
FiO ₂ increase >15%, % pts.	10	18	NS
Nasal trauma score first week, mean (SD)	3.1 (7.2)	11.8 (10.7)	<0.01
BPD at 36 weeks' gestation, % pts.	36	43	NS
Respiratory support, mean completed weeks (SD)	33.5 (2.88)	34.3 (3.51)	NS
Supplemental oxygen, mean completed weeks (SD)	36.9 (2.54)	38.0 (3.26)	0.06

NS, not significant; pts., patients; FiO₂, fraction of inspired oxygen; SD, standard deviation

CONCLUSION & KEY POINTS

NHF is associated with a similar rate of extubation failure at CPAP in premature infants \geq 28 weeks' GA

- Infants <28 weeks' GA experienced a higher rate of extubation failure
- NHF is associated with significantly less nasal trauma compared with CPAP, with no additional risk of adverse events

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom.



High-flow nasal cannulae in very preterm infants after extubation

AIM

- To compare the efficacy and safety of Nasal High Flow (NHF) therapy and nasal continuous positive airway pressure (nCPAP) for the noninvasive respiratory support of very preterm infants following extubation

METHOD

Patient group

- 303 infants <32 weeks' GA

Study design

- Multicenter, randomized, non-inferiority trial with a non-inferiority margin of 20%

Primary outcome

- Treatment failure within 7 days of extubation

Secondary outcomes

- Reintubation during the 7 day primary endpoint period defined as; requirement for supplemental oxygen at a GA of 36 weeks; pneumothorax following study entry; total number of days of any respiratory support following study entry; duration of oxygen supplementation following study entry; length of hospital admission
- The incidence, cause, and need for a change of treatment as a result of nasal trauma were also documented

Treatment regimen

- NHF (Fisher & Paykel Healthcare) at a starting flow rate of 5–6 L/min
 - Infants who met the failure criteria were “rescued” with nCPAP; OR
- nCPAP (Fisher & Paykel Healthcare or Hudson Respiratory Care) at a starting pressure of 7cm H₂O. A mechanical ventilator or a “bubble” system was used to generate nCPAP
 - Infants who met the failure criteria were reintubated

RESULTS

Primary outcome

- NHF was non-inferior to nCPAP with regard to the primary outcome, treatment failure within 7 days of extubation (risk difference 8.4%; 95% confidence interval [CI] -1.9, 18.7; see table)
- There were no significant differences between NHF and nCPAP groups in terms of the reasons for treatment failure
 - 48% of infants that met the failure criteria on NHF were “rescued” with nCPAP and did not require reintubation

Secondary outcomes

- There was no significant difference between the NHF and nCPAP groups for the secondary outcome of reintubation within 7 days of extubation (risk difference -7.4%; 95% CI -16.6, 1.8; P=0.12)
- NHF was associated with significantly lower rates of nasal trauma compared with nCPAP
- No significant differences were seen between the two treatment groups for other secondary outcomes and serious adverse events

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



	HFNC (n=152)	nCPAP (n=151)	P-value
Primary outcome (% pts)			
Treatment failure within 7 days of extubation	34.2	25.8	NS
Secondary outcomes			
Reintubation within 7 days of extubation (% infants)	17.8	25.2	NS
Oxygen supplementation at GA of 36 weeks (% pts)	30.9	34.4	NS
Pneumothorax following study entry (% pts)	0.7	2.6	NS
Nasal trauma (% infants)			
Any documented	39.5	54.3	0.01
Leading to change of treatment	5.3	17.9	0.001
Caused by randomized treatment	19.1	53.0	<0.001
Serious adverse events (% pts)			
Pneumothorax during randomized treatment	0.0	0.7	NS
Death before discharge	3.3	4.0	NS
pts, patients; IQR, interquartile range; NS, not significant.			

CONCLUSION & KEY POINTS

- NHF appears to have a similar efficacy and safety to CPAP when used post-extubation in infants ≥ 28 weeks' GA
- NHF is associated with significantly less nasal trauma compared with nCPAP therapy, with no additional risk of adverse events

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Heated, humidified high-flow nasal cannula versus nasal CPAP for respiratory support in neonates

AIM

- To determine the safety and efficacy of Nasal High Flow (NHF) therapy versus nasal Continuous Positive Airway Pressure (nCPAP) when used as noninvasive respiratory support in neonates with respiratory dysfunction

METHOD

Patient group

- 432 infants between 28-42 weeks' GA, with a birth weight ≥ 1000 g

Study design

- Multicenter RCT
- Infants were stratified according to birth weight (1000-1999g and ≥ 2000 g) and age at randomization (≤ 7 days vs > 7 days)

Primary outcome

- Failure of study support mode (defined as the need for intubation within the first 72 hours of treatment)

Secondary outcomes

- Rates of bronchopulmonary dysplasia (BPD), total ventilator days, days on supplemental oxygen, need for delayed intubation

Treatment regimen

- NHF (Fisher & Paykel Healthcare, Vapotherm or Hudson Respiratory Care devices) at a starting flow determined by patient weight; 1000-1999g infants received 3 L/min, 2000-2999g infants received 4 L/min, and ≥ 3000 g infants received 5 L/min OR
- nCPAP (bubble and Infant Flow nCPAP System, CareFusion) at a starting pressure of 5-6 cm H₂O
- Oxygenation and ventilation targets were 85-98% for oxygen saturation (SpO₂), and 40-65 mmHg for partial pressure of carbon dioxide (PCO₂)

RESULTS

Primary outcome

- There was no significant difference between groups in the primary outcome of failure of support within 72 hours (see table)
 - There were no significant differences in the reasons for early failure and intubation between groups: Increasing respiratory distress (83% in the NHF and nCPAP groups), increased FiO₂ (39% in the NHF group and 50% in the nCPAP group) and severe apnea (22% in the NHF group and 11% in the CPAP group)

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Secondary outcomes

- Infants managed with nCPAP had fewer days of any positive pressure support, and a shorter duration of study mode support compared to infants managed with heated humidified NHF (see table)
- At 7 days post-study entry, significantly more infants receiving NHF remained on the treatment than on nCPAP
 - There were no differences in the other endpoints measured
- Adverse events were similar between groups, and rates of failure were also similar between different devices used

Variable	Heated humidified NHF	nCPAP	P-value
Need for reintubation, % pts	15	11	NS
Remaining on therapy 7 days post-study entry, % pts	23	9	<0.001
BPD, % pts	20	16	NS
Median days on positive pressure support	6	4	<0.001
Median days on study mode	4	2	<0.001
Median days on supplemental O ₂	10	8	NS

BPD, bronchopulmonary dysplasia, NS, not significant; O₂, oxygen; pts, patients

CONCLUSIONS & KEY POINTS

- NHF appears to have a similar efficacy and safety to CPAP when used post-extubation in infants ≥ 28 weeks' GA
- There were no significant differences to the rate of adverse events between groups, and rates of failure were also similar between different devices used

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Heated, humidified high-flow nasal cannula vs nasal continuous positive airway pressure for respiratory distress syndrome of prematurity: a randomized clinical non-inferiority trial

AIM

- This trial investigated whether primary treatment with Nasal High Flow (NHF) therapy is non-inferior to Continuous Positive Airway Pressure (CPAP) in preterm infants with mild to moderate respiratory distress syndrome (RDS)

METHOD

Patient group

- Preterm infants 29 to <37 weeks' GA (GA) requiring noninvasive respiratory support for mild to moderate RDS

Study design

- Non-inferiority RCT with a margin of non-inferiority of 10% (determined by the absolute risk difference in the primary outcome)

Primary outcome

- Requirement for mechanical ventilation (MV) within 72 hours from commencing noninvasive respiratory support
 - Criteria for MV:
 - Persistent fraction of inspired oxygen (FiO_2) >0.4 to a target oxygen saturation (SpO_2) of 86-93% after surfactant
 - Severe apnea (>4 apnea episodes/hour or 2 apnea episodes/hour requiring positive pressure ventilation [PPV])
 - Acidosis (persistent partial pressure of carbon dioxide [$PaCO_2$] >70 mmHG, pH <7.2)

Treatment regimen

- NHF (Precision Flow, Vapotherm) starting flow 4-6 L/min OR
- nCPAP (SiPAP; Viasys Healthcare) starting pressure 4-6 cmH₂O
 - Switched from nCPAP to bilevel nCPAP in the following scenarios
- >4 episodes of apnea/hour
- >2 episodes/hour requiring PPV
- Evidence of increased work of breathing
- Bilevel CPAP started at a rate of 30 breaths/minute, inspiratory time 0.7-1 second, mean airway pressure of 6-8 cmH₂O

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



RESULTS

Intention to treat group

- 316 infants enrolled between January 2012 and June 2014 randomized to receive either
 - NHF (n = 158, mean GA 33.1 weeks, 52.5% female) OR
 - nCPAP (n = 158, mean GA 33.0 week, 47.5% female)

Primary outcome

- Requirement for MV within 72 hours from commencing respiratory support was similar in both treatment groups (see table)
 - 10.8% with NHF and 9.5% with nCPAP. NHF was non-inferior to nCPAP in this patient group

Secondary outcome

- No significant differences in secondary outcomes between treatment groups
 - Includes duration of respiratory support, need for surfactant, adverse events, the duration of hospitalization, full enteral feeding, weight or exclusive breast feeding at discharge

Outcome	NHF (n =158)	nCPAP (n = 158)	95% CI of risk difference or difference in medians	P value
Primary outcome				
MV within 72 hours, n (%)	17 (10.8)	15 (9.5)	-6.0 to 8.6	0.71
Median age at start of MV (IQR), h	27.0 (8.0-36.0)	7.0 (3.0-19.0)	-24.5-0.00	0.06
Median duration of MV (IQR), d	3.2 (1.2-5.0)	3.0 (1.2-6.0)	-1.25-2.25	0.72
Secondary outcomes				
Median duration of respiratory support (IQR), d	4.0 (2.0-6.0)	4.0 (2.0-7.0)	-1.0-0.5	0.45
Surfactant requirement, n (%)	70 (44.3)	73 (46.2)	-9.8-13.5	0.73
Any adverse event, n (%)	28 (17.7)	28 (17.7)	-9.0-9.0	>0.99
HFNC, heat humidified high flow nasal cannula; IQR, interquartile range; nCPAP, nasal continuous positive airway pressure (includes bilevel nCPAP)				

CONCLUSIONS & KEY POINTS

- NHF was non-inferior to CPAP as primary respiratory support for preterm infants 29 to <37 weeks' GA with mild to moderate RDS
- NHF and nCPAP showed similar efficacy and safety in this patient group

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Nasal high-flow therapy for primary respiratory support in preterm infants

AIM

- To compare the efficacy of nasal high flow (NHF) therapy and nasal continuous positive airway pressure (nCPAP) as primary respiratory support for preterm infants with respiratory distress

METHOD

Patient group

- Preterm infants 28 weeks to <37 weeks' GA, <24 hours old, with early respiratory distress
 - No prior endotracheal ventilation or surfactant replacement therapy
 - No urgent need for intubation or ventilation, no major congenital abnormality or pneumothorax

Study design

- Multi center non-inferiority RCT with a margin of non-inferiority of 10% (determined by the absolute risk difference in the primary outcome)
 - Nine NICUs in Australia and Norway
 - Infants stratified by GA (<32 weeks vs \geq 32 weeks) and study center

Primary outcome

- Treatment failure within 72 hours from randomization
 - Failure was determined if an infant receiving maximal support (NHF 8 L/min or nCPAP 8 cmH₂O) met one or more of the following criteria:
 - Fraction of inspired oxygen (FiO₂) \geq 0.4
 - pH \leq 7.2 plus partial pressure of carbon dioxide (PaCO₂) >60 mmHg (8.0 kPa)
 - \geq 2 episodes of apnea requiring positive pressure ventilation (PPV) within a 24-hour period or \geq 6 episodes requiring any intervention within a 6-hour period
 - Urgent need for intubation and mechanical ventilation

Secondary outcomes

- Reason(s) for treatment failure, use of mechanical ventilation within 72 hours from randomization or at any time during admission, nasal trauma, other complications including complications of prematurity and other measures of neonatal health and respiratory support use, cost of care

Adverse events

- Serious adverse events were defined as death before hospital discharge, pneumothorax or other air leak during treatment

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Treatment regimen

- NHF (Optiflow Junior; Fisher and Paykel Healthcare or Precision Flow; Vapotherm) at a starting flow rate of 6-8 L/min
 - Infants who met the criteria for treatment failure could receive nCPAP as rescue therapy (initiated at 7-8 cmH₂O)
 - Infants with continued treatment failure were intubated and ventilated OR
- nCPAP at a starting pressure 6-8 cmH₂O
 - Infants with treatment failure were intubated and ventilated

RESULTS

Study cessation

- Analysis of primary outcome data from the first 515 recruited infants (278 in the NHF group and 286 in the CPAP group) indicated a statistically significant difference between treatment groups
 - Continued recruitment was considered highly unlikely to show non-inferiority of NHF to nCPAP

Primary outcome

- Treatment failure within 72 hours of randomization was significantly higher with NHF than nCPAP
- 25.5% with NHF versus 13.3% with nCPAP (Risk difference 12.3% points, 95% confidence interval 5.8-18.7; P <0.001)

Outcome	NHF (n = 278)	nCPAP (n = 286)	Risk difference (95% CI)*	P value
	N (%)		Percentage points	
Primary intent-to-treat analysis				
Treatment failure within 72 hours	71 (25.5)	38 (13.3)	12.3 (5.8-18.7)	<0.001
GA <32 weeks'	46 (32.9)	27 (18.1)	14.7 (4.8-24.7)	0.004
GA ≥32 weeks'	25 (18.1)	11 (8.0)	10.1 (2.2-18.0)	0.01
Per-protocol analysis				
Treatment failure within 72 hours	64 (24.2)	36 (12.9)	11.3(4.8-17.8)	<0.001

*Positive values favor the nCPAP group, and negative values favor the NHF group. Apparent discrepancies in some of the risk differences are due to rounding. nCPAP, nasal continuous positive airway pressure; NHF, nasal high flow therapy; wk, weeks.

All clinical paper summaries are independently written by Biowrite Solutions. All rights reserved. No part of this publication may be reproduced by any process in any language without written consent of the copyright holder. Although great care has been taken to ensure that the information in this publication is accurate, none of Biowrite Solutions or Fisher & Paykel Healthcare shall be held responsible or in any way liable for the continued accuracy of the information, or for any errors, omissions or inaccuracies, or for any consequences arising therefrom



Secondary outcomes

- Intubation rate at 72 hours did not differ significantly between treatment groups
 - 15.5% with NHF versus 11.5% with nCPAP (use of nCPAP as rescue therapy may explain this finding)
- Median duration of respiratory support was 1 day longer with NHF (4 vs 3 days; P = 0.005)
- Supplemental oxygen requirement was more common with NHF (78.1% vs 69.6%; P = 0.02)
- The overall rate of adverse events, including rate of death, overall frequency of pneumothorax or other air leak from the lung, was similar between treatment groups
 - However, nasal trauma was significantly lower with NHF than with nCPAP (8.3% vs 18.5%; P < 0.001)

CONCLUSIONS

- Primary support with CPAP resulted in a significantly lower rate of treatment failure compared with NHF in preterm infants with respiratory distress
 - However there was no significant difference in the observed rate of intubation (secondary outcome) between treatment groups, suggesting that the use of NHF with “rescue CPAP” may result in similar rates of intubation





95% CONFIDENCE INTERVAL:

A statistical measure showing that 95% of the results for that parameter lie within the range quoted

100% RELATIVE HUMIDITY (RH):

The maximum amount of water a gas can hold at a given temperature.

BRONCHOPULMONARY DYSPLASIA (BPD):

A form of chronic lung disease that develops in premature neonates treated with oxygen and positive-pressure ventilation

BUBBLE CONTINUOUS POSITIVE AIRWAY PRESSURE (BCPAP):

Continuous positive airway pressure therapy delivered via a bubble generator

CONTINUOUS POSITIVE AIRWAY PRESSURE (CPAP):

A technique of respiratory therapy in which airway pressure is maintained above atmospheric pressure throughout the respiratory cycle by pressurization of the ventilatory circuit

DEAD SPACE:

A volume of gas that does not participate in gas exchange, is common to both the inspiratory and expiratory passages. There are different types of “dead space” including:

- **Alveolar dead space**
Volume of gas ventilating unperfused alveoli that has no blood perfusion (shunt or pulmonary embolism).
- **Anatomic dead space**
Volume of gas within the conducting zone of the lungs and upper airway. (Amount of volume that does not enter the alveoli.)
- **Mechanical dead space**
Expired air that is re-breathed through connecting tubing.
- **Physiological dead space**
Anatomic and alveolar dead space.

DISTENDING PRESSURE:

Pressure applied to the lungs to expand them. Can be applied using continuous positive or negative airway pressure to create a partial vacuum

ENDOTRACHEAL TUBE (ETT):

A tube inserted through the mouth or nose into the trachea to maintain an unobstructed airway

EXTUBATION:

Withdrawing an endotracheal tube (ETT) from a patient’s airway.

FRACTION OF INSPIRED OXYGEN (FIO₂):

The proportion of oxygen in the air that is inspired

FULL-TERM:

An infant born between 37 and 40 weeks gestation

FUNCTIONAL RESIDUAL VOLUME (FRC):

The volume in the lungs at the end-expiratory position

GA:

Gestational Age – Period of time between conception and birth.

HEATED, HUMIDIFIED GAS:

Air that has been heated and humidified prior to delivery by non-invasive ventilation, typically to 37°C and 100% relative humidity

HYPERCAPNIA:

The presence of an abnormally high level of carbon dioxide in the circulating blood

INFANT:

Children greater than 1 month to 2 years of age

INTUBATION:

The insertion of an ETT or tracheostomy tube into the trachea.

LOW BIRTH WEIGHT (LBW):

Birth weight less than 2500g

LUNG COMPLIANCE:

The ease of lung expansion

MECHANICAL VENTILATION (MV):

The use of an invasive artificial airway to mechanically assist or replace spontaneous breathing, when patients cannot do so on their own.

MINUTE VENTILATION (VE):

The volume of gas that moves in and out of the lungs in one minute; it is calculated by multiplying the exhaled tidal volume by the respiratory rate.

CONTINUOUS POSITIVE AIRWAY PRESSURE (CPAP):

A technique of respiratory therapy in which airway pressure is maintained above atmospheric pressure throughout the respiratory cycle by pressurization of the ventilator circuit.

NASAL INTERMITTENT POSITIVE PRESSURE VENTILATION (NIPPV):

A method of noninvasive ventilation that provides positive pressure to the back of the nose that is transferred to the lungs with intermittent breaths from a ventilator

NEONATE:

Premature infants or newborns less than 30 days old.

NICU:

Neonatal Intensive Care Unit

NON-INFERIORITY TRIAL:

A trial that is designed to determine whether the effect of a new treatment is not worse than a standard treatment

NONINVASIVE VENTILATION (NIV):

The delivery of ventilatory support without the need for an invasive artificial airway.

OXYGEN SATURATION (SPO₂):

Oxygen saturation as measured by pulse oximetry

PARTIAL PRESSURE OF CARBON DIOXIDE (PCO₂):

The part of total blood gas pressure exerted by carbon dioxide gas; a measure of how much carbon dioxide is dissolved in the blood and how well carbon dioxide is able to move out of the body.



PARTIAL PRESSURE OF OXYGEN (PAO₂):

The part of total blood gas pressure exerted by oxygen gas; a measure of how much oxygen is dissolved in the blood and how well oxygen is able to move from the airspace of the lungs into the blood

PEDIATRIC:

Referring to children up to 21 years of age; usually found in the PICU.

PICU:

Pediatric Intensive Care Unit

PNEUMOTHORAX:

Air or gas in the pleural space

POSITIVE END EXPIRATORY PRESSURE (PEEP):

The amount of pressure above atmospheric pressure present in the airway at the end of the expiratory cycle during mechanical ventilation

POSITIVE END INSPIRATORY PRESSURE (PIP):

The highest pressure applied to the lungs during inspiration

PRETERM:

An infant born before 37 weeks gestation regardless of their weight. They can be further divided into-

- **Moderate to late preterm**
32 to <37 weeks gestation
- **Very preterm**
28 to <32 weeks gestation
- **Extremely preterm**
<28 weeks gestation

RESPIRATORY DISTRESS SYNDROME (RDS):

A lung disease of the newborn, most frequently occurring in premature infants, that is caused by abnormally high alveolar surface tension as a result of a deficiency in lung surfactant; also called hyaline membrane disease

RESPIRATORY RATE:

The amount of breaths over a specified time period

SURFACTANT:

A substance produced in the lungs that tends to reduce the surface tension of the fluid in the lungs and helps make the small air sacs in the lung (alveoli) more stable

TIDAL VOLUME (VT):

Volume of air inspired or expired with each normal breath. The amount of gas delivered to a patient in one breath.

WORK OF BREATHING (WOB):

The force required to expand the lung against its elastic properties.



Notes



Notes

For more information please contact your local
Fisher & Paykel Healthcare representative