


F&P Neopuff CO₂ Detector

Instructions for Use



REF RD401

Setup (i)

During use (ii)

English	4
Français (French)	6
Deutsch (German)	8
Nederlands (Dutch)	10
Svenska (Swedish)	12
Italiano (Italian)	14
Norsk (Norwegian)	16
Suomi (Finnish)	18
Español (Spanish)	20
Português (Portuguese)	22
Português Brasil (Portuguese Brazilian)	24
Română (Romanian)	26
العربية (Arabic)	28
Ελληνικά (Greek)	30
Türkçe (Turkish)	32
Polski (Polish)	34
български (Bulgarian)	36
Dansk (Danish)	38
한국어 (Korean)	40
繁體中文版 (Traditional Chinese)	42
Íslenska (Icelandic)	44
Magyar (Hungarian)	46
Bahasa Indonesian (Indonesian)	48
Hrvatski (Croatian)	50
Česky (Czech)	52
Eesti keel (Estonian)	54
Latviešu (Latvian)	56
Lietuvių (Lithuanian)	58
Slovenčina (Slovak)	60
Slovenščina (Slovenian)	62

Intended purpose

The F&P Neopuff™ CO₂ Detector (RD401) provides colorimetric indication of exhaled CO₂ during the provision of invasive, or non-invasive ventilatory support to neonates and infants during T-piece resuscitation. This is a single use, prescription only device for use by healthcare professionals in a hospital environment. It is designed to be used in combination with approved Fisher & Paykel Healthcare products.

Instructions for use

The F&P Neopuff CO₂ Detector (RD401) is used with the F&P Neopuff T-piece resuscitator (RD900) when connected to the F&P ergonomic T-piece resuscitator circuit (RD1300). It is not to be used with a humidified and/or heated circuit.

Note: This device is an adjunct in patient assessment and is to be used in conjunction with other methods to determine clinical signs and symptoms.

Setup (i)

Do not use the detector if the device or packaging is damaged.
Do not use the detector if the device has expired.

1. Open the pouch at the tear tabs.
2. Ensure the detector sensor is colored blue-green.
Note: Discard if the detector sensor is permanently colored yellow.
3. Insert the detector into the duckbill port of the ergonomic T-piece circuit and ensure the detector is fully inserted.
4. Set up the Neopuff device as per the Neopuff user instructions.
5. Once the Neopuff T-piece resuscitator has been set up, do not remove the detector from the F&P ergonomic T-piece circuit.

During use (ii)

The F&P Neopuff CO₂ Detector **must not** come into contact with medication, including nebulized drugs.

6. Proceed to provide therapy to the patient as per the Neopuff user instructions.
The blue-green indicator circle will rapidly change from blue-green (inspiration) to yellow (expiration) in color when CO₂ is detected during therapy.
7. Dispose within 24 hrs of opening.

WARNING: The F&P Neopuff CO₂ Detector may not detect exhaled CO₂ if there is a significant level of leak between the patient and the T-piece.

Disposal

Dispose of this device according to hospital protocols.

Note: Reuse may result in transmission of infectious substances, interruption to treatment, serious harm and/or death.

Operating Range

The F&P Neopuff T-piece resuscitator should be used in the following operating range when using the F&P Neopuff CO₂ Detector:

- **Deliverable flow:** 8 to 15 L/min
- **PEEP:** 5 – 10 cmH₂O
- **PIP:** 20 – 40 cmH₂O
- **Operating temperature:** 18 to 40 °C

Technical specifications

- **Dead space:** 0 mL
- **Resistance to flow:** 0 cmH₂O @ 8 - 15 L/min
- **Patient weight range:** 400 to 5000 g.
- **Storage conditions:** -10 to +26 °C

The color indicators are as follows:

- **Blue-green:** <0.04% CO₂
- **Yellow:** >1.70% CO₂

Contraindications

There are no contraindications for using the Neopuff CO₂ Detector. Please refer to warnings.

General warnings

- The F&P Neopuff CO₂ Detector may not be interpreted properly by individuals with blue-green color blindness.
- Do not insert or remove the F&P Neopuff CO₂ Detector after PIP and PEEP have been set as this will change the pressures that are delivered. This may lead to patient harm.
- Failure to apply and use this product within the intended purpose, instructions for use, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient harm).
- California residents please be advised of the following pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: <http://www.fphcare.com/prop65>

Risks

- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.
- There are no known side effects associated with the use of the F&P Neopuff CO₂ Detector.

Symbol definitions

	Catalog number		Authorized European Representative
	Storage temperature range		Manufacturer
	Do not reuse		Consult instructions for use
	CE mark		Prescription only
	Lot number		Medical device
	Use-by date		Importer
	Country of origin and date of manufacture		Authorized Switzerland Representative
	Distributor		UK Responsible Person

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
Rx only

www.fphcare.com/neopuff-ifu

REF 626244 REV A 2023-10 © 2023 Fisher & Paykel Healthcare Limited.

F&P and Neopuff are trademarks of Fisher & Paykel Healthcare.

Additional copies of user instructions are available to the user on request.

For patent information, see www.fphcare.com/ip

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand
Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)**

Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务

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Japan (JP) Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)**

Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)**

Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań,

Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50

Switzerland (CH) CH REP Fisher & Paykel Healthcare GmbH,

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No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**

UK REP Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont

Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/**

Canada (CA) Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555

Denmark (DK) Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123

Ireland (IE) Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)**

Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)**

Tel: +46 8 564 76 680

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Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde