CE 0123 Rx only

fphcare.com/
neopuff-ifu

MD

F&P and Neopuff are trademarks of Fisher & Paykel Healthcare Limited.
REF 185045725 REV J 2023-05 © 2023 Fisher & Paykel Healthcare Limited

For patent information,
see www.fphcare.com/ip

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** **EC REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtabouef Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** **CH REP** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249, Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **UK REP** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

www.fphcare.com

Fisher & Paykel
HEALTHCARE






www.fphcare.com

Fisher & Paykel
HEALTHCARE

REF RD1300-10



F&P Neopuff™ Ergonomic T-piece Circuit

Instructions for use

1

2

English	6	Português (Brasil) (Brazilian Portuguese)	50
العربية (Arabic)	8	Română (Romanian)	52
Български (Bulgarian)	10	Русский (Russian)	54
Čeština (Czech)	12	Slovenčina (Slovak)	56
Dansk (Danish)	14	Slovenščina (Slovenian)	58
Deutsch (German)	16	Svenska (Swedish)	60
Ελληνικά (Greek)	18	ภาษาไทย (Thai)	62
Español (Spanish)	20	Türkçe (Turkish)	64
Eesti keel (Estonian)	22	اردو (Urdu)	66
Suomi (Finnish)	24	Tiếng Việt (Vietnamese)	68
Français (French)	26	繁體中文 (Traditional Chinese)	70
Hrvatski (Croatian)	28	简体中文 (Simplified Chinese)	72
Magyar (Hungarian)	30		
Bahasa Indonesia (Indonesian)	32		
Italiano (Italian)	34		
한국어 (Korean)	36		
Lietuvių (Lithuanian)	38		
Latviešu (Latvian)	40		
Nederlands (Dutch)	42		
Norsk (Norwegian)	44		
Polski (Polish)	46		
Português (Portuguese)	48		

RD1300-10 F&P Neopuff Ergonomic T-piece Circuit

English (en)

Intended purpose

The F&P Neopuff Ergonomic T-piece Circuit (RD1300-10) is a single use breathing circuit intended to provide ventilatory support to neonates and infants with respiratory insufficiency. The device is designed to connect to the F&P Neopuff T-piece Resuscitator or any other gas-powered resuscitator compliant to ISO-10651-5:2006. This device is designed for use in the hospital environment and must be prescribed by a physician. It is intended to be used by medical professionals including physicians, nurses, midwives and respiratory therapists. The intended population for use is neonates and infants.

Instructions for use

Consult the Labelling for your T-Piece Resuscitation Device for additional set-up & use instructions and additional warnings. An alternative non-gas powered means of resuscitation (such as a self-inflating bag) must be available.

Set-up

Do not use if package is damaged

- Connect the RD1300-10 to the T-piece Resuscitation Device.
 - Ensure the device is connected to the recommended flow source.
 - Adjust the flow rate to the prescribed input flow rate between 5 and 15 L/min.
 - DO NOT use outside this flow range**
 - Recommended flow rate is 8 L/min. For infants weighing ≥ 8 kg, a flow rate of 10 L/min or higher is recommended.**
- Occlude PEEP (Positive End Expiratory Pressure) cap and adjust blue PIP (Peak Inspiratory Pressure) knob on the Neopuff (or similar on other T-piece Resuscitators) to set PIP.
 - Test resuscitator function with the blue cap in place OR after connecting the test lung to the T-piece circuit. Test before using on the patient.**
- While the blue cap or test lung is still attached to the patient T-piece, adjust the PEEP cap to set PEEP

To resuscitate

Ensure you have removed blue cap or test lung before connecting to the mask or endotracheal tube.

- Resuscitate by placing and removing finger or thumb over the PEEP cap to allow inspiration and expiration

- Ensure pressures are monitored throughout resuscitation.**

Duckbill port

- The duckbill port may be used to deliver surfactant.
- The duckbill port may be used to perform suctioning.

Technical specifications

Input Flow Range	5 - 15 L/min
Circuit Length	1.6 m (63")
Patient Weight Range	Up to 10 kg (22 lb)
Dead Space	3.3 mL
Positive End-Expiratory Pressure (PEEP)	@5 L/min 1 to 5 cmH ₂ O [mbar] @8 L/min 1 to 9 cmH ₂ O [mbar] @10 L/min 2 to 15 cmH ₂ O [mbar] @15 L/min 3 to 25 cmH ₂ O [mbar]
T-piece Connection	15 mm female and 22 mm male
Intended T-piece Connections	Standard neonatal interface
Intended Duration of Use	Maximum of 14 days
Operating Temperature	-18 °C to 50 °C (-0.4 to +122 °F), up to 95% humidity
Storage Temperature	Storage: -10 to 50 °C (+14 to +122 °F), up to 95% humidity
Duckbill Port Catheter Sizing	5 Fr to 8 Fr (1.7 - 2.7 mm OD ; 0.06" - 0.11" OD)

Contraindications

There are no contraindications for using this device. Please refer to general warnings.

GENERAL WARNINGS

-  The F&P Neopuff Ergonomic T-piece Circuit is designed for SINGLE USE ONLY. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
 - Do not smoke, use naked flames and other sources of ignition near the unit when it is in use.

- The F&P Neopuff Ergonomic T-piece Circuit and F&P Neopuff T-piece Resuscitator are intended for use by medical professionals trained in infant resuscitation.
- DO NOT soak, wash, sterilize or re-use this product.
- DO NOT occlude the resuscitation tubing.
- DO NOT connect the F&P Neopuff Ergonomic T-piece circuit (RD1300-10) directly to a high pressure gas source. The F&P Neopuff Ergonomic T-piece circuit (RD1300-10) must be connected to the outlet port of the Neopuff or gas-powered resuscitator compliant to ISO 10651-5:2006. Bypassing the resuscitation device may result in serious patient harm (e.g. lung injury).
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65
- As this medical device uses an alternative small-bore connector design different from those specified in the ISO 80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need to be taken by the user to mitigate these reasonable foreseeable risks.

Risks

- There are no reported side effects associated with the use of this device.
- The use of this product is not without risk even when used as intended. Following all instructions and warnings provided, risks of hypoxia and barotrauma remain. These risks may result in serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

Disposal

Discard according to hospital protocols.

Symbol definitions

	Consult Instructions for use www.fphcare.com/neopuff-ifu		CE Mark
	Latex free		Do not reuse
	Authorized representative in the European Community		Not made with phthalates (DEHP, DBP, BBP)
	Date of manufacture and Country of manufacture		Use by date
	Storage temperature limit		Medical device
	Prescription only		Catalogue number
	Manufacturer		Batch code
	Distributor		Importer
	UKCA mark		Do not use if package is damaged

For patent information, see www.fphcare.com/ip

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