CE 0123 Rx only  fphcare.com/  
neopuff-ifu   

F&P and Neopuff are trademarks of Fisher & Paykel Healthcare Limited.  
REF 185043544 REV K 2022-11 © 2022 Fisher & Paykel Healthcare Limited

For patent information,  
see [www.fphcare.com/ip](http://www.fphcare.com/ip)

 Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland  
2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100  
Email: [info@fphcare.co.nz](mailto:info@fphcare.co.nz) Web: [www.fphcare.com](http://www.fphcare.com)

 **Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street,  
Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel  
do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil  
Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有  
限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486  
**France (FR)**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5,  
BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201  
Email: [c.s@fphcare.fr](mailto:c.s@fphcare.fr) **Germany (DE)** Fisher & Paykel Healthcare GmbH,  
Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)**  
Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110  
Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626  
**Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań,  
Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)**  
 Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland  
Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma:  
Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Östim Mahallesi  
1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK**  
**(GB)**  Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont  
Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)**  
Tel: 1800 446 3908 or +1 949 453 4000

**Austria (AT)** Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)**  
Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011  
**Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)**  
Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680


[www.fphcare.com](http://www.fphcare.com)

**Fisher & Paykel**  
HEALTHCARE

REF 900RD110


## F&P Neopuff™ Humidified Resuscitation Circuit

### Instructions for use


[www.fphcare.com](http://www.fphcare.com)



**Fisher & Paykel**  
HEALTHCARE

1

2

|                               |    |                                              |     |
|-------------------------------|----|----------------------------------------------|-----|
| English                       | 6  | Português (Brasil)<br>(Brazilian Portuguese) | 72  |
| العربية (Arabic)              | 9  | Română (Romanian)                            | 75  |
| Български (Bulgarian)         | 12 | Русский (Russian)                            | 78  |
| Čeština (Czech)               | 15 | Slovenčina (Slovak)                          | 81  |
| Dansk (Danish)                | 18 | Slovenščina (Slovenian)                      | 84  |
| Deutsch (German)              | 21 | Svenska (Swedish)                            | 87  |
| Ελληνικά (Greek)              | 24 | ภาษาไทย (Thai)                               | 90  |
| Español (Spanish)             | 27 | Türkçe (Turkish)                             | 93  |
| Eesti keel (Estonian)         | 30 | اردو (Urdu)                                  | 96  |
| Suomi (Finnish)               | 33 | Tiếng Việt (Vietnamese)                      | 99  |
| Français (French)             | 36 | 繁體中文 (Traditional Chinese)                   | 102 |
| Hrvatski (Croatian)           | 39 | 简体中文 (Simplified Chinese)                    | 105 |
| Magyar (Hungarian)            | 42 |                                              |     |
| Bahasa Indonesia (Indonesian) | 45 |                                              |     |
| Italiano (Italian)            | 48 |                                              |     |
| 한국어 (Korean)                  | 51 |                                              |     |
| Lietuvių (Lithuanian)         | 54 |                                              |     |
| Latviešu (Latvian)            | 57 |                                              |     |
| Nederlands (Dutch)            | 60 |                                              |     |
| Norsk (Norwegian)             | 63 |                                              |     |
| Polski (Polish)               | 66 |                                              |     |
| Português (Portuguese)        | 69 |                                              |     |

## 900RD110 F&P Neopuff Humidified Resuscitation Circuit

English (en)

### Intended purpose

The F&P Neopuff Humidified T-piece Circuit (900RD110) is a single use breathing circuit intended to provide heated, humidified ventilatory support to neonates and infants with respiratory insufficiency. The device is designed to connect to the F&P Neopuff Infant T-piece Resuscitator and F&P MR850 humidifier. This device is designed for use in the hospital environment and must be prescribed by a physician. It is intended to be used by medical professionals including physicians, nurses, midwives and respiratory therapists. The intended population for use is neonates and infants.

### Instructions for use

Consult the Labelling for your Neopuff T-piece Resuscitation Device, MR850/MR810\* Humidifier and MR290/MR225\* chambers for additional set-up & use instructions and additional warnings.

\*Not available in all markets

An alternative non-gas powered means of resuscitation (such as a self-inflating bag) must be available.

### Set-up

**Do not use if package is damaged**

#### Use Sterile water in the humidifier's Chamber

##### ① Fit a chamber to the humidifier.

Slide a MR290 chamber onto the humidifier base. Remove the blue caps. Hang the water bag, unwind the waterfeed set and spike water bag. The bag should be at least 50 cm above the chamber. Use USP sterile water for inhalation or equivalent.

**Or**

Slide a MR225 chamber onto the humidifier base. Remove the blue caps. Fill with 30 mL of sterile water or equivalent. This will allow for 60 minutes of running time (@ 8 L/min flow rate) and then refill as necessary during resuscitation.

##### ② Connect Circuit

- Connect the dry-line between the chamber inlet port and the Neopuff gas outlet.
- Connect the blue inspiratory limb to the remaining chamber port.
- Ensure the temperature probes and heater wire adaptors are fully inserted. If temperature probes are not used (i.e. using the MR810), ensure caps (BC410) are placed over the temperature probe ports.

##### ③ Connect to flow source

- Connect a blended oxygen/air supply or oxygen to the gas inlet port using gas supply line.
  - The Humidified Resuscitation System is to be used with medical grade oxygen/air**
- Adjust the flow rate to the prescribed input flow rate between 8 and 10 L/min.
  - DO NOT attempt to use product outside this flow range. Recommended flow rate 8 L/min. For infants weighing  $\geq 8$  kg, a flow rate of 10 L/min is recommended.**

##### ④ Set humidifier

Turn on the humidifier base.

MR850 – Check the humidifier is set to invasive mode. If not, hold the mode button until the invasive mode light is lit.

MR810 – Set the humidifier to high mode. Ensure the green heater wire light is lit.  Caution is required when delivering gases that are not body temperature saturated, to patients with bypassed airways.

##### ⑤ Set PIP

Ensure the test lung is inserted into the patient T-piece. Occlude PEEP (Positive End Expiratory Pressure) cap and adjust blue PIP (Peak Inspiratory Pressure) knob on the Neopuff (or similar on other T-piece Resuscitators) to set PIP.

**Connect the test lung to the T-piece circuit and test resuscitator function before connecting to the patient.**

##### ⑥ Set PEEP

While the test lung is still attached to the patient T-piece, adjust the PEEP cap to set PEEP.

Resuscitate by placing and removing finger or thumb over the PEEP cap to allow inspiration and expiration.

**Ensure pressures are monitored throughout resuscitation.**

## Technical specifications

|                                |                                                                                           |
|--------------------------------|-------------------------------------------------------------------------------------------|
| Input Flow Range               | 8–10 L/min                                                                                |
| Positive End-Expiratory (PEEP) | @8 L/min 1 to 10 cmH <sub>2</sub> O [mbar]<br>@10 L/min 2 to 13 cmH <sub>2</sub> O [mbar] |
| T-piece Connection             | 15 mm female (conical connector)                                                          |
| Circuit Length                 | 1.3 m                                                                                     |
| Patient Weight Range           | Up to 10 kg (22 lb)                                                                       |
| Dead Space                     | 6 mL                                                                                      |
| Intended T-piece Connections   | Standard neonatal interface                                                               |
| Intended Duration of Use       | Maximum of 7 days                                                                         |
| Operating Temperature          | 18 °C to 26 °C (+64 to +79 °F),<br>30 to 75% humidity                                     |
| Storage Temperature            | Storage: -10 to 50 °C (+14 to +122 °F), up to 95% humidity                                |
| Humidity                       | MR850 >33 mg/L<br>MR810 >10 mg/L                                                          |

## Contraindications

There are no contraindications for using this device. Please refer to general warnings.

## GENERAL WARNINGS

-  The humidified resuscitation system is designed for SINGLE USE ONLY. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- Do not smoke, use naked flames and other sources of ignition near the unit when it is in use.
  - The F&P Neopuff Humidified T-piece Circuit and F&P Neopuff T-piece Resuscitator are intended for use by medical professionals trained in infant resuscitation.
  - DO NOT soak, wash, sterilize or re-use this product.
  - DO NOT occlude the resuscitation tubing.
  - DO NOT operate the humidifier or the circuit with no gas flow.
  - Ensure the humidifier is mounted below the level of the patient.
  - DO NOT cover the circuit with materials such as blankets, towels or bed linen.
  - Avoid prolonged contact with patient's skin.
  - DO NOT stretch or milk tubing.
  - DO NOT fill the chamber above the maximum fill level line.
  - California residents please be advised of the following, pursuant to Proposition 65: This

product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)

- As this medical device uses an alternative small-bore connector design different from those specified in the ISO 80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need be taken by the user to mitigate these reasonable foreseeable risks.

## Risks

- There are no reported side effects associated with the use of this device.
- The use of this product is not without risk even when used as intended. Following all instructions and warnings provided, risks of hypoxia and barotrauma remain. These risks may result in serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority

## Disposal

Discard according to hospital protocols.

## Symbol definitions

|                |                                                                                                           |  |                                           |
|----------------|-----------------------------------------------------------------------------------------------------------|--|-------------------------------------------|
|                | Consult Instructions for use <a href="http://www.fphcare.com/neopuff-ifu">www.fphcare.com/neopuff-ifu</a> |  | CE Mark                                   |
|                | Latex free                                                                                                |  | Do not reuse                              |
|                | Authorized representative in the European Community                                                       |  | Not made with phthalates (DEHP, DBP, BBP) |
|                | Date of manufacture and country of manufacture                                                            |  | Use by date                               |
|                | Storage temperature limit                                                                                 |  | Medical device                            |
| <b>Rx only</b> | Prescription only                                                                                         |  | Catalogue number                          |
|                | Manufacturer                                                                                              |  | Batch code                                |
|                | Distributor                                                                                               |  | Importer                                  |
|                | UKCA mark                                                                                                 |  | Do not use if package is damaged          |

For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip)

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