CE 0123 Rx only

fphcare.com/
neopuff-ifu

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

REF 616858 REV B 2023-05 © 2023 Fisher & Paykel Healthcare Limited

www.fphcare.com

Fisher & Paykel
HEALTHCARE

Instructions for use

T-Piece Resuscitator


REF RD900AEU


I


1


2




3


4




8

9

English	6
العربية (Arabic)	9
Français (French)	12
Bahasa Indonesia (Indonesian)	15
한국어 (Korean)	18
Русский (Russian)	21
ภาษาไทย (Thai)	24
اردو (Urdu)	27
Tiếng Việt (Vietnamese)	30
繁體中文 (Traditional Chinese)	33

Intended Purpose

The F&P Neopuff T-piece Resuscitator (RD900AEU) is a reusable, manually operated, gas-powered resuscitator which provides ventilatory support to neonates and infants with respiratory insufficiency. This device is designed for use in the hospital environment and must be prescribed by a physician. It is intended to be used by medical professionals including physicians, nurses, midwives and respiratory therapists. The intended population for use is neonates and infants.

Clinical benefit

The F&P Neopuff T-piece Resuscitator provides life-saving ventilatory support to neonates and infants. It enables the delivery of controlled and consistent inflation pressures to a patient's airway, thereby helping to reduce the risk of lung damage and to initiate and maintain effective ventilation.

Instructions for use

1 Set up

The following procedure should be carried out prior to every use of the Neopuff to ensure that the device is functioning properly.

1. Check manometer reads zero with no gas flow. If not, the manometer requires calibration (refer to section 3.7) of the Technical Manual).
2. **Connect gas supply**
Connect a blended oxygen/air supply or oxygen to the gas inlet port using the Gas Supply Line.
3. **Connect T-Piece Circuit**
Connect the T-Piece Circuit to the gas outlet port. If using the Ergonomic T-piece  leave blue cap in place on the T-piece circuit or connect Test Lung to T-piece circuit. If using the Classic T-piece , connect Test Lung to T-piece circuit. Before use, inspect Test Lung for signs of damage such as discoloration.
Note: Replace test lung (RD020-01) if internal surfaces become contaminated or damaged.
4. **Check settings**
Adjust the flow rate to the prescribed input flow rate. Recommended flow rate is 8 L/min. For infants weighing ≥8 kg, a flow rate of 10 L/min or higher is recommended.
Note: Ensure the oxygen concentration of an oxygen/air supply is either monitored using an oxygen analyzer or preset using oxygen/air flow rate graphs.

5a. To check max pressure

Occlude PEEP (Positive End Expiratory Pressure) cap and turn PIP (Peak Inspiratory Pressure) control fully clockwise, until the

knob doesn't turn anymore.  or .

- 5b. Adjust max pressure control knob clockwise or counterclockwise to set desired max pressure.

Notes:

- The factory setting of the Max Pressure Relief is 40 cmH₂O [mbar].
- The Max Pressure Relief valve acts as an overall limit on the achievable circuit pressure. Resuscitation above 40 cmH₂O [mbar] cannot be achieved unless the Max Pressure Relief valve is adjusted.

Warning: The Max Pressure Relief can be adjusted up to a nominal 80 cmH₂O [mbar]. This should only be done in exceptional circumstances by a medical professional trained in infant resuscitation. Do not set the Max Pressure Relief above 80 cmH₂O [mbar].

6. To set PIP

While still occluding the PEEP cap, turn PIP control knob counterclockwise until the desired peak inspiratory pressure is set.  or .

7. To set PEEP

Adjust PEEP cap to the desired PEEP level.

Turn off gas supply and remove Test Lung from T-Piece. Ensure that the rigid plastic connector of the Test Lung is also removed from the T-Piece before attempting to connect a patient interface. Failing to do so may cause unacceptable delays during patient resuscitation.  or .

To resuscitate

8. Adjust gas supply to the flow rate used during set up.
9. Fit T-Piece to neonatal resuscitation mask and place over the baby's mouth and nose. Resuscitate by placing and removing finger or thumb over the PEEP cap to allow inspiration and expiration.  or .
10. OR Fit T-Piece to the endotracheal tube.
Resuscitate by placing and removing finger or thumb over the PEEP cap to allow inspiration and expiration.  or .

Ensure Pressures are monitored throughout resuscitation

Cleaning Instructions

Clean the exterior surfaces of the F&P Neopuff T-piece Resuscitator, using a cloth dampened with neutral detergent or an Ammonium Quat. disinfectant with a maximum alcohol content of 20% (such as Caviwipes, Sani-cloth HBm Ascepti-Wipes II, Spartan's TB-Cide Quat Wipes or Vernacare Tuffie Wipes Alcohol Free).

Clean the exterior surfaces of the Test Lung, Gas Supply Line and Gas Inlet Adapter using a cloth dampened with neutral detergent.

Do not clean the internal surfaces of the Test Lung (RD020-01),

Gas Supply Line (900RD008 & 900RD009) or Gas Inlet Adapter (900RD101).

Routine maintenance

A full technical description, including routine maintenance and service data, is contained in the Product Technical Manual available from your supplier, Fisher & Paykel Healthcare or fpkcare.com/neopuff-ifu

Disposal

Discard according to hospital protocols

Technical specifications

Patient Weight Range	Up to 10 kg (22 lb)
Manometer accuracy	± 2.0% of Full Scale Deflection
Manometer Range	-10 to 80 cmH ₂ O [mbar]
Peak Inspiratory Pressure (PIP)	@ 5 L/min approx. 2 to 70 cmH ₂ O [mbar]
	@ 8 L/min approx. 3 to 72 cmH ₂ O [mbar]
	@ 10 L/min approx. 4 to 73 cmH ₂ O [mbar]
	@ 15 L/min approx. 8 to 75 cmH ₂ O [mbar]
Positive End Expiratory Pressure (PEEP)	@ 5 L/min approx. 1 to 6 cmH ₂ O [mbar]
	@ 8 L/min approx. 1 to 10 cmH ₂ O [mbar]
	@ 10 L/min approx. 2 to 15 cmH ₂ O [mbar]
	@ 15 L/min approx. 4 to 17 cmH ₂ O [mbar]
Gas Inlet Flow Range	5 L/min (min) to 15 L/min (max)
Operating Time (400 L cylinder)	50 minutes (typical value based on a gas flow rate of 8 L/min)

Note: All performance figures listed above are representative only. PEEP values stated are based on typical clinical PIP settings. Higher PEEP values can be achieved if higher PIP values are set

Approved accessories and spare parts include.

- F&P Classic T-piece Circuit*
- F&P Ergonomic T-piece Circuit*
- F&P Humidified T-Piece Circuit*
- F&P Resuscitation Mask*
- F&P Gas Supply Lines
- F&P Gas Inlet Adapter
- F&P Test Lung
- F&P Mounting Brackets*

Note*: Not available in all markets, contact your local supplier or Fisher & Paykel Healthcare Representative

Storage Temperature: -10 to 50 °C (+14 to +122 °F), up to 95% humidity

Operating Conditions: -18°C to 50°C (-0.4 to +122°F), up to 95% humidity

The Neopuff™ T-piece Resuscitator has an expected service life of 10 years

Service Life of Test Lung (RD020-01): 12 Months

Service Life of Gas Supply Line (900RD008 & 900RD009): 12 Months

Service Life of Gas Inlet Adapter (900RD101): 12 Months

Contraindications

There are no contraindications for using this device. Please refer to general warnings

GENERAL WARNINGS

Please read and understand the instructions fully before using the Neopuff and related accessories

- The F&P Neopuff T-piece Resuscitator is to be used by medical professionals trained in infant resuscitation. It is recommended that users refer to local or international resuscitation guidelines to determine the suitability of different types of resuscitators for use, and resuscitation techniques.
- It is the responsibility of the purchaser to ensure all users of this device have been adequately trained in resuscitation techniques.
- The Neopuff T-piece Resuscitator must only be used after checking that correct pressures will be delivered to the patient.
- Do not smoke, use naked flames and other sources of ignition near the unit when it is in use.
- Connect to flow regulated oxygen or oxygen/air mixture only.
- The input gas flow rate is 5 to 15 L/min. The recommended operating gas flow rate is 8 L/min. Do not exceed the maximum flow rate of 15 L/min. Input flow ranges are circuit specific - refer to the circuit User Instructions.
- Use only recommended Fisher & Paykel Healthcare Neopuff T-piece Resuscitator accessories.
- The use of Fisher & Paykel Healthcare Gas Supply Lines is recommended
- Ensure all oxygen and air supplies are turned off and disconnected from the Neopuff before performing cleaning procedures. Explosion and fire hazards can exist when performing cleaning procedures in an oxygen-enriched environment.
- An alternative non-gas powered means of resuscitation (such as a self-inflating bag) must be available.
- For US Only: US Federal law restricts this device to sale by or

on the order of a licensed physician.

- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65
- Neopuff handle contains Nickel which may result in skin sensitization or an allergic skin reaction.
- Neopuff manometer contains Beryllium
- As this medical device uses an alternative small-bore connector design different from those specified in the ISO 80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need be taken by the user to mitigate these reasonable foreseeable risks.
- Do not use oil, grease or other substances that are incompatible with oxygen on any part of the F&P Neopuff T-piece Resuscitation system.
- Dropping the F&P Neopuff T-piece Resuscitator or other similar forms of impact may cause damage resulting in incorrect operation of the unit.
- Explosion and fire hazards can exist when performing cleaning procedures in an oxygen-enriched environment.

Risks

- There are no reported side effects associated with the use of this device.
- The use of this product is not without risk even when used as intended. Following all instructions and warnings provided, risks of hypoxia and barotrauma remain. These risks may result in serious injury or death
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

Symbol definitions

	Not made with phthalates (DEHP, DBP, BBP).		Storage temperature range
	Batch code		Catalogue number
	Consult instructions for use. www.fphcare.com/neopuff-ifu	Rx only	Prescription only

	Authorized representative in the European Community		Latex free
	Manufacturer		Medical device
	Date of manufacture and country of manufacture		Importer
	Distributor		UKCA mark
	Serial number		CE Mark
	Contains hazardous substances		

For patent information, see www.fphcare.com/ip

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