Infant T-Piece Resuscitator 900 Series









Instructions for Use




Neopuff is a trademark of Fisher & Paykel Healthcare Limited

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor
Australia (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/ 售后服务机构-费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 94 1590 355 Fax: +46 83 66 310 **France** [EC REP] Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH & Co. KG, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Russia** Tel: +7 495 782 2150 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

I Set-Up

II To Resuscitate

INFANT T-PIECE RESUSCITATOR 900 SERIES

ENGLISH en

Intended Use

The Fisher & Paykel Healthcare Neopuff™ Infant Resuscitator is a manually operated, gas-powered resuscitator which provides assisted respiratory breaths to newborn babies in delivery suites, nurseries and neonatal intensive care units.

Symbol	Definition
	Attention: Consult the Instructions for Use

General Warnings

- Please read and understand the instructions fully before using Neopuff and related accessories. The Neopuff is to be used only by persons trained in infant resuscitation. It is recommended that users refer to local or international resuscitation guidelines to determine the suitability of different types of resuscitator for use, and resuscitation technique.
- It is the responsibility of the purchaser to ensure that all users of this device have been adequately trained in resuscitation techniques.
- The Neopuff must only be used after checking that correct pressures will be delivered to the baby.
- Ensure no smoking, naked flames or sources of ignition are present while the unit is in use.
- Connect to flow-regulated oxygen or oxygen/air mixture only.
- Input gas flow rate 5 to 15 L/min. Recommended operating gas flow rate 8 L/min. Do not attempt to use a higher flow than 15 L/min. Input flow ranges are circuit specific, refer to circuit User Instructions.
- The Max Pressure Relief can be adjusted up to a nominal 80 cmH₂O [mbar], and should only be done in exceptional circumstances by persons trained in infant resuscitation. Do not attempt to set the Max Pressure Relief above 80 cmH₂O [mbar].
- Use only recommended Fisher & Paykel Healthcare Infant Resuscitator accessories.
- Use only a Fisher & Paykel Healthcare Gas Supply Line or approved equivalent.
- Ensure all oxygen and air supplies are turned off and disconnected from the Neopuff before performing cleaning procedures. Explosion and fire hazards can exist when performing cleaning procedures in an oxygen-enriched environment.
- An alternative means of resuscitation must be available.
- For US only: US Federal law restricts this device to sale by or on the order of a licensed physician.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

1 Set-Up

The following procedure should be carried out prior to every use of the Neopuff to ensure that the device is functioning correctly.

- 1 Check manometer reads zero with no gas flow. If not, the manometer requires calibration (refer to section 3.3.4 of the Technical Manual).

2 Connect Gas Supply

Connect an oxygen or blended oxygen/air supply to the gas inlet port using the Gas Supply Line.

3 Connect T-Piece Circuit

- Connect the T-Piece Circuit to the gas outlet port.
- If using Ergonomic T-piece  leave blue cap in place on the T-piece circuit or connect Test Lung to T-piece circuit. If using Classic T-piece , connect Test Lung to T-piece circuit. Before use, inspect Test Lung for signs of damage such as discoloration

4 Check Settings

Adjust the gas supply to desired flow rate between 5 and 15 L/min.

Note: Ensure the oxygen concentration of an oxygen/air supply is either monitored using an oxygen analyzer or preset using oxygen/air flow rate graphs.

5a To Check Max Pressure

Occlude PEEP¹ cap and turn PIP² control fully clockwise, until the knob doesn't turn anymore. ( or 

1. *Positive End Expiratory Pressure*
2. *Peak Inspiratory Pressure*

- ## 5b
- Adjust max pressure control knob clockwise or counterclockwise to set desired max pressure.

Note:

- The factory setting of the Max Pressure Relief is 40 cmH₂O [mbar].
- The Max Pressure Relief valve acts as an overall limit on the achievable circuit pressure. Resuscitation above 40 cmH₂O [mbar] cannot be achieved unless the Max Pressure Relief valve is adjusted.

6 To Set PIP

While still occluding the PEEP cap, turn PIP control knob counterclockwise until the desired peak inspiratory pressure is set. ( or 

7 To Set PEEP

Adjust PEEP cap to the desired PEEP level.

Turn off gas supply and remove Test Lung from T-Piece. Ensure that the rigid plastic connector of the Test Lung is also removed from the T-Piece before attempting to connect a mask or endotracheal tube. Failing to do so may cause unacceptable delays during patient resuscitation. (i or ii)

II To Resuscitate

8 Adjust gas supply to the desired flow rate.

9 Fit T-Piece to neonatal resuscitation mask and place over the baby's mouth and nose.

Resuscitate by placing and removing finger or thumb over the PEEP¹ cap to allow inspiration and expiration. (i or ii)

1. *Positive End Expiratory Pressure*

10 OR Fit T-Piece to the endotracheal tube.

Resuscitate by placing and removing finger or thumb over the PEEP¹ cap to allow inspiration and expiration. (i or ii)

1. *Positive End Expiratory Pressure*

Cleaning and Servicing

- Clean external surfaces of the Neopuff and Gas Supply Line using a damp cloth and mild soapy water or Isopropyl Alcohol.
- Dry all surfaces after cleaning with a clean soft cloth or paper towel.
- For more information on cleaning and maintenance of the Neopuff, Test Lung and other accessories, please refer to the Technical Manual (Part No. 185041597), available at fphcare.com

PERFORMANCE SPECIFICATIONS

Recommended Body Weight Range	Up to 10 kg
Manometer Range	-10 to 80 cmH ₂ O [mbar]
Peak Inspiratory Pressure (PIP)	@ 5 L/min approx. 2 to 70 cmH ₂ O [mbar]
	@ 8 L/min approx. 3 to 72 cmH ₂ O [mbar]
	@ 10 L/min approx. 4 to 73 cmH ₂ O [mbar]
	@ 15 L/min approx. 8 to 75 cmH ₂ O [mbar]
Positive End Expiratory Pressure (PEEP)	@ 5 L/min approx. 1 to 6 cmH ₂ O [mbar]
	@ 8 L/min approx. 1 to 10 cmH ₂ O [mbar]
	@ 10 L/min approx. 2 to 15 cmH ₂ O [mbar]
	@ 15 L/min approx. 4 to 17 cmH ₂ O [mbar]
Gas Inlet Flow Range	5 L/min (min) to 15 L/min (max)
Operating Time (400 L cylinder)	50 minutes (typical value based on a gas flow rate of 8 L/min)
<p>NOTE: All performance figures listed above are representative only. PEEP values stated are based on typical clinical PIP settings. Higher PEEP values can be achieved if higher PIP values are set.</p>	

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