CE 0123 Rx only  fphcare.com/
neopuff-ifu   

F&P and Neopuff are trademarks of Fisher & Paykel Healthcare Limited.
REF 185043891 REV I 2023-05 © 2023 Fisher & Paykel Healthcare Limited

For patent information,
see www.fphcare.com/ip

 Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box
14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz
Web: www.fphcare.com

 **Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street,
Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel
do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil
Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有
限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486
France (FR)  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5,
BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201
Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49,
73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032
India (IN) Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)**
Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel
Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464
Russia (RU) Tel: +7 495 782 21 50 **Switzerland (CH)**  Fisher & Paykel Healthcare
GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)**
Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret
Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara,
Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit
16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189
USA (US)/Canada (CA) Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)**
Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)**
Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346
Sweden (SE) Tel: +46 8 564 76 680


www.fphcare.com

Fisher & Paykel
HEALTHCARE

REF **RD803/RD804/RD805/RD806/
RD807**


F&P Neopuff™ Resuscitation Masks

Instructions for use


www.fphcare.com


Fisher & Paykel
HEALTHCARE

1

2

English	6	Português (Brasil)	50
العربية (Arabic)	8	(Brazilian Portuguese)	
български (Bulgarian)	10	Română (Romanian)	52
Čeština (Czech)	12	Русский (Russian)	54
Dansk (Danish)	14	Slovenčina (Slovak)	56
Deutsch (German)	16	Slovenščina (Slovenian)	58
Ελληνικά (Greek)	18	Svenska (Swedish)	60
Español (Spanish)	20	ภาษาไทย (Thai)	62
Eesti keel (Estonian)	22	Türkçe (Turkish)	64
Suomi (Finnish)	24	اردو (Urdu)	66
Français (French)	26	Tiếng Việt (Vietnamese)	68
Hrvatski (Croatian)	28	繁體中文 (Traditional Chinese)	70
Magyar (Hungarian)	30	简体中文 (Simplified Chinese)	72
Bahasa Indonesia (Indonesian)	32		
Italiano (Italian)	34		
한국어 (Korean)	36		
Lietuvių (Lithuanian)	38		
Latviešu (Latvian)	40		
Nederlands (Dutch)	42		
Norsk (Norwegian)	44		
Polski (Polish)	46		
Português (Portuguese)	48		

Intended Purpose

The F&P Neopuff Resuscitation Mask is a single use patient interface designed to deliver ventilatory support to neonates and infants with respiratory insufficiency. The device is designed to connect to a T-piece circuit or manual resuscitation device with an ISO 5356-1 15 mm conical female fitting. This device is designed for use in the hospital environment and must be prescribed by a physician. It is intended to be used by medical professionals including physicians, nurses, midwives and respiratory therapists. The intended population for use is neonates and infants.

Instructions for use

Set-up

Do not use if package is damaged

- 1 Fit an appropriately sized mask to the patient, covering the nose and mouth.
- 2 Insert the mask into the resuscitation device with a twisting motion.

Do not use this mask for therapies other than resuscitation

To Resuscitate

- 3 Create a seal over the patient's face to deliver therapy.

Ensure pressures are monitored throughout resuscitation

Disposal

Discard according to hospital protocols.

Contraindications

There are no contraindications for using this device. Please refer to general warnings.

Technical Specifications

Device connection	15 mm male	
Intended device connections	Common resuscitation devices	
Intended duration of use	Maximum of 14 days	
Operating temperature	-18°C to 50°C (-0.4 to +122°F), up to 95% humidity	
Storage temperature	-10°C to 50°C (+14 to +122°F), up to 95% humidity	
Mask size	Diameter	Dead space
RD803 (XS)	35 mm (1.38 in.)	13 mL
RD804 (S)	42 mm (1.65 in.)	21 mL
RD805 (M)	50 mm (1.97 in.)	34 mL
RD806 (L)	60 mm (2.36 in.)	58 mL
RD807 (XL)	72 mm (2.83 in.)	97 mL

GENERAL WARNINGS

The Resuscitation Mask is designed for **SINGLE USE ONLY**. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.

- The F&P Neopuff Resuscitation Mask is intended for use by medical professionals trained in infant resuscitation.
- DO NOT soak, wash, sterilize or re-use this product.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Risks

- There are no reported side effects associated with the use of this device.
- The use of this product is not without risk even when used as intended. Following all instructions and warnings provided, risks of hypoxia and barotrauma remain. These risks may result in serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

Symbol definitions

	Consult Instructions for use www.fphcare.com/neopuff-ifu		CE Mark
	Latex free		Do not reuse
	Authorized representative in the European Community		Not made with phthalates (DEHP, DBP, BBP)
	Date of manufacture and country of manufacture		Use by date
	Do not use if package is damaged		Medical device
	Storage temperature limit		Catalogue number
Rx only	Prescription only		Batch code
	Manufacturer		Importer
	Distributor		UKCA mark

For patent information, see www.fphcare.com/ip

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