









RT046 XS/S/M/L

Non-Vented Hospital Full Face Mask with Standard Elbow

Single Use 14 Rx only CE 0123

Fisher & Paykel HEALTHCARE

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

AUSTRALIA (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131 Tel: +61 3 9871 4900 **BRAZIL (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **CHINA (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路 31 号 G12 栋 301 号 电话: +86 20 32053486 **FRANCE (FR)** **EC REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât. F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 64465201 Email: c.s@fphcare.fr **GERMANY (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **HONG KONG (HK)** Tel: +852 2116 0032 **INDIA (IN)** Tel: +91 80 2309 6400 **JAPAN (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **KOREA (KR)** Tel: +82 2 6205 6900 **MEXICO (MX)** Tel: +52 55 9130 1626 **POLAND (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **RUSSIA (RU)** Tel: +7 495 782 21 50 **SWITZERLAND (CH)** **CH REP** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **TAIWAN (TW)** Tel: +886 2 8751 1739 **TURKEY (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **UK REP** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/CANADA (CA)** Tel: 1800 446 3908 or +1 949 453 4000 **AUSTRIA (AT)** Tel: 0800 29 31 23 **BENELUX (BE NL LU)** Tel: +31 40 216 3555 **DENMARK (DK)** Tel: +45 70 26 37 70 **FINLAND (FI)** Tel: +358 9 251 66 123 **IRELAND (IE)** Tel: 1800 409 011 **ITALY (IT)** Tel: +39 06 7839 2939 **NORWAY (NO)** Tel: +47 21 60 13 53 **SPAIN (ES)** Tel: +34 902 013 346 **SWEDEN (SE)** Tel: +46 8 564 76 680

English
Non-Vented Hospital Full Face Mask Standard Elbow Version

Intended Use
The Fisher & Paykel Healthcare single patient use masks are intended for use as an accessory to ventilators to enable non-invasive positive pressure ventilation (NPPV) therapy (CPAP or bi-level) to be delivered to spontaneously breathing adult patients (>30 kg) with respiratory insufficiency or respiratory failure who have been prescribed NPPV. The masks are to be fitted and therapy maintained by trained medical practitioners in a hospital/institutional environment with patient monitoring in place.

Technical Specifications	
Operating Pressure Range:	4 – 25 cmH ₂ O
Interface Connections:	ISO 5356-1 Conical Connectors
Mask Dead Space:	<325 cm ³
Resistance to flow through mask:	0.07 cmH ₂ O @ 50 L/min, 0.33 cmH ₂ O @ 100 L/min

- This product is intended for use for a maximum of 14 days.
- For use with pressurised gases provided by an external flow source or ventilator.
- For use with Ventilators/Devices delivering NIV therapy using dual limb circuit system.
- Can be used with a heated pass over humidifier to deliver humidified gases to the patient.
- When humidifying gases, use as part of Fisher & Paykel Healthcare non-invasive humidification systems F&P 850™ or F&P 950™.
- Note:** F&P 950™ System may not be available in all countries. F&P 950™ System is currently not cleared for sale in the USA.
- This product was not made with natural rubber latex.
- Dispose of mask according to hospital protocol.

Label – Mask Size Symbols

 = Extra Small  = Small  = Medium  = Large

 = Full Face  = Non Vented

Before Use

- Ensure you have read and understood the instructions completely.
- Remove all packaging prior to use.
- Do not use if package is damaged.
- Inspect the mask for damage. Discard the mask if any parts are broken or if the seal is torn.
- Ensure that the correct mask is being used for the therapy device. Refer to setup guide on the other side.
- Verify that the mask is the correct size for the patient.
- Verify that the therapy device (i.e. ventilator or flow source) including all alarms and safety systems are functioning correctly and that it is supplying the correct pressure(s).
- Clean the patient's face and then inspect for signs of redness, irritation or discomfort. Do not use if signs are present.
- Ensure adequate patient monitoring is in place.

Warnings

- Verify that the therapy device, including alarms and safety systems, are functioning correctly prior to use.
- This device is not fitted with an anti-asphyxiation valve and is not for use with single limb ventilation. Use with single limb ventilation may result in excessive carbon dioxide rebreathing in the event of flow disruption.
- The mask must only be worn when the therapy is being delivered.
- The mask must be fitted and therapy established by an appropriately trained medical practitioner or care provider.
- This mask may only be used in a hospital or clinical setting where the patient is adequately monitored by trained medical staff. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This mask is for single patient use. **Do not reuse.** Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death. Do not soak, wash, sterilise, or re-use this product. Avoid contact with chemicals, cleaning agents, or hand sanitisers.
- Facial hair, missing teeth/dentures and facial structure irregularities may compromise mask seal.
- California residents please be advised of the following, pursuant to Proposition 65:
This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Cautions

- Do not over tighten any of the headgear straps, overtightening may cause patient discomfort and/or leaks.
- Ensure the mask frame does not directly contact the forehead. If the mask frame touches the forehead the straps are overtightened.
- The RollFit™ Seal is designed to roll back and forth on the bridge of the nose to accommodate a range of face shapes. If the RollFit™ Seal is fully compressed, mask performance may be impaired.
- If the patient experiences skin redness, irritation or discomfort discontinue use and contact a physician.

Contraindications

Should not be used on patients who:

- Have cardiac or respiratory arrest, or severe hemodynamic instability.
- Are unconscious, unable to breathe spontaneously, uncooperative, unresponsive or unable to remove the mask.
- Have an upper airway obstruction, or an inability to clear secretions (impaired cough or swallow reflexes, excessive reflux, epistaxis, hiatal hernia).
- Have copious secretions, at risk of nausea/vomiting, or at high risk of aspiration of emesis.
- Have had head or facial surgery, trauma or burns.
- Have severe upper gastro-intestinal bleeding, or barotrauma (undrained pneumothorax).

If symptoms of these conditions occur discontinue treatment immediately.

Fitting Instructions

- Size the patient by aligning the blue line of the sizing guide between the eyes at eye level. The correct mask size is measured just below the lower lip on the chin.
- Connect the mask to the flow source. Ensure the flow source is turned on.
- Undo the lower clip(s) and slide the loose headgear over the patient's head.
- Connect the headgear clip(s) onto the mask frame.
- Adjust the upper headgear straps. If the mask frame touches the forehead this indicates that the straps are overtightened.
- If the patient has an orogastric or nasogastric tube, ensure the tube is placed under the integrated TubeFit zone, to achieve a good seal around the tube.
- Adjust the lower straps of the headgear.
- Gently pull the mask forward, allowing the seal to inflate so it can adjust to the patient's face to minimise leak.

Note: Readjust headgear/mask as required to achieve a good seal. Remember to gently pull the mask forward to inflate the seal.

For patent information, see www.fphcare.com/ip

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