










## INTENDED USE

The Fisher & Paykel Healthcare single patient use masks are intended for use as an accessory to ventilators to enable non-invasive positive pressure ventilation (NPPV) therapy (CPAP or bi-level) to be delivered to spontaneously breathing adult patients (>30 kg) with respiratory insufficiency or respiratory failure who have been prescribed NPPV. The masks are to be fitted and therapy maintained by trained medical practitioners in a hospital/institutional environment with patient monitoring in place.

## CONTRAINDICATIONS

This mask should not be used on patients who:

- are having cardiac or respiratory arrest, or severe hemodynamic instability
  - are unconscious, unable to breathe spontaneously, uncooperative, unresponsive, or unable to remove the mask
  - have an upper airway obstruction, or an inability to clear secretions (impaired cough or swallow reflexes, excessive reflux, epistaxis, hiatal hernia)
  - have copious secretions, at risk of nausea/vomiting, or at high risk of aspiration of emesis
  - have had head or facial surgery, trauma or burns
  - have severe upper gastro-intestinal bleeding, or barotrauma (undrained pneumothorax)
- If symptoms of these conditions occur discontinue treatment immediately.


## SIDE EFFECTS

After using this mask you may experience:

- redness, irritation, or pressure injury
- discomfort (eye irritation, claustrophobia, gastric insufflation, noise, dryness, nasal congestion, desynchrony)
- aspiration of fluid collected in the mask.

Please consult the appropriate medical practitioner if any of the above symptoms occur.

## TECHNICAL SPECIFICATIONS

|                                                                 |                                                                                    |
|-----------------------------------------------------------------|------------------------------------------------------------------------------------|
| Operating pressure range:                                       | 4 – 25 cmH <sub>2</sub> O                                                          |
| Interface connections:                                          | ISO 5356-1 conical connectors                                                      |
| Mask dead space:                                                | < 325 cm <sup>3</sup>                                                              |
| Anti-asphyxiation valve open to atmospheric pressure:           | 0.24 cmH <sub>2</sub> O                                                            |
| Anti-asphyxiation valve closed to atmospheric pressure:         | 0.80 cmH <sub>2</sub> O                                                            |
| Resistance to flow through mask:                                | 0.26 cmH <sub>2</sub> O @ 50 L/min,<br>0.63 cmH <sub>2</sub> O @ 100 L/min         |
| Pressure-flow curve for venting:                                |  |
| The venting flow rate may vary due to manufacturing variations. |                                                                                    |

- This mask was not made with natural rubber latex.
- This mask was not made with phthalates.

## SYMBOL DEFINITIONS

|                                                                                    |                                    |                                                                                     |                                              |                                                                                     |                                       |
|------------------------------------------------------------------------------------|------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------|
|  | Vented                             |  | Catalog Number                               |  | EU Authorised Representative          |
|  | Full Face                          |  | Batch Number                                 |  | European Conformity - TÜV SÜD         |
|  | Single-limb Compatible             |  | Date and Country of Manufacture<br>MX Mexico |  | UK Authorised Representative          |
|  | Size Extra Small                   |  | Manufacturer                                 |  | Packaging Materials                   |
|  | Size Small                         |  | Medical Device                               |  | Switzerland Authorised Representative |
|  | Size Medium                        |  | Expiry Date                                  |  | Warning                               |
|  | Size Large                         |  | Prescription Only                            |  | 14 Days' Maximum Use                  |
|  | Follow Instructions for Use        |  | Not Made with Phthalates (DEHP, BBP, DBP)    |  | Temperature Limits                    |
|  | Consult Instructions for Use       |  | Importer                                     |  | Do Not Use if Package is Damaged      |
|  | Not Made with Natural Rubber Latex |  | Keep Dry                                     |  | Distributor                           |
|  | Single Use                         |                                                                                     |                                              |                                                                                     |                                       |

## OPERATING INSTRUCTIONS

- This mask is intended for use for a maximum of 14 days.
- This mask is for single patient use.
- For use with ventilators and devices delivering NIV therapy using a single-limb circuit system.
- For use with pressurised gases provided by an external flow source or ventilator.
- Recommended for use as part of the Fisher & Paykel Healthcare noninvasive humidification systems, F&P 850 or F&P 950™, to deliver humidified gases to the patient.
- Note:** F&P 950 System may not be available in all countries.
- Dispose of mask according to hospital protocol.
- Notice to user: Any serious incident that has occurred in relation to the device should be reported to the competent authority in your country and your local Fisher & Paykel Healthcare representative.

## ⚠️ WARNINGS

- Appropriate patient monitoring, e.g., oxygen saturation, must be adhered to at all times. Failure to monitor the patient, e.g., in the event of flow disruption, may result in serious harm or death.
- Do not use this mask outside of a hospital/institutional environment in which the patient is adequately monitored by trained medical staff as this may result in serious harm or death.
- Verify that the therapy device, i.e., ventilator or flow source, including all alarms and safety systems are functioning correctly and that it is supplying the correct pressure(s). Failure to do so may result in excessive carbon dioxide rebreathing.
- This mask is fitted with an anti-asphyxiation valve **H** and is not for use with dual-limb ventilation. Use with dual-limb ventilation may result in compromised therapy.
- Do not block or seal the venting holes **I**. Doing so may result in excessive carbon dioxide rebreathing.
- Do not block the anti-asphyxiation valve openings **G**. Doing so may result in excessive carbon dioxide rebreathing in the event of flow disruption.
- Replace the mask every 14 days or sooner if noticeable damage occurs. Failure to do so may result in compromised therapy.
- Do not reuse, soak, wash or sterilise this mask. Avoid contact with chemicals, cleaning agents or hand sanitizers. Reusing or reprocessing the mask may result in transmission of infectious substances, inhalation of particulates or compromised therapy.
- There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.
- Smoking during oxygen therapy is dangerous and is likely to result in serious injury from fire.
- Use only lotions or salves that are labeled as oxygen compatible to avoid the risk of fire and burns.
- Do not lubricate fittings, connections, tubing or other accessories to avoid the risk of fire and burns.

## FITTING INSTRUCTIONS

- Size the patient by aligning the thin blue line of the sizing guide between the eyes at eye level. The correct mask size is measured just below the lower lip on the chin.
- Connect the mask to the flow source. Ensure the flow source is turned on.
- If using a flow source that requires pressure feedback, connect the pressure line to the pressure port **J** on the mask. Otherwise, ensure the pressure port is capped.
- Undo the headgear clip(s) **A** and slide the loose headgear over the patient's head.
- Connect the headgear clip(s) **A** onto the mask frame **B**.
- Adjust the upper headgear straps **C**. If the mask frame **B** touches the forehead this indicates that the straps are overtightened.
- If the patient has an orogastric or nasogastric tube, ensure the tube is placed under the integrated TubeFit™ zone **D** to achieve a good seal around the tube.
- Adjust the lower headgear straps **E**.
- Gently pull the mask forward, allowing the seal **F** to inflate so it can adjust to the patient's face to minimise leak.

**Note:** Readjust the headgear and mask as required to achieve a good seal. Remember to gently pull the mask forward to inflate the seal **F**.

## ⚠️ WARNINGS

- Remove all packaging prior to use. Packaging can occlude gas pathway resulting in compromised therapy.
- Do not use the mask if the mask is damaged. Using a damaged mask may result in compromised therapy.
- Do not use the mask if the packaging is damaged or unintentionally opened before use. Doing so may result in contamination leading to infection.
- Verify that the mask is the correct size for the patient. Incorrect sizing may result in compromised therapy, pressure injuries or eye infection.
- Only appropriately trained medical practitioners or care providers are to fit the mask and establish therapy. Fitment by an untrained person may result in a range of patient harms, some of which are serious.
- Do not overtighten any of the headgear straps **C E**. Overtightening may result in the mask frame **B** contacting the forehead, or the seal **F** being overly compressed. This may cause pressure injuries.

## NOTES

- Facial hair, missing teeth and dentures, and facial structure irregularities may compromise mask fit or seal.

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