CE 0123 Rx only  

REF 611734 REV F 2024-07 © 2024 Fisher & Paykel Healthcare Limited

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)**  Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

www.fphcare.com

Fisher & Paykel
HEALTHCARE

F&P Optiflow Junior 2

Next generation care

Instructions for use **Ventilator Transition Kit**

REF

OJR410VT

OJR412VT

OJR414VT

OJR416VT

OJR418VT


www.fphcare.com

Fisher & Paykel
HEALTHCARE

Sizing information

1. Prongs must not create a seal in the nares.
A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.

		WEIGHT GUIDE (KG)									
		500g	1	1.5	2	2.5	3	3.5	4	SPARES	
	OJR410VT									F&P Wigglepads™ 2 WJR110	
	XS										
	OJR412VT										
	S										

		WEIGHT GUIDE (KG)															
		1	2	3	4	5	6	7	8	9	10	15	20	25	30	SPARES	
	OJR414VT															F&P Wigglepads 2 WJR112	
	M																
	OJR416VT																
	L																
	OJR418VT																
	XL																

3.  Expected to fit patient  May fit patient







Technical specifications

Maximum flow rates (L/min)



Compatible Platform	 MR850		 F&P 950™*	
Compatible Circuit	RT265	RT266	950N80/J, 950N81/J, 950N82/J, 950N61/J	950P81/J
 OJR410VT XS	8	4	10	10
 OJR412VT S	9		10	10
 OJR414VT M	10		11	11
 OJR416VT L	15		28	34
 OJR418VT XL	15		31	36

- Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
- Maximum flow rates achieved are ventilator dependent. Refer to the manufacturer's instructions for correct usage.
- MR850 and F&P 950 flow rates are expressed in STPD.
- Refer to circuit user instructions for minimum flow rates.




*May not be available in all countries.

I**Connect
Cannula****1****2****3****II****Apply
Cannula****1****2****3**

III Remove Cannula

IV Replace F&P Wigglepads 2

F&P Optiflow Junior 2

Next generation care

English	6	Latviešu (Latvian)	71
العربية (Arabic)	9	Bahasa Melayu (Malay)	74
български (Bulgarian)	12	Nederlands (Dutch)	77
Česky (Czech)	15	Norsk (Norwegian)	80
Dansk (Danish)	18	Polski (Polish)	83
Deutsch (German)	21	Português (Portuguese)	87
Ελληνικά (Greek)	25	Português (Brasil) (Brazilian Portuguese)	90
Español (Spanish)	29	Română (Romanian)	93
Eesti keel (Estonian)	32	Русский (Russian)	96
Suomi (Finnish)	35	Slovenčina (Slovak)	100
Français (French)	38	Slovenščina (Slovenian)	103
Français Canadien (French Canadian)	41	Српски (Serbian)	106
עברית (Hebrew)	44	Svenska (Swedish)	109
Hrvatski (Croatian)	47	ไทย (Thai)	112
Magyar (Hungarian)	50	Türkçe (Turkish)	115
Bahasa Indonesia (Indonesian)	53	Українська (Ukrainian)	118
Íslenska (Icelandic)	56	اردو (Urdu)	121
Italiano (Italian)	59	Tiếng Việt (Vietnamese)	124
日本語 (Japanese)	62	简体中文版 (Simplified Chinese)	127
한국어 (Korean)	65	繁體中文版 (Traditional Chinese)	130
Lietuvių (Lithuanian)	68		

Optiflow, Wigglepads, and the Seahorse, Crab, Starfish, Octopus and Turtle images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

F&P Optiflow™ Junior 2 Ventilator Transition Kit

English (en)

Intended use

The Fisher & Paykel Healthcare Optiflow Junior 2 Ventilator Transition Kit allows the use of Optiflow Junior 2 cannula with approved Fisher & Paykel Healthcare breathing circuits when connected to a ventilator to deliver heated and humidified nasal high flow (NHF) therapy.

Optiflow Junior 2 nasal cannula is a single use nasal cannula intended for spontaneously breathing patients.

The intended pediatric subpopulation targeted for use of the F&P Optiflow Junior 2 nasal cannula range includes:

- Neonates, birth up to 1 month of age
- Infants, 1 month up to 2 years of age
- Children, 2 years up to 12 years of age

This product is designed for use in hospital environments and must be prescribed by a physician.

Contraindications

Suspected or confirmed abnormalities or trauma to the skull or airway involving any abnormal connection may allow pressure to be transmitted to unintended anatomical structures or tissues. Use of NHF in this clinical scenario can lead to serious injury or death.

Side effects

Nasal trauma and skin injury are side effects of using NHF interfaces.

GENERAL WARNINGS

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient may result in serious injury or death. For example, in the event of an interruption to gas flow.
- If using supplemental oxygen, keep ignition sources away from the patient.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate preexisting, untreated serious air leak syndrome which may lead to further serious injury or death.

- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of barotrauma, hypoxic injury and skin damage remain. These risks may result in serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

GENERAL CAUTIONS

- This cannula is intended to be used for a maximum of 14 days. Replace Wigglepads 2 as required.
- Using this product beyond 14 days may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Regularly monitor the patient to ensure skin integrity and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers. Secretions on the cannula and the prongs can be removed by gently wiping with a damp cloth.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious injury or death.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- DO NOT use in an MRI or similar magnetic-based scanner as the tubing contains stainless steel.
- Tubing may pose a risk of strangulation or airway restriction.
- DO NOT use if the product or its packaging has been tampered with.
- Ensure the patient does not lie on the tubing as this may apply pressure to the patient's ears or face.
- Product to only be used with medical grade gas supplies. The gas supply used with this device may unexpectedly fail to deliver oxygen or flow.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Only air and oxygen mixtures are intended to be used with the Optiflow Junior 2 Interface. The materials used may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.

Fitting Instructions

Connect Cannula 1

1. Select appropriate cannula size; A clear gap must be visible around each prong.
2. Connect the cannula to the adapter and the inspiratory limb of the breathing circuit. Remove the unheated extension if not required.
 - Connect the system to the ventilator and after following the manufacturer's instructions, ensure there is gas flow through the prongs.
3. Remove and discard the expiratory limb of the circuit. Disconnect the expiratory heater wire adapter if in use.

WARNINGS

- Failure to disconnect the heater wire adapter from the expiratory limb may increase the risk of fire or burns.
- A ventilator with Nasal High Flow mode and appropriate pressure limits must be used. Failure to do so may compromise therapy, leading to serious patient injury (e.g. hypoxia or barotrauma).

CAUTIONS

- The circuit temperature probe positioning and use of the unheated extension must be according to the circuit user instructions. Failure to remove the unheated extension when required will lead to excessive condensation.

Apply Cannula II

1. Prepare the patient's skin according to hospital protocol.
1. Remove the first backing tabs from the F&P Wigglepads 2 and avoid touching the adhesive.
2. Insert the cannula into the nares. Ensure the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.
3. Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.

WARNINGS

- DO NOT allow the prongs to seal in the nares. Occlusion may result in septal damage or barotrauma.

CAUTIONS

- DO NOT place the Wigglepads 2 on the patient's eyes, ears or injured skin.
- Ensure the cannula is placed directly on the Wigglepads 2. Direct skin contact caused by cannula misalignment may result in skin breakdown.

Remove Cannula III

Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2. Starting from the outside, peel towards the nose.

Replace F&P Wigglepads 2 IV

- Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.
1. Adhere the replacement Wigglepads 2 to the cannula, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.
 2. Remove the second backing tabs and stick onto the patient's cheeks.

Checks during operation

- Regular monitoring of the patient is necessary to ensure a slight gap between the cannula and septum is maintained, as well as the correct placement of the prongs in the nares. Reposition the cannula on the Wigglepads 2 if required.
- Regularly check skin integrity to prevent excessive pressure to the upper lip.
- To prevent nare occlusion, clean secretions from the cannula and the patient's face as required.
- Check the cannula remains secure. Replace the Wigglepads 2 if required.
- Ensure all connections are secure during use. Check the cannula is undamaged and that the flow path is maintained. Under excessive load, the cannula may disconnect to prevent forces being transferred to the patient.

CAUTIONS

- DO NOT wrap, insulate, stretch or crush the tubing as this may impair the performance of this product or compromise safety (including potentially causing patient injury).
- Monitor condensate to prevent occlusion or build-up of fluid. Drain away from the patient as required.

Technical specifications

DISPOSAL

There are no specific disposal requirements for this device. Disposal protocols appropriate to single use devices that may be contaminated with potentially infectious substances should be followed.

Operating conditions

- Ambient temperature range: 18 to 26 °C.

Note: This product may be used in an incubator.

Approved compatible equipment/accessories

- MR850 Humidifier in invasive mode with approved breathing circuits, accessories and chamber kits.
- F&P 950* Humidifier with approved breathing circuits, accessories and chamber kits.

Refer to the table for a complete list of approved breathing circuits and accessories.

*May not be available in all countries.

Operating flow rates

The operating flow rates for each F&P Optiflow Junior 2 cannula size are dependent on the breathing circuit and humidifier in use. Refer to the table for the operating flow rates for each cannula.

Sizing information notes

1. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.
3.  Expected to fit patient.  May fit patient.

Technical Specification Notes

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. Maximum flow rates achieved are ventilator dependent. Refer to the manufacturer's instructions for correct usage.
6. MR850 and F&P 950 flow rates are expressed in STPD.
7. Refer to circuit user instructions for minimum flow rates.

Symbol definitions

	Do not reuse	 0123	European Conformity
	Not made with phthalates (DEHP, DBP, BBP).		Storage temperature range
	Batch code		Medical device
	Consult instructions for use. fphcare.com/oj2-ifu		Prescription only
	Authorized representative in the European Community		This product is not made with natural rubber latex.
 YYYY-MM-DD	Use by date		Catalogue number
 YYYY-MM-DD	Date and Country of Manufacture NZ: New Zealand		Manufacturer
	Recyclable PET		Maximum 14 days use
	Importer		Distributor
	Switzerland authorized representative		UK responsible person
	Unique device identifier		

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