OPTIFLOW FILTERED
NASAL INTERFACE
WITH ANESTHESIA MASK
COMPATIBILITY

USER INSTRUCTIONS


REF **AA041**

CE 0123

Rx only

REF: 900942 REV D 2022-10 © 2022 Fisher & Paykel Healthcare Limited

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraisópolis, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât. F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

F&P 950, Optiflow Thrive, Optiflow Switch and Optiflow are trademarks of Fisher & Paykel Healthcare Ltd. For patent information, see www.fphcare.com/ip

Optiflow™ Switch Filtered Nasal Interface AA041 With Anesthesia mask compatibility Indications for Use

This product is a single-patient-use device that delivers respiratory gases to adult patients in hospitals and medical facilities.

Spontaneously Breathing Patients

This product is indicated for the delivery of Nasal High Flow (NHF) and Low Flow Oxygen to spontaneously breathing patients by appropriately qualified healthcare professionals.

Apneic Patients

This product is indicated for the delivery of Transnasal Humidified Rapid Insufflation Ventilatory Exchange ("THRIVE") to apneic patients in operating and procedure rooms, under the direction of an appropriately qualified anesthetist, anesthesiologist, or medical professional.

This product facilitates switching between the delivery of NHF/THRIVE and the delivery of anesthesia mask ventilation.

Contraindications

The following are contraindicated for this product. Failure to comply can lead to serious injury or death.

- Do not use this product with electrocautery or electrocautery devices on the head or neck. Refer to the FIRE DANGER information.
- Do not use where Continuous Positive Airway Pressure (CPAP) is contraindicated (e.g. pneumothorax, bullous lung disease, craniofacial trauma or surgery, airway foreign body, unstable hemodynamics).
- Do not use where nasal interfaces are contraindicated (e.g. nasal obstruction, nasal trauma).
- Do not administer inhalants, anesthetic gases or vapors through this product.
- Do not administer inhalants, anesthetic gases or vapors through an anesthesia mask while this product is fitted to the patient. The concentration of anesthetic agents administered to the patient may be altered by residual flow from the nasal interface. Gas composition measurements, such as concentration of anesthetic agent and Minimum Alveolar Concentration (MAC) may be inaccurate. This may result in delivery of an inadequate dose of anesthetic agent, or an excessive dose of anesthetic agent.

FIRE DANGER

This product is an open oxygen delivery system. Open oxygen delivery can increase the risk of a surgical fire occurring, causing serious injury or death. Extreme care must be taken. The following is advised.

- Contraindication:** Do not use this product with electrocautery or electrocautery devices on the head or neck.
- Warning:** Do not use this product where ignition sources and fuel are present. Use with ignition sources and fuel present completes the fire triangle, increasing the risk of fire. The following steps should be taken to reduce the risk of a surgical fire occurring:
 - Evaluate the oxygen needs of each individual patient and use the minimum supplementation required. Oxygen titration with a gas mixer should be considered.
 - Ensure the use of this product is communicated to all operating room personnel when conducting the procedure's fire risk assessment.
 - Ensure gas flow from the product is not in or near the surgical field AND the surgical site is free of all potential fuel sources, including alcohol skin preparation, gauze, sponges and drapes BEFORE potential ignition sources, such as electrocautery, electrocautery or laser devices are used.
 - Ensure the room is adequately ventilated as oxygen may accumulate over time.
 - Ensure gas flow from the product does not pool under drapes.
 - Follow the instructions for use of all surgical devices, including electrocautery, electrocautery and laser devices, regarding oxygen delivery.

Specifications

- Product is for single-patient-use only for a maximum period of 24 hours.
- System compatibilities and operating flow rates.

Compatible Systems	Anesthesia mask compatibility	Operating Flow Rate
F&P 950 Humidifier AA520 Flow Diverter AA452 Optiflow Oxygen Kit	✓	30-70 L/min
F&P Airvo™ 2 AA410 Airvo Tube and Chamber Kit	✓	10-50 L/min
AA020 Low Flow Tubing	✗	Up to 15 L/min

- Caution:** Flow rates above 15 L/min must be humidified. Failure to humidify can cause nasal and mucosal dryness and discomfort.
- Note:** When used with an Airvo 2 humidifier, the operating flow range is reduced:
 - All sizes: 10-50 L/min.
- Using an anesthesia mask with this product on Airvo 2 may result in a blockage alarm from the Airvo 2. This is part of normal operation. Flow from the Airvo 2 will resume automatically after removal of the anesthesia mask.
- Use only with compatible Fisher & Paykel Healthcare (F&P) humidifiers and accessories. **Warning:** Incompatible humidifiers and accessories which are used with this product may impair the

performance of this product or compromise safety (including potentially causing serious injury or death).

- Ambient operating temperature: 18-26 °C.
- Filtration efficiency: Viral: >99.996%, Bacterial: >99.999%.
- Filter is hydrophobic.
- This product is not intended to be sterile.
- This product was not made with natural rubber latex.

Setup and Operating Instructions

- Connect product to compatible system. Ensure all connections are secure.
- Ensure the product is sized correctly. **Warning:** Do not create a seal in the nares as this may present a risk of barotrauma or gastric insufflation. Set flow rate and ensure gas flow is exiting the prongs.
- The following practices must be adhered to when using an anesthesia mask with this product. **Warning:** Failure to do so can lead to excessive delivered pressure and cause serious injury or death:
 - Product must be used with a compatible flow diverter (e.g. AA520)
 - Ensure Adjustable Pressure Limiting (APL) valve on anesthesia machine is set to an appropriate relief pressure.
 - Monitor airway pressure and set maximum pressure alarm.
 - Remove product if ventilation with an anesthesia mask cannot be performed adequately.
- During use of this product, patients may desaturate rapidly with little warning and / or their arterial carbon dioxide may increase. **Warning:** Failure to follow the practices described below and appropriately address these situations can lead to serious injury or death.
 - Use appropriate patient monitoring, including oxygen saturation and carbon dioxide monitoring.
 - Ensure availability of appropriate measures to respond to a rapid desaturation and/or carbon dioxide increase.
 - In addition, during apnea, the following practices must also be adhered to:
 - Maintain a patent airway at all times (e.g. by applying jaw thrust).
 - Do not use product where a secure airway is indicated. This may increase the risk of pulmonary aspiration.

Warnings

- This product delivers oxygen. Oxygen titration with a gas mixer should be considered where hyperoxia is a concern.

Cautions

- When using this product with an anesthesia mask, residual flow of gas through the nasal interface may:
 - Alter the concentration of substances administered to the patient via the anesthesia mask (such as O₂ or bronchodilator).
 - Cause measurements to be inaccurate (such as tidal volume and gas composition).
- Do not use product where quantitative exhaled gas measurements are indicated. All oxygen delivery devices modify the patient's exhaled gas composition. Exhaled gas composition measurements may be significantly altered. Alternative monitoring methods such as arterial blood gas test should be considered.
- Do not reuse product on multiple patients. Reuse may result in transmission of infectious substances, interruption to treatment, serious injury or death.
- Do not use product if damaged. Product damage may result in impaired performance.
- Do not soak, wash or sterilize the product. Avoid contact with chemicals, cleaning agents and hand sanitizers as this may damage the product and/or cause leaks in the system and reduce the gas supply/humidity to the patient.
- Do not stretch or crush tubing, as this may occlude the gas flow or damage the tubing, resulting in impaired performance or gas leaks.
- Do not use product if packaging is not sealed or it is contaminated, as this presents a risk of infection.
- Do not use product in or near a magnetic resonance imaging (MRI) scanner.
- Dispose of product safely in accordance with standard hospital procedure. User may be exposed to respiratory tract fluids during disposal

Notes

- Refer to the instructions for use of compatible Fisher & Paykel Healthcare humidifiers and accessories for additional warnings, cautions, contraindications and system information.
- If a serious incident has occurred while using this product, please contact your local Fisher & Paykel Healthcare representative and Competent Authority.

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