F&P

Optiflow Junior 2+

Next generation care

Instructions for use
Nasal Cannula

REF OJR520


www.fphcare.com

Fisher & Paykel
HEALTHCARE

Sizing information

1. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.

OJR520 XXL	WEIGHT GUIDE (KG)						SPARES
	0	8	12.5	20	30	40	
			■	■	■	■	F&P Wigglepads™ 2 WJR114

3. ■ Expected to fit patient ■ May fit patient

Technical specifications

Maximum flow rates (L/min)

OJR520
XXL


MR850		AIRVO™ 2	F&P 950™*	
RT330	RT331	900PT561 900PT562*	950N40/J*	950P40/J*
36	45	50	36	60

- Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
- Maximum flow rates achieved with the RT331 are ventilator dependent. Refer to the manufacturer's instructions for correct use.
- MR850 and F&P 950 flow rates are expressed in STPD. Airvo 2 flow rates are expressed in BTPS.
- Refer to circuit user instructions for minimum flow rates.

*May not be available in all countries.





I

Apply
Cannula


II

Remove
Cannula

III

Replace F&P Wigglepads 2

F&P Optiflow Junior 2+

Next generation care

English	6	Latviešu (Latvian)	72
العربية (Arabic)	9	Bahasa Melayu (Malay)	75
български (Bulgarian)	12	Nederlands (Dutch)	78
Česky (Czech)	15	Norsk (Norwegian)	81
Dansk (Danish)	18	Polski (Polish)	84
Deutsch (German)	21	Português (Portuguese)	88
Ελληνικά (Greek)	25	Português (Brasil) (Brazilian Portuguese)	91
Español (Spanish)	29	Română (Romanian)	94
Eesti keel (Estonian)	32	Русский (Russian)	97
Suomi (Finnish)	35	Slovenčina (Slovak)	101
Français (French)	38	Slovenščina (Slovenian)	104
Français Canadien (French Canadian)	41	Српски (Serbian)	107
עברית (Hebrew)	44	Svenska (Swedish)	110
Hrvatski (Croatian)	47	ไทย (Thai)	113
Magyar (Hungarian)	50	Türkçe (Turkish)	116
Bahasa Indonesia (Indonesian)	53	Українська (Ukrainian)	119
Íslenska (Icelandic)	56	اردو (Urdu)	122
Italiano (Italian)	59	Tiếng Việt (Vietnamese)	125
日本語 (Japanese)	62	简体中文版 (Simplified Chinese)	128
한국어 (Korean)	65	繁體中文版 (Traditional Chinese)	131
Lietuvių (Lithuanian)	69		

Intended use

The Fisher & Paykel Healthcare Optiflow Junior 2+ nasal cannula is a single use nasal cannula intended for use with a nasal high flow therapy (NHF) system to deliver heated and humidified nasal high flow therapy to spontaneously breathing patients. This product is designed for use in hospital environments and must be prescribed by a physician.

The intended pediatric subpopulation targeted for use on the F&P Optiflow Junior 2+ Nasal Cannula includes:

- Infants, 1 month up to 2 years of age
- Children, 2 years up to 12 years of age

Contraindications

Suspected or confirmed abnormalities or trauma to the skull or airway involving any abnormal connection may allow pressure to be transmitted to unintended anatomical structures or tissues. Use of NHF in this clinical scenario can lead to serious injury or death.

Side effects

Nasal trauma and skin injury are side effects of using NHF interfaces.

GENERAL WARNINGS

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient may result in serious injury or death. For example, in the event of an interruption to gas flow.
- If using supplemental oxygen, keep ignition sources away from the patient.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate pre-existing, untreated serious air leak syndrome which may lead to further serious injury or death.
- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of barotrauma, hypoxic injury and skin damage remain. These risks may result in serious injury or death.

- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

GENERAL CAUTIONS

- This cannula is intended to be used for a maximum of 14 days. Replace Wigglepads 2 as required.
- Using this product beyond 14 days may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Regularly monitor the patient to ensure skin integrity and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers. Secretions on the cannula and the prongs can be removed by gently wiping with a damp cloth.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious injury or death.
- Only air and oxygen mixtures are intended to be used with the Optiflow Junior 2+ Interface. The materials used may not be compatible with anaesthetic or respirable gases, solutions/suspensions/ emulsions that have not been evaluated.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- Tubing may pose a risk of strangulation or airway restriction.
- DO NOT use if the product or its packaging has been tampered with.
- Ensure the patient does not lie on the tubing as this may apply pressure to the patient's ears or face.
- Product to only be used with medical grade gas supplies. The gas supply used with this device may unexpectedly fail to deliver oxygen or flow.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient injury).

Fitting instructions

Apply Cannula ①

1. Select appropriate cannula size; A clear gap must be visible around each prong.
 - Prepare the patient's skin according to hospital protocol.
 - Connect the system to the gas source and ensure there is gas flow through the prongs.
2. Remove the first backing tabs from the F&P Wigglepads 2 and avoid touching the adhesive.

- Insert the cannula into the nares. Ensure the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.
- Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.

WARNINGS

- DO NOT allow the prongs to seal in the nares. Occlusion may result in septal damage or barotrauma.
- When using the RT33I breathing circuit, a ventilator with Nasal High Flow mode and appropriate pressure limits must be used. Failure to do so may compromise therapy, leading to serious patient injury (e.g. hypoxia or barotrauma).

CAUTIONS

- DO NOT place the Wigglepads 2 on the patient's eyes, ears or injured skin.
- Ensure the cannula is placed directly on the Wigglepads 2. Direct skin contact caused by cannula misalignment may result in skin breakdown.

Remove Cannula

Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2. Starting from the outside, peel towards the nose.

Replace F&P Wigglepads 2

- Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.
- Adhere the replacement Wigglepads 2 to the cannula, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.
- Remove the second backing tabs and stick onto the patient's cheeks.

Checks during operation

- Regular monitoring of the patient is necessary to ensure a slight gap between the cannula and septum is maintained, as well as the correct placement of the prongs in the nares. Reposition the cannula on the Wigglepads 2 if required.
- Regularly check skin integrity to prevent excessive pressure to the upper lip.
- To prevent nare occlusion, clean secretions from the cannula and the patient's face as required.
- Check the cannula remains secure. Replace the Wigglepads 2 if required.
- Ensure all connections are secure during use. Check the cannula is undamaged and that the flow path is maintained. Under excessive load, the cannula may disconnect to prevent forces being transferred to the patient.

CAUTIONS

- DO NOT wrap, insulate, stretch or crush the tubing as this may impair the performance of this product or compromise safety (including potentially causing patient injury).
- Monitor condensate to prevent occlusion or build-up of fluid. Drain away from the patient as required.

Disposal

There are no specific disposal requirements for this device. Disposal protocols appropriate to single use devices that may be contaminated with potentially infectious substances should be followed.

Operating conditions

- Ambient temperature range: 18 to 26 °C.

Approved compatible equipment/accessories

- MR850 Humidifier in invasive mode with approved breathing circuits, accessories and chamber kits.
- AIRVO 2 Humidifier with approved breathing circuits, accessories and chamber kits.
- F&P 950 Humidifier* with approved breathing circuits, accessories and chamber kits.
- OPT014 Oxygen Tubing.
- Refer to the table for a complete list of approved breathing circuits and accessories.

Refer to the table for a complete list of approved breathing circuits and accessories.

* May not be available in all countries.

CAUTION

When using this product with an AIRVO 2 Humidifier, failure to use in Default mode may generate Check for Leaks alarms and may impair performance of this product or compromise safety (including potentially causing patient injury).

Operating flow rates

The operating flow rates are dependent on the breathing circuit and humidifier in use. Refer to the operating flow rates.

- OPT014 Oxygen Tubing flow rates 0.1-2 L/min for all cannula sizes.

Sizing notes

- Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
- Patient weight should only be used as a guide.
-  Expected to fit patient.  May fit patient.

Technical specification notes

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. Maximum flow rates achieved with the RT331 are ventilator dependant. Refer to the manufacturer's instructions for correct use.
6. MR850 and F&P 950 flow rates are expressed in STPD. Airvo 2 flow rates are expressed in BTPS.
7. Refer to circuit user instructions for minimum flow rates.

Symbol definitions

	Do not reuse
	Not made with phthalates (DEHP, DBP, BBP).
	Batch code
	Consult instructions for use. fphcare.com/oj2-ifu
	Authorized representative in the European Community
	Use by date YYYY-MM-DD
	Date and Country of Manufacture NZ: New Zealand
	Recyclable PET >PET<
	Importer
	Distributor
	European Conformity
	Storage temperature range -10°C (+14°F) / +50°C (+122°F)
	Catalogue number
	Prescription only

	This product is not made with natural rubber latex.
	Medical device
	Manufacturer
	Maximum 14 days use
	Switzerland authorized representative
	UK responsible person
	Unique device identifier

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英文名稱：“Fisher & Paykel” Nasal oxygen cannula (Non-Sterile)

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F&P, Optiflow, Wigglepads, and the Dolphin image shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

€ 0123 Rx only  

REF 612686 REV E 2024-08 © 2024 Fisher & Paykel Healthcare Limited

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900

Brazil (BR) Fisher & Paykel do Brasil, Rua Sampai Viana, 277 cj 21, Paraiso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002

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Germany (DE) Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0
Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63