F&P Optiflow Junior 2

Next generation care

Instructions for use **Blender Transition Kit**

REF

OJR410B

OJR412B

OJR414B

OJR416B

OJR418B


www.fphcare.com

Fisher & Paykel
HEALTHCARE

Sizing information

1. Prongs must not create a seal in the nares.
A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.

		WEIGHT GUIDE (KG)								
		500g	1	1.5	2	2.5	3	3.5	4	SPARES
	OJR410B XS									F&P Wigglepads[®] 2 WJR110
	OJR412B S									

		WEIGHT GUIDE (KG)														
		1	2	3	4	5	6	7	8	9	10	15	20	25	30	SPARES
	OJR414B M															F&P Wigglepads[®] 2 WJR112
	OJR416B L															
	OJR418B XL															

3.  Expected to fit patient  May fit patient

Technical specifications

Maximum flow rates (L/min)

Compatible Platform	 MR850	 F&P 950™*
Compatible Circuit	BC151 BC161 BC153 BC163 RT124* RT132 RT224* RT324	950N60/J* 950N62/J*
 OJR410B XS	8	10
 OJR412B S	9	10
 OJR414B M	10	11
 OJR416B L	15	28
 OJR418B XL	15	31






4. Flow rates above describe technical capability of the product when used at sea level.

Ensure clinical judgement is used when prescribing flow rates.



5. MR850 and F&P 950 flow rates are expressed in STPD.

6. Refer to circuit user instructions for minimum flow rates.



*May not be available in all countries.

I**Apply Pressure Relief****1****2****II****Connect Cannula****1****2****3****III****Apply Cannula****1****2****3**

IV Remove
Cannula

V Replace F&P
Wigglepads 2

English	6	Nederlands (Dutch)	87
العربية (Arabic)	9	Norsk (Norwegian)	91
български (Bulgarian)	12	Polski (Polish)	94
Česky (Czech)	16	Português (Portuguese)	98
Dansk (Danish)	19	Português (Brasil) (Brazilian Portuguese)	102
Deutsch (German)	22	Română (Romanian)	106
Ελληνικά (Greek)	26	Русский (Russian)	110
Español (Spanish)	30	Slovenčina (Slovak)	114
Eesti keel (Estonian)	34	Slovenščina (Slovenian)	118
Suomi (Finnish)	37	Српски (Serbian)	121
Français (French)	40	Svenska (Swedish)	125
Français Canadien (French Canadian)	44	ไทย (Thai)	128
עברית (Hebrew)	48	Türkçe (Turkish)	132
Hrvatski (Croatian)	51	Українська (Ukrainian)	135
Magyar (Hungarian)	54	اردو (Urdu)	139
Bahasa Indonesia (Indonesian)	58	Tiếng Việt (Vietnamese)	142
Íslenska (Icelandic)	62	简体中文版 (Simplified Chinese)	145
Italiano (Italian)	65	繁體中文版 (Traditional Chinese)	148
日本語 (Japanese)	69		
한국어 (Korean)	73		
Lietuvių (Lithuanian)	77		
Latviešu (Latvian)	80		
Bahasa Melayu (Malay)	83		

F&P Optiflow™ Junior 2 Blender Transition Kit

English 

Intended use

The Fisher & Paykel Healthcare Optiflow Junior 2 Blender Transition Kit allows the use of Optiflow Junior 2 cannula with approved Fisher & Paykel Healthcare breathing circuits when connected to a blender to deliver heated and humidified nasal high flow (NHF) therapy.

Optiflow Junior 2 nasal cannula is a single use nasal cannula intended for spontaneously breathing patients. This product is designed for use in hospital environments and must be prescribed by a physician.

The intended pediatric subpopulation targeted for use of the F&P Optiflow Junior 2 nasal cannula range includes:

- Neonates, birth up to 1 month of age
- Infants, 1 month up to 2 years of age
- Children, 2 up to 12 years of age

Contraindications

Suspected or confirmed abnormalities or trauma to the skull or airway involving any abnormal connection may allow pressure to be transmitted to unintended anatomical structures or tissues. Use of NHF in this clinical scenario can lead to serious injury or death.

Side effects

Nasal trauma and skin injury are side effects of using NHF interfaces.

GENERAL WARNINGS

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious injury or death.
- If using supplemental oxygen, keep ignition sources away from the patient.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate pre-existing, untreated serious air leak syndrome which may lead to further serious injury or death.
- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of barotrauma, hypoxic injury and skin damage remain. These risks may result in serious injury or death.

- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

GENERAL CAUTIONS

- This cannula is intended to be used for a maximum of 14 days. Replace Wigglepads 2 as required.
- Using this product beyond 14 days may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Regularly monitor the patient to ensure skin integrity and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers. Secretions on the cannula and the prongs can be removed by gently wiping with a damp cloth.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious injury or death.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- DO NOT use in an MRI or similar magnetic-based scanner as the tubing contains stainless steel.
- Tubing may pose a risk of strangulation or airway restriction.
- DO NOT use if the product or its packaging has been tampered with.
- Ensure the patient does not lie on the tubing as this may apply pressure to the patient's ears or face.
- Only air and oxygen mixtures are intended to be used with the Optiflow Junior 2 Interface. The materials used may not be compatible with anaesthetic or respirable gases, solutions/suspensions/emulsions that have not been evaluated.
- Product to only be used with medical grade gas supplies. The gas supply used with this device may unexpectedly fail to deliver oxygen or flow.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient injury).

Fitting instructions

Apply pressure relief 1

1. Connect the Optiflow Pressure Relief Manifold (blue cap) to the inlet port of the humidification chamber.
2. Connect the gas supply line to the Optiflow Pressure Relief Manifold.
 - Attach air/oxygen monitoring equipment as required. Refer to manufacturer's instructions for correct usage of monitoring equipment.

WARNINGS

- Ensure that the Optiflow Pressure Relief Manifold (blue cap) is used for nasal high flow therapy. Use of the incorrect Pressure Relief Manifold may affect the delivered flow and could result in serious patient injury (e.g. hypoxia).

CAUTIONS

- Do not use if the Optiflow Pressure Relief Manifold cap is missing or damaged.
- Ensure any unused ports have their caps in place before use.
- This device should be used only on the inlet (dry side) of the humidification chamber. Operation of the pressure relief manifold may be impaired if connected to the outlet port.
- Do not use an in-line nebulizer between the pressure relief manifold and the dry side of the humidification chamber.
- Do not use the Optiflow Pressure Relief Manifold (blue cap) with Bubble CPAP as excessive pressures may result.
- Do not bypass the Pressure Relief Manifold as this may cause injury.

Connect Cannula

1. Select appropriate cannula size; A clear gap must be visible around each prong
2. Connect the cannula to the adapter and the inspiratory limb of the breathing circuit.
 - Connect the system to the gas source and ensure there is gas flow through the prongs.
3. Remove and discard the expiratory limb of the circuit. Disconnect the expiratory heater wire adapter if in use.

WARNINGS

- Failure to disconnect the heater wire adapter from the expiratory limb may increase the risk of fire or burns.

Apply Cannula

- Prepare the patient's skin according to hospital protocol.
1. Remove the first backing tabs from the F&P Wigglepads 2 and avoid touching the adhesive.
 2. Insert the cannula into the nares. Ensure the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.
 3. Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.

WARNINGS

- DO NOT allow the prongs to seal in the nares. Occlusion may result in septal damage or barotrauma.

CAUTIONS

- DO NOT place the Wigglepads 2 on the patient's eyes, ears or injured skin.
- Ensure the cannula is placed directly on the Wigglepads 2. Direct skin contact caused by cannula misalignment may result in skin breakdown.

Remove Cannula

Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2. Starting from the outside, peel towards the nose.

Replace F&P Wigglepads 2

- Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.
1. Adhere the replacement Wigglepads 2 to the cannula, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.
 2. Remove the second backing tabs and stick onto the patient's cheeks.

Checks during operation

- Regular monitoring of the patient is necessary to ensure a slight gap between the cannula and septum is maintained, as well as the correct placement of the prongs in the nares. Reposition the cannula on the Wigglepads 2 if required.
- Regularly check skin integrity to prevent excessive pressure to the upper lip.
- To prevent nare occlusion, clean secretions from the cannula and the patient's face as required.
- Check the cannula remains secure. Replace the Wigglepads 2 if required.
- Ensure all connections are secure during use. Check the cannula is undamaged and that the flow path is maintained. Under excessive load, the cannula may disconnect to prevent forces being transferred to the patient.

CAUTIONS

- DO NOT wrap, insulate, stretch or crush the tubing as this may impair the performance of this product or compromise safety (including potentially causing patient injury).
- Monitor condensate to prevent occlusion or build-up of fluid. Drain away from the patient as required.

Disposal

There are no specific disposal requirements for this device. Disposal protocols appropriate to single use devices that may be contaminated with potentially infectious substances should be followed.

Technical specifications

Approved compatible equipment/accessories

- MR850 Humidifier in invasive mode with approved breathing circuits, accessories and chamber kits.
- F&P 950™ Humidifier with approved breathing circuits, accessories and chamber kits.

Refer to the table for a complete list of approved breathing circuits and accessories.

*May not be available in all countries

Operating flow rates

The operating flow rates for each F&P Optiflow Junior 2 cannula size are dependent on the breathing circuit and humidifier in use. Refer to the table for the operating flow rates for each cannula.

Sizing notes

1. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
2. Patient weight should only be used as a guide.
3.  Expected to fit patient.  May fit patient.

Technical specification notes

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. MR850 and F&P 950 flow rates are expressed in STPD.
6. Refer to circuit user instructions for minimum flow rates.

Pressure relief manifold

Nominal relief pressure 40 cmH₂O
(0.57 psi/3.9 kPa/0.039 bar)

1. Oxygen analyzer port: 15 mm female (behind cap)
2. Pressure relief manifold cap: blue
3. Inlet connection port: barbed air/oxygen
4. Outlet connection port: 22 mm female
5. Flow direction arrow

Operating conditions

Ambient temperature range: 18 to 26 °C.

NOTE: This product may be used in an incubator.

Symbol definitions

	Do not reuse	CE 0123	European Conformity
	Not made with phthalates (DEHP, DBP, BBP).	-10°C (+14°F)  +50°C (+122°F)	Storage temperature range
LOT	Batch code	REF	Catalogue number
MD	Medical device	Rx only	Prescription only
EC REP	Authorized representative in the European Community		This product is not made with natural rubber latex.
	Use by date YYYY-MM-DD		Manufacturer
	Date and Country of Manufacture NZ: New Zealand YYYY-MM-DD		Consult instructions for use fphcare.com/oj2-ifu
	Recyclable PET >PET<	14	Maximum 14 days use
	Importer		Distributor
CH REP	Switzerland authorized representative	UK REP	UK responsible person
UDI	Unique device identifier		

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型號：OJR410B, OJR412B, OJR414B, OJR416B, OJR418B

F&P, Optiflow, Wigglepads, and the Seahorse, Crab, Starfish, Octopus and Turtle images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

CE 0123 Rx only  

REF 611733 REV E 2024-08 © 2024 Fisher & Paykel Healthcare Limited

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900

Brazil (BR) Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002

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Germany (DE) Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0
Switzerland (CH)  Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63