


F&P 850™ System


F&P AIRVO™ 2


1 Apply Cannula

2 Remove Cannula

3 Replace Wigglepads™

This product is covered by or for use under one or more of the following patents or patent applications in the corresponding countries listed: Australia AU 345792, Europe EM 002100289-0001, EM 002100289-0002, EM 002100289-0003, EM 002132464-0001, EM 002132464-0002, EM 002132464-0003. The Bumblebee, Butterfly, Bird and Ladybird images above (collectively the Optiflow Critter Images) are trade marks of Fisher & Paykel Healthcare Limited, ©2019 Fisher & Paykel Healthcare Limited.

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English		
Nasal Cannula		
Maximum Patient Flow (L/min)		
Size MR850 AIRVO 2		
1. Apply Cannula		
2. Remove Cannula		
3. Replace Wigglepads™		
To Apply Cannula:		
• Select appropriate size. Patient weight is to be used as a guide only. Prongs should not fill the noses and a clear gap should be visible around each prong.		
• Ensure skin is dry/prepared.		
• Connect system to gas source. Place hand close to the nasal prongs to ensure that there is airflow exiting the prongs.		
• Remove first Wigglepad backing tabs without touching the adhesive.		
• Position cannula prongs high up into nares such that the cannula bridge rests just underneath the septum. Stick Wigglepads to the cheeks.		
• Remove second backing tabs and stick Wigglepads onto cheeks.		
Checks During Operation:		
• Check patient's septal integrity. Ensure cannula does not apply pressure to the septum. A slight gap should be visible between the cannula and septum. Check placement of prongs in nares. Reposition cannula on Wigglepads if required.		
• Check cannula remains secure on the face. Replace Wigglepads if required.		
• Ensure that the tubing is not applying excessive pressure to the ears and face.		
• Ensure that there are no secretions that could cause the prongs to seal in the nares.		
• Ensure that all connections are secure during use. Check cannula is undamaged and that the flow path is maintained. Under excessive load, the cannula may disconnect to prevent forces being transferred to the patient.		
System Specifications:		
• MR850 Humidifier with RT330 delivery tubing and chamber kit. MR850 Humidifier in invasive mode.		
OR		
AIRVO™ 2/myAIRVO™ 2 Humidifier with 900PT531 tubing and chamber kit.		
• Ambient Range: 18 – 26 °C.		
This product is not made with natural rubber latex, and does not contain PVC or Phthalates (DEHP, DBP, BBP).		
This product is intended to be used for a maximum of 7 days.		
Contraindications:		
This therapy must not be used where continuous positive airway pressure (CPAP) is contraindicated.		
Warnings:		
• Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.		
• Always ensure prongs are well positioned in nares. Failure to do so may result in hypoxia or septal damage.		
• Do not place Wigglepads on broken or damaged skin, patient's eyes or ears.		
• Tubing contains stainless steel. Do not use in Magnetic Resonance environment.		
• Tubing poses a potential strangulation hazard.		
• Do not use if packaging is not sealed.		
• MR850: This device must be used with the pressure relief valve supplied with the RT330 Optiflow Tubing Kit. Failure to use this safety device can allow unlimited pressure to be delivered to the patient.		
• Do not stretch or crush tube.		
• Do not soak, wash, sterilise, or re-use this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.		
• Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.		
• Use of a non-approved accessory could impair performance or compromise safety.		
• Do not use near a naked flame.		
• Check for condensate regularly. Drain as required.		
California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65		

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